



EDGEENDO®
Twice as Good. Half the Cost.

EdgeOne Platinum™ Heat Treated Fire-Wire™

EdgeOne Platinum™ Instruments

- VeriFile™
- Small
- Primary
- Medium
- Large

DIRECTIONS FOR USE

- EdgeOne Platinum™ Heat-Treated Files are not pre-sterilized and must be sterilized before use.
- EdgeOne Platinum™ Heat-Treated Files use the same motor, hand piece and settings as WaveOne® Gold Files.
- EdgeOne Platinum™ is designed to shape canals in a reverse-reciprocating motion and settings exactly like WaveOne® Gold.
- Use an in-and-out motion advancing the file 1-3 mm with each stroke.
- Fill the chamber with 17% EdgeLube™ EDTA Liquid.
- Coat each EdgeOne Platinum™ file with 17% EdgeGel™ EDTA Gel.

GlidePath:

- Fill the chamber with EdgeLube™ EDTA Liquid.
- Take #10 hand file to the estimated working length.
- Establish the working length with Apex Locator or X-ray.
- Established canal patency by taking the #10 hand file 1mm past the working length.
- Expand the GlidePath™ by taking a #15 hand file or EdgeGlidePath™ rotary or EdgeFind™ rotary files to working length.

Initial Shaping:

- Fill chamber with EdgeLube™ EDTA liquid.
- Always use the EdgeOne Platinum™ VeriFile™ as your first file for initial shaping of the canal by taking it to length with small in-and-out motions advancing the file apically 1-3mm per stroke.
- After the VeriFile™, rinse with EdgeLube™ EDTA liquid and recapitulate with a #10 hand file.

Finishing Canal:

- If the VeriFile™ goes down to working length without resistance, finish with either the Medium or Large file.
- If the VeriFile™ goes down to working length with moderate resistance, finish with the Primary.
- If the VeriFile™ goes down to working length with tight resistance, finish with the Small.
- If the VeriFile™ goes not go down to working length, alternate between the Small and VeriFile™ until the Small is to working length.

Disinfect Canal:

Use NaOCl or your current disinfecting technique.

Obturate Canal:

Use EdgeOne Platinum™ Gutta Percha Points.



1-855-985-3636 EdgeEndo.com