



EDGEENDO®

**EDGE TAPER
RETREAT™**

Product Type: D1, D2, D3

Specification: #20, #25, #30

Application: Use of root canal treatment, for cleaning and shaping root canal.

Product structure, usage:

Composed by file needle, file handle and rubber fortify parts. File needle is made with stainless steel. Fortify parts is made with silicon rubber. File needle have good elasticity, durability and corrosion resistance, can satisfy medical demand on root canal preparation and obturation for complicated root canal. File handle and rubber fortify parts can be sterilized by various liquor antiseptic or 145 degrees high temperature sterilizer.

Prohibition: Is not for use on people with nickel allergy

Cautions:

1. This product can only use by doctors with formal endodontic treatment training and knowledge
2. Beside hand-use instrument, root canal file also have mechanical instrument. Mechanical instrument for root canal treatment set with a bolt type connection, and can connect with a rotary motor. Rotary motor should have torque controller and rotation speed need to be restricted within an appropriate range, recommended to be within 150~600rpm.
3. Product is not sterilized, new products need to be sterilized before use;
4. Sterilization method include autoclave or use antiseptic;
5. If apply autoclave method, heating temperature should not over 145 degrees or it will result in file handle melting and fortify parts depreciate
6. For non-sterilized product, it can be used repeatedly for certain times, usually not recommended to use over 10 times. Number of use depends on root canal treatment condition, and usage. Instruments check should take place before every usage. If any instrument deformation is discovered, or file needle become blunt, or lost normal functioning, need to stop usage.
7. After usage this product need to be treat as medical discard.

Transport and storage:

Product transportation cannot damage product's original package.; store under indoor temperature; restricted for pressure

Manufacture date: See package

Period of Validity: 5 Years under recommended storage condition



EDGEENDO®
Twice as Good. Half the Cost.

1-855-985-3636

EDGEENDO.COM