

## Directions For Use

### EDGE C-FILES HAND INSTRUMENT

#### Composition

Edge C-Files™ are made with Stainless Steel.

#### Indications for Use

- Edge C-Files™ are used to form the initial GlidePath before EdgeEndo<sup>®</sup> NiTi rotary or reciprocating files.
- These instruments are for Dentist use only.

#### Contraindications

None Known.

#### Warnings

- Use of a rubber dam system before use is advised.
- Edge C-Files™ are intended for single use only to avoid file separation.

#### Precautions

- Exercise caution in canals that merge, divide, and/or have curvatures.
- For personal safety, wear protective equipment such as gloves, glasses, mask, etc.
- Lubricants such as EdgeLube™ and EdgeGel™ are recommended while using all hand instruments such as C-Files.
- To establish the working length and initial GlidePath take a #10 hand file to the working length.
- Clean flutes frequently and check for signs of distortion or wear.
- Edge C-File™ should be used in and in-out, watch-winding, or balanced-force motion.
- Edge C-File™ are used to form the initial GlidePath and to follow the canal until the working length is achieved.

#### Adverse Reactions

None known.

#### Step by step instructions

- Form straightline access to the canal orifice.
- Use a # 10 C-File to establish the working length with an apex locator and/or radiograph.
- Confirm patency with the #10 Edge C-File™
- Irrigate frequently with EdgeLube, 17% Liquid EDTA.
- If the #10 Edge C-File will not reach the working length, drop down and use a #6 Edge C-File™ to the working length, then #8, and then #10 Edge C-File™.
- Irrigate with EdgeLube™, 17% Liquid EDTA between each file.
- After the #10 is at the working length and patency is confirmed either 1) Complete the GlidePath with GlidePath systems EdgeGlidePath™ or EdgeFind™ 2) Start shaping the canal with EdgeEndo<sup>®</sup> rotary or reciprocating files



Rev. A 4/2017