EDGE  UTOPIA™

(

E
E
P
F
I
R
N
D
S
N
D
S
H

CHOOSE YOUR LANGUAGE

ENGLISH.....	5	SLOVENŠČINA	73
ESPAÑOL.....	10	SLOVENSKY	78
PORTUGUÊS	15	POLSKI	83
FRANÇAIS	21	EESTI KEEL	88
ITALIANO	26	LATVIEŠU VALODA	93
ROMÂNĂ	32	LIETUVIŠKAI	98
NEDERLANDS	37	ČESKY	104
DEUTSCH	42	MAGYAR	109
SVENSKA	48	TÜRKÇE	114
NORSK.....	53	БЪЛГАРСКИ	119
DANSKI	58	СРПСКИ	125
SUOMI	63	ΕΛΛΗΝΙΚΑ	130
HRVATSKI	68		

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ENGLISH

INTENDED PURPOSE

EdgeUtopia™ BioCeramic Sealer is a bioceramic root canal sealer, intended for exclusive use by dentists during endodontic obturation procedures. It is designed to fill and seal previously cleaned and shaped root canal systems, aiming to prevent reinfection and ensure long-term obturation. The device is supplied ready-to-use.

INDICATIONS

- Obturation of root canals of permanent teeth.
- Use during internal resorption treatment.

INTENDED POPULATION

Patients with permanent teeth with closed apices, diagnosed with pulpal or periapical disease.

INTENDED USERS

EdgeUtopia™ BioCeramic Sealer should be used by dentists trained to perform endodontic procedures.

INTENDED ENVIRONMENT

Dental offices.

CONTRAINDICATIONS

- Do not use the product after the expiration date.
- Do not use the product if the packaging is damaged.
- Do not use in patients with reported sensitivity to any of the formula components.
- Do not reuse the applicator tips.

COMPOSITION

Calcium silicate, calcium aluminate, calcium oxide, zirconium oxide, iron oxide, silicon dioxide, and dispersing agent.

TECHNICAL INFORMATION

- Biocompatible product.
- Ready-to-use presentation, no mixing required.
- Setting time: 30–240 minutes after insertion into the canal.

Important: Dry with paper points but avoid over-drying the canal. Moisture from dentinal tubules will initiate the material's setting reaction.

- High alkalinity.

- High radiopacity.
- No setting shrinkage.
- Can be used with conventional techniques, single cone, and thermoplastic obturation.

WARNINGS

- Avoid contact with eyes or skin. In case of contact, rinse immediately with water.
- Avoid contact with oral mucosa. In case of contact, rinse with water and avoid ingestion of the product. If any sensitivity persists, seek medical attention immediately.
- Safety and efficacy have not been established in vulnerable populations.
- Keep out of reach of children.

Note: Any serious incident related to the device must be reported to the manufacturer and the competent health authority in the country where the user and/or patient is located.

PRECAUTIONS

- Wear protective glasses, mask, and gloves when handling the product. In case of contact with eyes or skin, wash with water.
- Use absolute isolation in all indications.
- The product is moisture-sensitive, so close the packaging properly with adequate pressure to prevent drying of the product.
- Use the product only after infection control, as the material will not set in an acidic environment.
- Avoid excessive extrusion of the product into the periapical region.
- Disinfect the applicator tip prior to use. After use, the tip should be discarded.
- Clean and disinfect the syringe for subsequent uses.
- Return the syringe properly closed, along with the silica sachet, inside the packaging, which must be properly sealed and stored according to the "storage" section.
- For subsequent use, ensure the silica sachet inside the envelope remains blue.

TECHNIQUE OF USE

Use of Applicator Tips

1. A disposable plastic cover on the syringe is recommended to prevent cross-contamination.
2. Confirm the material flow from the syringe before direct application into the canal.
3. Position the tip and syringe assembly properly in the canal and press the plunger gently to avoid excessive extrusion of the product.
4. Before closing the syringe, pull back the plunger to stop the flow.
5. Clean excess material from the syringe tip and close with adequate pressure after use to prevent moisture entry that causes drying.

6. Clean and disinfect the syringe for subsequent uses to avoid cross-contamination; if contaminated with saliva or blood, discard it.

A. OBTURATION OF ROOT CANALS OF PERMANENT TEETH

Traditional Technique – Lateral Compaction

1. Administer anesthesia, place absolute isolation, then perform biomechanical preparation of the canal.

2. Remove the sealing tip and attach the applicator tip to the syringe.

3. Dry the canal only with paper points without over-drying.

4. Apply EdgeUtopia™ BioCeramic Sealer while gradually retracting the syringe to achieve complete canal filling.

Important: press the plunger gently to avoid excessive extrusion of the product.

5. Insert the main gutta-percha cone according to the last file used, followed by accessory cones, all coated with EdgeUtopia™ BioCeramic Sealer.

6. Take a radiograph to verify proper canal filling.

7. Cut the cone at the desired height with heated instruments, followed by vertical compaction.

8. Remove excess material from canal walls with water, perform coronal sealing, and restoration.

Important:

- EdgeUtopia™ BioCeramic Sealer does not undergo shrinkage setting.

- EdgeUtopia™ BioCeramic Sealer can be used with conventional, single cone, and thermoplastic obturation techniques.

- Use conventional gutta-percha removal techniques.

Single Cone Technique

1. Administer anesthesia, place absolute isolation, then perform biomechanical preparation of the canal.

2. Remove the sealing tip and attach the applicator tip to the syringe.

3. Dry the canal only with paper points without over-drying.

4. Apply EdgeUtopia™ BioCeramic Sealer while gradually retracting the syringe to achieve complete canal filling.

Important: press the plunger gently to avoid excessive extrusion of the product.

5. Insert the gutta-percha cone.

6. Take a radiograph to verify proper canal filling.

7. Cut the cone at the desired height with heated instruments, followed by vertical compaction.

8. Remove excess material from canal walls with appropriate irrigating solutions, perform coronal sealing, and restoration.

B. USE DURING INTERNAL RESORPTION TREATMENT

1. Administer anesthesia and place absolute isolation.

2. Remove granulation tissue from the resorption area with sharp curettes.
3. Neutralize the environment with calcium hydroxide paste or an intracanal dressing paste.
4. Remove calcium hydroxide or an intracanal dressing paste in the following session.
5. Dry the canal only with paper points without over-drying.
6. Insert EdgeUtopia™ BioCeramic Sealer with the applicator tip throughout the canal, prioritizing the resorption site.
7. Insert accessory gutta-percha cones according to the lateral compaction technique.
8. Cut cones at the top of the resorption area with heated instruments. Perform vertical compaction to improve cement flow at the resorption site.
9. Take a radiograph to verify correct filling of the resorption site and canal.
10. Obturate the remainder of the canal with gutta-percha.
11. Perform coronal sealing with glass ionomer or another preferred material and restore the tooth.
12. Take a radiograph and maintain follow-up for at least two years.

Obturation removal

1. Use conventional gutta-percha removal techniques.

STORAGE

- Keep the product in a dry, ventilated place between 15 and 30°C with relative humidity below 60%.
- Do not store in a refrigerator.
- Do not store the product near ammonia, ammonium nitrate, or chlorine-containing products. Avoid using disinfectant solutions containing any of these ingredients.

DISPOSAL METHOD

Dispose of the product in accordance with national laws and local authority recommendations for proper disposal. Improper disposal may present biological or physical risks; handle and dispose of as hazardous biological material if contaminated.

ATTENTION: This product must be used according to the instructions described in this manual. The manufacturer is not responsible for failure or damage caused by incorrect handling or use of this product or by its use in situations of non-compliance with this manual.

Symbology

	Consult instructions for use
	For care, see accompanying documents
	Maximum and minimum temperature limit
	Single use accessories
	Keep away from sunlight

	Keep dry
	Non-sterile
	Medical device
	CE Mark
	Manufacturer
	Date of manufacture

	Expiration date
	Batch number
	Catalogue number
	Unique device identifier (UDI)
	European Representative
	Swiss representative

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