INSTRUCTIONS FOR USE – EDGEFILE® X7

EN	INSTRUCTIONS FOR USE	2
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА	6
CS	NÁVOD K POUŽITÍ	10
DA	BRUGSANVISNING	14
DE	GEBRAUCHSANWEISUNG	18
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ.....	22
ES	INSTRUCCIONES DE USO	26
ET	KASUTUSJUHEND.....	30
FI	KÄYTTÖOHJEET	34
FR	NOTICE D’UTILISATION.....	38
HR	UPUTE ZA UPORABU	42
HU	HASZNALATI UTASITAS.....	46
IT	ISTRUZIONI PER L’USO	50
LT	NAUDOJIMO INSTRUKCIJA	54
LV	LIETOŠANAS INSTRUKCIJAS	58
NL	GEBRUIKSAANWIJZING	62
NO	BRUKSANVISNING.....	66
PL	INSTRUKCJA UŻYWANIA.....	70
PT	INSTRUÇÕES DE UTILIZAÇÃO	74
RO	INSTRUCȚIUNI DE UTILIZARE	78
SK	NÁVOD NA POUŽITIE.....	82
SL	NAVODILA ZA UPORABO	86
SV	BRUKSANVISNING.....	90
TR	KULLANIM TALIMATLARI	94


FOR DENTAL USE ONLY

INSTRUCTIONS FOR USE – EDGEFILE® X7


1) GENERAL DESCRIPTION

Composition:

The instrument is made of an Annealed Heat Treated (Fire-Wire™) nickel-titanium alloy and an endo stop made of silicone. All files have constant taper with a maximum flute diameter of 1mm and delivered in sterile state (sterilized by irradiation).

Tip size and taper are identified using two color marks: the colored ring on the handle indicates the tip diameter, and the stopper color indicates the taper.

The grooves in the shank type indicates the taper: .04 or .06. The handle length is 11 mm:

2 grooves = .04

3 grooves = .06

Parameters for use:

These instruments are compatible with endodontic motors compliant with ISO 1797.

Recommended speed 300-500 rpm continuous rotation. Torque limit : 300 g·cm (2.94 N·cm).

2) INTENDED USE

The EdgeFile® X7 endodontic instruments are intended for root canal treatment.

3) INDICATIONS FOR USE

EdgeFile® X7 are indicated in case of pulp or periapical diseases that require root canal treatment. All instruments are designed for preparation and shaping of root canals.

4) CONTRAINDICATIONS AND ADVERSE REACTIONS

None known.

5) WARNINGS AND PRECAUTIONS

- These instruments have not been tested on children, pregnant nor breastfeeding women.
- These endodontic instruments must not be used on individuals with a known allergic sensitivity to nickel-titanium and silicone.
- Single use instruments must not be retreated for subsequent use. This could compromise their integrity.
- Each instrument may be used to treat up to four (4) canals during a single procedure. Discard after use.
- For your own safety, wear personal protective equipment (masks, gloves, and goggles) during the treatment and reprocessing of the instrument.
- Do not use if sterile packaging seal is broken, damaged or wet.
- Do not use the instruments after expiration date.
- Reprocess the instrument appropriately before disposal.
- Carefully read the label, marking on the packaging and this instruction for use to ensure the correct identification and use of device.
- Use a dental dam when using endodontic instruments to avoid inhalation or ingestion by the patient.
- Take several radiographs from different angles to determine the anatomy of the root canals (length, width and curvature).
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- Cracks, deformations, signs of corrosion, loss in color or marking are signs that the device is no longer able to achieve the required performance level and should be discarded.
- Regularly clean flutes and inspect instrument during use, and discard if it shows any signs of distortion or wear.
- The packaging of the instruments should be carefully examined before opening (packaging integrity, no humidity, no visual degradation of instruments) to ensure that the packaging's integrity has not been compromised during transport and storage.
- If the packaging has been opened, damaged, become wet or the expiry date has passed, the sterile state of the instruments inside the packaging is not guaranteed. Discard the device(s) that are single use and consult the reprocessing instructions available on EdgeEndo's website for multiple use device(s).
- All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority according to local regulations.

6) PERFORMANCE CHARACTERISTICS

EdgeFile® X7 rotary files are compliant with all applicable testing requirements of the ISO 3630 series of standards.

7) PROTOCOL FOR USE

Simplified Technique Guide

- 1) Create straight line access to canal orifices.
- 2) Locate canals and explore using stainless steel hand instruments.
- 3) Irrigate before each hand or rotary file with liquid NaOCl or 17% EDTA.
- 4) Using a #10 K file or larger file to determine the working length with an apex locator and radiographs in each canal. This will also form a glide path in each canal.
- 5) Use #17/04 file in one or more passes, alternating with small-sized hand files if necessary, until working length is reached. If more coronal flare is desired it can be achieved by incorporating the EdgeGlidePath file used in a brushing motion.
- 6) Next use #25/04 to working length passively; if instrument has not reached working length use smaller additional shaping instrument #20/04 to working length alternating between #20/04 and #25/04.
- 7) After #25/04 reaches working length if clinician desires a larger apical or canal shape, take a #25/06 or other additional larger 04 or 06 instruments (#30, #35, #40, etc...).

Canal Preparation:

- 1) Perform steps 1 to 4 in the “Simplified Technique Guide” above.
- 2) Shape and clean the root canal using the 04 or 06 taper technique:

EdgeFile® X7 04 Taper Crown Down for All Canals:

- 1) Start with a 25/04 rotary file. Take the 25/04 to resistance or working length (whichever occurs first). If resistance is met before reaching the working length, then go to a 20/04. Take the 20/04 to resistance or working length (whichever occurs first).
- 2) Repeat going from 25/04 to 20/04 until one of the files goes to the working length. On occasion a 17/04 may be needed to be used to reach the working length. Then repeat going from 25/04 to 20/04 until one of the files go to the working length.
- 3) If this is the tip size you desire, then obturate. If not, take the next largest file to length. Keep taking the next largest size to length until you achieve the tip size you desire, then obturate. Between each rotary file recapitulate with a #10 or #15 tip hand file to maintain glide path and help lubricate to the canal terminus.

EdgeFile® X7 06 Taper Crown Down for Straight to Mildly Curved Canals:

If a .06 taper is desired, follow the same crown-down technique described for the .04 taper.

- 1) Start with a 25/06 rotary file. Take the 25/06 to resistance or working length (whichever occurs first). If resistance is met before reaching the working length, then go to a 20/06. Take the 20/06 to resistance or working length (whichever occurs first).
- 2) Repeat going from 25/06 to 20/06 until one of the files go to the working length. On occasion a 17/06 or 17/04 may be needed to be used to reach the working length. Then repeat going from 25/06 to 20/06 until one of the files go to the working length.
- 3) If this is the tip size you desire, then obturate. If not, take the next largest file to length. Keep taking the

INSTRUCTIONS FOR USE

EN

next largest size to length until you achieve the tip size you desire, then obturate. Between each rotary file recapitulate with a #10 or #15 tip hand file to maintain glide path and help lubricate to the canal terminus.

Obturation of Canal Systems:

- When using thermal carriers such as EdgeCore™ X7 or EdgeFill™ X7 use size verifiers to determine the proper sized carrier.
- When using a master gutta percha cone that matches the largest file taken to length, remember sometimes you may need to drop down in cone tip size if the corresponding gutta percha to your final rotary file does not go to length.

8) CONDITIONS AND TERM OF STORAGE, SHELF-LIFE

The shelf-life is 5 years for sterile instruments. The use-by date is indicated on the packaging. Sterility is preserved and guaranteed for the whole shelf life unless the packaging is broken, or the storage conditions are not respected.

9) DISPOSAL

When a device reaches the end of its life, make sure that it is discarded in accordance with the applicable laws and regulations.

10) SYMBOLS

For explanation of symbols for IFUs and labels, please consult the document “Glossary of the symbols used by EdgeEndo” available on EdgeEndo’s website.

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