



Safety Data Sheet

Edge Paper Points™

1. Identification of the substance/preparation and of the company/undertaking

- Product name	: Edge Paper Points™
- Supplier	: METABIOMED CO., LTD. Head Office & Osong Factory) 270, Osongsaengmyeong1-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, Korea Cheongju Factory) 136, Mochung-ro, Seowon-gu, Cheongju-si, Chungcheongbuk-do, Korea Tel : (82-43) 216-0433, 0434 Fax : (82-43) 216-0438

2. Hazards identification

- Physical/chemical hazards	: None
- Human health hazards	: None
- Environmental hazards	: None

3. Composition/information on ingredients

Information on ingredients
- Pure Pulp
Hazardous Ingredients
- None

4. First aids measures

- Inhalation	: N/A
- Ingestion	: N/A
- Skin contact	: N/A
- Eye contact	: N/A

5. Fire fighting measures

<u>Extinguishing media</u>	
- Suitable	: Use any means suitable for extinguishing surrounding fire.
- Not suitable	: N/A
- Special fire fighting procedures	: None
- Unusable fire/explosion	: None

hazards

- Hazardous thermal : None

(de)composition products

- Protection of fire fighters : None

6. Accidental release measures

- Personal precautions : None

- Environmental precautions : None

- Methods for cleaning up : N/A

7. Precautions for safe handling

- Handling : Keep container closed when not in use.

- Storage : Keep in tightly closed container, store in a cool and dry area. Protect against sunlight, thermal.

8. Exposure controls/personal protection

- Precautions for use : None

- Respiratory protection : Not necessary under normal conditions of use.

- Skin and body : N/A

- Hands : N/A

- Eyes : With sufficient extraction or closed system; safety goggles.

9. Physical and chemical properties

- Physical state : Solid

- Color : Whitish

- Odor : Odorless

- Boiling point : Decomposes

- Melting point : N/A

- Vapor pressure : N/A

- Specific gravity(H₂O=1) : 0.65~0.70

- Solubility in water : Insoluble

- pH : N/A

- Flash point : Not-regulated

- Autoignition temperature : N/A

- Lower explosion limit : N/A

- Upper explosion limit : N/A

- Viscosity : N/A

10. Stability and reactivity

- Stability : Stable

- Condition to avoid : N/A

- **Materials to avoid** : None
- **Hazardous reactions** : None
- **Hazardous decomposition products** : CO

11. Toxicological information

- **Acute toxicity** : N/A
- **Primary effects to skin** : Not a primary dermal irritant
- **Primary effects to eyes** : Not an ocular irritant
- **Primary effects to respiratory system** : Not an acute inhalation hazard
- **Sensitization** : N/A
- **Additional toxicological info** : N/A

12. Ecological information

- **General information** : N/A
- **Classification of water endangerment** : N/A

13. Disposal considerations

- **Product** : Not considered a RCRA hazardous waste as identified in Fed. Reg. 40CFR261. If this product cannot be used as directed by the label, contact the local Environmental Control Agency or Agent.

14. Transport information

- | | |
|--|--|
| - DOT shipping name : Not-regulated | - <u>Air transport ICAO-Ti and IATA-DGR class</u> : N/A |
| - DOT hazard class : Not-regulated | - <u>IATA-DGR class</u> : Not-regulated |
| - UN number : N/A | - UN number : N/A |
| - <u>Sea transport IMDG/GGV sea class</u> : N/A | |
| - UN number : N/A | |
| - EmS number : N/A | |
| - MFGA number : N/A | |

15. Regulatory information

- **TSCA Inventory Status** : This product and/or all of its components are included on the TSCA Inventory of Chemical Substances.

- **SARA 311/312 hazard Categories** : None

- **SARA 302 extremely hazardous substances** : None

- SARA 313 toxic chemicals	<u>Chemical name</u>	<u>CAS-No.</u>	<u>Weight %</u>
	None		

- **Classification according to EEC Guidelines** : Class a rule 6

- **Classification according to VBF** : N/A

- **Canadian regulations require the inclusion of the following statements under "Regulatory Information"; "This product has been classified in accordance with the hazard criteria of the controlled Products Regulations and the MSDS contains all the information required by the Controlled Products Regulations."**

16. Other information

The information contained in the Safety data sheet is correct to the best of our knowledge at the date of issue. It is intended as a guide for the safe use, handling, disposal, storage and transportation and is not intended as warranty or as a specification. The information relates only to the product specified and may not be suitable for combinations with other materials or in processes other than those specifically described herein.

HISTORY

- **Date of printing** : 3-8-2015

- **Date previous issue** : 3-8-2015

- **Date of issue** : 3-8-2015

- **Prepared by** : METABIOMED CO., LTD.