

OPERATIVE CARBIDES	ROUND											TAPER FISSURE PLAIN CUT			TAPER FISSURE PLAIN CUT LONG		TAPER FISSURE CROSS CUT					TAPER FISSURE CROSS CUT LONG		DOMED FISSURE PLAIN CUT		
																										
	1/4	1/2	1	2	3	4	5	6	7	8	10	169	170	169L	699	700	701	702	703	701L	1156	1157	1158			
	005	006	008	010	012	014	016	018	021	023	027	009	010	009	009	010	012	016	021	012	008	010	012			
	•Δ	•Δ	•	•Δ	•	•Δ	•	•Δ	•	•Δ	•	•	•	•	•	•	•Δ	•	•Δ	•Δ	•Δ	•	•Δ	•		
FG (19mm)																										
FG Surg (25mm)																										
RA (22mm)																										
RA Surg (26mm)																										
HP (44.5mm)																										



	ROUND									INVERTED CONE			DIABLO		BARREL		BARREL		
DIAMONDS																			
25 Per Pack																			
Head Diameter (mm)																			
Head Length (mm)																			
OAL (mm)																			
Extra Fine																			
Fine																			
Medium																			
Coarse																			
Extra Coarse																			

Head Diameter (mm)
Head Length (mm)
OAL (mm)
Extra Fine
Fine
Medium
Coarse
Extra Coarse

ROUND END TAPER								ROUND END TAPER LONG				ROUND END TAPER					
																	
012	014	015	016	018	019	021		016	018	020		010	012	014	016	018	
8.0	8.0	8.0	8.0	8.0	8.0	8.0		9.0	9.0	9.0		10.0	10.0	10.0	10.0	10.0	
22.0	22.0	22.0	22.0	22.0	22.0	22.0		22.0	22.0	22.0		24.0	24.0	24.0	24.0	24.0	
198012M	198014F		198016F	198018F				198L016F						199014F			
198012C	198014C	198015C	198016C	198018C		198021C		198L016C	198L018C	198L020C		199010C	199012C	199014C	199016C	199018M	
198012XC	198014XC	198015XC	198016XC	198018XC	198019XC	198021XC		198016XC	198018XC						199016XC	199018XC	

	ROUND END TAPER SHORT			INVERTED CONE		INVERTED CONE		INVERTED CONE		PEAR			FLAME						
																			
Head Diameter (mm)																			
Head Length (mm)																			
OAL (mm)																			
Extra Fine																			
Fine																			
Medium																			
Coarse																			
Extra Coarse																			
	199014SXC	199016SXC		225016M	233012M	234012M		237009M	237010C	237012C		249010C	249010F	249012C	249012XF	249012F	249014C	249014XF	249016C

	FLAME					BUD					BUD SLENDER		EGG						
																			
Head Diameter (mm)		010	012	014	016		016	018	020	023	025		025		016	018	019	023	025
Head Length (mm)		10.0	10.0	10.0	10.0		3.2	4.6	4.6	5.0	5.0		5.5		3.4	3.4		4.5	4.6
OAL (mm)		24.0	24.0	24.0	24.0		19.0	20.0	20.0	20.0	20.0		20.0		19.0	19.0		19.0	19.0
Extra Fine								257018XF		257023XF	257025XF							277025XF	
Fine		250010F	250012F	250014F		257016F	257018F	257020F	257023F	257025F								277023F	277025F
Medium			250012M				257018M		257023M					277018M				277023M	
Coarse			250012C				257018C	257020C	257023C				277016C	277018C				277023C	
Extra Coarse				250016C					257023XC			263025C				277019C		277023C	

	TORPEDO			GINGIVAL CURETTAGE				INTERDENTAL			GROSS REDUCTION CYLINDER SHORT		GROSS REDUCTION TAPER SHORT		FLAME NEEDLE SHORT HEAD		ROUND EDGE TAPER SHORT				
																					
Head Diameter (mm)	012	013		012	014	016	018		016	018		017		018		008	010		019	026	
Head Length (mm)	8.0			8.0	8.0	8.0	8.0		5.0	5.0		7.6		8.2		4.0	4.0		4.0	4.0	
OAL (mm)	22.0			22.0	22.0	22.0	22.0		19.0	19.0		20.0		20.0		23.0	17.0		18.0	18.0	
Extra Fine																					
Fine																					
Medium																					
Coarse																					
Extra Coarse																					
																					
	289012C	289013C		298012M	298014M	298016M	298018M		465016F	465018F		513017SC		519018SC		540008F	536010SF		544019SC		544026SF
				298012C	298014C	298016C	298018C														
						298016XC	298018XC														

	ROUND EDGE TAPER SHORT	ROUND EDGE TAPER SHORT	ROUND	ROUND END TAPER	ROUND END CYLINDER	TAPER CONE ROUNDED TIP	TAPER CONE UNCOATED TIP
Head Diameter (mm)	016	017	014	016	012	011	015
Head Length (mm)	8.0	8.0	1.4	8.0	8.0	10.0	10.0
OAL (mm)	22.0	20.0	23.0	22.0	22.0	21.0	21.0
Extra Fine							
Fine		546017SF					
Medium	546016C	546017SM	697014M	Z198016M	Z141012M	E164011C	E165014LC
Coarse		546017SC	697017M		Z141014M	E164015C	
Extra Coarse		546017SXC			Z139012M	E164018C	