



EDGE BURS®

Quick Guide

EdgeEndo® offers the strongest, sharpest, and most durable burs on the market, delivering unmatched concentricity. Crafted with cutting-edge Swiss-made precision machinery, EdgeBurs® include Carbides, Diamonds, Finishers, and Polishers—designed to meet virtually every dental procedure and need.

Better Carbides



Available in
10 & 100 pack

Better Diamonds



Available in
25 packs

Better Finishers



Available in
5 packs

Better Metal-Cutting



Available in
10 & 100 packs

Better Polishers



Available in
25 packs

Better Zirconia




Available in
25 packs



EDGEENDO®

Scan to
learn more



Rev A 02/2025



161
10.0
016

Type
Diameter
Head Length (mm)
Overall Length
Head Color

ROUND										INVERTED CONE			DIABLO		BARREL		BARREL					
DIAMONDS 25 Per Pack Head Diameter (mm) Head Length (mm) OAL (mm) Extra Fine Fine Medium Coarse Extra Coarse																						
	010	012	014	016	018	021	023	025		012	016		014	033	033	037						
	1.0	1.2	1.4	1.6	1.8	2.1	2.3	2.5		1.5	1.5		1.5	4.5	5.5	7.0						
	19.0	19.0	19.0	19.0	19.0	19.0	19.0	19.0		19.0	19.0		19.0	19.5	21.0	21.5						
	001010M	001012M	001014M	001016M	001018M		001023M	001025F		010012C	010016C		032014C	038033C	038033XC	039033C	039037C					
	001010C	001012C	001014C	001016C	001018C	001021C	001023C										039037XC					
					001018XC		001023XC															
WHEEL - ROUNDED RIM				FLAT END CYLINDER		FLAT END CYLINDER			FLAT END CYLINDER		BEVELED CYLINDER		BEVELED CYLINDER		ROUND END CYLINDER							
 Head Diameter (mm) Head Length (mm) OAL (mm) Extra Fine Fine Medium Coarse Extra Coarse																						
	036	040	042		010		010	012	014	012	012	012	016	012								
	1.2	1.3	1.4		3.0		4.0	4.0	4.0	8.0	8.0	10.0	10.0	6.0								
	19.0	17.5	17.5		19.0		19.0	19.0	19.0	22.0	22.0	24.0	24.0	18.5								
	067036C	067040SC	067042SC		108010F	108010M	109010M	109012M	109014C	111012C	130012C	131012C	131016C	140012SXC								
		067040SXC					109010C	109012C														
ROUND END CYLINDER		ROUND END CYLINDER- SHORT		KS CYLINDER SHORT		END CUTTING CYLINDER SHORT			POINTED CONE		POINTED CONE											
 Head Diameter (mm) Head Length (mm) OAL (mm) Extra Fine Fine Medium Coarse Extra Coarse																						
	012	014	012	014	014	012	012	014	016	014	016		014	016	017	010	014	014	016	016	017	
	8.0	8.0	8.0	8.0	8.0	7.0	21.0	21.0	21.0	8.0	8.0		10.0	10.0	10.0	10.0	10.0	10.0	10.0	10.0	10.0	
	22.0	22.0	22.0	22.0	19.5	20.0				22.0	22.0		24.0	24.0	24.0	24.0	24.0	24.0	24.0	24.0	24.0	
	141012C	141014C	141012XC	141014XC	141014SXC	146012SC	150012SM	150014SM	150016SM	165014XF	165014F	165016F	165014M	165014C	165014XC	166010F	166010M	166014F	166014M	166016XF	166016F	166017F
POINTED CONE		FLAT END TAPER		FLAT END TAPER		ROUND END TAPER SHORT		ROUND END TAPER		ROUND END TAPER		ROUND END TAPER SHORT										
 Head Diameter (mm) Head Length (mm) OAL (mm) Extra Fine Fine Medium Coarse Extra Coarse																						
	012	014	016	018	014	016	018	016	018	017	024	016	022	025			023					
	11.5	8.0	8.0	8.0	8.0	10.0	10.0	10.0	10.0	6.0	7.0	6.0	8.0	8.0			8.5					
	24.5	22.0	22.5	22.5	22.0	22.5	22.5	24.0	24.0	18.5	18.5	16.5	22.0	22.0			20.0					
	167012XF	172014M	172016M	172018C	172014XC	172016XC	172018XC	173016C	173018C	197017SXC	197024SXC	197016SHC	198022C	198025C	198022XC	198025XC	198023SC					
	167012C	172014C	172016C	172018C	172014XC	172016XC	172018XC	173016C	173018C													

Head Diameter (mm)
Head Length (mm)
OAL (mm)
Extra Fine
Fine
Medium
Coarse
Extra Coarse

ROUND END TAPER SHORT			INVERTED CONE		INVERTED CONE		INVERTED CONE		PEAR			FLAME						
																		
014	016		016	012	012	009	010	012	010	012	010	012	014	016				
10.0	10.0		4.0	2.8	4.0	1.5	3.0	3.0	8.0	8.0	8.0	8.0	8.0	8.0				
22.0	22.0		19.0	19.0	19.0	19.0	19.0	19.0	22.0	22.0	22.0	22.0	22.0	22.0				
199014SXC	199016SXC		225016M	233012M	234012M	237009M	237010C	237012M	249010F	249010C	249012XF	249012F	249012M	249012C	249012XC	249014XF	249014C	249016C

TORPEDO	GINGIVAL CURETTAGE						INTERDENTAL	GROSS REDUCTION CYLINDER SHORT	GROSS REDUCTION TAPER SHORT	FLAME NEEDLE SHORT HEAD	ROUND EDGE TAPER SHORT		
 012 8.0 22.0	 013	 012 8.0 22.0	 014 8.0 22.0	 016 8.0 22.0	 018 8.0 22.0	 016 5.0 19.0	 018 5.0 19.0	 017 7.6 20.0	 018 8.2 20.0	 008 4.0 23.0	 010 4.0 17.0	 019 4.0 18.0	 026 4.0 18.0
289012C 289013C	298012M 298012C	298014M 298014C	298016F 298016C 298016XC	298018M 298018C 298018XC		465016F 465018F		513017SC	519018SC	540008F	536010SF	544019SC	544026SF

ROUND EDGE TAPER SHORT		ROUND EDGE TAPER SHORT		ROUND		ROUND END TAPER		ROUND END CYLINDER			TAPER CONE ROUNDED TIP			TAPER CONE UNCOATED TIP										
	016 8.0 22.0		017 8.0 20.0	ZIRCONIA CUTTING DIAMONDS		014 1.4 23.0		017 1.7 23.0		016 8.0 22.0		012 8.0 22.0		014 8.0 22.0		012 4.0 19.0		011 10.0 21.0		015 10.0 21.0		018 10.0 21.0		015 10.0 21.0
546016C		546017SF 546017SM 546017SC 546017SXC			697014M	697017M	Z198016M	Z141012M	Z141014M	Z139012M	E164011C	E164015C	E164018C	E165014LC										