



EDGEBURS® QUICK GUIDE

BETTER CARBIDES



\$15.88/10 pack
\$109.88/100 pack

BETTER DIAMONDS



\$39.88/25 pack

BETTER FINISHERS



\$19.88/5 pack

BETTER METAL-CUTTING



EdgeOrca


Replaces Your Current Metal-Cutter
\$19.88/10 pack • \$158.88/100 pack

BETTER POLISHERS



\$39.88/25 pack

BETTER ZIRCONIA



\$39.88/25 pack

BETTER BURS, ^{up to} 1/2 THE PRICE

EdgeEndo® is selling the strongest, sharpest and most durable burs available with unparalleled concentricity. They are manufactured using the latest state-of-the-art Swiss-made precision machine tools. EdgeBurs® include Carbides, Diamonds, Finishers and Polishers. We have a bur for essentially every procedure and dental need.



EDGE DENTAL®
SPECIALTIES

1-855-985-3636

EDGEENDO.COM

	ROUND										TAPER FISSURE PLAIN CUT			TAPER FISSURE PLAIN CUT LONG		TAPER FISSURE CROSS CUT					TAPER FISSURE CROSS CUT LONG		DOMED FISSURE PLAIN CUT		
OPERATIVE CARBIDES																									
Head Dia. (mm)	1/4	1/2	1	2	3	4	5	6	7	8	10	169	170	169L	699	700	701	702	703	701L	1156	1157	1158		
FG (19mm)	Δ	Δ		Δ		Δ		Δ		Δ							Δ			Δ		Δ			
FG Surg (25mm)				Δ		Δ		Δ		Δ							Δ					Δ			
RA (22mm)				Δ		Δ		Δ		Δ															
RA Surg (26mm)																									
HP (44.5mm)						β		β		$\beta\Delta$							β	$\beta\Delta$	$\beta\Delta$						

	DOMED FISSURE CROSS CUT		FLAT FISSURE PLAIN CUT			FLAT FISSURE CROSS CUT			FLAT FISSURE CROSS CUT LONG		PEAR			PEAR LONG		INVERTED CONE				8 FLUTES		AMALGAM PREP	
Head Dia. (mm)	1557	1558	56	57	58	556	557	558	557L	330	331	332	331L	33 1/2	34	35	37	1171 1/2	1172	245			
FG (19mm)	Δ	Δ				Δ	Δ	Δ		Δ	Δ					Δ			Δ	Δ			
FG Surg (25mm)																							
FG Short (16mm)																							
RA (22mm)																		Δ	Δ				

	ROUND		BULLET			EGG			UF EGG		TAPER POINTED				TAPER ROUND END			FLAME		LONG FLAME	
FINISHING CARBIDES																					
5 Per Pack	7006	7008	7801	7802	7803	7404	7406	7408	9406	7610	7611	7612	7613	7642	7664	7675	7104	7106	48L	48L	
US #	41	41	247	247	247	379	379	379	379UF	132	133	134	135	375R	375R	375R	46	46	48L	48L	
Fig #	018	023	009	010	012	014	018	023	081	008	010	014	014	010	014	016	014	018	010	012	
Diameter	1.8	2.3	4.0	3.4	3.4	3.3	3.7	3.9	3.4	3.0	4.2	6.0	9.0	7.0	8.2	8.0	3.5	3.7	8.0	8.0	
Head Length (mm)																					

	NEEDLE				UF NEEDLE		10 BLADE SAFE END POINTED TAPER		12 BLADE SAFE END TAPER		ORCA CARBIDE TURBO						ORCA DX	
US #	7901	7902	7903	9903	EF9-UF	OF1-UF	1558T	1557T	331T	330T	245T	557T	PR-DX					
Fig #	246	246	246	246-UF	135SE-UF	379-UF	012	010	010	008	009	010	012					
Diameter	009	010	012	012	014	023	4.0	3.9	1.9	2.0	2.7	3.7	4.0					
Head Length (mm)	3.8	3.8	3.9	5.5	9.0	4.3	Δ	Δ		Δ		Δ						

	SHINE - FINE GRIT			REVEAL - MEDIUM GRIT			SPECIALTY CARBIDES			151-Z		ENDO ACCESS		LINDEMANN CROSS CUT		D LINDEMANN CROSS CUT		LINDEMANN	
POLISHERS																			
25 Per Roll	Point	Cup	Knife	Point	Cup	Knife	151-Z	151L-Z	152	152	162	162A	161						
Type																			
Diameter	5.5	6.0	10.0	5.5	6.0	10.0	016	016	016	016	10.0	9.2	10.0						
Head Length (mm)	10.7	7.4	5.1	10.7	7.4	5.1	11	11	9.0	10.0	016	016	016						
Overall Length	28.0	24.8	22.6	28.0	24.8	22.6	23.0	28.0	23										
Head Color	\bullet	\bullet	\bullet	\circ	\circ	\circ													

	ROUND									INVERTED CONE			DIABLO		BARREL		BARREL		
DIAMONDS																			
25 Per Pack																			
Head Diameter (mm)	010 012 014 016 018 021 023 025									012 016			014		033		033 037		
Head Length (mm)	1.0 1.2 1.4 1.6 1.8									1.5 1.5			1.5		4.5		5.5 7.0		
OAL (mm)	19.0 19.0 19.0 19.0 19.0 19.0 19.0 19.0									19.0 19.0			19.0		19.5		21.0 21.5		
Extra Fine																			
Fine																			
Medium	001010M 001012M 001014M 001016M 001018M									001023M			001025F						
Coarse	001010C 001012C 001014C 001016C 001018C									001021C 001023C			001023M		038033C 038033XC		039033C 039037C 039037XC		
Extra Coarse	001018XC 001018XC									001023XC									
	WHEEL - ROUNDED RIM				FLAT END CYLINDER		FLAT END CYLINDER			FLAT END CYLINDER		BEVELED CYLINDER		BEVELED CYLINDER		ROUND END CYLINDER			
Head Diameter (mm)	036 040 042				010		010 012 014			012		012		012 016		012			
Head Length (mm)	1.2 1.3 1.4				3.0		4.0 4.0 4.0			8.0		10.0		10.0 24.0		6.0			
OAL (mm)	19.0 17.5 17.5				19.0		19.0 19.0 19.0			22.0		22.0		24.0 24.0		18.5			
Extra Fine	067036C 067040SC 067042SC				108010F		109010M 109010C 109012M 109012C 109014C			111012C		130012C		131012C 131016C		140012SXC			
Fine																			
Medium					108010M														
Coarse	067040SXC																		
Extra Coarse																			
	ROUND END CYLINDER		ROUND END CYLINDER- SHORT		KS CYLINDER SHORT		END CUTTING CYLINDER SHORT			POINTED CONE		POINTED CONE							
Head Diameter (mm)	012 014		014		012		012 014 016			014 016		010 014 016 017							
Head Length (mm)	8.0 8.0		8.0		7.0		21.0 21.0 21.0			8.0 8.0		10.0 10.0 10.0 24.0							
OAL (mm)	22.0 22.0		19.5		20.0		21.0 21.0 21.0			22.0 22.0		24.0 24.0 24.0 24.0							
Extra Fine	141012C 141014C		141014SXC		146012SC		150012SM 150014SM 150016SM			165014XF 165014F 165014M 165014C 165014XC		166014XF 166014F 166016XF 166016F 166017F							
Fine	141012XC 141014XC									165016F		166010F 166010M 166014M 166014C 166016XC							
Medium																			
Coarse																			
Extra Coarse																			
	POINTED CONE		FLAT END TAPER			FLAT END TAPER			ROUND END TAPER SHORT		ROUND END TAPER		ROUND END TAPER SHORT						
Head Diameter (mm)	012		014 016 018			016 018			017 024		016		022 025						
Head Length (mm)	11.5		8.0 8.0 8.0			10.0 10.0			6.0 7.0		6.0		8.0 8.0						
OAL (mm)	24.5		22.0 22.5 22.5			24.0 24.0			18.5 18.5		16.5		22.0 22.0						
Extra Fine	167012XF										197016SHC		198022C 198025C						
Fine													198022XC 198025XC						
Medium	167012C		172014M 172016M			173016C 173018C			197017SXC 197024SXC				198023XC						
Coarse			172014C 172016C 172018C			173018XC													
Extra Coarse			172014XC 172016XC 172018XC																

DIAMONDS	ROUND END TAPER								ROUND END TAPER LONG			ROUND END TAPER				
	Head Diameter (mm)	012	014	015	016	018	019	021	016	018	020	010	012	014	016	018
Head Length (mm)	8.0	8.0	8.0	8.0	8.0	8.0	8.0	9.0	9.0	9.0	10.0	10.0	10.0	10.0	10.0	
OAL (mm)	22.0	22.0	22.0	22.0	22.0	22.0	22.0	22.0	22.0	22.0	24.0	24.0	24.0	24.0	24.0	
Extra Fine																
Fine		198014F		198016F	198018F			198L016F					199014F			
Medium	198012M	198014M		198016M	198018M										199018M	
Coarse	198012C	198014C	198015C	198016C	198018C		198021C	198L016C	198L018C	198L020C	199010C	199012C	199014C	199016C	199018C	
Extra Coarse	198012XC	198014XC	198015XC	198016XC	198018XC	198019XC	198021XC	198016XC	198018XC					199016XC	199018XC	

	ROUND END TAPER SHORT		INVERTED CONE		INVERTED CONE		INVERTED CONE		PEAR			FLAME			
	Head Diameter (mm)	014	016	016	016	012	012	012	012	009	010	012	010	012	014
Head Length (mm)	10.0	10.0	4.0	19.0	2.8	19.0	4.0	19.0	1.5	3.0	3.0	8.0	8.0	8.0	8.0
OAL (mm)	22.0	22.0	19.0		19.0		19.0		19.0	19.0	19.0	22.0	22.0	22.0	22.0
Extra Fine															
Fine												249010F	249012F	249014F	
Medium			225016M		233012M		234012M		237009M		237012M	249010C	249012C	249014C	249016C
Coarse										237010C	237012C	249010C	249012C	249014C	249016C
Extra Coarse	199014SXC	199016SXC										249010C	249012C	249014C	249016C

	FLAME				BUD					BUD SLENDER		EGG				
	Head Diameter (mm)	010	012	014	016	016	018	020	023	025	025	016	018	019	023	025
Head Length (mm)	10.0	10.0	10.0	10.0	3.2	4.6	4.6	5.0	5.0	5.5	3.4	3.4		4.5	4.6	
OAL (mm)	24.0	24.0	24.0	24.0	19.0	20.0	20.0	20.0	20.0	20.0	19.0	19.0		19.0	19.0	
Extra Fine						257018XF		257023XF	257025XF					277023XF	277025XF	
Fine	250010F	250012F	250014F		257016F	257018F	257020F	257023F	257025F					277018F	277023F	
Medium		250012M				257018M		257023M						277018M	277023M	
Coarse		250012C		250016C		257018C	257020C	257023C		263025C				277016C	277023C	
Extra Coarse							257020C	257023C	257023XC					277016C	277023C	

	TORPEDO		GINGIVAL CURETTAGE				INTERDENTAL			GROSS REDUCTION CYLINDER SHORT		GROSS REDUCTION TAPER SHORT		FLAME NEEDLE SHORT HEAD		ROUND EDGE TAPER SHORT	
	Head Diameter (mm)	012	013	012	014	016	018	016	018	017	018	008	010	019	026		
Head Length (mm)	8.0		8.0	8.0	8.0	8.0	5.0	5.0	7.6	8.2	4.0	4.0	4.0	4.0			
OAL (mm)	22.0		22.0	22.0	22.0	22.0	19.0	19.0	20.0	20.0	23.0	17.0	18.0	18.0			
Extra Fine																	
Fine					298016F		465016F	465018F			540008F			544026SF			
Medium			298012M	298014M	298016M	298018M											
Coarse	289012C	289013C	298012C	298014C	298016C	298018C			513017SC	519018SC			536010SF	544019SC			
Extra Coarse					298016XC	298018XC											

	ROUND EDGE TAPER SHORT	ROUND EDGE TAPER SHORT	ROUND		ROUND END TAPER	ROUND END CYLINDER			TAPER CONE ROUNDED TIP			TAPER CONE UNCOATED TIP	
	Head Diameter (mm)	016	017	014	017	016	012	014	012	011	015	018	015
Head Length (mm)	8.0	8.0	1.4	1.7	8.0	8.0	8.0	4.0	10.0	10.0	10.0	10.0	
OAL (mm)	22.0	20.0	23.0	23.0	22.0	22.0	22.0	19.0	21.0	21.0	21.0	21.0	
Extra Fine													
Fine		546017SF											
Medium		546017SM			697014M	697017M	Z198016M	Z141012M	Z141014M	Z139012M			
Coarse	546016C	546017SC									E164011C	E164015C	E164018C
Extra Coarse		546017SXC											E165014LC

ZIRCONIA CUTTING DIAMONDS