

**EDGEENDO<sup>®</sup>**

# **PRODUCT CATALOG**

**2024**

**FREE GROUND SHIPPING ON ALL ORDERS!**

**100% MONEY BACK GUARANTEE\***

EdgeEndo.com • 1-855-985-3636

Click to follow us on:   

\*See website for details

# TABLE OF CONTENTS

## Shaping

EdgeUtopia Files . . . . .	5
EdgeX7 Utopia . . . . .	5
EdgeSequel Utopia . . . . .	7
EdgeTaper Blaze Utopia . . . . .	9
EdgeOne Blaze Utopia . . . . .	12

## Files

EdgeCoil . . . . .	14
EdgeCoil Glidepath . . . . .	14
EdgeEvolve . . . . .	15
EdgeFile . . . . .	16
EdgeFind. . . . .	22
EdgeGlidepath . . . . .	23
EdgeMorphology . . . . .	24
EdgeOne Fire . . . . .	25
EdgeOne Glidepath . . . . .	27
EdgeOne Platinum . . . . .	29
EdgeGlidepath . . . . .	29
EdgeSequel Sapphire. . . . .	30
EdgeTaper . . . . .	33
EdgeTaper Platinum. . . . .	34
EdgeV-Taper HT . . . . .	36
EdgeTaper Retreat. . . . .	39
EdgeGates . . . . .	40

HandFiles . . . . .	41
---------------------	----

## Power

EdgeApex HP . . . . .	46
EdgeApex . . . . .	48
EdgePack . . . . .	49
EdgeFlow . . . . .	52

## Irrigation

EdgePRO. . . . .	55
Irrigation Syringes . . . . .	56

Cleansers. . . . .	57
Ultrasonics. . . . .	59

## Obturation

Obturators. . . . .	62
Points . . . . .	71
Sealers . . . . .	84
Temp . . . . .	86

## Supplies

EndoRing . . . . .	87
--------------------	----

(Click Category Above To Jump To Section)





# EDGE X7 UTOPIA™

**Enhanced Technology. Remarkably Beautiful. Ultra Secure Handle.**

The New EdgeX7 Utopia™ delivers a Technology-Enhanced EdgeFile X7® with a Remarkably-Beautiful gold blade featuring our FireWire Blaze heat treatment. Plus, we added our new ultra secure, Ultra-Fit handle and shank.

## Pre-Sterilized



Parabolic Cross Section

Available in 3 lengths: 21, 25 & 29 mm, 04 & 06 File Size

Maximum flute diameter 1 mm



High Cutting Efficiency



Coated Finish



No Bounce Back



High Strength


## FireWire Blaze

FireWire Blaze—Years to develop and worth every minute. We wanted to outdo our original FireWire, and with FireWire Blaze, we believe we have! Our FireWire Blaze proprietary heat treatment provides a super sharp cutting edge and accurate blade surface. The quality is visible!

Utopia Alternative for EdgeFile® X7

**QTY 21 mm \$35.95 /Pre-Sterilized 6 pk**

	17.04		PURPLE	X7U170421SE
	17.06		PURPLE	X7U170621SE
	20.04		YELLOW	X7U200421SE
	20.06		YELLOW	X7U200621SE
	25.04		RED	X7U250421SE
	25.06		RED	X7U250621SE
	30.04		BLUE	X7U300421SE
	30.06		BLUE	X7U300621SE
	35.04		GREEN	X7U350421SE
	35.06		GREEN	X7U350621SE
	40.04		BLACK	X7U400421SE
	40.06		BLACK	X7U400621SE
	45.04		WHITE	X7U450421SE
	45.06		WHITE	X7U450621SE



**QTY 25 mm \$35.95 /Pre-Sterilized 6 pk**

	17.04		PURPLE	X7U170425SE
	17.06		PURPLE	X7U170625SE
	20.04		YELLOW	X7U200425SE
	20.06		YELLOW	X7U200625SE
	25.04		RED	X7U250425SE
	25.06		RED	X7U250625SE
	30.04		BLUE	X7U300425SE
	30.06		BLUE	X7U300625SE
	35.04		GREEN	X7U350425SE
	35.06		GREEN	X7U350625SE
	40.04		BLACK	X7U400425SE
	40.06		BLACK	X7U400625SE
	45.04		WHITE	X7U450425SE
	45.06		WHITE	X7U450625SE

**QTY 29 mm \$35.95 /Pre-Sterilized 6 pk**

	17.04		PURPLE	X7U170429SE
	17.06		PURPLE	X7U170629SE
	20.04		YELLOW	X7U200429SE
	20.06		YELLOW	X7U200629SE
	25.04		RED	X7U250429SE
	25.06		RED	X7U250629SE
	30.04		BLUE	X7U30049SE
	30.06		BLUE	X7U300629SE
	35.04		GREEN	X7U350429SE
	35.06		GREEN	X7U350629SE
	40.04		BLACK	X7U400429SE
	40.06		BLACK	X7U40069SE
	45.04		WHITE	X7U450429SE
	45.06		WHITE	X7U450629SE



[Click Here to Order](#)

**Scan or Click For More Information**

Utopia Alternative for EdgeFile® X7

### Assorted Packs .04 \$35.95 /6 pk

QTY (17, 20, 25, 25, 30, 35)


	21MM	X7UASST0421SE
	25MM	X7UASST0425SE

### Assorted Packs .06 \$35.95 /6 pk

QTY (17, 20, 25, 25, 30, 35)

	21MM	X7UASST0621SE
	25MM	X7UASST0625SE

Glide Path      Initial Shaping      Optional Shaping



[Click Here to Order](#)


Scan or Click For  
More Information

# EDGESEQUEL UTOPIA™

## Enhanced Technology. Remarkably Beautiful. Ultra Secure Handle.

EdgeSequel Utopia™ delivers a Technology-Enhanced Constant 02, 04, 06 Taper System with a Remarkably-Beautiful blue blade featuring our FireWire Blaze Blue heat treatment. Plus, we added our ultra secure, Ultra-Fit handle and shank.

### Pre-Sterilized



Triangular Cross Section  
Available in 3 lengths: 21, 25 & 31 mm



Highly Flexible



No Bounce Back



High Strength

---

### FireWire Blaze Blue


















FireWire Blaze Blue is our newest proprietary heat treated NiTi. You will find that its flexibility, and feel provides unparalleled performance and shapes. It Is Truly, Amazing!

Utopia Alternative for EndoSequence® and Vortex®





**QTY 19 mm \$35.95 /Pre-Sterilized 6 pk**


	25.08	YELLOW	ES250819SE
--	-------	--------	------------

**QTY 21 mm \$35.95 /Pre-Sterilized 6 pk**














	10.02		PURPLE	ES100221SE
	10.04		PURPLE	ES100421SE
	15.02		WHITE	ES150221SE
	15.04		WHITE	ES150421SE
	15.06		WHITE	ES150621SE
	20.02		YELLOW	ES200221SE
	20.04		YELLOW	ES200421SE
	25.04		RED	ES250421SE
	25.06		RED	ES250621SE
	30.02		BLUE	ES300221SE
	30.04		BLUE	ES300421SE
	30.06		BLUE	ES300621SE
	35.04		GREEN	ES350421SE
	35.06		GREEN	ES350621SE
	40.04		BLACK	ES400421SE
	40.06		BLACK	ES400621SE
	50.04		YELLOW	ES500421SE

**QTY 25 mm \$35.95 /Pre-Sterilized 6 pk**

	10.02		PURPLE	ES100225SE
	10.04		PURPLE	ES100425SE
	15.02		WHITE	ES150225SE
	15.04		WHITE	ES150425SE



**QTY 25 mm \$35.95 /Pre-Sterilized 6 pk**

	15.06		WHITE	ES150625SE
	20.02		YELLOW	ES200225SE
	20.04		YELLOW	ES200425SE
	25.04		RED	ES250425SE
	25.06		RED	ES250625SE
	30.02		BLUE	ES300225SE
	30.04		BLUE	ES300425SE
	30.06		BLUE	ES300625SE
	35.04		GREEN	ES350425SE
	35.06		GREEN	ES350625SE
	40.04		BLACK	ES400425SE
	40.06		BLACK	ES400625SE
	50.04		YELLOW	ES500425SE











**Scan or Click For More Information**

[Click Here to Order](#)

Vortex® is a registered trademark of Dentsply Tulsa Dental Specialties. EndoSequence® is a registered trademark of Brassler USA

Utopia Alternative for EndoSequence® and Vortex®

**QTY 31 mm \$35.95 /Pre-Sterilized 6 pk**

	15.04		WHITE	ES150431SE
	25.04		RED	ES250431SE
	25.06		RED	ES250631SE
	30.04		BLUE	ES300431SE
	30.06		BLUE	ES300631SE
	35.06		GREEN	ES350631SE
	40.04		BLACK	ES400431SE
	50.04		YELLOW	ES500431SE

**Assorted Packs \$25.88 /3 pk**

QTY (15/04, 25/04, 30/04)

	21MM	ASST15300421SE
	25MM	ASST15300425SE
	31MM	ASST15300431SE

**Assorted Packs \$35.95 /6 pk**

04 taper;21mm; Size 15, 20, 25, 30, 35

QTY + 25/.08 19mm

	21MM	V2ASST0421SE
	25MM	V2ASST0425SE
	31MM	V2ASST0431SE

**Assorted Packs \$25.88 /3 pk**

QTY (15/04, 25/04, 25/06)

	21MM	ASST15250621SE
	25MM	ASST15250625SE
	31MM	ASST15250631SE

[Click Here to Order](#)



**Scan or Click For  
More Information**

Vortex® is a registered trademark of Dentsply Tulsa Dental Specialties. EndoSequence® is a registered trademark of Brassler USA




# EDGETAPER BLAZE UTOPIA™

## Enhanced Technology. Remarkably Beautiful. Ultra Secure Handle.

The NEW EdgeTaper Blaze Utopia™ is an alternative to ProTaper Gold®, with more efficient, resilient blade, than EdgeTaper Platinum™. It can be used with the same technique as ProTaper Gold® and EdgeTaper Platinum™. It features our FireWire Blaze Heat-Treatment making it efficient, safe, and flexible. Plus, we added our new ultra secure, Ultra-Fit Handle and all EdgeUtopia™ files are provided in pre-sterilized blister packs.

### Pre-Sterilized



Tri-Shaped Cross Section  
Available in 3 lengths: 21, 25 & 31 mm




### FireWire Blaze






FireWire Blaze—Years to develop and worth every minute. We wanted to outdo our original FireWire, and with FireWire Blaze, we believe we have! Our FireWire Blaze proprietary heat treatment provides a super sharp cutting edge and accurate blade surface. The quality is visible!



**QTY 19 mm \$35.95 /Pre-Sterilized 6 pk**

	SX		YELLOW	ETBUSX19SU
--	----	---	--------	------------






**QTY 21 mm \$35.95 /Pre-Sterilized 6 pk**

	S1		PURPLE	ETBUS121SU
	S2		WHITE	ETBUS221SU
	F1		YELLOW	ETBUF121SU
	F2		RED	ETBUF221SU
	F3		BLUE	ETBUF321SU






**Assorted Packs \$35.95 /6 pk**

	21MM	ETBUA21SU
	25MM	ETBUA25SU
	31MM	ETBUA31SU

**QTY 25 mm \$35.95 /Pre-Sterilized 6 pk**

	S1		PURPLE	ETBUS125SU
	S2		WHITE	ETBUS225SU
	F1		YELLOW	ETBUF125SU
	F2		RED	ETBUF225SU
	F3		BLUE	ETBUF325SU

**QTY 31 mm \$35.95 /Pre-Sterilized 6 pk**

	S1		PURPLE	ETBUS131SU
	S2		WHITE	ETBUS231SU
	F1		YELLOW	ETBUF131SU
	F2		RED	ETBUF231SU
	F3		BLUE	ETBUF331SU



**Scan or Click For More Information**

[Click Here to Order](#)

ProTaper Gold® is a registered trademark of Dentsply Tulsa Dental Specialties.




# EDGEONE BLAZE UTOPIA™

## Enhanced. Beautiful. Utopia.

EdgeOne Blaze Utopia™ delivers enhanced reciprocating technology. EdgeOne Blaze Utopia™ is our premier single shaping file reciprocating system. EdgeOne Blaze Utopia™ is an alternative to WaveOne Gold®. It features our FireWire Blaze Heat-Treatment making it efficient, safe, and flexible. Plus, we added our new ultra secure, Ultra-Fit Handle and it's pre-sterilized.

### Pre-Sterilized Single File Shaping System



S-Shaped Cross Section  
For Increased Debris Removal

Available in 3 lengths: 21, 25 & 31 mm



Highly  
Flexible



No Bounce  
Back



High  
Strength

---

### FireWire Blaze

FireWire Blaze—Years to develop and worth every minute. We wanted to outdo our original FireWire, and with FireWire Blaze, we believe we have! Our FireWire Blaze proprietary heat treatment provides a super sharp cutting edge and accurate blade surface. The quality is visible!



**QTY 21 mm \$25.88 /Pre-Sterilized 3 pk**

	R-GLIDE	<input type="checkbox"/>	WHITE	EOBURG1421SU
	R20	<input checked="" type="checkbox"/>	YELLOW	EOBU2021SU
	R25	<input checked="" type="checkbox"/>	RED	EOBU2521SU
	R35	<input checked="" type="checkbox"/>	GREEN	EOBU3521SU
	R45	<input type="checkbox"/>	WHITE	EOBU4521SU

**QTY 25 mm \$25.88 /Pre-Sterilized 3 pk**

	R-GLIDE	<input type="checkbox"/>	WHITE	EOBURG1425SU
	R20	<input checked="" type="checkbox"/>	YELLOW	EOBU2025SU
	R25	<input checked="" type="checkbox"/>	RED	EOBU2525SU
	R35	<input checked="" type="checkbox"/>	GREEN	EOBU3525SU
	R45	<input type="checkbox"/>	WHITE	EOBU4525SU

**QTY 31 mm \$25.88 /Pre-Sterilized 3 pk**

	R-GLIDE	<input type="checkbox"/>	WHITE	EOBURG1431SU
	R20	<input checked="" type="checkbox"/>	YELLOW	EOBU2031SU
	R25	<input checked="" type="checkbox"/>	RED	EOBU2531SU
	R35	<input checked="" type="checkbox"/>	GREEN	EOBU3531SU
	R45	<input type="checkbox"/>	WHITE	EOBU4531SU

**Assorted Packs \$34.88 /4 pk**

QTY R20, R25, R35, R45

	21MM	EOBUA21SU
	25MM	EOBUA25SU
	31MM	EOBUA31SU



**Scan or Click For More Information**

[Click Here to Order](#)

WaveOne® Gold is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGE COIL™

## Heat-Treated FireWire™ NiTi

(Unique one file system - Reverse Reciprocating)

QTY **\$35.88 /3 pk**

<input type="checkbox"/>	21mm	ECR21
<input type="checkbox"/>	25mm	ECR25
<input type="checkbox"/>	31mm	ECR31



[Click Here to Order](#)

# EDGE COIL GLIDEPATH™

Reciprocating GlidePath

QTY **21 mm \$25.88 /3 pk**

<input type="checkbox"/>	16.02	<input type="checkbox"/>	WHITE	ECGP21
--------------------------	-------	--------------------------	-------	--------

QTY **25 mm \$25.88 /3 pk**

<input type="checkbox"/>	16.02	<input type="checkbox"/>	WHITE	ECGP25
--------------------------	-------	--------------------------	-------	--------

QTY **31 mm \$25.88 /3 pk**

<input type="checkbox"/>	16.02	<input type="checkbox"/>	WHITE	ECGP31
--------------------------	-------	--------------------------	-------	--------



[Click Here to Order](#)



Scan or Click For  
More Information



**QTY 21 mm \$35.95 /6 pk**

	17.04	☐	WHITE	EE170421
	17.06	☐	WHITE	EE170621
	17.08	☐	WHITE	EE170821
	17.10	☐	WHITE	EE171021
	20.04	■	YELLOW	EE200421
	20.06	■	YELLOW	EE200621
	20.08	■	YELLOW	EE200821
	20.10	■	YELLOW	EE201021
	25.04	■	RED	EE250421
	25.06	■	RED	EE250621
	25.08	■	RED	EE250821
	25.10	■	RED	EE251021
	30.04	■	BLUE	EE300421
	30.06	■	BLUE	EE300621
	30.08	■	BLUE	EE300821
	30.10	■	BLUE	EE301021
	35.04	■	GREEN	EE350421
	35.06	■	GREEN	EE350621
	35.08	■	GREEN	EE350821
	35.10	■	GREEN	EE351021
	40.04	■	BLACK	EE400421
	40.06	■	BLACK	EE400621
	40.08	■	BLACK	EE400821
	40.10	■	BLACK	EE401021

**QTY 25 mm \$35.95 /6 pk**

	17.04	☐	WHITE	EE170425
	17.06	☐	WHITE	EE170625
	17.08	☐	WHITE	EE170825
	17.10	☐	WHITE	EE171025
	20.04	■	YELLOW	EE200425
	20.06	■	YELLOW	EE200625
	20.08	■	YELLOW	EE200825
	20.10	■	YELLOW	EE201025
	25.04	■	RED	EE250425
	25.06	■	RED	EE250625
	25.08	■	RED	EE250825
	25.10	■	RED	EE251025
	30.04	■	BLUE	EE300425
	30.06	■	BLUE	EE300625
	30.08	■	BLUE	EE300825
	30.10	■	BLUE	EE301025
	35.04	■	GREEN	EE350425
	35.06	■	GREEN	EE350625
	35.08	■	GREEN	EE350825
	35.10	■	GREEN	EE351025
	40.04	■	BLACK	EE400425
	40.06	■	BLACK	EE400625
	40.08	■	BLACK	EE400825
	40.10	■	BLACK	EE401025

[Click Here to Order](#)

GT® and GTX® are registered trademark of Dentsply Tulsa Dental Specialties.



**Scan or Click For  
More Information**




# EDGEFILE® X1

## Heat-Treated FireWire™ NiTi




(WaveOne® Users)






QTY **21 mm \$35.95 /3 pk**

	20.06		YELLOW	X1200621RF
	25.06		RED	X1250621RF
	40.06		BLACK	X1400621RF

QTY **25 mm \$35.95 /3 pk**

	20.06		YELLOW	X1200625RF
	25.06		RED	X1250625RF
	40.06		BLACK	X1400625RF

QTY **29 mm \$35.95 /3 pk**

	20.06		YELLOW	X1200629RF
	25.06		RED	X1250629RF
	40.06		BLACK	X1400629RF

### Assorted Pack .06 \$35.95 /3 pk

QTY Sizes 20, 25, 40

	21MM	ASSTX121
	25MM	ASSTX125
	29MM	ASSTX129



Scan or Click For  
More Information

[Click Here to Order](#)

WaveOne® is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGEFILE® X3







## Heat-Treated FireWire™ NiTi

Rotary Files







(ProTaper® Users)



**QTY 21 mm \$35.95 /6 pk**

	17.06 (N1)		PURPLE	X3170621RF
	17.04 (N2)		WHITE	X3170421RF
	20.06 (C1)		YELLOW	X3200621RF
	25.06 (C2)		RED	X3250621RF
	30.06 (C3)		BLUE	X3300621RF
	40.06 (C4)		BLACK	X3400621RF

**25 mm \$35.95 /6 pk**

	17.06 (N1)		PURPLE	X3170625RF
	17.04 (N2)		WHITE	X3170425RF
	20.06 (C1)		YELLOW	X3200625RF
	25.06 (C2)		RED	X3250625RF
	30.06 (C3)		BLUE	X3300625RF
	40.06 (C4)		BLACK	X3400625RF

**QTY 29 mm \$35.95 /6 pk**

	17.06 (N1)		PURPLE	X3170629RF
	17.04 (N2)		WHITE	X3170429RF
	20.06 (C1)		YELLOW	X3200629RF
	25.06 (C2)		RED	X3250629RF
	30.06 (C3)		BLUE	X3300629RF
	40.06 (C4)		BLACK	X3400629RF



**Scan or Click For  
More Information**


[Click Here to Order](#)

ProTaper® is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGEFILE® X7

## Excellent Strength and Flexibility

EdgeEndo®'s Leading Proprietary Global Platform:  
Most recognized and proven solution!



## Heat Treated FireWire™ NiTi

FireWire™ NiTi alloy is an innovation in file metallurgy. It makes our files more flexible and significantly increases resistance to cyclic fatigue. FireWire™ NiTi enables EdgeEndo® files to not “bounce back”, preserving canal anatomy, and carefully follows the canal as they shape.\*

\*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph:<http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit [www.edgeendo.com/matter-of-temperature/](http://www.edgeendo.com/matter-of-temperature/)

# EDGEFILE® X7

Heat-Treated FireWire™ NiTi

Rotary Files  
(Universal)

**QTY 21 mm \$35.95 /6 pk**

	17.04	□	WHITE	X7170421RF
	17.06	□	WHITE	X7170621RF
	20.04	■	YELLOW	X7200421RF
	20.06	■	YELLOW	X7200621RF
	25.04	■	RED	X7250421RF
	25.06	■	RED	X7250621RF
	30.04	■	BLUE	X7300421RF
	30.06	■	BLUE	X7300621RF
	35.04	■	GREEN	X7350421RF
	35.06	■	GREEN	X7350621RF
	40.04	■	BLACK	X7400421RF
	40.06	■	BLACK	X7400621RF
	45.04	□	WHITE	X7450421RF
	45.06	□	WHITE	X7450621RF



**QTY 29 mm \$35.95 /6 pk**

	17.04	□	WHITE	X7170429RF
	17.06	□	WHITE	X7170629RF
	20.04	■	YELLOW	X7200429RF
	20.06	■	YELLOW	X7200629RF
	25.04	■	RED	X7250429RF
	25.06	■	RED	X7250629RF
	30.04	■	BLUE	X7300429RF
	30.06	■	BLUE	X7300629RF
	35.04	■	GREEN	X7350429RF
	35.06	■	GREEN	X7350629RF
	40.04	■	BLACK	X7400429RF
	40.06	■	BLACK	X7400629RF
	45.04	□	WHITE	X7450429RF
	45.06	□	WHITE	X7450629RF

**QTY 25 mm \$35.95 /6 pk**

	17.04	□	WHITE	X7170425RF
	17.06	□	WHITE	X7170625RF
	20.04	■	YELLOW	X7200425RF
	20.06	■	YELLOW	X7200625RF
	25.04	■	RED	X7250425RF
	25.06	■	RED	X7250625RF
	30.04	■	BLUE	X7300425RF
	30.06	■	BLUE	X7300625RF
	35.04	■	GREEN	X7350425RF
	35.06	■	GREEN	X7350625RF
	40.04	■	BLACK	X7400425RF
	40.06	■	BLACK	X7400625RF
	45.04	□	WHITE	X7450425RF
	45.06	□	WHITE	X7450625RF



[Click Here to Order](#)

Scan or Click For  
More Information



### Assorted Pack .04 \$35.95 /6 pk

QTY (17, 20, 25, 30, 35, 40)

	21MM	ASSTX70421
	25MM	ASSTX70425
	29MM	ASSTX70429



### Assorted Pack .06 \$35.95 /6 pk

QTY (17, 20, 25, 30, 35, 40)

	21MM	ASSTX70621
	25MM	ASSTX70625
	29MM	ASSTX70629

[Click Here to Order](#)



Scan or Click For  
More Information

# EDGEFILE® XR

Non Heat-Treated NiTi

Retreatment Files



QTY **15 mm \$35.95 /4 pk**

	25.12 (R1)		XRR1251215
	25.08 (R2)	■ RED	XRR2250815

QTY **23 mm \$35.95 /4 pk**

	25.04 (R4)	■ BLACK	XRR4250423
--	------------	---------	------------

QTY **19 mm \$35.95 /4 pk**

	25.06 (R3)	■ GREY	XRR3250619
--	------------	--------	------------




[Click Here to Order](#)






Scan or Click For  
More Information






**QTY 21 mm \$31.95 /6 pk**

	P1		PURPLE	EFNDP121
	P2		WHITE	EFNDP221
	P3		YELLOW	EFNDP321

**QTY 31 mm \$31.95 /6 pk**

	P1		PURPLE	EFNDP131
	P2		WHITE	EFNDP231
	P3		YELLOW	EFNDP331

**QTY 25 mm \$31.95 /6 pk**

	P1		PURPLE	EFNDP125
	P2		WHITE	EFNDP225
	P3		YELLOW	EFNDP325

**Assorted Packs \$31.95 /3 pk**

QTY (P1, P2, P3)

	21MM	EFNDASST21
	25MM	EFNDASST25
	31MM	EFNDASST31



**Scan or Click For  
More Information**

[Click Here to Order](#)

PathFile® is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGEGLIDEPATH™

## Heat-Treated FireWire™ NiTi

Rotary Files  
(ProGlider® Users)

**21 mm \$25.88 /4 pk**

16.02  WHITE EGP0221

**25 mm \$25.88 /4 pk**

16.02  WHITE EGP0225

**31 mm \$25.88 /4 pk**

16.02  WHITE EGP0231



**Scan or Click For  
More Information**

[Click Here to Order](#)

ProGlider® is a registered trademark of Dentsply Tulsa Dental Specialties.



## Heat-Treated FireWire™ NiTi

(TruNatomy® Users)



QTY **19 mm \$25.88 /3 pk**

<input type="checkbox"/>	ORIFICE	<input type="checkbox"/>	PURPLE	MOROR19HT
--------------------------	---------	--------------------------	--------	-----------

QTY **21 mm \$25.88 /3 pk**

<input type="checkbox"/>	GLIDER	<input type="checkbox"/>	WHITE	MORGL21HT
<input type="checkbox"/>	SMALL	<input type="checkbox"/>	YELLOW	MORSM21HT
<input type="checkbox"/>	PRIMARY	<input type="checkbox"/>	RED	MORPRI21HT
<input type="checkbox"/>	MEDIUM	<input type="checkbox"/>	BLUE	MORMED21HT
<input type="checkbox"/>	LARGE	<input type="checkbox"/>	BLACK	MORLG21HT

**Assorted Packs \$45.45/ 5pk**

(Orifice, Glider, Small, Primary, Medium)

<input type="checkbox"/>	21MM	ASSTMORP21
<input type="checkbox"/>	25MM	ASSTMORP25
<input type="checkbox"/>	31MM	ASSTMORP31

QTY **25 mm \$25.88 /3 pk**

<input type="checkbox"/>	GLIDER	<input type="checkbox"/>	WHITE	MORGL25HT
<input type="checkbox"/>	SMALL	<input type="checkbox"/>	YELLOW	MORSM25HT
<input type="checkbox"/>	PRIMARY	<input type="checkbox"/>	RED	MORPRI25HT
<input type="checkbox"/>	MEDIUM	<input type="checkbox"/>	BLUE	MORMED31HT
<input type="checkbox"/>	LARGE	<input type="checkbox"/>	BLACK	MORLG31HT

QTY **31 mm \$25.88 /3 pk**

<input type="checkbox"/>	GLIDER	<input type="checkbox"/>	WHITE	MORGL31HT
<input type="checkbox"/>	SMALL	<input type="checkbox"/>	YELLOW	MORSM31HT
<input type="checkbox"/>	PRIMARY	<input type="checkbox"/>	RED	MORPRI31HT
<input type="checkbox"/>	MEDIUM	<input type="checkbox"/>	BLUE	MORMED31HT
<input type="checkbox"/>	LARGE	<input type="checkbox"/>	BLACK	MORLG31HT



[Click Here to Order](#)

TruNatomy® is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGEONE FIRE™

## Excellent Strength and Flexibility

The Standard in Reciprocating Technology



Parallelogram shaped cross-section

Available in 3 lengths: 21, 25, 31mm



Highly flexible



No bounce back



High strength

## Heat Treated FireWire™ NiTi

FireWire™ NiTi alloy is an innovation in file metallurgy. It makes our files more flexible and significantly increases resistance to cyclic fatigue. FireWire™ NiTi enables EdgeEndo® files to not “bounce back”, preserving canal anatomy, and carefully follows the canal as they shape.\*

\*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph:<http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit [www.edgeendo.com/matter-of-temperature/](http://www.edgeendo.com/matter-of-temperature/) WaveOne® Gold is a registered trademark of Dentsply Tulsa Dental Specialties.





# EDGEONE FIRE™

## Heat-Treated FireWire™ NiTi





Reciprocating File  
(WaveOne® Gold Users)







QTY **21 mm \$25.88 /3 pk**

	SMALL		YELLOW	EOFSM21
	PRIMARY		RED	EOFPRI21
	MEDIUM		GREEN	EOFMED21
	LARGE		WHITE	EOFLG21

QTY **25 mm \$25.88 /3 pk**

	SMALL		YELLOW	EOFSM25
	PRIMARY		RED	EOFPRI25
	MEDIUM		GREEN	EOFMED25
	LARGE		WHITE	EOFLG25

QTY **31 mm \$25.88 /3 pk**

	SMALL		YELLOW	EOFSM31
	PRIMARY		RED	EOFPRI31
	MEDIUM		GREEN	EOFMED31
	LARGE		WHITE	EOFLG31

**Assorted Packs \$34.16 /4 pk**

QTY (20, 25, 35, 45))

	21MM	ASSTEOF21
	25MM	ASSTEOF25
	31MM	ASSTEOF31



Scan or Click For  
More Information

[Click Here to Order](#)

WaveOne® Gold is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGE GLIDEPATH™

for EdgeOne Fire™

**21 mm \$22.88 /3 pk**

<input type="checkbox"/>	16.02	<input type="checkbox"/>	WHITE	EOFGP21
--------------------------	-------	--------------------------	-------	---------

**25 mm \$22.88 /3 pk**

<input type="checkbox"/>	16.02	<input type="checkbox"/>	WHITE	EOFGP25
--------------------------	-------	--------------------------	-------	---------

**31 mm \$22.88 /3 pk**

<input type="checkbox"/>	16.02	<input type="checkbox"/>	WHITE	EOFGP31
--------------------------	-------	--------------------------	-------	---------

Reciprocating GlidePath



[Click Here to Order](#)



# EDGEONE PLATINUM™

## Heat-Treated FireWire™ NiTi

(WaveOne® Gold Users)



QTY **21 mm \$25.88 /3 pk**

	VERIFILE™	<input type="checkbox"/>	WHITE	EOPVER21
	SMALL	<input type="checkbox"/>	YELLOW	EOPSM21
	PRIMARY	<input type="checkbox"/>	RED	EOPPRI21
	MEDIUM	<input type="checkbox"/>	GREEN	EOPMED21
	LARGE	<input type="checkbox"/>	WHITE	EOPLG21

QTY **25 mm \$25.88 /3 pk**

	VERIFILE™	<input type="checkbox"/>	WHITE	EOPVER25
	SMALL	<input type="checkbox"/>	YELLOW	EOPSM25
	PRIMARY	<input type="checkbox"/>	RED	EOPPRI25
	MEDIUM	<input type="checkbox"/>	GREEN	EOPMED25
	LARGE	<input type="checkbox"/>	WHITE	EOPLG25

QTY **31 mm \$25.88 /3 pk**

	VERIFILE™	<input type="checkbox"/>	WHITE	EOPVER31
	SMALL	<input type="checkbox"/>	YELLOW	EOPSM31
	PRIMARY	<input type="checkbox"/>	RED	EOPPRI31
	MEDIUM	<input type="checkbox"/>	GREEN	EOPMED31
	LARGE	<input type="checkbox"/>	WHITE	EOPLG31

### Assorted Packs \$42.45 /5 pk

QTY (15, 20, 25, 35, 45)

	21MM	ASSTEOP21
	25MM	ASSTEOP25
	31MM	ASSTEOP31



Scan or Click For  
More Information

[Click Here to Order](#)

WaveOne® Gold is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGE GLIDEPATH™

for EdgeOne Platinum™

Reciprocating GlidePath



QTY **21 mm \$22.88 /3 pk**

<input type="text"/>	16.02	<input type="checkbox"/>	WHITE	EOGP21
----------------------	-------	--------------------------	-------	--------

QTY **25 mm \$22.88 /3 pk**

<input type="text"/>	16.02	<input type="checkbox"/>	WHITE	EOGP25
----------------------	-------	--------------------------	-------	--------


QTY **31 mm \$22.88 /3 pk**

<input type="text"/>	16.02	<input type="checkbox"/>	WHITE	EOGP31
----------------------	-------	--------------------------	-------	--------

[Click Here to Order](#)

# EDGESEQUEL SAPPHIRE™

**Bodybuilder Strength. Gymnast Flexibility.**



Triangular Cross Section  
Available in 3 lengths: 21, 25 & 31 mm



Highly Flexible



No Bounce Back



High Strength

## Heat Treated FireWire™ NiTi

FireWire™ NiTi alloy is an innovation in file metallurgy. It makes our files more flexible and significantly increases resistance to cyclic fatigue. FireWire™ NiTi enables EdgeEndo® files to not “bounce back”, preserving canal anatomy, and carefully follows the canal as they shape.\*

## The Strength of EdgeSequel Sapphire™

Cyclic Fatigue Testing - Average Time to Failure in Seconds\*

EDGESEQUEL  
SAPPHIRE™

633 seconds

Vortex Blue®

299.9 seconds

EndoSequence®

80.9 seconds

\*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph: <http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance>  
Vortex® and Vortex® Blue are registered trademarks of Dentsply Tulsa Dental Specialties. EndoSequence® is a registered trademark of Brasseler USA

# EDGESEQUEL

## SAPPHIRE™

Heat-Treated FireWire™ NiTi

Rotary Files

(Vortex®, Vortex® Blue & EndoSequence® Users)



QTY **21 mm \$25.88 /4 pk**

	15.04	□	WHITE	ES150421
	15.06	□	WHITE	ES150621
	20.04	■	YELLOW	ES200421
	20.06	■	YELLOW	ES200621
	25.04	■	RED	ES250421
	25.06	■	RED	ES250621
	30.04	■	BLUE	ES300421
	30.06	■	BLUE	ES300621
	35.04	■	GREEN	ES350421
	35.06	■	GREEN	ES350621
	40.04	■	BLACK	ES400421
	40.06	■	BLACK	ES400621
	45.04	□	WHITE	ES450421
	45.06	□	WHITE	ES450621
	50.04	■	YELLOW	ES500421
	50.06	■	YELLOW	ES500621
	55.04	■	RED	ES550421
	60.04	■	BLUE	ES600421
	70.04	■	GREEN	ES700421
	80.04	■	BLACK	ES800421
	90.04	□	WHITE	ES900421

QTY **25 mm \$25.88 /4 pk**

	15.04	□	WHITE	ES150425
	15.06	□	WHITE	ES150625
	20.04	■	YELLOW	ES200425
	20.06	■	YELLOW	ES200625
	25.04	■	RED	ES250425
	25.06	■	RED	ES250625
	30.04	■	BLUE	ES300425
	30.06	■	BLUE	ES300625
	35.04	■	GREEN	ES350425
	35.06	■	GREEN	ES350625
	40.04	■	BLACK	ES400425
	40.06	■	BLACK	ES400625
	45.04	□	WHITE	ES450425
	45.06	□	WHITE	ES450625
	50.04	■	YELLOW	ES500425
	50.06	■	YELLOW	ES500625
	55.04	■	RED	ES550425
	60.04	■	BLUE	ES600425
	70.04	■	GREEN	ES700425
	80.04	■	BLACK	ES800425
	90.04	□	WHITE	ES900425



[Click Here to Order](#)

Scan or Click For  
More Information

Vortex® is a registered trademark of Dentsply Tulsa Dental Specialties. EndoSequence® is a registered trademark of Brassler USA

# EDGESEQUEL

## SAPPHIRE™

### Heat-Treated FireWire™ NiTi

Rotary Files

(Vortex®, Vortex® Blue & EndoSequence® Users)



QTY **31 mm \$25.88 /4 pk**

	15.04	□	WHITE	ES150431
	15.06	□	WHITE	ES150631
	20.04	■	YELLOW	ES200431
	20.06	■	YELLOW	ES200631
	25.04	■	RED	ES250431
	25.06	■	RED	ES250631
	30.04	■	BLUE	ES300431
	30.06	■	BLUE	ES300631
	35.04	■	GREEN	ES350431
	35.06	■	GREEN	ES350631
	40.04	■	BLACK	ES400431
	40.06	■	BLACK	ES400631
	45.04	□	WHITE	ES450431
	45.06	□	WHITE	ES450631
	50.04	■	YELLOW	ES500431
	50.06	■	YELLOW	ES500631
	55.04	■	RED	ES550431
	60.04	■	BLUE	ES600431
	70.04	■	GREEN	ES700431
	80.04	■	BLACK	ES800431
	90.04	□	WHITE	ES900431

### Assorted Packs .04 \$25.88 /4 pk

QTY (15, 20, 25, 30)

	21MM	ASST15300421
	25MM	ASST15300425
	31MM	ASST15300431

### Assorted Packs .04 \$25.88 /4 pk

QTY (35, 40, 45, 50)

	21MM	ASST35500421
	25MM	ASST35500425
	31MM	ASST35500431

### Assorted Packs .06 \$25.88 /4 pk

QTY (15, 20, 25, 30)

	21MM	ASST15300621
	25MM	ASST15300625
	31MM	ASST15300631

### Assorted Packs .06 \$25.88 /4 pk

QTY (35, 40, 45, 50)

	21MM	ASST35500621
	25MM	ASST35500625
	31MM	ASST35500631




[Click Here to Order](#)

Scan or Click For  
More Information

Vortex® is a registered trademark of Dentsply Tulsa Dental Specialties. EndoSequence® is a registered trademark of Brassler USA

# EDGETAPER™

## The Taper Has Been Edge-ified



Highly flexible



No bounce back



High strength

## The Amazing Strength of EdgeTaper™

Cyclic Fatigue Testing - Average Time to Failure in Seconds\*

EDGETAPER  
PLATINUM™

677 seconds

ProTaper Gold®

312 seconds

EDGETAPER™

127 seconds

ProTaper®

119 seconds








\*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph: <http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance>  
ProTaper® Gold and ProTaper are registered trademarks of Dentsply Tulsa Dental Specialties










QTY **19 mm \$35.95 /6 pk**

	SX	ETSX19
--	----	--------








QTY **21 mm \$35.95 /6 pk**

	S1	 PURPLE	ETS121
	S2	 WHITE	ETS221
	F1	 YELLOW	ETF121
	F2	 RED	ETF221
	F3	 BLUE	ETF321
	F4	 BLACK	ETF421
	F5	 YELLOW	ETF521

QTY **31 mm \$35.95 /6 pk**

	S1	 PURPLE	ETS131
	S2	 WHITE	ETS231
	F1	 YELLOW	ETF131
	F2	 RED	ETF231
	F3	 BLUE	ETF331
	F4	 BLACK	ETF431
	F5	 YELLOW	ETF531

QTY **25 mm \$35.95 /6 pk**

	S1	 PURPLE	ETS125
	S2	 WHITE	ETS225
	F1	 YELLOW	ETF125
	F2	 RED	ETF225
	F3	 BLUE	ETF325
	F4	 BLACK	ETF425
	F5	 YELLOW	ETF525

**Assorted Packs \$35.95 /6 pk**

QTY (S1, S2, F1, F2, F3, SX)

	21MM	ETASST21
	25MM	ETASST25
	31MM	ETASST31

[Click Here to Order](#)




**Scan or Click For  
More Information**

ProTaper® is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGE TAPER

## PLATINUM™

### Edge-ified Sequence with Enhanced Strength and Flexibility



Convex Triangular cross-section

Available in 3 lengths: 21, 25, 31 mm



Highly  
flexible



No bounce  
back



High  
strength

### Heat Treated FireWire™ NiTi

FireWire™ NiTi alloy is an innovation in file metallurgy. It makes our files more flexible and significantly increases resistance to cyclic fatigue. FireWire™ NiTi enables EdgeEndo® files to not “bounce back”, preserving canal anatomy, and carefully follows the canal as they shape.\*

\*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph: <http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit [www.edgeendo.com/matter-of-temperature/](http://www.edgeendo.com/matter-of-temperature/) ProTaper Gold® is a registered trademark of Dentsply Tulsa Dental Specialties.



# EDGE TAPER

## PLATINUM™

Heat-Treated FireWire™ NiTi

(ProTaper® Gold Users)



QTY **19 mm \$35.95 /6 pk**

	SX	ETSX19HT
--	----	----------

QTY **21 mm \$35.95 /6 pk**

	S1	PURPLE	ETS121HT
	S2	WHITE	ETS221HT
	F1	YELLOW	ETF121HT
	F2	RED	ETF221HT
	F3	BLUE	ETF321HT
	F4	BLACK	ETF421HT
	F5	YELLOW	ETF521HT

QTY **31 mm \$35.95 /6 pk**

	S1	PURPLE	ETS131HT
	S2	WHITE	ETS231HT
	F1	YELLOW	ETF131HT
	F2	RED	ETF231HT
	F3	BLUE	ETF331HT
	F4	BLACK	ETF431HT
	F5	YELLOW	ETF531HT

QTY **25 mm \$35.95 /6 pk**

	S1	PURPLE	ETS125HT
	S2	WHITE	ETS225HT
	F1	YELLOW	ETF125HT
	F2	RED	ETF225HT
	F3	BLUE	ETF325HT
	F4	BLACK	ETF425HT
	F5	YELLOW	ETF525HT

**Assorted Packs \$35.95 /6 pk**

	QTY	(S1, S2, F1, F2, F3, SX)
	21MM	ETP21MM
	25MM	ETP25MM
	31MM	ETP31MM

[Click Here to Order](#)





Scan or Click For  
More Information


ProTaper® Gold is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGE V-TAPER HT™

**Strong, Advanced, Flexible, Efficient**



Parabolic Cross Section  
Available in 3 lengths: 21, 25 & 3 mm













EdgeV-Taper HT™ features our Heat-Treated™ Fire Wire NiTi making them super strong and flexible. They are designed with a Safe-Core™ Parabolic Cross-Section. The Safe-Core Parabolic Cross-Section is excellent for deep-hole drilling and reaming as done in root canal procedures. The Safe-Core™ Parabolic Cross-Section combines the attributes of being a highly efficient and flexible instrument while being extremely safe and resistant to fracture, alleviating your stress and resulting in a more enjoyable procedure.

EdgeV-Taper HT™ system includes tip sizes 20, 22, 25, 30, 35, 40, 45 and 50. They have a changing variable taper along the length of the file. Lengths available are 21mm, 25mm and 31mm.











\*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph:<http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit [www.edgeendo.com/matter-of-temperature/](http://www.edgeendo.com/matter-of-temperature/) ProTaper Gold® is a registered trademark of Dentsply Tulsa Dental Specialties.













**QTY 21 mm \$35.95 /6 pk**

	20.06		YELLOW	VF200621
	25.05		RED	VF250521
	22.07		YELLOW	VF220721
	25.08		RED	VF250821
	30.06		BLUE	VF300621
	30.08		BLUE	VF300821
	35.06		GREEN	VF350621
	40.06		BLACK	VF400621
	45.06		WHITE	VF450621
	50.06		YELLOW	VF500621

**QTY 31 mm \$35.95 /6 pk**

	20.06		YELLOW	VF200631
	25.05		RED	VF250531
	22.07		YELLOW	VF220731
	25.08		RED	VF250831
	30.06		BLUE	VF300631
	30.08		BLUE	VF300831
	35.06		GREEN	VF350631
	40.06		BLACK	VF400631
	45.06		WHITE	VF450631
	50.06		YELLOW	VF500631

**QTY 25 mm \$35.95 /6 pk**

	20.06		YELLOW	VF200625
	25.05		RED	VF250525
	22.07		YELLOW	VF220725
	25.08		RED	VF250825
	30.06		BLUE	VF300625
	30.08		BLUE	VF300825
	35.06		GREEN	VF350625
	40.06		BLACK	VF400625
	45.06		WHITE	VF450625
	50.06		YELLOW	VF500625

[Click Here to Order](#)



**Scan or Click For  
More Information**

# EDGE V-TAPER<sup>HT</sup> GLIDEPATH



QTY **17 mm \$35.95 /6 pk**

	25.09		VFGP250917
--	-------	--	------------

QTY **21 mm \$35.95 /6 pk**

	13.03	<input checked="" type="checkbox"/>	PINK	VFGP130321
	14.03	<input type="checkbox"/>	WHITE	VFGP140321
	17.04	<input type="checkbox"/>	WHITE	VFGP170421

QTY **25 mm \$35.95 /6 pk**

	13.03	<input checked="" type="checkbox"/>	PINK	VFGP130325
	14.03	<input type="checkbox"/>	WHITE	VFGP140325
	17.04	<input type="checkbox"/>	WHITE	VFGP170425

QTY **29 mm \$35.95 /6 pk**

	14.03	<input type="checkbox"/>	WHITE	VFGP140329
--	-------	--------------------------	-------	------------

QTY **31 mm \$35.95 /6 pk**

	13.03	<input checked="" type="checkbox"/>	PINK	VFGP130331
	17.04	<input type="checkbox"/>	WHITE	VFGP170431

[Click Here to Order](#)



Scan or Click For  
More Information

# EDGE TAPER RETREAT™



QTY **16 mm \$35.95 /6 pk**

<input type="text"/>	D1	<input type="checkbox"/>	WHITE	ETRD116
----------------------	----	--------------------------	-------	---------

QTY **18 mm \$35.95 /6 pk**

<input type="text"/>	D2	<input type="checkbox"/>	WHITE	ETRD218
----------------------	----	--------------------------	-------	---------

QTY **22 mm \$35.95 /6 pk**

<input type="text"/>	D6	<input type="checkbox"/>	WHITE	ETRD322
----------------------	----	--------------------------	-------	---------

QTY **Assorted \$35.95 /6 pk**

<input type="text"/>	ASST D1-D3	<input type="checkbox"/>	WHITE	ETRASSTD123
----------------------	------------	--------------------------	-------	-------------

[Click Here to Order](#)



**QTY 28 mm \$11.88 /6 pk**

	#1	EG0128
	#2	EG0228
	#3	EG0328
	#4	EG0428
	#5	EG0528
	#6	EG0628

**QTY 32 mm \$11.88 /6 pk**

	#1	EG0132
	#2	EG0232
	#3	EG0332
	#4	EG0432
	#5	EG0532
	#6	EG0632
















**Assorted Packs \$12.88 /6 pk**

QTY (#1, #2, #3, #4, #5, #6)

	28MM	EGASST28
	32MM	EGASST32

[Click Here to Order](#)






**QTY 21 mm \$3.88 /6 pk**

6		PINK	KF621
8		GREY	KF821
10		PURPLE	KF1021
15		WHITE	KF1521
20		YELLOW	KF2021
25		RED	KF2521
30		BLUE	KF3021
35		GREEN	KF3521
40		BLACK	KF4021
45		WHITE	KF4521
50		YELLOW	KF5021
55		RED	KF5521
60		BLUE	KF6021
70		GREEN	KF7021
80		BLACK	KF8021
















**QTY 25 mm \$3.88 /6 pk**

6		PINK	KF625
8		GREY	KF825
10		PURPLE	KF1025
15		WHITE	KF1525
20		YELLOW	KF2025
25		RED	KF2525
30		BLUE	KF3025
35		GREEN	KF3525
40		BLACK	KF4025
45		WHITE	KF4525



50		YELLOW	KF5025
55		RED	KF5525
60		BLUE	KF6025
70		GREEN	KF7025
80		BLACK	KF8025

**QTY 31 mm \$3.88 /6 pk**

6		PINK	KF631
8		GREY	KF831
10		PURPLE	KF1031
15		WHITE	KF1531
20		YELLOW	KF2031
25		RED	KF2531
30		BLUE	KF3031
35		GREEN	KF3531
40		BLACK	KF4031
45		WHITE	KF4531
50		YELLOW	KF5031
55		RED	KF5531
60		BLUE	KF6031
70		GREEN	KF7031
80		BLACK	KF8031

[Click Here to Order](#)

# EDGEK-FILE™ NITi

## Non Heat-Treated NiTi

(Tulsa Dental K-File® Users)



### QTY 21 mm \$5.88 /6 pk

15	□	WHITE	NITIKF1521
20	■	YELLOW	NITIKF2021
25	■	RED	NITIKF2521
30	■	BLUE	NITIKF3021
35	■	GREEN	NITIKF3521
40	■	BLACK	NITIKF4021
45	□	WHITE	NITIKF4521
50	■	YELLOW	NITIKF5021
55	■	RED	NITIKF5521
60	■	BLUE	NITIKF6021

### QTY 31 mm \$5.88 /6 pk

15	□	WHITE	NITIKF1531
20	■	YELLOW	NITIKF2031
25	■	RED	NITIKF2531
30	■	BLUE	NITIKF3031
35	■	GREEN	NITIKF3531
40	■	BLACK	NITIKF4031
45	□	WHITE	NITIKF4531
50	■	YELLOW	NITIKF5031
55	■	RED	NITIKF5531
60	■	BLUE	NITIKF6031


### QTY 25 mm \$5.88 /6 pk

15	□	WHITE	NITIKF1525
20	■	YELLOW	NITIKF2025
25	■	RED	NITIKF2525
30	■	BLUE	NITIKF3025
35	■	GREEN	NITIKF3525
40	■	BLACK	NITIKF4025
45	□	WHITE	NITIKF4525
50	■	YELLOW	NITIKF5025
55	■	RED	NITIKF5525
60	■	BLUE	NITIKF6025








[Click Here to Order](#)










# EDGE BROACHES™



QTY **21 mm \$6.88 /10 pk**

#0		PURPLE	EB4XF021
#1		WHITE	EB3XF121
#2		YELLOW	EB2XF221
#3		RED	EB1XF321
#4		BLUE	EB0XF421
#5		GREEN	EBM0521
#6		BLACK	EBC0621
ASST #0-6		EBASST0621	
ASST #2-4		EBASST2421	

QTY **25 mm \$6.88 /10 pk**

#0		PURPLE	EB4XF025
#1		WHITE	EB3XF125
#2		YELLOW	EB2XF225
#3		RED	EB1XF325
#4		BLUE	EB0XF425
#5		GREEN	EBM0525
#6		BLACK	EBC0625
ASST #0-6		EBASST0625	
ASST #2-4		EBASST2425	

[Click Here to Order](#)

# EDGE C-FILE™

## Stainless Steel



QTY **21 mm \$5.88 /6 pk**

	6		PINK	CF621
	8		GREY	CF821
	10		PURPLE	CF1021
	12		BROWN	CF1221
	15		WHITE	CF1521

QTY **25 mm \$5.88 /6 pk**

	6		PINK	CF625
	8		GREY	CF825
	10		PURPLE	CF1025
	12		BROWN	CF1225
	15		WHITE	CF1525

[Click Here to Order](#)



**QTY 21 mm \$3.88 /6 pk**

	15	<input type="checkbox"/>	WHITE	HED1521
	20	<input type="checkbox"/>	YELLOW	HED2021
	25	<input type="checkbox"/>	RED	HED2521
	30	<input type="checkbox"/>	BLUE	HED3021
	35	<input type="checkbox"/>	GREEN	HED3521
	40	<input type="checkbox"/>	BLACK	HED4021
	45	<input type="checkbox"/>	WHITE	HED4521
	50	<input type="checkbox"/>	YELLOW	HED5021

**QTY 25 mm \$3.88 /6 pk**

	15	<input type="checkbox"/>	WHITE	HED1525
	20	<input type="checkbox"/>	YELLOW	HED2025
	25	<input type="checkbox"/>	RED	HED2525
	30	<input type="checkbox"/>	BLUE	HED3025
	35	<input type="checkbox"/>	GREEN	HED3525
	40	<input type="checkbox"/>	BLACK	HED4025
	45	<input type="checkbox"/>	WHITE	HED4525
	50	<input type="checkbox"/>	YELLOW	HED5025

**QTY 31 mm \$3.88 /6 pk**

	15	<input type="checkbox"/>	WHITE	HED1531
	20	<input type="checkbox"/>	YELLOW	HED2031
	25	<input type="checkbox"/>	RED	HED2531
	30	<input type="checkbox"/>	BLUE	HED3031
	35	<input type="checkbox"/>	GREEN	HED3531
	40	<input type="checkbox"/>	BLACK	HED4031
	45	<input type="checkbox"/>	WHITE	HED4531
	50	<input type="checkbox"/>	YELLOW	HED5031

[Click Here to Order](#)

# EDGEAPEX HP™



- Pre-programmed with top selling EdgeEndo files PLUS 10 fully customizable Memory Program slots
- Integrated Apex Locator option to use as a standalone Apex locator with handfile
- Efficient brushless motor, low noise, long service life
- Cordless portable endo motor with combined length determination
- 360° rotation of contra angle
- Adopt real-time feedback technology and dynamic torque control, effectively preventing file separation.
- Wireless Charging
- Auto On/Off, Torque Reverse, Apical Reverse, Apical Slow Down
- Max Torque 5.0 Ncm
- Visual torque on file in use
- Power 5.0 Ncm
- Change Dominant Hand
- Range: 100 - 2500

# EDGEAPEX HP™

Qty **\$995.00/ UNIT**

	1 UNIT	EAPEXHP
--	--------	---------

## Additional Accessories included in package:

- Measuring Wire
- Lip Hook (x2)
- File Clip (x4)
- Touch Probe (x2)
- Power Adapter
- Protective Silicon Cover (x2)
- Spray Nozzle
- O-Rings (x2)
- Disposable Sleeves (1 pack)



## Accessories

### MEASURING WIRE

QTY **\$17.88/ 1 PER PACK**

	EAPEXMW
--	---------



### CONTRA ANGLE

QTY **\$199.88/ 1 PER PACK**

	EAPEXCA
--	---------



### LIP HOOK

QTY **\$14.88/ 5 PER PACK**

	EAPEXLH
--	---------



### PROTECTIVE SILICONE COVER

QTY **\$10.88/ 3 PER PACK**

	EAPEXSC
--	---------



### FILE CLIP

QTY **\$59.88/ 5 PER PACK**

	EAPEXFC
--	---------



### CHARGING BASE

QTY **\$74.88/ 1 PER PACK**

	EAPEXCB
--	---------



### TOUCH PROBE

QTY **\$8.88/ 1 PER PACK**

	EAPEXTP
--	---------



[Click Here to Order](#)

# EDGEAPEX™

QTY **\$495.00/ UNIT**

	1 UNIT	EEALU
--	--------	-------

### Additional Accessories included in package:

- Base
- Measuring Wire
- File Clip (4)
- Touch Probe (2)
- Pulp Tester Probe (2)
- Adapter
- Tester
- Enlarged display of apical foramen
- Instruction card for apical stop
- Apex Locator troubleshooting card



## Apex Locator Accessories

### MEASURING WIRE

QTY **\$17.88/ 1 PER PACK**

	EEALUMW
--	---------



### LITHIUM BATTERY

QTY **\$17.88/ 1 PER PACK**

	EEALULB
--	---------



### LIP HOOK

QTY **\$14.88/ 5 PER PACK**

	EEALULH5
--	----------



### TESTER

QTY **\$14.88/ 1 PER PACK**

	EEALUT
--	--------



### FILE CLIP

QTY **\$59.88/ 5 PER PACK**

	EEALUFC5
--	----------



### ENDO RULER

QTY **\$14.88/ 1 PER PACK**

	EEALUER
--	---------



### TOUCH PROBE

QTY **\$8.88/ 1 PER PACK**

	EEALUMF
--	---------



[Click Here to Order](#)

# EDGEPACK™



- Wireless Gutta Percha Obturation Device
- Ergonomic design and ambidextrous accessible
- Easy to use
- 0.35mm pre-bend tip
- Improved penetrability into narrow root canals
- High precision temperature control
- Smooth cutting of GP from high-efficiency conversion system
- Compact size for easy cleaning and storage
- Visual and auditory prompts for safe operation
- Preset temperatures: 150°C, 180°C, 200°C, 230°C
- Time related features: Warm up time within 0.2 seconds, Auto shut off after 10 seconds,
- 4 hour battery life with each charge

# EDGEPACK™

QTY **445.00/ UNIT**

<input type="text"/>	1 UNIT	EEDPU
----------------------	--------	-------

## Additional Accessories included in package:

- Main unit
- Charging Stand
- Power Adapter
- Heat Insulation Sleeve of Tips (4)
- Instruction Manual
- Qualified Certification
- Packing List
- EdgePack Tips
- 4004 (1)
- 4504 (1)
- 5506 (1)
- 5508 (1)



## Accessories

### EDGEPACK CHARGING STAND

QTY **\$59.88/ UNIT**

<input type="text"/>	1 UNIT	EEDPUCS
----------------------	--------	---------



[Click Here to Order](#)



# EdgePack Tips

## EDGEPACK 3504 TIP

**\$49.88/ 1 PER PACK**

QTY

EEDPUT3504



## EDGEPACK 5508 TIP

**\$49.88/ 1 PER PACK**

QTY

EEDPUT5508



## EDGEPACK 4004 TIP

**\$49.88/ 1 PER PACK**

QTY

EEDPUT4004



## EDGEPACK 5508L TIP

**\$49.88/ 1 PER PACK**

QTY

EEDPUT5508L



## EDGEPACK 4504 TIP

**\$49.88/ 1 PER PACK**

QTY

EEDPUT4504



## EDGEPACK 5510 TIP

**\$49.88/ 1 PER PACK**

QTY

EEDPUT5510



## EDGEPACK 5004 TIP

**\$49.88/ 1 PER PACK**

QTY

EEDPUT5004



## EDGEPACK 6012 TIP

**\$49.88/ 1 PER PACK**

QTY

EEDPUT6012



## EDGEPACK 5506 TIP

**\$49.88/ 1 PER PACK**

QTY

EEDPUT5506



[Click Here to Order](#)

# EDGEFLOW™



- Wireless and electric
- Excellent gutta-percha fluidity, uniform and smooth gutta-percha output without bubbles.
- Precise and intuitive display of remaining gutta-percha
- Auto heating stop after timeout; Auto shutdown after timeout
- Charging base with battery
- Temperature monitoring design
- Dual batteries deliver longer endurance, enabling filling about 100 root canals.
- Re-usable Injection needles available in various sizes and lengths
- Preset temperatures: 100°C, 120°C, 150°C, 180°C, 200°C
- Injection Speeds: Slow, Medium, Fast

# EDGEFLOW™

QTY **\$595.00/ UNIT**

<input type="text"/>	1 UNIT	EEPTO
----------------------	--------	-------

## Additional Accessories included in package:

- Main unit
- Charging Stand
- Power Adapter
- Battery (2)
- Heat Insulation Hood
- Wrench of Needle
- Gutta-Percha Ring
- Instruction Manual
- Qualified Certification
- Packing List
- Injection Needle (4)
  - E23G-M (1)
  - E23G-L (1)
  - E25G-L (2)



## Accessories

### OBTURATOR WRENCH

QTY **\$429.88/ 1 PER PACK**

<input type="text"/>	EEPTOW
----------------------	--------



### EDGEFLOW™ BATTERY

QTY **\$29.88/ 1 PER PACK**

<input type="text"/>	EEPTOB
----------------------	--------



### EDGEFLOW™ HEAT INSULATION HOOD

QTY **\$17.88/ 3 PER PACK**

<input type="text"/>	EEPTOIH3
----------------------	----------



[Click Here to Order](#)

# EdgeFlow Re-Useable Needle Tips

**NEEDLE TIP E20G-S**  
**\$33.88/ 1 PER PACK**  
QTY   
EE20GSIN



**NEEDLE TIP E23G-L**  
**\$33.88/ 1 PER PACK**  
QTY   
EE23GLIN



**NEEDLE TIP E20G-M**  
**\$33.88/ 1 PER PACK**  
QTY   
EE20GMIN



**NEEDLE TIP E25G-M**  
**\$33.88/ 1 PER PACK**  
QTY   
EE25GMIN



**NEEDLE TIP E20G-L**  
**\$33.88/ 1 PER PACK**  
QTY   
EE20GLIN



**NEEDLE TIP E25G-L**  
**\$33.88/ 1 PER PACK**  
QTY   
EE25GLIN



**NEEDLE TIP E23G-M**  
**\$33.88/ 1 PER PACK**  
QTY   
EE23GMIN



**NOTE: Gutta Percha Bars/Pellets not currently available for purchase through EdgeEndo®.**

[Click Here to Order](#)

**EDGEPRO**<sup>®</sup>  
Supreme Clean. Light Speed.

LAUNCHING  
**ENDO**  
FORWARD



### **Cleaning. Enhanced.**

EdgePRO<sup>®</sup> provides outstanding Cleaning, Debridement and Disinfection by removing infected tissue, biofilm and smear layer killing up to 99% of bacteria commonly found in the root canal.(1,2,3,4) This should give you the confidence you need to treat even the most complex anatomies.

### **Endodontics. Streamlined.**

EdgePRO<sup>®</sup> lets "You be You" – shape and obturate the canals the way you want without the need to change your technique. EdgePRO<sup>™</sup> offers the potential for streamlined workflow and one-visit treatment resulting in reduced treatment time and patient visits.

### **Technology. Affordable.**

EdgePRO<sup>®</sup> has state-of-the-art technology at more affordable pricing that is as low as 1/3 the cost of other irrigation devices. EdgePRO<sup>™</sup> gives your practice access to the future of Endodontics and an edge in winning more referral business.

**Call 1-855-TO-CLEAN to order**  
**Scan or visit [EdgeProEndo.com](http://EdgeProEndo.com)**  
**for more information**






1) Christo JE, Zilm PS, Sullivan T, Cathro PR. Efficacy of low concentrations of sodium hypochlorite and low-powered Er,Cr:YSGG laser activated irrigation against an Enterococcus faecalis biofilm. International Endodontic Journal. 2016;49(3):279-286. doi:10.1111/iej.12447 2) Gordon W, Atabakhsh VA, Meza F, et al. The antimicrobial efficacy of the ErCrYSGG laser with radial emitting tips on root canal dentin walls infected with E faecalis. JADA. 2007;138(July):992-1002. doi:10.14219/jada.archive.2007.0297 3) Internal data, FDA 510(K) 4a) 1. Schoop, U et al. "The Impact of an ErCrYSGG Laser with Radial-Firing Tips on Endodontic Treatment." Lasers in medical science 24.1 (2009): 59-65. Web. 5 Dec. 2011. 4b) Arnabat, Josep et al. "Bactericidal Activity of Erbium, Chromium: Yttrium-scandium-gallium-garnet Laser in Root Canals." Lasers Medical Science 25 (2010): 805-810. Web

# Irrigation Tips & Syringes

## Edge® Irrigation Needle Tips




QTY **\$9.88/ 100PK**

	23ga. Blue		10818365020260
	27ga. Yellow		10818365020277
	30ga. Purple		10818365020284



## Edge® Irrigation Probe Tips

QTY **\$15.88/ 100PK**

	23ga. Blue		10818365020291
	27ga. Yellow		10818365020307
	30ga. Purple		10818365020314



## Edge® Pre-Assembled Irrigation Syringe with 1" Notched Tip

QTY **\$20.88/ 100PK**

	23ga. Blue	00818365021062
	27ga. Yellow	00818365021055
	30ga. Purple	00818365021079



## Edge® Luer-Lock Syringe

QTY **\$11.88/ 100PK**

	3CC	10818365020383
--	-----	----------------



[Click Here to Order](#)

# EDGE ICE™

QTY **\$15.88**

6 OZ. SPRAY BOTTLE 818365020065

[Click Here to Order](#)



# EDGEENDO®

## 2% CHLORHEXIDINE GLUCONATE

QTY **\$29.88**

480ML SOLUTION BOTTLE 818365020041

[Click Here to Order](#)



# EDGEENDO®

## 17% EDTA

QTY **\$29.88**

480ML SOLUTION BOTTLE 818365020058

[Click Here to Order](#)



# EDGEENDO®

## 5% SODIUM HYPOCHLORITE

QTY **\$32.88**

480ML SOLUTION BOTTLE 818365020034

[Click Here to Order](#)





# EDGEMIX™

QTY **\$59.88**

16OZ BOTTLE 18365020416

[Click Here to Order](#)



# EDGELUBE™

QTY **\$10.99**

100ML BOTTLE (17% EDTA SOLUTION) 8806387892083

[Click Here to Order](#)



# EDGE GEL™

QTY **\$12.99**

5.4ML BOTTLE (19% EDTA GEL) 8806387892076

[Click Here to Order](#)





# EDGETIP™

## EdgeTip #1

Diamond-coated ball tip, length 20mm, taper 5%.

QTY **\$49.88**

<input type="text"/>	EDGETIP #1	ETBD
----------------------	------------	------



## EdgeTip #2

Diamond-coated steel tip 76µm, length 18mm, taper 5%

QTY **\$49.88**

<input type="text"/>	EDGETIP #2	ET18D
----------------------	------------	-------



## EdgeTip #3

Diamond-coated retreatment tip, 30µm, length 20mm, taper 5%

QTY **\$49.88**

<input type="text"/>	EDGETIP #3	ET20D
----------------------	------------	-------



## EdgeTip #4 & #5

Long Diamond-coated retreatment tip, 30µm, 40mm, taper 4%

QTY **\$49.88**

<input type="text"/>	EDGETIP #4	ET40D4
<input type="text"/>	EDGETIP #5	ET40D5



## EdgeTip #6

Short Titanium- Niobium tip, length 15mm taper 4%

QTY **\$49.88**

<input type="text"/>	EDGETIP #6	ED25S
----------------------	------------	-------



[Click Here to Order](#)

# EDGETIP™

## EdgeTip #7

Titanium- Niobium tip, length 20mm, taper 3%

QTY **\$49.88**

<input type="text"/>	EDGETIP #7	ET25
----------------------	------------	------



## EdgeTip #8

Titanium- Niobium tip, 25mm, taper 3%

QTY **\$49.88**

<input type="text"/>	EDGETIP #8	ET25L
----------------------	------------	-------



[Click Here to Order](#)

# EDGE SURGICAL TIP™

## EdgeSurgical Tip #1

Diamond-coated universal tip 30µm, length 3mm, taper 9%

QTY **\$49.88**

<input type="text"/>	EDGESURGICAL TIP #1	AS3D1
----------------------	---------------------	-------



## EdgeSurgical Tip #2

Diamond-coated universal tip 30µm, taper 9%

QTY **\$49.88**

<input type="text"/>	EDGESURGICAL TIP #2	AS3D2
----------------------	---------------------	-------



## EdgeSurgical Tip #3

Left-oriented retro surgery tip, diamond-coated 30µm, length 5mm, taper 7%

QTY **\$49.88**

<input type="text"/>	EDGESURGICAL TIP #3	P15LD
----------------------	---------------------	-------



## EdgeSurgical Tip #4

Left-oriented tip, diamond-coated 30µm, length 3mm, taper 10%

QTY **\$49.88**

<input type="text"/>	EDGESURGICAL TIP #4	ASLD
----------------------	---------------------	------



## EdgeSurgical Tip #5

Right-oriented tip, diamond-coated 30µm, length 6mm, taper 7%

QTY **\$49.88**

<input type="text"/>	EDGESURGICAL TIP #5	P15RD
----------------------	---------------------	-------



## EdgeSurgical Tip #6

Right-oriented tip, diamond-coated 30µm, length 3mm, taper 10%

QTY **\$49.88**

<input type="text"/>	EDGESURGICAL TIP #6	ASRD
----------------------	---------------------	------










[Click Here to Order](#)

# EDGECORE™

20/pk

## EDGECORE 20/pk

QTY **100.88**

20		YELLOW	EE2020C
25		RED	EE2520C
30		BLUE	EE3020C
35		GREEN	EE3520C
40		BLACK	EE4020C
45		WHITE	EE4520C
50		YELLOW	EE5020C










[Click Here to Order](#)

# EDGEFILL™

20/pk

## EDGEFILL 20/pk

QTY **\$100.88**

20		YELLOW	EE2020F
25		RED	EE2520F
30		BLUE	EE3020F
35		GREEN	EE3520F
40		BLACK	EE4020F
45		WHITE	EE4520F
50		YELLOW	EE5020F














[Click Here to Order](#)



## EDGECORE for EdgeCoil™












QTY **\$35.95/5 pk + 1 Verifier**

	20		YELLOW	EC20C
	25		RED	EC25C
	30		BLUE	EC30C
	35		GREEN	EC35C
	40		BLACK	EC40C
	45		WHITE	EC45C
	50		YELLOW	EC50C
	55		RED	EC55C
	60		BLUE	EC60C
	70		GREEN	EC70C
	80		BLACK	EC80C

[Click Here to Order](#)

## EDGEFILL for EdgeCoil™






QTY **\$35.95 /6PK**

	20		YELLOW	EC20C
	25		RED	EC25C
	30		BLUE	EC30C
	35		GREEN	EC35C
	40		BLACK	EC40C
	45		WHITE	EC45C
	50		YELLOW	EC50C
	55		RED	EC55C
	60		BLUE	EC60C
	70		GREEN	EC70C
	80		BLACK	EC80C

[Click Here to Order](#)








**EDGECORE for EdgeEvolve™**  
**\$35.95 /5 pk + 1 Verifier**

	20		YELLOW	EE2006C
	25		RED	EE2506C
	30		BLUE	EE3006C
	35		GREEN	EE3506C
	40		BLACK	EE4006C

[Click Here to Order](#)

**EDGEFILL for EdgeEvolve™**  
**\$35.95 /6 pk**

	20		YELLOW	EE2006O
	25		RED	EE2506O
	30		BLUE	EE3006O
	35		GREEN	EE3506O
	40		BLACK	EE4006O

[Click Here to Order](#)



## EDGECORE for EdgeFile® X1

**\$35.95 /5 pk + 1 Verifier**

QTY (for WaveOne® files only not X1 files)

	20.06		YELLOW	X12006C
	25.08		RED	X12508C
	40.08		BLACK	X14008C

[Click Here to Order](#)

## EDGECORE for EdgeFile® X3

**\$35.95 /5 pk + 1 Verifier**

QTY

	20		YELLOW	X32006C
	25		RED	X32506C
	30		BLUE	X33006C
	40		BLACK	X34006C
	50		YELLOW	X35006C

[Click Here to Order](#)

## EDGECORE for EdgeFile® X5

**\$35.95 /5 pk + 1 Verifier**

QTY

	20.04		YELLOW	X52004C
	20.06		YELLOW	X52006C
	20.08		YELLOW	X52008C
	30.04		BLUE	X53004C
	30.06		BLUE	X53006C
	30.08		BLUE	X53008C
	40.04		BLACK	X54004C
	40.06		BLACK	X54006C
	40.08		BLACK	X54008C

[Click Here to Order](#)

## EDGECORE for EdgeFile® X7

**\$35.95 /5 pk + 1 Verifier**

QTY




	20		YELLOW	X72006C
	25		RED	X72506C
	30		BLUE	X73006C
	35		GREEN	X73506C
	40		BLACK	X74006C
	45		WHITE	X74506C
	50		YELLOW	X75006C

[Click Here to Order](#)



## EDGEFILL for EdgeFile® X1






QTY **\$35.95 /6 pk**

	20		YELLOW	X12006O
	25		RED	X12508O
	40		BLACK	X12008O

[Click Here to Order](#)

## EDGEFILL for EdgeFile® X3










QTY **\$35.95 /6 pk**

	20		YELLOW	X32006O
	25		RED	X32506O
	30		BLUE	X33006O
	40		BLACK	X34006O
	50		YELLOW	X35006O

[Click Here to Order](#)

## EDGEFILL for EdgeFile® X5








QTY **\$35.95 /6 pk**

	20.04		YELLOW	X52004O
	20.06		YELLOW	X52006O
	20.08		YELLOW	X52008O
	30.04		BLUE	X53004O
	30.06		BLUE	X53006O
	30.08		BLUE	X53008O
	40.04		BLACK	X54004O
	40.06		BLACK	X54006O
	40.08		BLACK	X54008O

[Click Here to Order](#)

## EDGEFILL for EdgeFile® X7

QTY **\$35.95 /6 pk**

	20		YELLOW	X72006O
	25		RED	X72506O
	30		BLUE	X73006O
	35		GREEN	X73506O
	40		BLACK	X74006O
	45		WHITE	X74506O
	50		YELLOW	X75006O





[Click Here to Order](#)





**EDGECORE for EdgeMorphology™**





QTY **\$35.95 /5 pk + 1 Verifier**

	SMALL		YELLOW	MORSMC
	PRIMARY		RED	MORPRIC
	MEDIUM		BLUE	MORMEDC
	LARGE		BLACK	MORLGC

[Click Here to Order](#)

**EDGEFILL for EdgeMorphology™**

QTY **\$35.95 /6 pk**

	SMALL		YELLOW	MORSMO
	PRIMARY		RED	MORPRIO
	MEDIUM		BLUE	MORMEDO
	LARGE		BLACK	MORLGO

[Click Here to Order](#)





# EDGEONE

For EdgeOne Fire™ and EdgeOne Platinum™



**EDGECORE for EdgeOne**





QTY **\$35.95 /5 pk + 1 Verifier**

	SMALL		YELLOW	EOPSM20C
	PRIMARY		RED	EOPPRI25C
	MEDIUM		GREEN	EOPMED35C
	LARGE		WHITE	EOPLG45C

[Click Here to Order](#)

**EDGEFILL for EdgeOne**

QTY **\$35.95 /6 pk**

	SMALL		YELLOW	EOPSM20O
	PRIMARY		RED	EOPPRI25O
	MEDIUM		GREEN	EOPMED35O
	LARGE		WHITE	EOPLG45O

[Click Here to Order](#)

# EDGESEQUEL

## SAPPHIRE™



### EDGECORE for EdgeSequel Sapphire™

\$35.95 /5 pk + 1 Verifier

20		YELLOW	ES2006C
25		RED	ES2506C
30		BLUE	ES3006C
35		GREEN	ES3506C
40		BLACK	ES4006C
45		WHITE	ES4506C
50		YELLOW	ES5006C
55		RED	ES5506C
60		BLUE	ES6006C
70		GREEN	ES7006C
80		BLACK	ES8006C
90		WHITE	ES9006C

[Click Here to Order](#)

### EDGEFILL for EdgeSequel Sapphire™

\$35.95 /6 pk

20		YELLOW	ES2006O
25		RED	ES2506O
30		BLUE	ES3006O
35		GREEN	ES3506O
40		BLACK	ES4006O
45		WHITE	ES4506O
50		YELLOW	ES5006O
55		RED	ES5506O
60		BLUE	ES6006O
70		GREEN	ES7006O
80		BLACK	ES8006O
90		WHITE	ES9006O

[Click Here to Order](#)

# EDGETAPER™

For EdgeTaper™ and EdgeTaper Platinum™



### EDGECORE for EdgeTaper™

\$35.95 /5 pk + 1 Verifier

F1		YELLOW	ET2006C
F2		RED	ET2506C
F3		BLUE	ET3006C
F4		BLACK	ET4006C
F5		YELLOW	ET5006C

[Click Here to Order](#)

### EDGEFILL for EdgeTaper™

\$35.95 /6 pk

F1		YELLOW	ET2006O
F2		RED	ET2506O
F3		BLUE	ET3006O
F4		BLACK	ET4006O
F5		YELLOW	ET5006O

[Click Here to Order](#)



**EDGECORE for EdgeV-File HT™**  
**\$35.95 /5 pk + 1 Verifier**

20		YELLOW	VF2006C
25		RED	VF2506C
30		BLUE	VF3006C
35		GREEN	VF3506C
40		BLACK	VF4006C
45		WHITE	VF4506C
50		YELLOW	VF5006C

[Click Here to Order](#)

**EDGEFILL for EdgeV-File HT™**  
**\$35.95 /6PK**

20		YELLOW	VF2006O
25		RED	VF2506O
30		BLUE	VF3006O
35		GREEN	VF3506O
40		BLACK	VF4006O
45		WHITE	VF4506O
50		YELLOW	VF5006O




[Click Here to Order](#)



## EDGEVERIFY X1






(For WaveOne® Users)

**\$10.95 /6 pk**

20.06		YELLOW	X12006V
25.08		RED	X12508V
40.08		BLACK	X14008V









## EDGEVERIFY X3

**\$10.95 /6 pk**

20		YELLOW	X32006V
25		RED	X32506V
30		BLUE	X33006V
40		BLACK	X34006V
50		YELLOW	X35006V








## EDGEVERIFY X5

**\$10.95 /6 pk**

20.04		YELLOW	X52004V
20.06		YELLOW	X52006V
20.08		YELLOW	X52008V
30.04		BLUE	X53004V
30.06		BLUE	X53006V
30.08		BLUE	X53008V
40.04		BLACK	X54004V
40.06		BLACK	X54006V
40.08		BLACK	X54008V

## EDGEVERIFY X7

**\$10.95 /6 pk**

20		YELLOW	X72006V
25		RED	X72506V
30		BLUE	X73006V
35		GREEN	X73506V
40		BLACK	X74006V
45		WHITE	X74506V
50		YELLOW	X75006V

[Click Here to Order](#)

WaveOne® is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGEFILE® X7

## Bioceramic Gutta Percha



### .04 Gutta Percha Points

QTY **\$39.88 60/pk**

	17	<input type="checkbox"/>	X71704BCGP
	20	<input type="checkbox"/>	X72004BCGP
	25	<input type="checkbox"/>	X72504BCGP
	30	<input type="checkbox"/>	X73004BCGP
	35	<input type="checkbox"/>	X73504BCGP
	40	<input type="checkbox"/>	X74004BCGP
	45	<input type="checkbox"/>	X74004BCGP

[Click Here to Order](#)

### .06 Gutta Percha Points

QTY **\$39.88 60/pk**

	17	<input type="checkbox"/>	X71706BCGP
	20	<input type="checkbox"/>	X72006BCGP
	25	<input type="checkbox"/>	X72506BCGP
	30	<input type="checkbox"/>	X73006BCGP
	35	<input type="checkbox"/>	X73506BCGP
	40	<input type="checkbox"/>	X74006BCGP
	45	<input type="checkbox"/>	X74506BCGP

[Click Here to Order](#)

# EDGETAPER

## PLATINUM™

## Bioceramic Gutta Percha



### Gutta Percha Points

QTY **\$39.88 60/pk**

	F1	<input type="checkbox"/>	ETF1BCGP
	F2	<input type="checkbox"/>	ETF2BCGP
	F3	<input type="checkbox"/>	ETF3BCGP
	F4	<input type="checkbox"/>	ETF4BCGP
	F5	<input type="checkbox"/>	ETF5BCGP

[Click Here to Order](#)

# EDGEENDO



## Paper Points

QTY **\$7.99 /pk 100 Points**

	15.04	☐	WHITE	8806387821878
	15.06	☐	WHITE	8806387821922
	20.04	■	YELLOW	8806387821885
	20.06	■	YELLOW	8806387821861
	25.04	■	RED	8806387821731
	25.06	■	RED	8806387821663
	30.04	■	BLUE	8806387821748
	30.06	■	BLUE	8806387821670
	35.04	■	GREEN	8806387821755
	35.06	■	GREEN	8806387821687
	40.04	■	BLACK	8806387821892
	40.06	■	BLACK	8806387821939
	45.04	☐	WHITE	8806387821908
	45.06	☐	WHITE	8806387821946
	50.04	■	YELLOW	8806387821915
	50.06	■	YELLOW	8806387821953
	55.04	■	RED	8806387822028
	60.04	■	BLUE	8806387822035
	70.04	■	GREEN	8806387822042
	80.04	■	BLACK	8806387822059

## Assorted Paper Points

QTY **\$6.99 /60 pk**

	15-40.04	ASST	8806387821717
	15-40.06	ASST	8806387821724

## Gutta Percha Points

QTY **\$12.99 /pk 60 Points**

	15.04	☐	WHITE	8806387800996
	15.06	☐	WHITE	8806387801146
	20.04	■	YELLOW	8806387801009
	20.06	■	YELLOW	8806387801153
	25.04	■	RED	8806387801016
	25.06	■	RED	8806387801160
	30.04	■	BLUE	8806387801023
	30.06	■	BLUE	8806387801177
	35.04	■	GREEN	8806387801030
	35.06	■	GREEN	8806387801184
	40.04	■	BLACK	8806387801047
	40.06	■	BLACK	8806387801191
	45.04	☐	WHITE	8806387801054
	45.06	☐	WHITE	8806387801207
	50.04	■	YELLOW	8806387801061
	50.06	■	YELLOW	8806387801214
	55.04	■	RED	8806387801078
	60.04	■	BLUE	8806387801085
	70.04	■	GREEN	8806387801092
	80.04	■	BLACK	8806387801108

## Assorted Gutta Percha Points

QTY **\$11.99 /60 pk**

	15-40.04	ASST	8806387801115
	15-40.06	ASST	8806387801269

[Click Here to Order](#)

[Click Here to Order](#)

# EDGEACCESSORY™

## Paper Points

**\$7.99 200 points/pk**

	XXF, XF, F, M, C, XC	ASSORTED	EASSTPP
--	----------------------	----------	---------

[Click Here to Order](#)



## Gutta Percha Points

**\$12.99 100 points/pk**

	XXF, XF, F, M, C, XC	ASSORTED	EASSTGP
--	----------------------	----------	---------

[Click Here to Order](#)




# EDGE COIL™

## Paper Points

**\$7.99 /pk 100 Points**

	15 - 40	ASSORTED	ECASSTPP1540
	45 - 80	ASSORTED	ECASSTPP4580

[Click Here to Order](#)



## Gutta Percha Points

**\$12.99 /pk 60 Points**

	15 - 40	ASSORTED	ECASSTGP1540
	45 - 80	ASSORTED	ECASSTGP4580

[Click Here to Order](#)







## DIA-SYSTEM Paper Points

QTY **\$6.99 / 100 pk**

	20.04		YELLOW	GTP2004PP
	20.06		YELLOW	GTP2006PP
	20.08		YELLOW	GTP2008PP
	20.10		YELLOW	GTP2010PP
	25.04		RED	GTP2504PP
	25.06		RED	GTP2506PP
	25.08		RED	GTP2508PP
	25.10		RED	GTP2510PP
	30.04		BLUE	GTP3004PP
	30.06		BLUE	GTP3006PP
	30.08		BLUE	GTP3008PP
	30.10		BLUE	GTP3010PP
	35.04		GREEN	GTP3504PP
	35.06		GREEN	GTP3506PP
	35.08		GREEN	GTP3508PP
	35.10		GREEN	GTP3510PP
	40.04		BLACK	GTP4004PP
	40.06		BLACK	GTP4006PP
	40.08		BLACK	GTP4008PP
	40.10		BLACK	GTP4010PP

## DIA-SYSTEM Gutta Percha Points

QTY **\$12.99 / 60 pk**

	20.04		YELLOW	GTP2004GP
	20.06		YELLOW	GTP2006GP
	20.08		YELLOW	GTP2008GP
	20.10		YELLOW	GTP2010GP
	25.04		RED	GTP2504GP
	25.06		RED	GTP2506GP
	25.08		RED	GTP2508GP
	25.10		RED	GTP2510GP
	30.04		BLUE	GTP3004GP
	30.06		BLUE	GTP3006GP
	30.08		BLUE	GTP3008GP
	30.10		BLUE	GTP3010GP
	35.04		GREEN	GTP3504GP
	35.06		GREEN	GTP3506GP
	35.08		GREEN	GTP3508GP
	35.10		GREEN	GTP3510GP
	40.04		BLACK	GTP4004GP
	40.06		BLACK	GTP4006GP
	40.08		BLACK	GTP4008GP
	40.10		BLACK	GTP4010GP

[Click Here to Order](#)



# EDGEONE FIRE™

(Use Exactly Like WaveOne™ Gold)











## Paper Points

## Gutta Percha Points

QTY **\$7.99 /100 pk**

QTY **\$12.99 /60 pk**

	SMALL		YELLOW	EOFSMPP
	PRIMARY		RED	EOFPRIPP
	MEDIUM		GREEN	EOFMEDPP
	LARGE		WHITE	EOFLGPP

	SMALL		YELLOW	EOFSMGP
	PRIMARY		RED	EOFPRIGP
	MEDIUM		GREEN	EOFMEDGP
	LARGE		WHITE	EOFLGPP

[Click Here to Order](#)

[Click Here to Order](#)

WaveOne® Gold is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGEONE

## PLATINUM™



### Paper Slim Points

QTY **\$11.99 / 100 pk**

	VERIFILE™	<input type="checkbox"/>	WHITE	EOPVERPP
	SMALL	<input type="checkbox"/>	YELLOW	EOPSMPP
	PRIMARY	<input type="checkbox"/>	RED	EOPPRIPP
	MEDIUM	<input type="checkbox"/>	GREEN	EOPMEDPP
	LARGE	<input type="checkbox"/>	WHITE	EOPLGPP

[Click Here to Order](#)

### Gutta Percha Slim Points

QTY **\$11.99 / 60 pk**

	VERIFILE™	<input type="checkbox"/>	WHITE	EOPVERGP
	SMALL	<input type="checkbox"/>	YELLOW	EOPSMGP
	PRIMARY	<input type="checkbox"/>	RED	EOPPRIGP
	MEDIUM	<input type="checkbox"/>	GREEN	EOPMEDGP
	LARGE	<input type="checkbox"/>	WHITE	EOPLGGP

[Click Here to Order](#)

# EDGESEQUEL

## SAPPHIRE™



### Paper Points

**\$7.99 100/pk**

	15.04		WHITE	ESS1504PP
	15.06		WHITE	ESS1506PP
	20.04		YELLOW	ESS2004PP
	20.06		YELLOW	ESS2006PP
	25.04		RED	ESS2504PP
	25.06		RED	ESS2506PP
	30.04		BLUE	ESS3004PP
	30.06		BLUE	ESS3006PP
	35.04		GREEN	ESS3504PP
	35.06		GREEN	ESS3506PP
	40.04		BLACK	ESS4004PP
	40.06		BLACK	ESS4006PP
	45.04		WHITE	ESS4504PP
	45.06		WHITE	ESS4506PP
	50.04		YELLOW	ESS5004PP
	50.06		YELLOW	ESS5006PP
	55.04		RED	ESS5504PP
	60.04		BLUE	ESS6004PP
	70.04		GREEN	ESS7004PP
	80.04		BLACK	ESS8004PP

[Click Here to Order](#)



### Gutta Percha Points

**\$12.99 60/pk**

	15.04		WHITE	ESS1504GP
	15.06		WHITE	ESS1506GP
	20.04		YELLOW	ESS2004GP
	20.06		YELLOW	ESS2006GP
	25.04		RED	ESS2504GP
	25.06		RED	ESS2506GP
	30.04		BLUE	ESS3004GP
	30.06		BLUE	ESS3006GP
	35.04		GREEN	ESS3504GP
	35.06		GREEN	ESS3506GP
	40.04		BLACK	ESS4004GP
	40.06		BLACK	ESS4006GP
	45.04		WHITE	ESS4504GP
	45.06		WHITE	ESS4506GP
	50.04		YELLOW	ESS5004GP
	50.06		YELLOW	ESS5006GP
	55.04		RED	ESS5504GP
	60.04		BLUE	ESS6004GP
	70.04		GREEN	ESS7004GP
	80.04		BLACK	ESS8004GP






[Click Here to Order](#)

# EDGETAPER ENCOREPLATINUM™



## Paper Slim Points

QTY **\$7.99 100/pk**






	X1		ETEX1PP
	X2		ETEX2PP
	X3		ETEX3PP
	X4		ETEX4PP
	X5		ETEX5PP

[Click Here to Order](#)





## Gutta Percha Slim Points

QTY **\$12.99 60/pk**





	X1		ETEX1GP
	X2		ETEX2GP
	X3		ETEX3GP
	X4		ETEX4GP
	X5		ETEX5GP

[Click Here to Order](#)



**Paper Points**





QTY **\$7.99 /60 pk**

	SMALL		YELLOW	EMORSMP
	PRIMARY		YELLOW	EMORPRIP
	MEDIUM		YELLOW	EMORMEDP
	LARGE		YELLOW	EMORLGPP

[Click Here to Order](#)

**Gutta Percha Points**

QTY **\$12.99 /60 pk**

	SMALL		YELLOW	EMORSMP
	PRIMARY		YELLOW	EMORPRIP
	MEDIUM		YELLOW	EMORMEDP
	LARGE		YELLOW	EMORLGPP






[Click Here to Order](#)

# EDGETAPER and EDGETAPER™ PLATINUM™



## Paper Slim Points






QTY **\$6.99 100/pk**

	F1		ETF1GP
	F2		ETF2GP
	F3		ETF3GP
	F4		ETF4GP
	F5		ETF5GP

[Click Here to Order](#)

## Gutta Percha Slim Points

QTY **\$11.99 60/pk**











	F1		ETF1PP
	F2		ETF2PP
	F3		ETF3PP
	F4		ETF4PP
	F5		ETF5PP

[Click Here to Order](#)

# EDGE V-FILE<sup>HT</sup>™













## Paper Points \$7.99 /100 pk

QTY				
	20.06		YELLOW	EVHTPP2006
	25.06		RED	EVHTPP2506
	22.07		YELLOW	EVHTPP2207
	25.08		RED	EVHTPP2508
	30.06		BLUE	EVHTPP3006
	30.08		BLUE	EVHTPP3008
	35.06		GREEN	EVHTPP3506
	40.06		BLACK	EVHTPP4006
	45.06		WHITE	EVHTPP4506
	50.06		YELLOW	EVHTPP5006

[Click Here to Order](#)



## Gutta Percha Points \$12.99 /60 pk

QTY				
	20.06		YELLOW	EVHTGP2006
	25.06		RED	EVHTGP2506
	22.07		YELLOW	EVHTGP2207
	25.08		RED	EVHTGP2508
	30.06		BLUE	EVHTGP3006
	30.08		BLUE	EVHTGP3008
	35.06		GREEN	EVHTGP3506
	40.06		BLACK	EVHTGP4006
	45.06		WHITE	EVHTGP4506
	50.06		YELLOW	EVHTGP5006

[Click Here to Order](#)

# GTP Greater Taper Points



## DIA-SYSTEM Paper Points \$9.99 100/pk

20.04	YELLOW	GTP2004GP
20.06	YELLOW	GTP2006GP
20.08	YELLOW	GTP2008GP
20.10	YELLOW	GTP2010GP
25.04	RED	GTP2504GP
25.06	RED	GTP2506GP
25.08	RED	GTP2508GP
25.10	RED	GTP2510GP
30.04	BLUE	GTP3004GP
30.06	BLUE	GTP3006GP
30.08	BLUE	GTP3008GP
30.10	BLUE	GTP3010GP
35.04	GREEN	GTP3504GP
35.06	GREEN	GTP3506GP
35.08	GREEN	GTP3508GP
35.10	GREEN	GTP3510GP
40.04	BLACK	GTP4004GP
40.06	BLACK	GTP4006GP
40.08	BLACK	GTP4008GP
40.10	BLACK	GTP4010GP

[Click Here to Order](#)

## DIA-SYSTEM Gutta Percha Points \$15.99 60/pk

20.04	YELLOW	GTP2004GP
20.06	YELLOW	GTP2006GP
20.08	YELLOW	GTP2008GP
20.10	YELLOW	GTP2010GP
25.04	RED	GTP2504GP
25.06	RED	GTP2506GP
25.08	RED	GTP2508GP
25.10	RED	GTP2510GP
30.04	BLUE	GTP3004GP
30.06	BLUE	GTP3006GP
30.08	BLUE	GTP3008GP
30.10	BLUE	GTP3010GP
35.04	GREEN	GTP3504GP
35.06	GREEN	GTP3506GP
35.08	GREEN	GTP3508GP
35.10	GREEN	GTP3510GP
40.04	BLACK	GTP4004GP
40.06	BLACK	GTP4006GP
40.08	BLACK	GTP4008GP
40.10	BLACK	GTP4010GP

[Click Here to Order](#)



# Cell Packs













**Sterile**

Made from specially formulated highly absorbent papers, our Cell Packs are millimeter marked to provide preliminary guidance in measuring the root canal depth before proceeding clinically with the application of Gutta Percha Points



## DiaDent Paper Points

QTY **\$11.40 120/pk**

	15.04		WHITE	MP258-703
	15.06		WHITE	MP259-703
	20.04		YELLOW	MP258-704
	20.06		YELLOW	MP259-704
	25.04		RED	MP258-705
	25.06		RED	MP259-705
	30.04		BLUE	MP258-706
	30.06		BLUE	MP259-706
	35.04		GREEN	MP258-707
	35.06		GREEN	MP259-707
	40.04		BLACK	MP258-708
	40.06		BLACK	MP259-708

[Click Here to Order](#)

# EDGESEAL™

QTY **\$29.95**

<input type="text"/>	13.5G DUAL SYRINGE	818365020010
----------------------	--------------------	--------------



Resin-based root canal sealer designed to work with various obturation methods. Provides excellent biocompatibility. Easy to mix paste. Hermetic sealing ability. Non-staining to teeth. Insoluble in tissue fluids. Good radiopacity.

[Click Here to Order](#)

# EDGECAOH™

QTY **\$19.95**

<input type="text"/>	2.2G SYRINGE	00818365020003
----------------------	--------------	----------------



Temporary root-canal filler for weeping canals that forms a hard-tissue barrier (apexification). Internal and external root resorption. Pre-mixed paste. Syringe tip for easy delivery of paste into canal. 10 Disposable tips included to prevent cross-contamination.

[Click Here to Order](#)

# EDGE UTOPIA™

## BioCeramic Sealer 2g

QTY **\$69.88**

<input type="text"/>	2G SYRINGE + 20 DISPOSABLE TIPS	EUBCSEAL2G
----------------------	---------------------------------	------------



Excellent flow characteristics can be used in both warm vertical and single cone obturation techniques. Hydrophilic properties are ideal for setting in a moist environment. No Mixing Required. Has a setting time of 120-240 minutes after insertion into the canal and does not shrink when set.

[Click Here to Order](#)

# NeoPRO™

## BioCeramic Sealer

QTY **\$79.88**

<input type="text"/>	2G SYRINGE + 15 Disposable Tips	NEOPRO2G
----------------------	---------------------------------	----------



NeoPRO™ is a bioactive paste consisting of an extremely fine, inorganic powder of tricalcium dicalcium silicate in an organic medium. The product is packaged ready-to-use. No mixing is required. NeoPRO™ is designed to set in vivo in the presence of moisture provided by the surrounding tissues.

- **Bioactive bioceramic, Does not discolor teeth, Radiopaque, Resin-free**

Working Time: > 30 min. No mixing is required. NeoPRO™ BioCeramic Sealer can be applied immediately and introduced directly into the root canal.

Initial Setting Time: at 37°C, in vivo (or moist environment): 11 hrs ± 1 hr.

[Click Here to Order](#)

# EDGE BIO CERAMIC™ SEALER

QTY **\$129.88**

<input type="text"/>	2G SYRINGE + 15 Disposable Tips	818365020072
----------------------	---------------------------------	--------------

### EdgeBIOCERAMIC™ Tips

QTY **\$17.88**

<input type="text"/>	15 TIPS	EEBCT15
----------------------	---------	---------



EdgeBioCeramic™ Sealer features a unique resin-free formula making it highly biocompatible and hydrophilic. Unlike conventional sealers, EdgeBioCeramic™ Sealer does not shrink and therefore it is not necessary to compact the gutta percha.

- **Biocompatible, Premixed, Hydrophilic, Highly antibacterial, Radiopaque, Re-treatable, Zero shrinkage, Cost effective.**

Working Time: No mixing is required. EdgeBioCeramic™ Sealer can be applied immediately and introduced directly into the root canal.

Setting Time: 4 hours, however in very dry canals the setting time can be more than 10 hours.

[Click Here to Order](#)

# EDGE BIO CERAMIC™ THERMALFLOW™



QTY \$129.88

<input type="text"/>	1.5G SYRINGE + 15 Disposable Tips	818365020232
----------------------	-----------------------------------	--------------

EdgeBioCeramic™ ThermalFlow™ Root Canal Sealer features a lower viscosity when heated and is more radiopaque. Ideal for use in warm obturation techniques!

- **Biocompatible, Premixed, Hydrophilic, Highly antibacterial, Radiopaque, Re-treatable, Zero shrinkage, Cost effective.**

Working Time: No mixing is required. EdgeBioCeramic™ ThermalFlow™ Root Canal Sealer can be applied immediately and introduced directly into the root canal.

[Click Here to Order](#)

# EDGE BIO CERAMIC™ RETROFILL™ AND PERFORATION REPAIR



QTY \$189.88

<input type="text"/>	0.5G SYRINGE	818365020089
----------------------	--------------	--------------

EdgeBioCeramic™ RetroFill™ is a unique bioceramic repair material that exhibits superior healing and handling properties. EdgeBioCeramic™ RetroFill™ is pre-mixed and ready to use, unlike MTA. Simply syringe out a small amount of material and deliver to the repair site.

Working Time: No mixing is required. The setting reaction begins as soon as the material is placed in contact with a moist environment. The working time is no more than 30 minutes.

Setting Time: Minimum of 2 hours in normal conditions but can take longer to set in extremely dry canals.

[Click Here to Order](#)

# EDGE TEMP™ PLUS

QTY \$10.95

<input type="text"/>	40G	WHITE PASTE	ETEMPPLUS
----------------------	-----	-------------	-----------



Hydraulic temporary restorative material. Designed for fabricating temporary fillings and other temporary prosthetics until the permanent restoration can be delivered. Setting time: 30 minutes.

[Click Here to Order](#)


# ENDORING®

## EndoRing® II Kit with REDUCED GLARE Metal Ruler and e-Foam® Rotary HD Inserts

QTY **\$49.88**

	1 KIT	ER2GLEE-S
--	-------	-----------

The Kit comes with 1 EndoRing II Assembly – 1 EndoRing Cup, 3 Finger Rings (S,M,L), 1 Reusable GelWell Cup, 1 GelWell Adaptor Ring PLUS 4 Single Use GelWell Cups, 6 Silver e-Foam HD Inserts, 1 Reduced Glare Black Metal Ruler



## EndoRing® II Kit with Metal Ruler

QTY **\$41.88**

	1 KIT	ER2BEE-S
--	-------	----------

The Kit comes with 1 EndoRing II Assembly – 1 EndoRing Cup, 3 Finger Rings (S,M,L), 1 Reusable GelWell Cup, 1 GelWell Adaptor Ring PLUS 4 Single Use GelWell Cups, 6 Silver e-Foam HD Inserts, 1 Reduced Glare Black Metal Ruler



## Foam Refills

QTY **\$39.88 48/pk**

	Silver Rotary HD e-foam® Insert Refill	ERFSHDEE-S
	Blue e-foam® Insert Refill	ERFBEE-S
	Yellow e-foam® Insert Refill	ERFYEE-S
	Combo e-foam® Insert Refill (Yellow & blue)	ERFEE-S

[Click Here to Order](#)



## EDGEBURS®

Scan or Click For EdgeBurs® Catalog



Our mission is to deliver dental products and solutions, at a substantially lower cost, which in turn benefits practitioners and patients everywhere. Innovation is the heart of EdgeEndo®. We are dedicated to the pursuit of bringing leading-edge products to the industry. This thought and modern business model is revolutionizing dental practices across the U.S. and the world.


We believe premium technology shouldn't have to come with a premium price tag. We sell direct to dental professionals, without unnecessary costs, which characterizes other dental companies and are reflected in inflated prices.



**Scan or Click  
For Clinical Cases**



**Scan or Click  
For Clinical Studies**



**Scan or Click  
For Articles**

**EdgeEndo.com** • 1-855-985-3636

Click to follow us on:   