


EDGE X7

UTOPIA™

Utopia Alternative to EdgeFile® X7



\$35.95
6/ Pre-Sterilized Pack

- Premium version of EdgeFile® X7
- More Efficient Cutting
- Constant .04 or .06 Taper
- 21, 25, 29mm Lengths
- Assorted Pack Available

EDGE DIAMOND

UTOPIA™

Utopia Alternative to ProTaper Ultimate™




\$35.95
6/ Pre-Sterilized Pack

- Rotary or Forward Reciprocation
- Continuously Changing Diamond Shaped Cross-section
- Minimally Invasive Design
- Progressive and Regressive Variable Taper
- 21mm & 25mm lengths
- Assorted Pack Available

EDGE TAPER BLAZE

UTOPIA™

Utopia Alternative to ProTaper Gold® & EdgeTaper Platinum™




\$35.95
6/ Pre-Sterilized Pack

- Premium Version Of Edgetaper Platinum™
- Electropolished Blade For Smoother Navigation
- Improved Cutting Efficiency
- Variable Taper System
- 21, 25, 31mm Lengths
- Assorted Pack Available

EDGE SEQUEL

UTOPIA™

Utopia Alternative to Vortex® and EndoSequence® & EdgeSequel Sapphire



\$35.95
6/ Pre-Sterilized Pack

- Premium Version Of EdgeSequel Sapphire™
- Shorter Handle For Improved Posterior Access
- Constant .02, .04 Or .06 Taper
- #25/.08 19mm Also Available
- 21, 25, 31mm Lengths
- Assorted Pack Available

EDGE ONE BLAZE

UTOPIA™

Utopia Alternative to WaveOne® Gold & EdgeOne Fire



\$25.88
3/ Pre-Sterilized Pack

- Premium Version Of Edgeone Fire™
- Our Greatest Single File Shaping System Yet
- Resilient To Unwinding
- Variable Taper System
- 21, 25, 31mm Lengths
- Assorted Pack Available

Our Next Big....Little Thing!

EDGE FILE® X7



#1 Fan Favorite File System

\$35.95
6/ PACK

HEAT-TREATED FIREWIRE™ NiTi

- Universal file system
- 2x more cyclic fatigue resistant than other file systems
- Tips 17-45, 04 & 06 constant taper
- 21, 25, & 29mm
- Available Assorted Pack

EDGE V-TAPER HT™



VTaper® 2H & Endosequence® VT Scout™ Alternative

\$35.95
6 PER PACK

HEAT-TREATED FIREWIRE™ NiTi

- Identical to SS White V-Taper® 2H
- Minimally Invasive
- Expansive size portfolio available
- 21, 25, 31mm

EDGE TAPER PLATINUM™



ProTaper® Gold Alternative

\$35.95
6 PER PACK

HEAT-TREATED FIREWIRE™ NiTi

- 2x more resistant to cyclic fatigue* compared to ProTaper® Gold.
- SX, S1, S2, F1, F2, F3, F4, F5
- 21, 25, & 31mm
- Available Assorted Pack

EDGE TAPER™




ProTaper® Alternative

\$35.95
6/ PACK

- 6x more resistant than ProTaper®
- SX, S1, S2, F1, F2, F3, F4, F5
- 21, 25, & 31mm
- Available Assorted Pack

EDGE SEQUEL SAPPHIRE™



Vortex® Blue and EndoSequence® CM Alternative

\$25.88
4 PER PACK

HEAT-TREATED FIREWIRE™ NiTi

- 2x more resistant to cyclic fatigue* compared to Vortex® Blue, 8x more resistant than Vortex® and EndoSequence®
- Constant taper
- Tips 04 15-90 taper, 06 15-50 taper
- 21, 25, & 31mm
- Available Assorted Pack

EDGE ONE PLATINUM™



WaveOne® Gold Alternative

\$25.88
3 PER PACK

HEAT-TREATED FIREWIRE™ NiTi

- Reverse Reciprocating
- Triangular Cross-section
- Verifile, Small, Primary, Medium & Large
- 21, 25, & 31mm
- Available Assorted Pack

EDGE ONE FIRE™



WaveOne® Gold Alternative

\$25.88
3 PER PACK

HEAT-TREATED FIREWIRE™ NiTi

- 5X more resistant to cyclic fatigue* compared to WaveOne® Gold
- Reverse Reciprocating
- Parallelogram Cross-section
- Small, Primary, Medium & Large
- 21, 25, & 31mm
- Available Assorted Pack

*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph: <https://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> ProTaper®, ProTaper Gold®, Vortex®, Vortex Blue®, ProTaper Ultimate™ and WaveOne® Gold are registered trademarks of Dentsply Tulsa Dental Specialties. EndoSequence®, EndoSequence® CM, EndoSequence VT Scout™ are registered trademarks of Brasseler USA. V-Taper® 2H is a registered trademark of SS White.

ACCESS • SHAPING • IRRIGATION • OBTURATION • SEALING

EDGE SURGICAL TIP™



\$49.88
Per Tip

EDGE TIP™



\$49.88
Per Tip



Take Advantage of
Our Newly Lowered
Prices on Select Burs!



EDGE K-FILE™

\$5.15
6 pack

- Sizes: 6-80
- Alternative to Tulsa Dental K-file®
- 21, 25, & 31mm Lengths



EDGE C-FILE™

\$5.95
6 pack

- Sizes: 6, 8, 10, 12, 15
- 21 & 25mm

Also Available: EdgeHeadstoms - \$5.15 • EdgeK-File NiTi - \$7.88



EDGEAPEX HP™

\$1095.00

- Pre-programmed with top selling EdgeEndo files
- 10 fully customizable Memory Program slots
- Integrated Apex Locator during instrumentation
- Option to use Apex locator with handfile
- 360° rotation of contra angle



EDGEICE™

Pulp Vitality
Refrigerant Spray

\$15.88



EDGEMIX™

EDTA-based formula
enhanced with
Chlorhexidine for added
bacterial properties

\$59.88

Also Available: EdgeLube, EDTA, Sodium Hypochlorite, Chlorhexidine



EDGECORE™

Pink Gutta Percha Core
Thermal Obturator

\$35.95

5 pack & 1 verifier



EDGEFILL™


Plastic Carrier
Thermal Obturator

\$35.95

6 Pack

Also Available: Universal sized EdgeCore and EdgeFill 20pks Available - \$100.88

NEW!



EdgeEndo® Gutta Percha Points - \$12.99

Optimized Apical Fit with Coronal Venting

- Tight Apical Fit
- Extra Space In The Coronal Third To Facilitate The Upward Venting Of Sealer
- Minimizes The Risk Of Apical Extrusion



EdgeEndo® Paper Points

Universal Constant Taper

\$7.99



EdgeEndo® Gutta Points

Universal Constant Taper

\$12.99

EDGE BIO CERAMIC™
SEALER



\$129.88

Per 2g syringe +
15 reduced waste tips

#1 Fan Favorite sealer with
zero shrinkage

EDGE BIO CERAMIC™
THERMALFLOW™



\$129.88

Per 1.5g syringe +
15 reduced waste tips

Low viscosity ideal for
warm obturation

EDGE UTOPIA®
BioCeramic Sealer



\$72.88

Per 2g syringe + 20 tips

High pH of 12.5 for optimal
antibacterial properties

EDGE BIO CERAMIC RETROFILL™
AND PERFORATION REPAIR MATERIAL



\$189.88

Per 0.5g syringe

For all repair procedures and
pediatric pulpotomy

EDGE TEMP™ PLUS



Hydraulic Temporary Restorative
Alternative to Cavit™ by 3M

\$10.95



EDGEFLOW™

\$695.00

- Excellent gutta-percha fluidity
- Uniform and smooth gutta-percha output without bubbles.
- Precise and intuitive display of remaining guttapercha
- Simple preset heating with the press of a button.



EDGEPACK™

\$545.00

- High precision temperature control
- 0.35mm pre-bend tip
- Variety of Tapers Available in Tip Sizes 35-60
- Improved penetrability into narrow root canals
- Ensures smooth cutting of gutta percha

Cavit™ is a registered trademark of 3M. Tulsa Dental K-File® is a registered trademark of DENTSPLY Tulsa Dental Specialties



Product
Evaluation
Request

Fast & Free Shipping On All Orders
Hassle Free Returns • Straightforward Pricing
• Exceptional Customer Service



Full
Product
Portfolio



EDGEENDO®

Rev. A 02/2025