



**EDGEENDO®**

# 2025 RESIDENT PRICING PRODUCT CATALOG

The Endodontic Graduate Program from EdgeEndo® is a special purchase program for first and second year residents. We're excited to provide you this opportunity to work with our products at a special discounted price. You may use this opportunity during residency and until December 31<sup>st</sup> of your graduation year.

Participants must create an account if they do not already have one. All products must be ordered by phone 1-855-985-3636. Limits are shown for each product category and are available to Mix/Match.

**[Click For More Information](#)**

## Created by an Endodontist for Endodontists

Established in 2012 by American Endodontist Charles Goodis, EdgeEndo® was founded with the vision to make root canal treatments safe, easy, and cost-effective overcoming the shortfalls of manufacturers offering outdated technology with a high price tag.

**The solution** – Innovation combined with high-quality and cost-effectiveness to support the growth of endodontic practices by saving chair time and delivering outstanding patient outcomes. In just ten years, EdgeEndo® revolutionized the endodontic industry in the USA becoming a leader in endodontic solutions!

**Safe** – Healthy Patients - Happy Dentists. EdgeEndo® solutions are developed with the patient at the center. Our FireWire heat treatment processes provides excellent resistance to cyclic fatigue and flexibility to carefully follow the canal, preserving canal anatomy and dentin and allowing for faster procedures. All for increased safety and productivity.

**Easy** – Changing habits is always a challenge for any dental professional. Our solutions are designed for ease of use and are compatible with most endodontic techniques. EdgeEndo® continues expanding its product offering to cover the entire endodontic procedure; diagnosis, instrumentation, obturation and restoration. No need for multiple sources and complexity. We make it easy for the dentists.


**Cost-Effective** – EdgeEndo® has proven that quality does not need to come with a premium price tag. Our goal is to make root canals accessible to anyone, without compromises.

**EdgeEndo® remains steadfast in its mission to provide game-changing products and solutions**

**Innovation** - EdgeEndo® is the inventor of FireWire and FireWire Blaze, a unique combination of heat treatment and cryogenic applications allowing NiTi instruments to deliver incredible flexibility and resistance to cyclic fatigue.

**Quality** – Thousands of clinicians trust EdgeEndo® and our mission to deliver a seamless and consistent experience with all our products. The control of manufacturing processes is our obsession, and our instruments are all manufactured in the USA, Switzerland, and France under strict quality controls.

**Revolutionize endodontic procedures with EdgeEndo®:  
the partner dental professionals have been searching for!**



# Table of Contents

## EdgePower Suite

|                       |    |
|-----------------------|----|
| EdgeApex HP . . . . . | 4  |
| EdgeFlow . . . . .    | 6  |
| EdgeApex . . . . .    | 9  |
| EdgePack . . . . .    | 10 |

## Diagnosis & Isolation

|                   |    |
|-------------------|----|
| EdgeIce. . . . .  | 13 |
| EdgeDam . . . . . | 13 |

## Supplies

|                    |    |
|--------------------|----|
| EndoRing . . . . . | 14 |
|--------------------|----|

## Initial Access

|                              |    |
|------------------------------|----|
| EdgeTip. . . . .             | 15 |
| EdgeBurs - Carbide . . . . . | 16 |
| EdgeGates. . . . .           | 19 |
| EdgeBroaches . . . . .       | 20 |
| EdgeHedstroms . . . . .      | 21 |
| EdgeTaper Retreat. . . . .   | 21 |
| BC Solvent. . . . .          | 22 |

## Glidepath

|                                        |    |
|----------------------------------------|----|
| EdgeFind. . . . .                      | 23 |
| EdgeGlidepath . . . . .                | 23 |
| EdgeOneGlidepath . . . . .             | 24 |
| EdgeOne Blaze Utopia - RGlide. . . . . | 24 |

## Access Refinement

|                              |    |
|------------------------------|----|
| EdgeSurgical Tip . . . . .   | 25 |
| EdgeBurs - Diamond . . . . . | 26 |

## Working Length Determination

|                           |    |
|---------------------------|----|
| EdgeK-File NiTi . . . . . | 30 |
| EdgeK-File . . . . .      | 31 |
| EdgeC-File. . . . .       | 32 |

## File Systems

### Rotary • Constant Taper

|                             |    |
|-----------------------------|----|
| EdgeX7 Utopia . . . . .     | 33 |
| EdgeFile X7 . . . . .       | 35 |
| EdgeFile X3 . . . . .       | 37 |
| EdgeSequel Utopia . . . . . | 39 |

|                              |    |
|------------------------------|----|
| EdgeSequel Sapphire. . . . . | 42 |
| EdgeEvolve . . . . .         | 44 |

### Rotary • Variable Taper

|                                  |    |
|----------------------------------|----|
| EdgeDiamond Utopia. . . . .      | 47 |
| EdgeV-Taper HT . . . . .         | 49 |
| EdgeTaper Blaze Utopia . . . . . | 51 |
| EdgeTaper Platinum. . . . .      | 53 |
| EdgeTaper . . . . .              | 55 |

### Reciprocating • Variable Taper

|                                |    |
|--------------------------------|----|
| EdgeOne Blaze Utopia . . . . . | 57 |
| EdgeOne Fire . . . . .         | 59 |

### Reciprocating • Constant Taper

|                            |    |
|----------------------------|----|
| EdgeOne Platinum . . . . . | 61 |
| EdgeFile X1 . . . . .      | 63 |

## Irrigation

|                                     |    |
|-------------------------------------|----|
| EdgeLube . . . . .                  | 65 |
| EdgeGel . . . . .                   | 65 |
| Irrigation Tips & Syringes. . . . . | 66 |
| Clensers . . . . .                  | 67 |

## Sealers

|                                       |    |
|---------------------------------------|----|
| EdgeBioCeramic Sealer. . . . .        | 68 |
| EdgeUtopia BioCeramic Sealer. . . . . | 68 |
| EdgeBioCeramic Thermal Flow . . . . . | 68 |
| EdgeSeal. . . . .                     | 69 |

## Obturation

|                     |    |
|---------------------|----|
| EdgeCore . . . . .  | 69 |
| Edge Fill . . . . . | 69 |
| Points. . . . .     | 70 |

## Intracanal Medicament & Repair

|                                            |    |
|--------------------------------------------|----|
| EdgeCaOH. . . . .                          | 75 |
| EdgeTemp Plus. . . . .                     | 75 |
| EdgeBioCeramic Retrofill & Repair. . . . . | 75 |
| EdgeUtopia RRM. . . . .                    | 76 |

## Starter Bundles . . . . . 78

## Category Limit Quantities. . . . . 79

# EDGEAPEX HP™



- Built in Apex Locator with option to use as a standalone Apex locator with handfile
- Pre-programmed with top selling EdgeEndo files PLUS 10 fully customizable Memory Program slots
- Efficient brushless motor, low noise, long service life
- Cordless portable endo motor with combined length determination
- 360° rotation of contra angle
- Adopt real-time feedback technology and dynamic torque control, effectively preventing file separation
- Wireless Charging
- Torque Control
- Auto On/Off, Torque Reverse, Apical Reverse, Apical Slow Down
- Visual torque on file in use
- Max Torque 5.0 Ncm
- Change Dominant Hand
- Range: 100 - 2500 RPM

**\$799.88**

# EDGEAPEX HP™

**\$799.88**

|        |        |
|--------|--------|
| 1 UNIT | EAPXHP |
|--------|--------|

## Accessories included in package:

- Measuring Wire
- Lip Hook (x2)
- File Clip (x4)
- Touch Probe (x2)
- Power Adapter
- Protective Silicon Cover (x2)
- Spray Nozzle
- O-Rings (x2)
- Disposable Sleeves (1 pack)

**Warranty:** 2 years for handpiece motor, charging base, and power adapter. 1 year for contra angle attachment. 6 months for accessories.



## Accessories

|                                 |
|---------------------------------|
| <b>MEASURING WIRE - \$17.88</b> |
| 1 PACK                          |
| EAPXMW                          |



|                                |
|--------------------------------|
| <b>CONTRA ANGLE - \$199.88</b> |
| 1 PACK                         |
| EAPXCA                         |



|                             |
|-----------------------------|
| <b>TOUCH PROBE - \$8.88</b> |
| 1 PACK                      |
| EAPXTP                      |



|                            |
|----------------------------|
| <b>FILE CLIP - \$59.88</b> |
| 5 PACK                     |
| EAPXFC                     |



|                           |
|---------------------------|
| <b>LIP HOOK - \$14.88</b> |
| 5 PACK                    |
| EAPXLH                    |



|                                |
|--------------------------------|
| <b>CHARGING BASE - \$74.88</b> |
| 1 PACK                         |
| EAPXCB                         |



# EDGEFLOW™



- Wireless and electric
- Excellent gutta-percha fluidity, uniform and smooth gutta-percha output without bubbles.
- Precise and intuitive display of remaining gutta-percha
- Auto heating stop after timeout; Auto shutdown after timeout
- Charging base with battery
- Temperature monitoring design
- Dual batteries deliver longer endurance, enabling filling about 100 root canals.
- Re-usable Injection needles available in various sizes and lengths
- Preset temperatures: 100°C, 120°C, 150°C, 180°C, 200°C
- Injection Speeds: Slow, Medium, Fast

**\$499.88**

# EDGEFLOW™

\$499.88

|        |       |
|--------|-------|
| 1 UNIT | EEPTO |
|--------|-------|

**Accessories included in package:**

- Main Unit
- Charging Stand
- Power Adapter
- Battery (x2)
- Heat Insulation Hood
- Wrench of Needle
- Gutta-Percha Ring
- Instruction Manual
- Qualified Certification
- Packing List
- Injection Needles (x4)
  - E23G-M (x1)
  - E23G-L (x1)
  - E25G-L (x2)

**Warranty:** 1 year for device.



## Accessories

|                         |                |
|-------------------------|----------------|
| <b>OBTURATOR WRENCH</b> | <b>\$13.88</b> |
| 1 PACK                  |                |
| EEPTOW                  |                |



|                          |                |
|--------------------------|----------------|
| <b>EDGEFLOW™ BATTERY</b> | <b>\$28.88</b> |
| 1 PACK                   |                |
| EEPTOB                   |                |



|                                       |                |
|---------------------------------------|----------------|
| <b>EDGEFLOW™ HEAT INSULATION HOOD</b> | <b>\$16.88</b> |
| 3 PACK                                |                |
| EEPTOHIH3                             |                |



# EdgeFlow™ Re-Useable Needle Tips

|                          |                |
|--------------------------|----------------|
| <b>NEEDLE TIP E20G-S</b> | <b>\$31.88</b> |
| <b>1 PACK</b>            |                |
| EE20GSIN                 |                |



|                          |                |
|--------------------------|----------------|
| <b>NEEDLE TIP E23G-L</b> | <b>\$31.88</b> |
| <b>1 PACK</b>            |                |
| EE23GLIN                 |                |



|                          |                |
|--------------------------|----------------|
| <b>NEEDLE TIP E20G-M</b> | <b>\$31.88</b> |
| <b>1 PACK</b>            |                |
| EE20GMIN                 |                |



|                          |                |
|--------------------------|----------------|
| <b>NEEDLE TIP E23G-M</b> | <b>\$31.88</b> |
| <b>1 PACK</b>            |                |
| EE25GMIN                 |                |



|                          |                |
|--------------------------|----------------|
| <b>NEEDLE TIP E20G-L</b> | <b>\$31.88</b> |
| <b>1 PACK</b>            |                |
| EE20GLIN                 |                |



|                          |                |
|--------------------------|----------------|
| <b>NEEDLE TIP E23G-L</b> | <b>\$31.88</b> |
| <b>1 PACK</b>            |                |
| EE25GLIN                 |                |



|                          |                |
|--------------------------|----------------|
| <b>NEEDLE TIP E23G-M</b> | <b>\$31.88</b> |
| <b>1 PACK</b>            |                |
| EE23GMIN                 |                |



## EDGEPELLETS™

|                 |                |
|-----------------|----------------|
| <b>100 PACK</b> | <b>\$22.49</b> |
| 00810184450092  |                |



# EDGEAPEX™

**\$549.88**

|        |       |
|--------|-------|
| 1 UNIT | EEALU |
|--------|-------|



## Accessories included in package:

- Base
- Measuring Wire
- File Clip (x4)
- Touch Probe (x2)
- Pulp Tester Probe (x2)
- Adapter
- Tester
- Enlarged display of apical foramen
- Instruction card for apical stop
- Apex Locator troubleshooting card

**Warranty:** 1 year for device. 6 months for accessories.

## Apex Locator Accessories

|                       |                |
|-----------------------|----------------|
| <b>MEASURING WIRE</b> | <b>\$18.73</b> |
| 1 PACK                |                |
| EEALUMW               |                |



|               |                |
|---------------|----------------|
| <b>TESTER</b> | <b>\$15.59</b> |
| 1 PACK        |                |
| EEALUT        |                |



|                 |                |
|-----------------|----------------|
| <b>LIP HOOK</b> | <b>\$15.59</b> |
| 5 PACK          |                |
| EEALULH5        |                |



|                   |                |
|-------------------|----------------|
| <b>ENDO RULER</b> | <b>\$15.59</b> |
| 1 PACK            |                |
| EEALUER           |                |



|                    |                |
|--------------------|----------------|
| <b>TOUCH PROBE</b> | <b>\$11.40</b> |
| 1 PACK             |                |
| EEALUMF            |                |



|                  |                |
|------------------|----------------|
| <b>FILE CLIP</b> | <b>\$62.72</b> |
| 5 PACK           |                |
| EEALUFC5         |                |



# EDGE PACK™



- Wireless Gutta Percha Obturation Device
- Ergonomic design and ambidextrous accessible
- Easy to use
- 0.35mm pre-bend tip
- Improved penetrability into narrow root canals
- High precision temperature control
- Smooth cutting of GP from high-efficiency conversion system
- Compact size for easy cleaning and storage
- Visual and auditory prompts for safe operation
- Preset temperatures: 150°C, 180°C, 200°C, 230°C
- Time related features: Warm up time within 0.2 seconds, Auto shut off after 10 seconds,
- 4 hour battery life with each charge

**\$429.88**

# EDGEPACK™

**\$429.88**

|        |       |
|--------|-------|
| 1 UNIT | EEDPU |
|--------|-------|

**Accessories included in package:**

- Charging Stand
- Power Adapter
- WP4004 Tip
- WP4504 Tip
- WP5506 Tip
- WP5508 Tip
- Heat Insulation Sleeve of Tips (x4)
- Instruction Manual and Qualified Certification

**Warranty:** 1 year for device and tips.



## Accessories

|                                |         |                |
|--------------------------------|---------|----------------|
| <b>EDGEPACK CHARGING STAND</b> |         | <b>\$58.88</b> |
| 1 UNIT                         | EEDPUCS |                |



# EdgePack Tips

|                          |                |
|--------------------------|----------------|
| <b>EDGEPACK 3504 TIP</b> | <b>\$33.96</b> |
| 1 PACK                   |                |
| EEDPUT3504               |                |



|                          |                |
|--------------------------|----------------|
| <b>EDGEPACK 5508 TIP</b> | <b>\$33.96</b> |
| 1 PACK                   |                |
| EEDPUT5508               |                |



|                          |                |
|--------------------------|----------------|
| <b>EDGEPACK 4004 TIP</b> | <b>\$33.96</b> |
| 1 PACK                   |                |
| EEDPUT4004               |                |



|                           |                |
|---------------------------|----------------|
| <b>EDGEPACK 5508L TIP</b> | <b>\$33.96</b> |
| 1 PACK                    |                |
| EEDPUT5508L               |                |




|                          |                |
|--------------------------|----------------|
| <b>EDGEPACK 4504 TIP</b> | <b>\$33.96</b> |
| 1 PACK                   |                |
| EEDPUT4504               |                |



|                          |                |
|--------------------------|----------------|
| <b>EDGEPACK 5510 TIP</b> | <b>\$33.96</b> |
| 1 PACK                   |                |
| EEDPUT5510               |                |



|                          |                |
|--------------------------|----------------|
| <b>EDGEPACK 5004 TIP</b> | <b>\$33.96</b> |
| 1 PACK                   |                |
| EEDPUT5004               |                |



|                          |                |
|--------------------------|----------------|
| <b>EDGEPACK 6012 TIP</b> | <b>\$33.96</b> |
| 1 PACK                   |                |
| EEDPUT6012               |                |



|                          |                |
|--------------------------|----------------|
| <b>EDGEPACK 5506 TIP</b> | <b>\$33.96</b> |
| 1 PACK                   |                |
| EEDPUT5506               |                |




|                    |              |                |
|--------------------|--------------|----------------|
| <b>EdgeIce™</b>    |              | <b>\$13.88</b> |
| 6 OZ. SPRAY BOTTLE | 818365020065 |                |

A quick way to check cold sensitivity and pulp vitality. Significantly colder than ethyl chloride or ice. EdgeIce™ features a directional nozzle for easy placement on micro applicator or cotton pellet.



## EDGE DAM™

|                  |                          |
|------------------|--------------------------|
| <b>Latex 5x5</b> | <b>52 Pack - \$10.69</b> |
| Medium           | EDBLM5552                |
| Heavy            | EDBLH5552                |
| <b>Latex 6x6</b> | <b>36 Pack- \$10.69</b>  |
| Medium           | EDBLM6636                |
| Heavy            | EDBLH6636                |



|                       |                          |
|-----------------------|--------------------------|
| <b>Latex Free 5x5</b> | <b>50 Pack - \$63.79</b> |
| Medium                | EDVLFM5550               |
| Heavy                 | EDVLFH5550               |
| <b>Latex Free 6x6</b> | <b>50 Pack - \$63.79</b> |
| Medium                | EDVLFM6650               |
| Heavy                 | EDVLFH6650               |



# ENDORING®

## EndoRing® II Kit with REDUCED GLARE Metal Ruler and e-Foam® Rotary HD Inserts \$42.88

|       |           |
|-------|-----------|
| 1 KIT | ER2GSEE-S |
|-------|-----------|

1 ring assembly & 6 silver HD foam inserts



## EndoRing® II Kit with Metal Ruler \$35.88

|       |          |
|-------|----------|
| 1 KIT | ER2BEE-S |
|-------|----------|

1 ring assembly, 3 blue foam inserts & 3 yellow foam inserts



## Foam Refills 48 Pack

**\$31.88**

|                                               |            |
|-----------------------------------------------|------------|
| Silver Rotary HD e-foam® Insert Refill        | ERFSHDEE-S |
| Blue e-foam® Insert Refill                    | ERFBEE-S   |
| Yellow e-foam® Insert Refill                  | ERFYEE-S   |
| Combo e-foam® Insert Refill (Yellow & blue) - | ERFEE-S    |



EdgeTips™ 1-5 are diamond-coated for more efficiency. We don't use zirconium nitride. EdgeTips™ 6-8 are a special titanium alloy and their special geometries allow precision work in areas of limited access.

- Find calcified or hidden canals.
- Remove intracanal obstructions.
- Remove dental restorations and cores.
- Dislodge posts and broken instruments.
- Remove remove pulp stones.
- Unique shaft designs for better access and visibility.
- Length and diameters designed to progress down the root.

**EdgeTip #1** - Diamond-coated ball tip, length 20mm, taper 5%.

ETBD



**EdgeTip #2** - Diamond-coated steel tip 76µm, length 18mm, taper 5%

ET18D



**EdgeTip #3** - Diamond-coated retreatment tip, 30µm, length 20mm, taper 5%

ET20D



**EdgeTip #4** - Long Diamond-coated retreatment tip, 30µm, 40mm, taper 4%

ET40D4



**EdgeTip #6** - Short Titanium- Niobium tip, length 15mm taper 4%

ED25S



**EdgeTip #7** - Titanium- Niobium tip, length 20mm, taper 3%

ET25



**EdgeTip #8** - Titanium- Niobium tip, 25mm, taper 3%


ET25L



## OPERATIVE CARBIDES




ROUND



| US NO.                 | ISO NO.                | 1/4                      | 1/2                 | 1                    | 2                | 3                    | 4                |
|------------------------|------------------------|--------------------------|---------------------|----------------------|------------------|----------------------|------------------|
| <b>FG (19MM)</b>       | <b>500 314 001 001</b> | <b>005</b>               | <b>006</b>          | <b>008</b>           | <b>010</b>       | <b>012</b>           | <b>014</b>       |
|                        | 10 PACK                | \$9.30<br>FG1QTR10       | \$16.63<br>FG1HLF10 | \$16.63<br>FG110     | \$16.63<br>FG210 | \$9.30<br>FG310      | \$16.63<br>FG410 |
| <b>FG SURG. (25MM)</b> | <b>500 316 001 001</b> |                          | <b>006</b>          |                      | <b>010</b>       |                      | <b>014</b>       |
|                        | 10 PACK                | \$16.63<br>FG1HLF-SURG10 |                     | \$16.63<br>FG2SURG10 |                  | \$16.63<br>FG4SURG10 |                  |
| <b>RA (22MM)</b>       | <b>500 204 001 001</b> |                          | <b>006</b>          |                      | <b>010</b>       | <b>012</b>           | <b>014</b>       |
|                        | 10 PACK                |                          | \$16.63<br>RA1HLF10 |                      | \$16.63<br>RA210 | \$10.35<br>RA310     | \$16.63<br>RA410 |

ROUND



| US NO.                 | ISO NO.                | 5                | 6                    | 7                | 8                    | 10                |
|------------------------|------------------------|------------------|----------------------|------------------|----------------------|-------------------|
| <b>FG (19MM)</b>       | <b>500 314 001 001</b> | <b>016</b>       | <b>018</b>           |                  | <b>023</b>           |                   |
|                        | 10 PACK                | \$16.63<br>FG510 | \$16.63<br>FG610     |                  | \$16.63<br>FG810     |                   |
| <b>FG SURG. (25MM)</b> | <b>500 316 001 001</b> |                  | <b>018</b>           |                  | <b>023</b>           |                   |
|                        | 10 PACK                |                  | \$16.63<br>FG6SURG10 |                  | \$16.63<br>FG8SURG10 |                   |
| <b>RA (22MM)</b>       | <b>500 204 001 001</b> |                  | <b>018</b>           | <b>021</b>       | <b>023</b>           | <b>027</b>        |
|                        | 10 PACK                |                  | \$16.63<br>RA610     | \$10.35<br>RA710 | \$16.63<br>RA810     | \$16.63<br>RA1010 |



INVERTED CONE



| US NO.           | ISO NO.                | 34                | 35                |
|------------------|------------------------|-------------------|-------------------|
| <b>FG (19MM)</b> | <b>500 314 010 001</b> | <b>008</b>        | <b>010</b>        |
|                  | 10 PACK                | \$16.63<br>FG3410 | \$16.63<br>FG3510 |



AMALGAM PREP



| US NO.                 | ISO NO.                | 245                  |
|------------------------|------------------------|----------------------|
| <b>FG (19MM)</b>       | <b>500 314 234 006</b> | <b>009</b>           |
|                        | 10 PACK                | \$16.63<br>FG24510   |
| <b>FG SHORT (16MM)</b> | <b>500 313 234 006</b> | <b>009</b>           |
|                        | 10 PACK                | \$16.63<br>FG245SS10 |



PEAR

| US NO.                 | ISO NO.                | 330                  | 331               | 332               |
|------------------------|------------------------|----------------------|-------------------|-------------------|
| <b>FG (19MM)</b>       | <b>500 314 237 001</b> | <b>008</b>           | <b>010</b>        | <b>012</b>        |
|                        | 10 PACK                | \$16.63<br>FG33010   | \$9.30<br>FG33110 | \$9.30<br>FG33210 |
| <b>FG SHORT (16MM)</b> | <b>500 313 237 001</b> | <b>008</b>           |                   |                   |
|                        | 10 PACK                | \$16.63<br>FG330SS10 |                   |                   |



PEAR LONG HEAD

| US NO.           | ISO NO.                | 331L                |
|------------------|------------------------|---------------------|
| <b>FG (19MM)</b> | <b>500 314 238 006</b> | <b>010</b>          |
|                  | 10 PACK                | \$16.63<br>FG331L10 |



FLAT FISSURE CROSS CUT

| US NO.                 | ISO NO.                | 556                | 557                    | 558                    |
|------------------------|------------------------|--------------------|------------------------|------------------------|
| <b>FG (19MM)</b>       | <b>500 314 107 007</b> | <b>008</b>         | <b>010</b>             | <b>012</b>             |
|                        | 10 PACK                | \$16.63<br>FG55610 | \$16.63<br>FG55710     | \$89.30<br>FG55810     |
| <b>FG SURG (25MM)</b>  | <b>500 316 107 007</b> |                    | <b>010</b>             | <b>012</b>             |
|                        | 10 PACK                |                    | \$16.63<br>FG557SURG10 | \$12.44<br>FG558SURG10 |
| <b>FG SHORT (16MM)</b> | <b>500 313 107 007</b> |                    | <b>010</b>             |                        |
|                        | 10 PACK                |                    | \$ 16.63<br>FG557SS10  |                        |



FLAT FISSURE  
CROSS CUT LONG HEAD

| US NO.           | ISO NO.                | 557L                |
|------------------|------------------------|---------------------|
| <b>FG (19MM)</b> | <b>500 314 110 007</b> | <b>010</b>          |
|                  | 10 PACK                | \$16.63<br>FG557L10 |



DOMED FISSURE PLAIN  
CUT

| US NO.           | ISO NO.                | 1157               | 1158               |
|------------------|------------------------|--------------------|--------------------|
| <b>FG (19MM)</b> | <b>500 314 137 006</b> | <b>010</b>         | <b>012</b>         |
|                  | 10 PACK                | \$9.30<br>FG115710 | \$9.30<br>FG115810 |



DOMED FISSURE CROSS  
CUT

| US NO.           | ISO NO.                | 1557               | 1558               |
|------------------|------------------------|--------------------|--------------------|
| <b>FG (19MM)</b> | <b>500 314 137 007</b> | <b>010</b>         | <b>012</b>         |
|                  | 10 PACK                | \$9.30<br>FG155710 | \$9.30<br>FG155810 |

## SPECIALTY



151-Z

| FIG NO.             | 151-Z             |
|---------------------|-------------------|
| DIAMETER (1/10MM)   | 016               |
| H/LENGTH (MM)       | 11.0              |
| OVERALL LENGTH (MM) | 23.0              |
| 5 PACK              | \$23.97<br>151ZB5 |



ENDO ACCESS

| FIG NO.             | 152              |
|---------------------|------------------|
| DIAMETER (1/10MM)   | 016              |
| H/LENGTH (MM)       | 9.0              |
| OVERALL LENGTH (MM) | 23               |
| 5 PACK              | \$23.97<br>152B5 |



LINDEMANN  
CROSS CUT

| US NO.          | 162                    |
|-----------------|------------------------|
| H/LENGTH (MM)   | 10.0                   |
| FG SURG. (25MM) | 016                    |
| 5 PACK          | \$23.97<br>FGSURG162B5 |



D LINDEMANN  
CROSS CUT

| US NO.          | 162A                    |
|-----------------|-------------------------|
| H/LENGTH (MM)   | 9.2                     |
| FG SURG. (25MM) | 016                     |
| 5 PACK          | \$23.97<br>FGSURG162AB5 |



LINDEMANN

| US NO.          | 161                    |
|-----------------|------------------------|
| H/LENGTH (MM)   | 10.0                   |
| FG SURG. (25MM) | 016                    |
| 5 PACK          | \$23.97<br>FGSURG161B5 |

| 6 Pack | \$7.88 | 28 mm  | 32 mm  |
|--------|--------|--------|--------|
| #1     |        | EG0128 | EG0132 |
| #2     |        | EG0228 | EG0232 |
| #3     |        | EG0328 | EG0332 |
| #4     |        | EG0428 | EG0432 |
| #5     |        | EG0528 | EG0532 |
| #6     |        | EG0628 | EG0632 |



| ASSORTED PACK |          | \$7.88 |
|---------------|----------|--------|
| 28MM          | EGASST28 |        |
| 32MM          | EGASST32 |        |



| 6 Pack <span style="float: right;">\$7.21</span> |                                                                                          |          |          |
|--------------------------------------------------|------------------------------------------------------------------------------------------|----------|----------|
| Size                                             | Color                                                                                    | 21mm     | 25mm     |
| 0                                                |  PURPLE | EB4XF021 | EB4XF025 |
| 1                                                |  WHITE  | EB3XF121 | EB3XF125 |
| 2                                                |  YELLOW | EB2XF221 | EB2XF225 |
| 3                                                |  RED    | EB1XF321 | EB1XF325 |
| 4                                                |  BLUE   | EB0XF421 | EB0XF425 |
| 5                                                |  GREEN  | EBM0521  | EBM0525  |
| 6                                                |  BLACK  | EBC0621  | EBC0625  |

| Assorted Pack 10 Pack <span style="float: right;">\$7.21</span>                |            |            |
|--------------------------------------------------------------------------------|------------|------------|
| Size                                                                           | 21mm       | 25mm       |
| 0-6                                                                            | EBASST0621 | EBASST0625 |
| <b>INCLUDED:</b> #0 (X1), #1 (X1), #2 (X2), #3 (X2), #4 (X2), #5 (X1), #6 (X1) |            |            |
| 2-4                                                                            | EBASST2421 | EBASST2425 |
| <b>INCLUDED:</b> #2 (X3), #4 (X4), #5 (X3)                                     |            |            |

# EDGEHEDSTROMS™

Stainless Steel



| 6 Pack |                                 |         |         | \$5.39  |
|--------|---------------------------------|---------|---------|---------|
| Size   | Color                           | 21mm    | 25mm    | 31mm    |
| 15     | <input type="checkbox"/> WHITE  | HED1521 | HED1525 | HED1531 |
| 20     | <input type="checkbox"/> YELLOW | HED2021 | HED2025 | HED2031 |
| 25     | <input type="checkbox"/> RED    | HED2521 | HED2525 | HED2531 |
| 30     | <input type="checkbox"/> BLUE   | HED3021 | HED3025 | HED3031 |
| 35     | <input type="checkbox"/> GREEN  | HED3521 | HED3525 | HED3531 |
| 40     | <input type="checkbox"/> BLACK  | HED4021 | HED4025 | HED4031 |
| 45     | <input type="checkbox"/> WHITE  | HED4521 | HED4525 | HED4531 |
| 50     | <input type="checkbox"/> YELLOW | HED5021 | HED5025 | HED5031 |

## EDGETAPER RETREAT™

Rotary Retreatment Files



EdgeTaper Retreat™ files are designed for the different needs of unfilling the coronal third, the middle third and the apical third before canal reshaping.

| 16 mm 6 Pack |                                |         | \$19.99 |
|--------------|--------------------------------|---------|---------|
| D1           | <input type="checkbox"/> WHITE | ETRD116 |         |

| 18 mm 6 Pack |                                |         | \$19.99 |
|--------------|--------------------------------|---------|---------|
| D2           | <input type="checkbox"/> WHITE | ETRD218 |         |

| 22 mm 6 Pack |                                |         | \$19.99 |
|--------------|--------------------------------|---------|---------|
| D3           | <input type="checkbox"/> WHITE | ETRD322 |         |

| Assorted 6 Pack |                                |             | \$19.99 |
|-----------------|--------------------------------|-------------|---------|
| D1-D3           | <input type="checkbox"/> WHITE | ETRASSTD123 |         |

# EDGE U UTOPIA™ BC SOLVENT

**EdgeUtopia™ BC Solvent is not just another solution,  
it's *THE* solution for dissolving Bioceramic Sealers**

## **Dissolves Bioceramic Sealers and Calcium Hydroxide.**


EdgeUtopia™ BC Solvent is engineered to dissolve the toughest bioceramic sealers and calcium hydroxide, outperforming other traditional solutions on the market. No more struggle with difficult to remove materials.

## **Tough on Sealer. Gentle on Teeth.**

While other products fall short on dissolving tough materials without damaging sensitive structures, EdgeUtopia BC Solvent works precisely where it's needed most, while maintaining the integrity of the surrounding tooth tissues.

## **EdgeUtopia™ BC Solvent - Preserving What Matters.**

|                |                |
|----------------|----------------|
| <b>\$68.06</b> |                |
| 240mL Bottle   | 10818207024937 |



Rev A 03/2025

# EDGE FIND™

## Non Heat-Treated Endonol™ NiTi

| 6 Pack |                                                                                          | \$19.99  |          |          |
|--------|------------------------------------------------------------------------------------------|----------|----------|----------|
| Size   | Color                                                                                    | 21mm     | 25mm     | 31mm     |
| P1     |  PURPLE | EFNDP121 | EFNDP131 | EFNDP125 |
| P2     |  WHITE  | EFNDP221 | EFNDP231 | EFNDP225 |
| P3     |  YELLOW | EFNDP321 | EFNDP331 | EFNDP325 |



| Assorted 3 Pack | \$34.52 | 21mm       | 25mm       | 31mm       |
|-----------------|---------|------------|------------|------------|
| P1, P2, P3      |         | EFNDASST21 | EFNDASST25 | EFNDASST31 |

# EDGE GLIDE PATH™

## Heat-Treated FireWire™ NiTi

Rotary GlidePath



| 3 pack |                                                                                           | \$12.88 |
|--------|-------------------------------------------------------------------------------------------|---------|
| Size   | Color                                                                                     |         |
| 21mm   |  WHITE | EGP0221 |
| 25mm   |  WHITE | EGP0225 |
| 31mm   |  WHITE | EGP0231 |

The EdgeGlidePath™ heat treated FireWire NiTi file features the flexibility and strength that EdgeEndo® is known for, replacing as many as 8 instruments to create a glidepath faster, with far less effort and with a 1mm maximum flute diameter. Replaces #15, 20, 25, 30, 35 hand files in a step-back technique and a #1, 2 & 3 Gates Glidden. Used exactly like Tulsa Dental ProGlider®\* at substantially lower cost.

\*ProGlider® is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGEONE GLIDEPATH™

For EdgeOne Fire™



Reciprocating GlidePath

| 3 pack |                                | \$12.88 |
|--------|--------------------------------|---------|
| Size   | Color                          |         |
| 21mm   | <input type="checkbox"/> WHITE | EOFGP21 |
| 25mm   | <input type="checkbox"/> WHITE | EOFGP25 |
| 31mm   | <input type="checkbox"/> WHITE | EOFGP31 |

# EDGEONE GLIDEPATH™

For EdgeOne Platinum™



Reciprocating GlidePath

| 3 pack |                                | \$12.88 |
|--------|--------------------------------|---------|
| Size   | Color                          |         |
| 21mm   | <input type="checkbox"/> WHITE | EOGP21  |
| 25mm   | <input type="checkbox"/> WHITE | EOGP25  |
| 31mm   | <input type="checkbox"/> WHITE | EOGP31  |

# EDGEONE BLAZE UTOPIA™

R-Glide



Reciprocating GlidePath

| 3 pack |                                | \$12.88      |
|--------|--------------------------------|--------------|
| Size   | Color                          |              |
| 21mm   | <input type="checkbox"/> WHITE | EOBURG1421SU |
| 25mm   | <input type="checkbox"/> WHITE | EOBURG1425SU |
| 31mm   | <input type="checkbox"/> WHITE | EOBURG1431SU |

\*ProGlider® is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGE SURGICAL TIP™

\$32.88

EdgeSurgical Tips™ are designed for use in endodontic microsurgery. Made with a diamond coating, these instruments cut dentin faster without need for repeated ultrasonics. The dentin surface is left smooth yet microscopically rough which results in better adaptation of filling materials.

- Diamond coating for efficient cutting of dentin.
- Irrigation strategically placed at the tip for control.
- Tip angles and lengths designed for access to difficult areas.
- We don't use zirconium nitride.

## EdgeSurgical Tip #1

Diamond-coated universal tip 30µm, length 3mm, taper 9%

AS3D1



## EdgeSurgical Tip #2

Diamond-coated universal tip 30µm, taper 9%

AS3D2



## EdgeSurgical Tip #3

Left-oriented retro surgery tip, diamond-coated 30µm, length 5mm, taper 7%

P15LD



## EdgeSurgical Tip #4

Left-oriented tip, diamond-coated 30µm, length 3mm, taper 10%

ASLD



## EdgeSurgical Tip #5

Right-oriented tip, diamond-coated 30µm, length 6mm, taper 7%

P15RD



## EdgeSurgical Tip #6

Right-oriented tip, diamond-coated 30µm, length 3mm, taper 10%

ASRD



## DIAMONDS



**001 | 801**

**Round**

| <b>Head Diameter 1/10mm</b> | 010                       | 012                       | 014                       | 016                       | 018                        | 021                       | 023                        | 025                       |
|-----------------------------|---------------------------|---------------------------|---------------------------|---------------------------|----------------------------|---------------------------|----------------------------|---------------------------|
| <b>L1: Head Length L1mm</b> | 1.0                       | 1.2                       | 1.4                       | 1.6                       | 1.8                        | 2.1                       | 2.3                        | 2.5                       |
| <b>L2: OAL mm</b>           | 19.0                      | 19.0                      | 19.0                      | 19.0                      | 19.0                       | 19.0                      | 19.0                       | 19.0                      |
| <b>ISO</b>                  |                           |                           |                           |                           |                            |                           |                            |                           |
| <b>806 314 001 514...</b>   |                           |                           |                           |                           |                            |                           |                            | 025<br>\$33.88<br>001025F |
| <b>806 314 001 524...</b>   | 010<br>\$33.88<br>001010M | 012<br>\$33.88<br>001012M | 014<br>\$33.88<br>001014M | 016<br>\$33.88<br>001016M | 018<br>\$33.88<br>001018M  |                           | 023<br>\$33.88<br>001023M  |                           |
| <b>806 314 001 534...</b>   | 010<br>\$33.88<br>001010C | 012<br>\$33.88<br>001012C | 014<br>\$33.88<br>001014C | 016<br>\$33.88<br>001016C | 018<br>\$33.88<br>001018C  | 021<br>\$33.88<br>001021C | 023<br>\$33.88<br>001023C  |                           |
| <b>806 314 001 544...</b>   |                           |                           |                           |                           | 018<br>\$33.88<br>001018XC |                           | 023<br>\$33.88<br>001023XC |                           |



**032 | 813**

**Diabolo**

| <b>Head Diameter 1/10mm</b> | 014                       |
|-----------------------------|---------------------------|
| <b>L1: Head Length L1mm</b> | 1.5                       |
| <b>L2: OAL mm</b>           | 19.0                      |
| <b>ISO</b>                  |                           |
| <b>806 314 032 534...</b>   | 014<br>\$33.88<br>032014C |



**544 | 845KR Short**

**Round Edge Taper**

| <b>Head Diameter 1/10mm</b> | 019                        | 026                        |
|-----------------------------|----------------------------|----------------------------|
| <b>L1: Head Length L1mm</b> | 4.0                        | 4.0                        |
| <b>L2: OAL mm</b>           | 18.0                       | 18.0                       |
| <b>ISO</b>                  | Short                      | Short                      |
| <b>806 313 544 514...</b>   |                            | 026<br>\$33.88<br>544026SF |
| <b>806 313 544 534...</b>   | 019<br>\$33.88<br>544019SC |                            |



**546 | 847KR**

**Round Edge Taper**



|                      |                             |
|----------------------|-----------------------------|
| Head Diameter 1/10mm | 017                         |
| L1: Head Length L1mm | 8.0                         |
| L2: OAL mm           | 20.0                        |
| ISO                  | Short                       |
| 806 314 546 514...   | 017<br>\$33.88<br>546017SF  |
| 806 314 546 524...   | 017<br>\$33.88<br>546017SM  |
| 806 314 546 534...   | 017<br>\$33.88<br>546017SC  |
| 806 314 546 544...   | 017<br>\$33.88<br>546017SXC |



**546 | 847KR**

**Round Edge Taper**



|                      |                           |
|----------------------|---------------------------|
| Head Diameter 1/10mm | 016                       |
| L1: Head Length L1mm | 8.0                       |
| L2: OAL mm           | 22.0                      |
| ISO                  |                           |
| 806 313 546 534...   | 016<br>\$33.88<br>546016C |



**197 | 855 Short**

**Round End Taper**



|                      |                             |                             |
|----------------------|-----------------------------|-----------------------------|
| Head Diameter 1/10mm | 017                         | 024                         |
| L1: Head Length L1mm | 6.0                         | 7.0                         |
| L2: OAL mm           | 18.5                        | 18.5                        |
| ISO                  | Short                       | Short                       |
| 806 314 197 544...   | 017<br>\$33.88<br>197017SXC | 024<br>\$33.88<br>197024SXC |



**197 | 855 XShort**

**Round End Taper**



|                      |                             |
|----------------------|-----------------------------|
| Head Diameter 1/10mm | 016                         |
| L1: Head Length L1mm | 6.0                         |
| L2: OAL mm           | 16.5                        |
| ISO                  | XShort                      |
| 806 313 197 534...   | 016<br>\$33.88<br>197016SHC |



**198** | 856

**Round End Taper**

|                      |                             |                             |                             |                           |                             |                             |                             |
|----------------------|--------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|
| Head Diameter 1/10mm | 012                                                                                                          | 014                                                                                                          | 015                                                                                                          | 016                                                                                                          | 018                                                                                                            | 019                                                                                                            | 021                                                                                                            |
| L1: Head Length L1mm | 8.0                                                                                                          | 8.0                                                                                                          | 8.0                                                                                                          | 8.0                                                                                                          | 8.0                                                                                                            | 8.0                                                                                                            | 8.0                                                                                                            |
| L2: OAL mm           | 22.0                                                                                                         | 22.0                                                                                                         | 22.0                                                                                                         | 22.0                                                                                                         | 22.0                                                                                                           | 22.0                                                                                                           | 22.0                                                                                                           |
| ISO                  |                                                                                                              |                                                                                                              |                                                                                                              |                                                                                                              |                                                                                                                |                                                                                                                |                                                                                                                |
| 806 316 198 514...   |                                                                                                              |  014<br>\$33.88<br>198014F  |                                                                                                              |  016<br>\$33.88<br>198016F  |  018<br>\$33.88<br>198018F  |                                                                                                                |                                                                                                                |
| 806 316 198 524...   |  012<br>\$23.97<br>198012M  |  014<br>\$33.88<br>198014M  |                                                                                                              |  016<br>\$23.97<br>198016M  |  018<br>\$23.97<br>198018M  |                                                                                                                |                                                                                                                |
| 806 316 198 534...   |  012<br>\$33.88<br>198012C  |  014<br>\$33.88<br>198014C  |  015<br>\$33.88<br>198015C  |  016<br>\$33.88<br>198016C  |  018<br>\$33.88<br>198018C  |                                                                                                                |  021<br>\$33.88<br>198021C  |
| 806 316 198 544...   |  012<br>\$33.88<br>198012XC |  014<br>\$33.88<br>198014XC |  015<br>\$33.88<br>198015XC |  016<br>\$33.88<br>198016XC |  018<br>\$33.88<br>198018XC |  019<br>\$23.97<br>198019XC |  021<br>\$33.88<br>198021XC |



**198** | 856

**Round End Taper**

|                      |                               |                               |
|----------------------|------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|
| Head Diameter 1/10mm | 022                                                                                                              | 025                                                                                                              |
| L1: Head Length L1mm | 8.0                                                                                                              | 8.0                                                                                                              |
| L2: OAL mm           | 22.0                                                                                                             | 22.0                                                                                                             |
| ISO                  |                                                                                                                  |                                                                                                                  |
| 806 315 198 534...   |  022<br>\$33.88<br>(468-0698) |  025<br>\$33.88<br>(468-0700) |
| 806 315 198 544...   |  022<br>\$33.88<br>(468-0699) |  025<br>\$33.88<br>(468-0701) |



**198** | 856 Short

**Round End Taper**

|                      |                             |
|----------------------|------------------------------------------------------------------------------------------------------------------|
| Head Diameter 1/10mm | 023                                                                                                              |
| L1: Head Length L1mm | 8.5                                                                                                              |
| L2: OAL mm           | 20.0                                                                                                             |
| ISO                  | Short                                                                                                            |
| 806 314 198 534...   |  023<br>\$23.97<br>198023SC |



## 198 | 856L Long Round End Taper

|                      |                               |                               |                              |
|----------------------|---------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|
| Head Diameter 1/10mm | 016                                                                                                           | 018                                                                                                           | 020                                                                                                          |
| L1: Head Length L1mm | 9.0                                                                                                           | 9.0                                                                                                           | 9.0                                                                                                          |
| L2: OAL mm           | 22.0                                                                                                          | 22.0                                                                                                          | 22.0                                                                                                         |
| ISO                  |                                                                                                               |                                                                                                               |                                                                                                              |
| 806 315 198L 514...  |  016<br>\$33.88<br>198L016F  |                                                                                                               |                                                                                                              |
| 806 315 198L 534...  |  016<br>\$33.88<br>198L016C  |  018<br>\$33.88<br>198L018C  |  020<br>\$33.88<br>198L020C |
| 806 315 198L 544...  |  016<br>\$33.88<br>198L016XC |  018<br>\$33.88<br>198L018XC |                                                                                                              |



## 199 | 850 Round End Taper

|                      |                              |                              |                              |                               |                               |
|----------------------|---------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|
| Head Diameter 1/10mm | 010                                                                                                           | 012                                                                                                           | 014                                                                                                           | 016                                                                                                              | 018                                                                                                              |
| L1: Head Length L1mm | 10.0                                                                                                          | 10.0                                                                                                          | 10.0                                                                                                          | 10.0                                                                                                             | 10.0                                                                                                             |
| L2: OAL mm           | 24.0                                                                                                          | 24.0                                                                                                          | 24.0                                                                                                          | 24.0                                                                                                             | 24.0                                                                                                             |
| ISO                  |                                                                                                               |                                                                                                               |                                                                                                               |                                                                                                                  |                                                                                                                  |
| 806 316 199 524...   |                                                                                                               |                                                                                                               |                                                                                                               |                                                                                                                  |  018<br>\$33.88<br>199018M  |
| 806 316 199 534...   |  010<br>\$33.88<br>199010C |  012<br>\$33.88<br>199012C |  014<br>\$33.88<br>199014C |  016<br>\$33.88<br>199016C  | 018<br>\$25.01<br>199018C                                                                                        |
| 806 316 199 544...   |                                                                                                               |                                                                                                               |                                                                                                               |  016<br>\$33.88<br>199016XC |  018<br>\$33.88<br>199018XC |



## 199 | 850 Short Round End Taper

|                      |                              |                              |
|----------------------|-----------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
| Head Diameter 1/10mm | 014                                                                                                             | 016                                                                                                             |
| L1: Head Length L1mm | 10.0                                                                                                            | 10.0                                                                                                            |
| L2: OAL mm           | 22.0                                                                                                            | 22.0                                                                                                            |
| ISO                  | Short                                                                                                           | Short                                                                                                           |
| 806 315 199 514...   |  014<br>\$33.88<br>199014SXC |  016<br>\$33.88<br>199016SXC |



## 263 | 369 Bud Slender

|                      |                            |
|----------------------|-----------------------------------------------------------------------------------------------------------------|
| Head Diameter 1/10mm | 025                                                                                                             |
| L1: Head Length L1mm | 5.5                                                                                                             |
| L2: OAL mm           | 20.0                                                                                                            |
| ISO                  |                                                                                                                 |
| 806 314 263 534...   |  025<br>\$33.88<br>263025C |

# EDGEK-FILE™ NITI



Non Heat-Treated NiTi

Hand File

Alternative to Tulsa Dental K-File®

| 6 Pack |                                                                                           |            |            | \$8.25     |
|--------|-------------------------------------------------------------------------------------------|------------|------------|------------|
| Size   | Color                                                                                     | 21mm       | 25mm       | 31mm       |
| 15     |  WHITE   | NITIKF1521 | NITIKF1525 | NITIKF1531 |
| 20     |  YELLOW  | NITIKF2021 | NITIKF2025 | NITIKF2031 |
| 25     |  RED     | NITIKF2521 | NITIKF2525 | NITIKF2531 |
| 30     |  BLUE    | NITIKF3021 | NITIKF3025 | NITIKF3031 |
| 35     |  GREEN   | NITIKF3521 | NITIKF3525 | NITIKF3531 |
| 40     |  BLACK   | NITIKF4021 | NITIKF4025 | NITIKF4031 |
| 45     |  WHITE   | NITIKF4521 | NITIKF4525 | NITIKF4531 |
| 50     |  YELLOW | NITIKF5021 | NITIKF5025 | NITIKF5031 |
| 55     |  RED   | NITIKF5521 | NITIKF5525 | NITIKF5531 |
| 60     |  BLUE  | NITIKF6021 | NITIKF6025 | NITIKF6031 |

Tulsa Dental K-File® is a registered trademark of Dentsply Tulsa Dental Specialties



| <b>6 Pack</b> |                                                                                            |             |             |             | <b>\$5.39</b> |
|---------------|--------------------------------------------------------------------------------------------|-------------|-------------|-------------|---------------|
| <b>Size</b>   | <b>Color</b>                                                                               | <b>21mm</b> | <b>25mm</b> | <b>31mm</b> |               |
| 6             |  PINK     | KF621       | KF625       | KF631       |               |
| 8             |  GREY     | KF821       | KF825       | KF831       |               |
| 10            |  PURPLE   | KF1021      | KF1025      | KF1031      |               |
| 15            |  WHITE    | KF1521      | KF1525      | KF1531      |               |
| 20            |  YELLOW   | KF2021      | KF2025      | KF2031      |               |
| 25            |  RED      | KF2521      | KF2525      | KF2531      |               |
| 30            |  BLUE     | KF3021      | KF3025      | KF3031      |               |
| 35            |  GREEN    | KF3521      | KF3525      | KF3531      |               |
| 40            |  BLACK  | KF4021      | KF4025      | KF4031      |               |
| 45            |  WHITE  | KF4521      | KF4525      | KF4531      |               |
| 50            |  YELLOW | KF5021      | KF5025      | KF5031      |               |
| 55            |  RED    | KF5521      | KF5525      | KF5531      |               |
| 60            |  BLUE   | KF6021      | KF6025      | KF6031      |               |
| 70            |  GREEN  | KF7021      | KF7025      | KF7031      |               |
| 80            |  BLACK  | KF8021      | KF8025      | KF8031      |               |

# EDGE C-FILE™



**Stainless Steel**

Hand File

| 6 Pack |                                                                                          | \$6.23 |        |
|--------|------------------------------------------------------------------------------------------|--------|--------|
| Size   | Color                                                                                    | 21mm   | 25mm   |
| 6      |  PINK   | CF621  | CF625  |
| 8      |  GREY   | CF821  | CF825  |
| 10     |  PURPLE | CF1021 | CF1025 |
| 12     |  BROWN  | CF1221 | CF1225 |
| 15     |  WHITE  | CF1521 | CF1525 |



# EDGE X7 UTOPIA™



**\$19.99**  
6 Pack

## Enhanced Technology. Remarkably Beautiful.

- Enhanced version of EdgeFile® X7
- Features proprietary FireWire Blaze Heat Treatment for better cutting efficiency
- Constant Taper File
- Pre-Sterilized



- Parabolic Shapped Cross-Section
- Available in 3 lengths: 21, 25 & 29 mm
- File Size: 04 & 06
- Maximum flute diameter 1 mm
- Speed: 500 RPM
- Torque: 4.0 Ncm



High Cutting  
Efficiency



No Bounce  
Back



High  
Strength

# EDGE X7



## UTOPIA™



| Pre-Sterilized 6 Pack |        |             |             |             | Points       |              | Obturators |          |
|-----------------------|--------|-------------|-------------|-------------|--------------|--------------|------------|----------|
| Size                  | Color  | 21mm        | 25mm        | 29mm        | Paper Points | Gutta Percha | EdgeCore   | EdgeFill |
| 17.04                 | PURPLE | X7U170421SE | X7U170425SE | X7U170429SE | SPP1504-1    | X71704GP-1   | N/A        | N/A      |
| 20.04                 | YELLOW | X7U170621SE | X7U170625SE | X7U170629SE | SPP2004-1    | X72004GP-1   | X72006C    | X72006O  |
| 25.04                 | RED    | X7U200421SE | X7U200425SE | X7U200429SE | SPP2504-1    | X72504GP-1   | X72506C    | X72506O  |
| 30.04                 | BLUE   | X7U200621SE | X7U200625SE | X7U200629SE | SPP3004-1    | X73004GP-1   | X73006C    | X73006O  |
| 35.04                 | GREEN  | X7U250421SE | X7U250425SE | X7U250429SE | SPP3504-1    | X73504GP-1   | X73506C    | X73506O  |
| 40.04                 | BLACK  | X7U250621SE | X7U250625SE | X7U250629SE | SPP4004-1    | X74004GP-1   | X74006C    | X74006O  |
| 45.04                 | WHITE  | X7U300421SE | X7U300425SE | X7U300429SE | SPP4504-1    | X74504GP-1   | X74506C    | X74506O  |
| 17.06                 | PURPLE | X7U300621SE | X7U300625SE | X7U300629SE | SPP1506-1    | X71706GP-1   | N/A        | N/A      |
| 20.06                 | YELLOW | X7U350421SE | X7U350425SE | X7U350429SE | SPP2006-1    | X72006GP-1   | X72006C    | X72006O  |
| 25.06                 | RED    | X7U350621SE | X7U350625SE | X7U350629SE | SPP2506-1    | X72506GP-1   | X72506C    | X72506O  |
| 30.06                 | BLUE   | X7U400421SE | X7U400425SE | X7U400429SE | SPP3006-1    | X73006GP-1   | X73006C    | X73006O  |
| 35.06                 | GREEN  | X7U400621SE | X7U400625SE | X7U400629SE | SPP3506-1    | X73506GP-1   | X73506C    | X73506O  |
| 40.06                 | BLACK  | X7U450421SE | X7U450425SE | X7U450429SE | SPP4006-1    | X74006GP-1   | X74006C    | X74006O  |
| 45.06                 | WHITE  | X7U450621SE | X7U450625SE | X7U450629SE | SPP4506-1    | X74506GP-1   | X74506C    | X74506O  |

Pair EdgeDiamond Utopia™ with:

# EDGE UTOPIA™

## BioCeramic Sealer 2g

818365020072



### Assorted Packs .04 6 Pack

(17, 20, 25, 25, 30, 35)



| Length | Files         | Paper Points | Gutta Percha |
|--------|---------------|--------------|--------------|
| 21MM   | X7UASST0421SE | SPP04AST-1   | X7GP04AST-1  |
| 25MM   | X7UASST0425SE |              |              |



# EDGEFILE® X7



**\$19.99**  
6 Pack

## Unmatched Strength and Flexibility

- Universal File System
- #1 Fan Favorite
- FireWire NiTi Heat Treatment provides super elasticity with no shape memory
- Constant Tapered Rotary Files
- Non-Sterile
- Tighter pitch and flutes which deflect more torsional forces along the length of the file
- Non-cutting tip with flutes extending to 0.25mm—offering precision and safety in one design

Constant Taper .04 or .06 until wire diameter reaches 1.0mm

Non-cutting tip

FireWire NiTi Heat Treatment


Calibration Rings

ISO Tip Size Indicator Band

Short handle 11mm

Silicone Stopper:  
Red .04 Taper  
Green .06 Taper

Taper identification bands:  
2 stripes = .04 Taper  
3 stripes = .06 Taper




- Parabolic Shaped Cross-Section
- Available in 3 lengths: 21, 25 & 29 mm
- Constant Taper: .04 and .06
- Rotary files
- Speed: 450 RPM
- Torque: 2.9 NCm



Capable of  
90° Curves



No Bounce  
Back



High  
Strength

# EDGEFILE® X7



| 6 Pack |        |            |            |            | Points       |              | Obturators |          |
|--------|--------|------------|------------|------------|--------------|--------------|------------|----------|
| Size   | Color  | 21mm       | 25mm       | 29mm       | Gutta Percha | Paper Points | EdgeCore   | EdgeFill |
| 17.04  | WHITE  | X7170421RF | X7170425RF | X7170429RF | SPP1504-1    | X71704GP-1   | N/A        | N/A      |
| 20.04  | YELLOW | X7170621RF | X7170625RF | X7170629RF | SPP2004-1    | X72004GP-1   | X72006C    | X72006O  |
| 25.04  | RED    | X7200421RF | X7200425RF | X7200429RF | SPP2504-1    | X72504GP-1   | X72506C    | X72506O  |
| 30.04  | BLUE   | X7200621RF | X7200625RF | X7200629RF | SPP3004-1    | X73004GP-1   | X73006C    | X73006O  |
| 35.04  | GREEN  | X7250421RF | X7250425RF | X7250429RF | SPP3504-1    | X73504GP-1   | X73506C    | X73506O  |
| 40.04  | BLACK  | X7250621RF | X7250625RF | X7250629RF | SPP4004-1    | X74004GP-1   | X74006C    | X74006O  |
| 45.04  | WHITE  | X7300421RF | X7300425RF | X7300429RF | SPP4504-1    | X74504GP-1   | X74506C    | X74506O  |
| 17.06  | WHITE  | X7300621RF | X7300625RF | X7300629RF | SPP1506-1    | X71706GP-1   | N/A        | N/A      |
| 20.06  | YELLOW | X7350421RF | X7350425RF | X7350429RF | SPP2006-1    | X72006GP-1   | X72006C    | X72006O  |
| 25.06  | RED    | X7350621RF | X7350625RF | X7350629RF | SPP2506-1    | X72506GP-1   | X72506C    | X72506O  |
| 30.06  | BLUE   | X7400421RF | X7400425RF | X7400429RF | SPP3006-1    | X73006GP-1   | X73006C    | X73006O  |
| 35.06  | GREEN  | X7400621RF | X7400625RF | X7400629RF | SPP3506-1    | X73506GP-1   | X73506C    | X73506O  |
| 40.06  | BLACK  | X7450421RF | X7450425RF | X7450429RF | SPP4006-1    | X74006GP-1   | X74006C    | X74006O  |
| 45.06  | WHITE  | X7450621RF | X7450625RF | X7450629RF | SPP4506-1    | X74506GP-1   | X74506C    | X74506O  |

Pair EdgeFile® X7 with:

**EDGEBIOCERAMIC™**  
SEALER

818365020072



| Assorted Pack .04 6 Pack<br>(17, 20, 25, 30, 35, 40) |            |              |              |
|------------------------------------------------------|------------|--------------|--------------|
| Length                                               | Files      | Paper Points | Gutta Percha |
| 21mm                                                 | ASSTX70421 | SPP04AST-1   | X7GP04AST-1  |
| 25mm                                                 | ASSTX70425 |              |              |
| 29mm                                                 | ASSTX70429 |              |              |

| Assorted Pack .06 6 Pack<br>(17, 20, 25, 30, 35, 40) |            |              |              |
|------------------------------------------------------|------------|--------------|--------------|
| Length                                               | Files      | Paper Points | Gutta Percha |
| 21mm                                                 | ASSTX70621 | SPP06AST-1   | X7GP06AST-1  |
| 25mm                                                 | ASSTX70625 |              |              |
| 29mm                                                 | ASSTX70629 |              |              |



# EDGEFILE® X3



**\$19.99**  
6 Pack

## Engineered for Strength with Flexible Precision

- Designed with a Safe-Core™ Parabolic Cross- Section
- .04 or .06 constant taper file until flute diameter reaches 1.0 then taper is 0
- Safe unwinding feature
- Super elasticity. No Shape Memory
- FireWire NiTi Heat Treated
- Variable Tapered Rotary Files
- Features its own dedicated Glidepath system
- Non-Sterile
- Alternative to ProTaper® Next<sup>1</sup> systems



- Parabolic shaped cross-section
- Available in 3 lengths: 21, 25 & 29 mm
- Speed 300-500 RPM
- Torque 2.94 Ncm



Highly  
Flexible



No Bounce  
Back



High  
Strength

1) ProTaper® Next is a registered trademark of Dentsply Sirona.

# EDGEFILE® X3



| 6 Pack     |        |            |            |            | Points       |              | Obturators |          |
|------------|--------|------------|------------|------------|--------------|--------------|------------|----------|
| Size/Taper | Color  | 21mm       | 25mm       | 29mm       | Paper Points | Gutta Points | EdgeCore   | EdgeFill |
| N1 17.06   | PURPLE | X3170621RF | X3170625RF | X3170629RF | SPP1506-1    | X71706GP-1   | N/A        | N/A      |
| N2 17.04   | WHITE  | X3170421RF | X3170425RF | X3170429RF | SPP1504-1    | X71704GP-1   | N/A        | N/A      |
| C1 20.06   | YELLOW | X3200621RF | X3200625RF | X3200629RF | SPP2006-1    | X72006GP-1   | X72006C    | X72006O  |
| C2 25.06   | RED    | X3250621RF | X3250625RF | X3250629RF | SPP2506-1    | X72506GP-1   | X72506C    | X72506O  |
| C3 30.06   | BLUE   | X3300621RF | X3300625RF | X3300629RF | SPP3006-1    | X73006GP-1   | X73006C    | X73006O  |
| C4 40.06   | BLACK  | X3400621RF | X3400625RF | X3400629RF | SPP4006-1    | X74006GP-1   | X74006C    | X74006O  |

Pair EdgeFile® X3 with:

**EDGE BIO CERAMIC™**  
SEALER  
818365020072



| EdgeVerify™ for<br>EdgeFile® X3 6 Pack |           |
|----------------------------------------|-----------|
| Size                                   | Item Code |
| 20                                     | X32006V   |
| 25                                     | X32506V   |
| 30                                     | X33006V   |
| 40                                     | X34006V   |
| 50                                     | X35006V   |

# EDGESEQUEL



# UTOPIA™



**\$19.99**  
6 Pack

## Enhanced Technology. Remarkably Beautiful.

- Features proprietary FireWire Blaze Blue Heat Treatment
- Improved cutting efficiency compared to EdgeSequel Sapphire™
- 12mm handle length for better posterior access
- Constant Taper File
- Pre-Sterilized
- Alternative to EndoSequence® and Vortex®




\*19mm orifice opener available in #25/.08

Vortex® is a registered trademark of Dentsply Tulsa Dental Specialties. EndoSequence® is a registered trademark of Brasseler USA

# EDGESEQUEL



## UTOPIA™



| 19mm Pre-Sterilized 6 Pack |        |            |            |            |              |              |            |          |
|----------------------------|--------|------------|------------|------------|--------------|--------------|------------|----------|
| Size/Taper                 | Color  |            |            |            |              |              |            |          |
| 25.08                      | YELLOW | ES250819SE |            |            |              |              |            |          |
| Pre-Sterilized 6 Pack      |        |            |            |            | Points       |              | Obturators |          |
| Size/Taper                 | Color  | 21mm       | 25mm       | 31mm       | Paper Points | Gutta Percha | EdgeCore   | EdgeFill |
| 10.02                      | PURPLE | ES100221SE | ES100225SE | N/A        | N/A          | N/A          | N/A        | N/A      |
| 15.02                      | WHITE  | ES150221SE | ES150225SE | N/A        | N/A          | N/A          | N/A        | N/A      |
| 20.02                      | YELLOW | ES200221SE | ES200225SE | N/A        | N/A          | N/A          | N/A        | N/A      |
| 30.02                      | BLUE   | ES300221SE | ES300425SE | N/A        | N/A          | N/A          | N/A        | N/A      |
| 10.04                      | PURPLE | ES100421SE | ES100425SE | N/A        | N/A          | N/A          | N/A        | N/A      |
| 15.04                      | WHITE  | ES150421SE | ES150425SE | ES150431SE | SPP1504-1    | SGP1504-1    | N/A        | N/A      |
| 20.04                      | YELLOW | ES200421SE | ES200421SE | N/A        | SPP2004-1    | SGP2004-1    | ES2006C    | ES2006O  |
| 25.04                      | RED    | ES250421SE | ES250425SE | ES250431SE | SPP2504-1    | SGP2504-1    | ES2506C    | ES2506O  |
| 30.04                      | BLUE   | ES300421SE | ES300425SE | ES300431SE | SPP3004-1    | SGP3004-1    | ES3006C    | ES3006O  |
| 35.04                      | GREEN  | ES350421SE | ES350425SE | N/A        | SPP3504-1    | SGP3504-1    | ES3506C    | ES3506O  |
| 40.04                      | BLACK  | ES400421SE | ES400425SE | ES400431SE | SPP4004-1    | SGP4004-1    | ES4006C    | ES4006O  |
| 50.04                      | YELLOW | ES500421SE | ES500425SE | ES500431SE | SPP5004-1    | SGP5004-1    | ES5006C    | ES5006O  |
| 15.06                      | WHITE  | ES150621SE | ES150625SE | N/A        | SPP1506-1    | SGP1506-1    | N/A        | N/A      |
| 25.06                      | RED    | ES250621SE | ES250625SE | ES250631SE | SPP2506-1    | SGP2506-1    | ES2506C    | ES2506O  |
| 30.06                      | BLUE   | ES300621SE | ES300625SE | ES300631SE | SPP3006-1    | SGP3006-1    | ES3006C    | ES3006O  |
| 35.06                      | GREEN  | ES350621SE | ES350625SE | ES350631SE | SPP3506-1    | SGP3506-1    | ES3506C    | ES3506O  |
| 40.06                      | BLACK  | ES400621SE | ES400625SE | N/A        | SPP4006-1    | SGP4006-1    | ES4006C    | ES4006O  |

| Assorted Pack 3 Pack<br>(15.04, 25.04, 30.04)    |                |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------|
| 21mm                                                                                                                                                                                                                                                                                                | ASST15300421SE |
| 25mm                                                                                                                                                                                                                                                                                                | ASST15300425SE |
| 31mm                                                                                                                                                                                                                                                                                                | ASST15300431SE |

| Assorted Pack 3 Pack<br>(15.04, 25.04, 25.06)    |                |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------|
| 21mm                                                                                                                                                                                                                                                                                                | ASST15250621SE |
| 25mm                                                                                                                                                                                                                                                                                                | ASST15250625SE |
| 31mm                                                                                                                                                                                                                                                                                                | ASST15250631SE |

| Assorted Pack 6 Pack<br>(15.04, 20.04, 25.04, 30.04, 35.04 & 25.08 19mm)       |              |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| 21mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | V2ASST0421SE |
| 25mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | V2ASST0425SE |
| 31mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | V2ASST0431SE |

Pair EdgeSequel Utopia™ with:

EDGE



UTOPIA™

BioCeramic Sealer 2g

818365020072




# EDGESEQUEL SAPPHIRE™



**\$12.88**  
**4 Pack**

## Bodybuilder Strength. Gymnast Flexibility.

- Excellent Cyclic Fatigue Compared To Other Leading Files\*
- FireWire NiTi Heat Treated
- Constant Tapered Files
- Non-Sterile
- Alternative to Vortex®, VortexBlue®, and EndoSequence®



## Cyclic Fatigue Testing - Average Time to Failure in Seconds

**EDGESEQUEL**  
SAPPHIRE™

633Seconds

Vortex Blue®

299.9 Seconds

EndoSequence®

80.9 Seconds

\*Vortex®, Vortex® Blue, are registered trademarks of Dentsply Tulsa Dental Specialties. EndoSequence® is a registered trademark of Brassler USA \*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph:<http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit [www.edgeendo.com/matter-of-temperature](http://www.edgeendo.com/matter-of-temperature)

# EDGESEQUEL<sup>™</sup>

## SAPPHIRE<sup>™</sup>



| 4 Pack |        |          |          |          | Points       |              | Obturers |          |
|--------|--------|----------|----------|----------|--------------|--------------|----------|----------|
| Size   | Color  | 21mm     | 25mm     | 31mm     | Paper Points | Gutta Points | EdgeCore | EdgeFill |
| 25.08  | RED    | ES250821 | ES250825 | ES250831 | SPP1504-1    | SGP1504-1    | N/A      | N/A      |
| 15.04  | WHITE  | ES150421 | ES150425 | ES150431 | SPP2004-1    | SGP2004-1    | ES2006C  | ES2006O  |
| 20.04  | YELLOW | ES150621 | ES150625 | ES150631 | SPP2504-1    | SGP2504-1    | ES2506C  | ES2506O  |
| 25.04  | RED    | ES200421 | ES200425 | ES200431 | SPP3004-1    | SGP3004-1    | ES3006C  | ES3006O  |
| 30.04  | BLUE   | ES200621 | ES200625 | ES200631 | SPP3504-1    | SGP3504-1    | ES3506C  | ES3506O  |
| 35.04  | GREEN  | ES250421 | ES250425 | ES250431 | SPP4004-1    | SGP4004-1    | ES4006C  | ES4006O  |
| 40.04  | BLACK  | ES250621 | ES250625 | ES250631 | SPP4504-1    | SGP4504-1    | ES4506C  | ES4506O  |
| 45.04  | WHITE  | ES300421 | ES300425 | ES300431 | SPP5004-1    | SGP5004-1    | ES5006C  | ES5006O  |
| 50.04  | YELLOW | ES300621 | ES300625 | ES300631 | ESS5504GP    | ESS5504PP    | ES5506C  | ES5506O  |
| 55.04  | RED    | ES350421 | ES350425 | ES350431 | ESS6004GP    | ESS6004PP    | ES6006C  | ES6006O  |
| 60.04  | BLUE   | ES350621 | ES350625 | ES350631 | ESS7004GP    | ESS7004PP    | ES7006C  | ES7006O  |
| 70.04  | GREEN  | ES400421 | ES400425 | ES400431 | ESS8004GP    | ESS8004PP    | ES8006C  | ES8006O  |
| 80.04  | BLACK  | ES400621 | ES400625 | ES400631 | N/A          | N/A          | ES9006C  | ES9006O  |
| 90.04  | WHITE  | ES450421 | ES450425 | ES450431 | SPP1506-1    | SGP1506-1    | N/A      | N/A      |
| 15.06  | WHITE  | ES450621 | ES450625 | ES450631 | SPP2006-1    | SGP2006-1    | ES2006C  | ES2006O  |
| 20.06  | YELLOW | ES500421 | ES500425 | ES500431 | SPP2506-1    | SGP2506-1    | ES2506C  | ES2506O  |
| 25.06  | RED    | ES500621 | ES500625 | ES500631 | SPP3006-1    | SGP3006-1    | ES3006C  | ES3006O  |
| 30.06  | BLUE   | ES550421 | ES550425 | ES550431 | SPP3506-1    | SGP3506-1    | ES3506C  | ES3506O  |
| 35.06  | GREEN  | ES600421 | ES600425 | ES600431 | SPP4006-1    | SGP4006-1    | ES4006C  | ES4006O  |
| 40.06  | BLACK  | ES700421 | ES700425 | ES700431 | SPP4506-1    | SGP4506-1    | ES4506C  | ES4506O  |
| 45.06  | WHITE  | ES800421 | ES800425 | ES800431 | SPP5006-1    | SGP5006-1    | ES5006C  | ES5006O  |
| 50.06  | YELLOW | ES900421 | ES900425 | ES900431 | SPP5006-1    | SGP5006-1    | ES5006C  | ES5006O  |

Pair EdgeSequel Sapphire<sup>™</sup> with:

**EDGE BIO CERAMIC<sup>™</sup>**  
SEALER

818365020072



| Assorted Packs .04 4 Pack |              |
|---------------------------|--------------|
| Sizes: (15, 20, 25, 30)   |              |
| 21mm                      | ASST15300421 |
| 25mm                      | ASST15300425 |
| 31mm                      | ASST15300431 |
| Sizes: (35, 40, 45, 50)   |              |
| 21mm                      | ASST35500421 |
| 25mm                      | ASST35500425 |
| 31mm                      | ASST35500431 |

| Assorted Packs .06 4 Pack |              |
|---------------------------|--------------|
| Sizes: (15, 20, 25, 30)   |              |
| 21mm                      | ASST15300621 |
| 25mm                      | ASST15300625 |
| 31mm                      | ASST15300631 |
| Sizes: (35, 40, 45, 50)   |              |
| 21mm                      | ASST35500621 |
| 25mm                      | ASST35500625 |
| 31mm                      | ASST35500631 |

# EDGE EVOLVE™




**\$18.99**  
**6 Pack**

## The Evolution of the Revolution

- Minimal dentin with a conservative flute and taper. All EdgeEvolve files, even our largest 40.10, only have a max diameter of 1mm. FireWire NiTi Heat Treated
- Extremely Flexible Constant Taper File with a non-cutting tip
- Non-Sterile
- Available in .04, .06, .08, .10 Tapers
- Alternative to GTX®<sup>1</sup>

Taper Identification Bands:  
2 stripes = .04 Taper  
3 stripes = .06 Taper  
4 stripes = .08 Taper  
5 stripes = .10 Taper



- Triangular shaped cross-section
- Available in 3 lengths: 21, 25 & 29mm
- Speed: 300-500 RPM
- Torque: 2.94 Ncm



Extremely  
Flexible




No Bounce  
Back



High  
Strength

1) GTX® is a registered trademark of Dentsply Tulsa Dental Specialties.

# EDGE EVOLVE™



| 6 Pack     |        |          |          |          | Points       |              | Obturators |          |
|------------|--------|----------|----------|----------|--------------|--------------|------------|----------|
| Size/Taper | Color  | 21mm     | 25mm     | 29mm     | Paper Points | Gutta Percha | EdgeCore   | EdgeFill |
| 17.04      | WHITE  | EE170421 | EE170425 | EE170429 | SPP1504-1    | SGP1504-1    | N/A        | N/A      |
| 17.06      | WHITE  | EE170621 | EE170625 | EE170629 | SPP1506-1    | SGP1506-1    | N/A        | N/A      |
| 17.08      | WHITE  | EE170821 | EE170825 | EE170829 | N/A          | N/A          | N/A        | N/A      |
| 17.10      | WHITE  | EE171021 | EE171025 | EE171029 | N/A          | N/A          | N/A        | N/A      |
| 20.04      | YELLOW | EE200421 | EE200425 | EE200429 | SPP2004-1    | SGP2004-1    | EE2006C    | EE2006O  |
| 20.06      | YELLOW | EE200621 | EE200625 | EE200629 | SPP2006-1    | SGP2006-1    | EE2006C    | EE2006O  |
| 20.08      | YELLOW | EE200821 | EE200825 | EE200829 | N/A          | N/A          | EE2006C    | EE2006O  |
| 20.10      | YELLOW | EE201021 | EE201025 | EE201029 | N/A          | N/A          | EE2006C    | EE2006O  |
| 25.04      | RED    | EE250421 | EE250425 | EE250429 | SPP2504-1    | SGP2504-1    | EE2506C    | EE2506O  |
| 25.06      | RED    | EE250621 | EE250625 | EE250629 | SPP2506-1    | SGP2506-1    | EE2506C    | EE2506O  |
| 25.08      | RED    | EE250821 | EE250825 | EE250829 | N/A          | N/A          | EE2506C    | EE2506O  |
| 25.10      | RED    | EE251021 | EE251025 | EE251029 | N/A          | N/A          | EE2506C    | EE2506O  |
| 30.04      | BLUE   | EE300421 | EE300425 | EE300429 | SPP3004-1    | SGP3004-1    | EE3006C    | EE3006O  |
| 30.06      | BLUE   | EE300621 | EE300625 | EE300629 | SPP3006-1    | SGP3006-1    | EE3006C    | EE3006O  |
| 30.08      | BLUE   | EE300821 | EE300825 | EE300829 | N/A          | N/A          | EE3006C    | EE3006O  |
| 30.10      | BLUE   | EE301021 | EE301025 | EE301029 | N/A          | N/A          | EE3006C    | EE3006O  |
| 35.04      | GREEN  | EE350421 | EE350425 | EE350429 | SPP3504-1    | SGP3504-1    | EE3506C    | EE3506O  |
| 35.06      | GREEN  | EE350621 | EE350625 | EE350629 | SPP3506-1    | SGP3506-1    | EE3506C    | EE3506O  |
| 35.08      | GREEN  | EE350821 | EE350825 | EE350829 | N/A          | N/A          | EE3506C    | EE3506O  |
| 35.10      | GREEN  | EE351021 | EE351025 | EE351029 | N/A          | N/A          | EE3506C    | EE3506O  |
| 40.04      | BLACK  | EE400421 | EE400425 | EE400429 | SPP4004-1    | SGP4004-1    | EE4006C    | EE4006O  |
| 40.06      | BLACK  | EE400621 | EE400625 | EE400629 | SPP4006-1    | SGP4006-1    | EE4006C    | EE4006O  |
| 40.08      | BLACK  | EE400821 | EE400825 | EE400829 | N/A          | N/A          | EE4006C    | EE4006O  |
| 40.10      | BLACK  | EE401021 | EE401025 | EE401029 | N/A          | N/A          | EE4006C    | EE4006O  |

|                                                                                                                                                                      |            |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| <b>Assorted Pack 20</b> 4 Per Pack<br>(20.04, 20.06, 20.08, 20.10) <span style="background-color: yellow; display: inline-block; width: 50px; height: 15px;"></span> |            |
| 21mm                                                                                                                                                                 | EEASST2021 |
| 25mm                                                                                                                                                                 | EEASST2025 |
| 29mm                                                                                                                                                                 | EEASST2029 |

|                                                                                                                                                                   |            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| <b>Assorted Pack 25</b> 4 Per Pack<br>(25.04, 25.06, 25.08, 25.10) <span style="background-color: red; display: inline-block; width: 50px; height: 15px;"></span> |            |
| 21mm                                                                                                                                                              | EEASST2521 |
| 25mm                                                                                                                                                              | EEASST2525 |
| 29mm                                                                                                                                                              | EEASST2529 |

|                                                                                                                                                                    |            |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| <b>Assorted Pack 30</b> 4 Per Pack<br>(30.04, 30.06, 30.08, 30.10) <span style="background-color: blue; display: inline-block; width: 50px; height: 15px;"></span> |            |
| 21mm                                                                                                                                                               | EEASST3021 |
| 25mm                                                                                                                                                               | EEASST3025 |
| 29mm                                                                                                                                                               | EEASST3029 |

|                                                                                                                                                                     |            |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| <b>Assorted Pack 40</b> 4 Per Pack<br>(40.04, 40.06, 40.08, 40.10) <span style="background-color: black; display: inline-block; width: 50px; height: 15px;"></span> |            |
| 21mm                                                                                                                                                                | EEASST4021 |
| 25mm                                                                                                                                                                | EEASST4025 |
| 29mm                                                                                                                                                                | EEASST4029 |

Pair EdgeEvolve™ with:

EDGE BIO CERAMIC™  
SEALER


818365020072



**\$19.99**  
6 Pack

## Elegance in Innovation, Strength in Versatility

- Minimally Invasive Progressive/Regressive Variable Tapered File System
- Ability to perform in **both rotary or forward reciprocation** settings
- Rotary Motion: Delivers maximum cutting efficiency and debris removal
- Forward Reciprocation Motion: Maximizes safety while still allowing for efficient cutting
- Unique geometry, innovative cutting angles and optimized FireWire™ Blaze Heat-Treatment\*\* allows files to navigate curvatures with ease while maintaining dentin preservation in the coronal third
- Alternative to ProTaper Ultimate™ & Minimally Invasive Alternative to ProTaper™\*



- Diamond Cross-Section
- Available in 21mm and 25mm lengths
- Progressive/Regressive Variable Taper File
- Pre-Sterilized



Incredibly Strong



Very Versatile



Minimally Invasive

\*ProTaper Ultimate™ and ProTaper™ are registered trademarks of Dentsply Sirona. \*\*Slider is not heat-treated Best Product Rating Dental Product Shopper Vol 19 No 2

# EDGE DIAMOND



## UTOPIA™



Best Product Rating Dental Product Shopper Vol 19 No 2

| 21mm Pre-Sterilized 6 Pack |      |       |        |            | EdgeTaper™ Points |              | EdgeTaper™ Obturators |          |
|----------------------------|------|-------|--------|------------|-------------------|--------------|-----------------------|----------|
|                            | Size | Taper | Color  | Item Code  | Paper Points      | Gutta Points | EdgeCore              | EdgeFill |
| SLIDER*                    | 16   | 02v   | PURPLE | EDU1602V21 | N/A               | N/A          | N/A                   | N/A      |
| SHAPER                     | 20   | 04v   | WHITE  | EDU2004V21 | N/A               | N/A          | N/A                   | N/A      |
| F1                         | 20   | 07v   | YELLOW | EDU2007V21 | ETPPF1-1          | ETGPF1-1     | ET2006C               | ET2006O  |
| F2                         | 25   | 07v   | RED    | EDU2507V21 | ETPPF2-1          | ETGPF2-1     | ET2506C               | ET2506O  |
| F3                         | 30   | 07v   | BLUE   | EDU3007V21 | ETPPF3-1          | ETGPF3-1     | ET3006C               | ET3006O  |
| FX                         | 35   | 07v   | GREEN  | EDU3507V21 | ETPPF4-1          | ETGPF4-1     | ET4006C               | ET4006O  |

| 25mm Pre-Sterilized 6 Pack |      |       |        |            | EdgeTaper™ Points |              | EdgeTaper™ Obturators |          |
|----------------------------|------|-------|--------|------------|-------------------|--------------|-----------------------|----------|
|                            | Size | Taper | Color  | Item Code  | Paper Points      | Gutta Points | EdgeCore              | EdgeFill |
| SLIDER*                    | 16   | 02v   | PURPLE | EDU1602V25 | N/A               | N/A          | N/A                   | N/A      |
| SHAPER                     | 20   | 04v   | WHITE  | EDU2004V25 | N/A               | N/A          | N/A                   | N/A      |
| F1                         | 20   | 07v   | YELLOW | EDU2007V25 | ETPPF1-1          | ETGPF1-1     | ET2006C               | ET2006O  |
| F2                         | 25   | 07v   | RED    | EDU2507V25 | ETPPF2-1          | ETGPF2-1     | ET2506C               | ET2506O  |
| F3                         | 30   | 07v   | BLUE   | EDU3007V25 | ETPPF3-1          | ETGPF3-1     | ET3006C               | ET3006O  |
| FX                         | 35   | 07v   | GREEN  | EDU3507V25 | ETPPF4-1          | ETGPF4-1     | ET4006C               | ET4006O  |

\*Slider is not Heat Treated

| Assorted Pack 6 Pack |            |
|----------------------|------------|
|                      |            |
| 21mm                 | EDUASST21V |
| 25mm                 | EDUASST25V |

Pair EdgeDiamond Utopia™ with:

# EDGE UTOPIA™

## BioCeramic Sealer 2g

818365020072




# EDGE V-TAPER<sup>HT</sup>



**\$19.99**  
6 Pack

## Strong. Advanced. Flexible. Efficient.

- Designed with a Safe-Core™ Parabolic Cross- Section
- The Safe-Core™ Parabolic Cross- Section combines the attributes of being a highly efficient and flexible instrument while being resistant to fracture, minimizing your stress and resulting in a more reliable procedure.
- FireWire NiTi Heat Treated
- Variable Tapered Rotary Files
- Features its own dedicated Glidepath system
- Non-Sterile
- Alternative to DC Taper®<sup>1</sup>



- Parabolic shaped cross-section
- Available in 3 lengths: 21, 25 & 31 mm
- Speed 300-400 RPM
- Torque - 4.0- 5.2 Ncm



Highly  
Flexible



No Bounce  
Back



High  
Strength

1) DC Taper® is a registered trademark of SS White Dental®.

**Rotary • Variable Taper • File System**

# EDGE V-TAPER<sup>HT</sup>



| Orifice Opener - 6 Pack |       |            |
|-------------------------|-------|------------|
| Size/Taper              | Color | 17mm       |
| 25.09                   | GOLD  | VFGP250917 |

| Glidepaths 6 Pack |       |            |            |            |            |
|-------------------|-------|------------|------------|------------|------------|
| Size/Taper        | Color | 21mm       | 25mm       | 29mm       | 31mm       |
| 13.03 GP          | PINK  | VFGP130321 | VFGP130325 | N/A        | VFGP130331 |
| 14.03 GP          | WHITE | VFGP140321 | VFGP140325 | VFGP140329 | N/A        |
| 17.04 GP          | WHITE | VFGP170421 | VFGP170425 | N/A        | VFGP170431 |

| 6 Pack     |        |          |          |          | Points       |              | Obturers |          |
|------------|--------|----------|----------|----------|--------------|--------------|----------|----------|
| Size/Taper | Color  | 21mm     | 25mm     | 31mm     | Paper Points | Gutta Points | EdgeCore | EdgeFill |
| 20.06      | YELLOW | VF200621 | VF200625 | VF200631 | EVHTPP2006   | EVHTGP2006   | VF2006C  | VF2006O  |
| 25.06      | RED    | VF250521 | VF250525 | VF250531 | EVHTPP2506   | EVHTGP2506   | VF2506C  | VF2506O  |
| 30.06      | BLUE   | VF220721 | VF220725 | VF220731 | EVHTPP3006   | EVHTGP3006   | VF3006C  | VF3006O  |
| 35.06      | GREEN  | VF250821 | VF250825 | VF250831 | EVHTPP3506   | EVHTGP3506   | VF3506C  | VF3506O  |
| 40.06      | BLACK  | VF300621 | VF300625 | VF300631 | EVHTPP4006   | EVHTGP4006   | VF4006C  | VF4006O  |
| 45.06      | WHITE  | VF300821 | VF300825 | VF300831 | EVHTPP4506   | EVHTGP4506   | VF4506C  | VF4506O  |
| 50.06      | YELLOW | VF350621 | VF350625 | VF350631 | EVHTPP5006   | EVHTGP5006   | VF5006C  | VF5006O  |
| 22.07      | YELLOW | VF400621 | VF400625 | VF400631 | EVHTPP2207   | EVHTGP2207   | VF2006C  | VF2006O  |
| 25.08      | RED    | VF450621 | VF450625 | VF450631 | EVHTPP2508   | EVHTGP2508   | VF2506C  | VF2506O  |
| 30.08      | BLUE   | VF500621 | VF500625 | VF500631 | EVHTPP3008   | EVHTGP3008   | VF3006C  | VF3006O  |

Pair EdgeV-Taper HT™ with:

**EDGE BIO CERAMIC™**  
SEALER

818365020072



**EdgeVerify™ for  
EdgeV-Taper HT™ 6 Pack**

| Size | Item Code |
|------|-----------|
| 20   | VF2006V   |
| 25   | VF2506V   |
| 30   | VF3006V   |
| 35   | VF3506V   |
| 40   | VF4006V   |
| 45   | VF4506V   |
| 50   | VF5006V   |


# EDGETAPER BLAZE UTOPIA™



**\$19.99**  
6 Pack

## Enhanced Technology. Remarkably Beautiful.

- Features proprietary FireWire Blaze Heat Treatment
- Electropolished for better cutting efficiency
- Shorter handle length for better posterior access
- Non-cutting tip easily glides to the working length
- Pre-Sterilized
- Alternative to ProTaper Gold®
- Optimized version of EdgeTaper Platinum™



- Triangular Shaped Cross-Section
- Available in 3 lengths: 21, 25 & 31 mm
- Speed: 300 RPM
- Torque: 5.1 Ncm (SX&S1), 2.5 Ncm (S2&F1), 3.1 Ncm (F2&F3)

\*ProTaper® Gold is a registered trademark of Dentsply Tulsa Dental Specialties



High  
Cutting  
Efficiency



Electropolished  
Finish



High  
Strength

**Rotary • Variable Taper • File System**

# EDGETAPER BLAZE



## UTOPIA™



| Pre-Sterilized 6 Pack |                                                                                          |            |
|-----------------------|------------------------------------------------------------------------------------------|------------|
| Size                  | Color                                                                                    | 19mm       |
| SX                    |  YELLOW | ETBUSX19SU |

| Pre-Sterilized 6 Pack |                                                                                           |            |            |            | EdgeTaper™ Points |              | EdgeTaper™ Obturators |          |
|-----------------------|-------------------------------------------------------------------------------------------|------------|------------|------------|-------------------|--------------|-----------------------|----------|
| Size                  | Color                                                                                     | 21mm       | 25mm       | 31mm       | Paper Points      | Gutta Percha | EdgeCore              | EdgeFill |
| S1                    |  PURPLE  | ETBUS121SU | ETBUS125SU | ETBUS131SU | N/A               | N/A          | N/A                   | N/A      |
| S2                    |  WHITE   | ETBUS221SU | ETBUS225SU | ETBUS231SU | N/A               | N/A          | N/A                   | N/A      |
| F1                    |  YELLOW | ETBUF121SU | ETBUF125SU | ETBUF131SU | ETPPF1-1          | ETGPF1-1     | ET2006C               | ET2006O  |
| F2                    |  RED   | ETBUF221SU | ETBUF225SU | ETBUF231SU | ETPPF2-1          | ETGPF2-1     | ET2506C               | ET2506O  |
| F3                    |  BLUE  | ETBUF321SU | ETBUF325SU | ETBUF331SU | ETPPF3-1          | ETGPF3-1     | ET3006C               | ET3006O  |

| Assorted Packs 6 Pack<br>(SX, S1, S2, F1, F2, F3)      |           |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| 21mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | ETBUA21SU |
| 25mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | ETBUA25SU |
| 31mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | ETBUA31SU |

Pair EdgeTaper Blaze Utopia™ with:

# EDGE UTOPIA™

## BioCeramic Sealer 2g

818365020072



# EDGETAPER


## PLATINUM™



**\$19.99**  
6 Pack

### Edge-ified Sequence with Enhanced Strength and Flexibility

- 2x The Cyclic Fatigue as ProTaper Gold®
- 6x The Cyclic Fatigue as ProTaper®
- FireWire NiTi Heat Treated
- Variable Tapered Files
- Non-Sterile



- Convex Triangular Cross- Section
- Available in 3 lengths: 21, 25 & 31 mm
- Speed: 300-400 RPM
- Torque: 4.0-5.2 Ncm



Highly  
Flexible



No Bounce  
Back




High  
Strength

**Rotary • Variable Taper • File System**

# EDGE TAPER

## PLATINUM™



| 6 Pack |                                                                                          |          |
|--------|------------------------------------------------------------------------------------------|----------|
| Size   | Color                                                                                    | 19mm     |
| SX     |  YELLOW | ETSX19HT |

| 6 Pack |                                                                                            |          |          |          | Points       |              | Obturators |          |
|--------|--------------------------------------------------------------------------------------------|----------|----------|----------|--------------|--------------|------------|----------|
| Size   | Color                                                                                      | 21mm     | 25mm     | 31mm     | Paper Points | Gutta Percha | EdgeCore   | EdgeFill |
| S1     |  PURPLE   | ETS121HT | ETS125HT | ETS131HT | N/A          |              |            |          |
| S2     |  WHITE    | ETS221HT | ETS225HT | ETS231HT | N/A          |              |            |          |
| F1     |  YELLOW   | ETS121HT | ETF125HT | ETF131HT | ETPPF1-1     | ETGPF1-1     | ET2006C    | ET2006O  |
| F2     |  RED      | ETF221HT | ETF225HT | ETF231HT | ETPPF2-1     | ETGPF2-1     | ET2506C    | ET2506O  |
| F3     |  BLUE   | ETF321HT | ETF325HT | ETF331HT | ETPPF3-1     | ETGPF3-1     | ET3006C    | ET3006O  |
| F4     |  BLACK  | ETF421HT | ETF425HT | ETF431HT | ETPPF4-1     | ETGPF4-1     | ET4006C    | ET4006O  |
| F5     |  YELLOW | ETF521HT | ETF525HT | ETF531HT | ETPPF5-1     | ETGPF5-1     | ET5006C    | ET5006O  |

| Assorted Packs                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                               |              |
|----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| Files                            |       | Paper Points                                                                                                                                                                                                                                                                                                                                                                                                                                  | Gutta Percha |
| 6 Pack: (SX, S1, S2, F1, F2, F3) |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |      |              |
| 21mm                             | ETP21MM                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | ETPPAST                                                                                                                                                                                                                                                                                                                                                                                                                                       | ETGPAST      |
| 25mm                             | ETP25MM                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                                                                                                                                               |              |
| 31mm                             | ETP31MM                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                                                                                                                                               |              |

Pair EdgeTaper Platinum™ with:

# EDGE BIO CERAMIC™

## SEALER

818365020072




# EDGE TAPER™



**\$19.99**  
**6 Pack**

## Edge-ified Sequence with Enhanced Strength

- 6x The Cyclic Fatigue as ProTaper®<sup>1</sup>
- Variable Tapered Files
- Non-Sterile
- Non-Heat Treated



- Convex Triangular Cross- Section
- Available in 3 lengths: 21, 25 & 31 mm
- Speed: 300-400 RPM
- Torque: 4.0-5.2 Ncm



No Bounce  
Back



High  
Strength

## Cyclic Fatigue Testing - Average Time to Failure in Seconds

**EDGE TAPER™**  
PLATINUM™

677 Seconds

ProTaper Gold®

312 Seconds

**EDGE TAPER™**

127 Seconds

ProTaper®

119 Seconds

1) ProTaper® is a registered trademark of Dentsply Tulsa Dental Specialties. \*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph:<http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit [www.edgeendo.com/matter-of-temperature](http://www.edgeendo.com/matter-of-temperature)

## Rotary • Variable Taper • File System

# EDGE TAPER™

Non-Heat Treated



| 6 Pack |                                                                                          |        |
|--------|------------------------------------------------------------------------------------------|--------|
| Size   | Color                                                                                    | 19mm   |
| SX     |  YELLOW | ETSX19 |

| 6 Pack |                                                                                            |        |        |        | Points       |              | Obturators |          |
|--------|--------------------------------------------------------------------------------------------|--------|--------|--------|--------------|--------------|------------|----------|
| Size   | Color                                                                                      | 21mm   | 25mm   | 31mm   | Paper Points | Gutta Percha | EdgeCore   | EdgeFill |
| S1     |  PURPLE   | ETS121 | ETS125 | ETS131 | N/A          |              |            |          |
| S2     |  WHITE    | ETS221 | ETS225 | ETS231 | N/A          |              |            |          |
| F1     |  YELLOW   | ETF121 | ETF125 | ETF131 | ETPPF1-1     | ETGPF1-1     | ET2006C    | ET2006O  |
| F2     |  RED    | ETF221 | ETF225 | ETF231 | ETPPF2-1     | ETGPF2-1     | ET2506C    | ET2506O  |
| F3     |  BLUE   | ETF321 | ETF325 | ETF331 | ETPPF3-1     | ETGPF3-1     | ET3006C    | ET3006O  |
| F4     |  BLACK  | ETF421 | ETF425 | ETF431 | ETPPF4-1     | ETGPF4-1     | ET4006C    | ET4006O  |
| F5     |  YELLOW | ETF521 | ETF525 | ETF531 | ETPPF5-1     | ETGPF5-1     | ET5006C    | ET5006O  |

EDGE BIO CERAMIC™  
SEALER



| Assorted Packs                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |          |                                                                                                                                                                                                                                                                                                                                                                                                                                               |              |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| Files       |          | Paper Points                                                                                                                                                                                                                                                                                                                                                                                                                                  | Gutta Percha |
| <b>6 Pack:</b> (SX, S1, S2, F1, F2, F3)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |          |      |              |
| 21mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | ETASST21 | ETPPAST                                                                                                                                                                                                                                                                                                                                                                                                                                       | ETGPAST      |
| 25mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | ETASST25 |                                                                                                                                                                                                                                                                                                                                                                                                                                               |              |
| 31mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | ETASST31 |                                                                                                                                                                                                                                                                                                                                                                                                                                               |              |


# EDGEONE BLAZE UTOPIA™



## Enhanced Technology. Remarkably Beautiful.

**\$12.88**  
3 Pack

- Premier Single File Shaping Reciprocating System
- Features proprietary FireWire Blaze Heat Treatment
- Electropolished for better cutting efficiency
- Shorter handle length for better posterior access
- Variable Taper File
- Pre-Sterilized
- Alternative to WaveOne Gold® and EdgeOne Fire™



- S-Shaped Cross-Section - For Increased Debris Removal
- Available in 3 lengths: 21, 25 & 31 mm
- CCW 170°/ CW 60°



Highly  
Flexible



No Bounce  
Back



High  
Strength

VOL. 18 NO. 4 Dental Product Shopper. "WaveOne® Gold is a registered trademarks of Dentsply Tulsa Dental Specialties.

EDGEONE BLAZE



UTOPIA™



VOL. 18 NO. 4 Dental Product Shopper.

| Pre-Sterilized 3 Pack |                                            |              |              |              | Points       |              | Obturators |           |
|-----------------------|--------------------------------------------|--------------|--------------|--------------|--------------|--------------|------------|-----------|
| Size                  | Color                                      | 21mm         | 25mm         | 31mm         | Paper Points | Gutta Points | EdgeCore   | EdgeFill  |
| R-GLIDE*              | <input type="checkbox"/> WHITE             | EOBURG1421SU | EOBURG1425SU | EOBURG1431SU | N/A          |              | N/A        |           |
| SM R20                | <input checked="" type="checkbox"/> YELLOW | EOBU2021SU   | EOBU2025SU   | EOBU2031SU   | SPP2006-1    | EORGP20-1    | EOPSM20C   | EOPSM20O  |
| PRI R25               | <input checked="" type="checkbox"/> RED    | EOBU2521SU   | EOBU2525SU   | EOBU2531SU   | SPP2506-1    | EORGP25-1    | EOPPRI25C  | EOPPRI25O |
| MED R35               | <input checked="" type="checkbox"/> GREEN  | EOBU3521SU   | EOBU3525SU   | EOBU3531SU   | SPP3506-1    | EORGP35-1    | EOPMED35C  | EOPMED35O |
| LG R45                | <input type="checkbox"/> WHITE             | EOBU4521SU   | EOBU4525SU   | EOBU4531SU   | SPP4506-1    | EORGP45-1    | EOPLG45C   | EOPLG45O  |

| Assorted Pack* 4 Pack<br>(R20, R25, R35, R45) |           |
|-----------------------------------------------|-----------|
| 21MM                                          | EOBUA21SU |
| 25MM                                          | EOBUA25SU |
| 31MM                                          | EOBUA31SU |

Pair EdgeOne Blaze Utopia™ with:

\*R-Glide Sold Separately



BioCeramic Sealer 2g

818365020072





# EDGEONE FIRE™



**\$12.88**  
3 Pack

## The Standard in Reciprocating Technology

- 5x the Cyclic Fatigue as WaveOne® Gold
- FireWire NiTi Heat Treated
- Variable Tapered Reverse Reciprocating Files
- Non-Sterile
- Alternative to WaveOne® Gold



- Parallelogram shaped cross-section
- Available in 3 lengths: 21, 25 & 31 mm
- Designed to shape canals in reverse-reciprocating motion CCW 150°/CW 60° or WaveOne® Gold<sup>1</sup> settings



Highly Flexible



No Bounce Back



High Strength

## Cyclic Fatigue Testing - Average Time to Failure in Seconds

**EDGEONE**  
FIRE™

2844 Seconds

WaveOne® Gold

486 Seconds

WaveOne® Gold is a registered trademarks of Dentsply Tulsa Dental Specialties. \*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph:<http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit [www.edgeendo.com/matter-of-temperature](http://www.edgeendo.com/matter-of-temperature)

# EDGEONE

## FIRE™



| 3 Pack       |                                |         |         |         |
|--------------|--------------------------------|---------|---------|---------|
| Size         | Color                          | 21mm    | 25mm    | 31mm    |
| GLIDEPATH 19 | <input type="checkbox"/> WHITE | EOFGP21 | EOFGP25 | EOFGP31 |

| 3 Pack |                                            |          |          |          | Points       |              | Obturators |           |
|--------|--------------------------------------------|----------|----------|----------|--------------|--------------|------------|-----------|
| Size   | Color                                      | 21mm     | 25mm     | 31mm     | Paper Points | Gutta Percha | EdgeCore   | EdgeFill  |
| SM 20  | <input checked="" type="checkbox"/> YELLOW | EOFSM21  | EOFSM25  | EOFSM31  | EOFSMPP-1    | EOFSMGP-1    | EOPSM20C   | EOPSM20O  |
| PRI 25 | <input checked="" type="checkbox"/> RED    | EOFPRI21 | EOFPRI25 | EOFPRI31 | EOFPRIPP-1   | EOFPRIGP-1   | EOPPRI25C  | EOPPRI25O |
| MED 35 | <input checked="" type="checkbox"/> GREEN  | EOFMED21 | EOFMED25 | EOFMED31 | EOFMEDPP-1   | EOFMEDGP-1   | EOPMED35C  | EOPMED35O |
| LG 45  | <input type="checkbox"/> WHITE             | EOFLG21  | EOFLG25  | EOFLG31  | EOFLGGPP-1   | EOFLGGP-1    | EOPLG45C   | EOPLG45O  |

Pair EdgeOne Fire™ with:

## EDGE BIO CERAMIC™

### SEALER

818365020072



| Assorted Packs                                                                                                                             |           |                                                                                                                                      |              |
|--------------------------------------------------------------------------------------------------------------------------------------------|-----------|--------------------------------------------------------------------------------------------------------------------------------------|--------------|
| Files <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/> |           | Paper Points                                                                                                                         | Gutta Percha |
| 4 Pack (SM, PRI, MED, LG)                                                                                                                  |           | <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/> |              |
| 21mm                                                                                                                                       | ASSTEOF21 | EOFPPAST-1                                                                                                                           | EOFGPAST-1   |
| 25mm                                                                                                                                       | ASSTEOF25 |                                                                                                                                      |              |
| 31mm                                                                                                                                       | ASSTEOF31 |                                                                                                                                      |              |

## Constant Taper. Reliable Reciprocation.

- Constant taper reciprocating system
- Verifile dictates sequence
- Verifile immediately following 10 and 15 handfiles
- Non-Sterile
- Alternative to WaveOne® Gold<sup>1</sup>



- Triangular shaped cross-section
- Available in 3 lengths: 21, 25 & 31 mm
- Designed to shape canals in reverse-reciprocating motion
- CCW 150°/CW 60° or WaveOne® Gold<sup>1</sup> settings



Highly Flexible



No Bounce Back



High Strength

<sup>1</sup>WaveOne® Gold is a registered trademark of Dentsply Tulsa Dental Specialties.



| 3 Pack    |                                |        |        |        |
|-----------|--------------------------------|--------|--------|--------|
| Size      | Color                          | 21mm   | 25mm   | 31mm   |
| GLIDEPATH | <input type="checkbox"/> WHITE | EOGP21 | EOGP25 | EOGP31 |

| 3 Pack   |                                            |          |          |          | Points       |              | Obturators |           |
|----------|--------------------------------------------|----------|----------|----------|--------------|--------------|------------|-----------|
| Size     | Color                                      | 21mm     | 25mm     | 31mm     | Paper Points | Gutta Percha | EdgeCore   | EdgeFill  |
| VERIFILE | <input type="checkbox"/> WHITE             | EOPVER21 | EOPVER25 | EOPVER31 | N/A          | N/A          | N/A        | N/A       |
| SMALL    | <input checked="" type="checkbox"/> YELLOW | EOPSM21  | EOPSM25  | EOPSM31  | SPP2006-1    | X72006GP-1   | EOPSM20C   | EOPSM20O  |
| PRIMARY  | <input checked="" type="checkbox"/> RED    | EOPPRI21 | EOPPRI25 | EOPPRI31 | SPP2506-1    | X72506GP-1   | EOPPRI25C  | EOPPRI25O |
| MEDIUM   | <input checked="" type="checkbox"/> GREEN  | EOPMED21 | EOPMED25 | EOPMED31 | SPP3506-1    | X73506GP-1   | EOPMED35C  | EOPMED35O |
| LARGE    | <input type="checkbox"/> WHITE             | EOPLG21  | EOPLG25  | EOPLG31  | SPP4506-1    | X74506GP-1   | EOPLG45C   | EOPLG45O  |

Pair EdgeOne Platinum<sup>™</sup> with:

**EDGE**BIOCERAMIC<sup>™</sup>  
SEALER

818365020072



| Assorted Pack                                                                                                                                                       |           |                                                                                                                                      |              |  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|--------------------------------------------------------------------------------------------------------------------------------------|--------------|--|
| Files <input type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/> |           | Paper Points                                                                                                                         | Gutta Percha |  |
| 5 Pack (VER, SM, PRI, MED, LG)                                                                                                                                      |           | <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/> |              |  |
| 21mm                                                                                                                                                                | ASSTEOP21 | EOFPPAST-1                                                                                                                           | EOFGPAST-1   |  |
| 25mm                                                                                                                                                                | ASSTEOP25 |                                                                                                                                      |              |  |
| 31mm                                                                                                                                                                | ASSTEOP31 |                                                                                                                                      |              |  |



# EDGEFILE® X1



**\$18.99**  
3 Pack

## Constant Taper. Consistent Reciprocation.

- Tighter pitch and flutes which deflect more torsional forces along the length of file due to the parabolic cross-section
- All X1 files (20.06 small, 25.06 primary, 40.06 large) have a constant 6% taper until the flute diameter is 1.0 mm then it goes to 0 taper
- FireWire NiTi Heat Treated
- Short handle 11mm
- Non-Sterile
- Alternative to WaveOne®<sup>1</sup>



- Parabolic Shaped Cross-Section
- Non - Cutting Tip
- Available in 3 lengths: 21, 25 & 31 mm
- Constant Tapered Reciprocating file system should use the same motor settings as WaveOne®<sup>1</sup>



Highly Flexible



No Bounce Back



High Strength

<sup>1</sup>) WaveOne® is a registered trademark of Dentsply Sirona.

# EDGEFILE® X1



| 3 Pack     |                                                                                          |            |            |            | Points       |              | Obturators |          |
|------------|------------------------------------------------------------------------------------------|------------|------------|------------|--------------|--------------|------------|----------|
| Size/Taper | Color                                                                                    | 21mm       | 25mm       | 29mm       | Gutta Percha | Paper Points | EdgeCore   | EdgeFill |
| 20.06      |  YELLOW | X1200621RF | X1200625RF | X1200629RF | SPP2006-1    | X72006GP-1   | X72006C    | X72006O  |
| 25.06      |  RED    | X1250621RF | X1250625RF | X1250629RF | SPP2506-1    | X72506GP-1   | X72506C    | X72506O  |
| 40.06      |  BLACK  | X1400621RF | X1400625RF | X1400629RF | SPP4006-1    | X74006GP-1   | X74006C    | X74006O  |

Pair EdgeFile® X1 with:

**EDGEBIOCERAMIC™**  
SEALER

818365020072



**Assorted Pack 3 Pack**  
(20, 25, 40)



|      |          |
|------|----------|
| 21mm | ASSTX121 |
| 25mm | ASSTX125 |
| 29mm | ASSTX129 |

|                                                  |               |
|--------------------------------------------------|---------------|
| <b>EdgeLube™ - \$7.99</b><br>(17% EDTA SOLUTION) |               |
| 100ML BOTTLE                                     | 8806387892083 |



Presence of calcification in the canal. Narrow canals on posterior teeth. Smear layer removal.

Removal of inorganic substances by Chelation (EDTA-Ca<sup>2+</sup>). Easier instrumentation by softening dentin canal wall. Produces a canal wall free of smear layer.

|                                             |               |
|---------------------------------------------|---------------|
| <b>EdgeGel™ - \$10.99</b><br>(19% EDTA GEL) |               |
| 5.4ML BOTTLE                                | 8806387892076 |



Root canal preparation with files and reamers. Presence of calcification in the canal. Narrow canals in posterior teeth. Smear layer removal.

Removal of inorganic substances by chelation (EDTA-Ca<sup>2+</sup>). Easier filling with lubricating action. Produces a canal wall free of smear layer.

# Irrigation Tips & Syringes

|                                                          |                                                                                   |                |
|----------------------------------------------------------|-----------------------------------------------------------------------------------|----------------|
| <b>Edge® Irrigation Needle Tips - \$8.88</b><br>100 Pack |                                                                                   |                |
| 23ga. Blue                                               |  | 10818365020260 |
| 27ga. Yellow                                             |  | 10818365020277 |
| 30ga. Purple                                             |  | 10818365020284 |



|                                                          |                                                                                    |                |
|----------------------------------------------------------|------------------------------------------------------------------------------------|----------------|
| <b>Edge® Irrigation Probe Tips - \$14.88</b><br>100 Pack |                                                                                    |                |
| 23ga. Blue                                               |   | 10818365020291 |
| 27ga. Yellow                                             |   | 10818365020307 |
| 30ga. Purple                                             |  | 10818365020314 |



|                                                                                         |                                                                                     |                |
|-----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|----------------|
| <b>Edge® Pre-Assembled Irrigation Syringe with 1" Notched Tip - \$19.88</b><br>100 Pack |                                                                                     |                |
| 23ga. Blue                                                                              |  | 00818365021062 |
| 27ga. Yellow                                                                            |  | 00818365021055 |
| 30ga. Purple                                                                            |  | 00818365021079 |

|                                                     |                                                                                     |                |
|-----------------------------------------------------|-------------------------------------------------------------------------------------|----------------|
| <b>Edge® Luer-Lock Syringe - \$9.88</b><br>100 Pack |                                                                                     |                |
| 3CC                                                 |  | 10818365020383 |

|                                             |              |
|---------------------------------------------|--------------|
| <b>2% Chlorhexidine Gluconate - \$29.88</b> |              |
| 480ML BOTTLE                                | 818365020041 |



|                           |              |
|---------------------------|--------------|
| <b>17% EDTA - \$29.88</b> |              |
| 480ML BOTTLE              | 818365020058 |



|                                        |              |
|----------------------------------------|--------------|
| <b>5% Sodium Hypochlorite- \$32.88</b> |              |
| 480ML BOTTLE                           | 818365020034 |



**EDGEMIX™**

|                |             |
|----------------|-------------|
| <b>\$59.88</b> |             |
| 16OZ BOTTLE    | 18365020416 |



EDTA-based formula enhanced with Chlorhexidine for added anti-bacterial properties. Easy to use. 2-3 mL per canal. In instances where calcifications or a lot of debris is present, the higher end dosage should be utilized.

- Clears 27% more dentin tubules than EDTA
- Kills 99.99% bacteria in 10 seconds
- Does not form a precipitant when mixed with NaOCl

# EDGE BIO CERAMIC™ SEALERS

## EDGE BIO CERAMIC™ SEALER

2g syringe + 15 disposable tips

### Cold Obturation

EdgeBioCeramic Sealer™ trusted in over **50 million cases worldwide** since 2008, EdgeBioCeramic Sealer™ offers proven, reliable performance. Easy-to-use, premixed, and resin-free, it's hydrophilic with zero shrinkage. The reduced waste tip preserves material by 62%, maximizing efficiency and cost savings.

#### Features:

- Radiopacity: High
- Backed by over 17 years of clinical use and extensive studies
- Ideal viscosity for single cone hydraulic condensation
- Composition: Tricalcium Silicate, Dicalcium Silicate, Zirconium Oxide, Tantalum Oxide

**\$101.00**

Sealer Sku: 818365020072  
Disposable Tips Sku: EEBCT15

**\$18.73/ 15 Tips**



## EDGE UTOPIA™ BioCeramic Sealer

2g syringe + 20 disposable tips


### Warm & Cold Obturation

EdgeUtopia™ Sealer with its superior flow makes it **optimal for both warm or cold obturation techniques**. Its high pH of 12.5 offers powerful antibacterial properties, keeping the canal infection-free. Designed for moist environments, its hydrophilic nature ensures a secure seal in wet conditions.

#### Features:

- Radiopacity: Higher
- Our most cost effective Bioceramic sealer
- Ideal for both cold and warm obturation techniques
- Composition: Tricalcium silicate, dicalcium silicate, tricalcium aluminate, calcium oxide, zirconium oxide, silicon dioxide, polyethylene glycol and iron oxide

Sealer Sku: EUBCSEAL2G  
**\$59.88**



## EDGE BIO CERAMIC™ THERMALFLOW™

1.5g syringe + 15 disposable tips

### Warm Obturation

EdgeBioCeramic™ ThermalFlow™ is an **optimized version** of EdgeBioCeramic Sealer™ with the same resin-free, hydrophilic, and zero shrinkage benefits, but with **lower viscosity for enhanced warm obturation techniques**. The reduced waste tip preserves material by 62%, maximizing efficiency and cost savings.

#### Features:

- Radiopacity: Highest
- Ideal for Warm Obturation with enhanced flowability
- Lowest viscosity and highest radiopacity
- Composition: Tricalcium Silicate, Dicalcium Silicate, Zirconium Oxide, Tantalum Oxide

**\$157.88**

Sealer Sku: 818365020232  
Disposable Tips Sku: EEBCT15

**\$18.73/ 15 Tips**



# EDGESEAL™



13.5g Dual Syringe

Resin-based root canal sealer designed to work with various obturation methods. Provides excellent biocompatibility. Easy to mix paste. Hermetic sealing ability. Non-staining to teeth. Insoluble in tissue fluids. Good radiopacity.

Working Time: 35 minutes at 23°C (73°F)

Setting Time: 45 minutes at 37°C (99°F)

Film Thickness: 3.3μ

**\$22.99**

818365020010

## EDGECORE™

Pink Gutta Percha Core  
Thermal Obturator



5 Pack & 1 Verifier

## EDGEFILL™

Plastic Carrier  
Thermal Obturator



6 Pack

## 20 Pack Obturators - \$64.88

| Size | EdgeCore | EdgeFill |
|------|----------|----------|
| 20   | EE2020C  | EE2020O  |
| 25   | EE2520C  | EE2520O  |
| 30   | EE3020C  | EE3020O  |
| 35   | EE3520C  | EE3520O  |
| 40   | EE4020C  | EE4020O  |
| 45   | EE4520C  | EE4520O  |
| 50   | EE5020C  | EE5020O  |

## NEW EdgeEndo® Gutta Percha: Optimized Performance, Seamless Results Optimized For Use with Our Most Popular NiTi File Systems!



### Seamless Integration

Tailored to work effortlessly with EdgeEndo file systems for streamlined endodontic procedures.



### Proprietary Stiffening Process

Provides improved handling, less tip crumpling, and improved displacement of sealer.



### Optimized Apical Fit with Coronal Venting

Designed for a tight apical fit with extra space in the coronal third to facilitate the upward venting of sealer thus minimizing the risk of apical extrusion. Additionally, this design helps prevent gutta percha from “hanging up” in longer canals promoting more efficient obturation.

## Gutta Percha Quick Reference


These are the suggested gutta percha points for EdgeOne Fire, all EdgeTaper Families and EdgeOne Blaze Utopia NiTi Systems.

| EDGEONE<br>FIRE™<br>Gutta Percha                                                                                                                                                  |         |            |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|------------|
| 29.5mm- 60 Pack                                                                                                                                                                   |         | Item Codes |
| YELLOW                                                                                                                                                                            | SMALL   | EOFSGMP    |
| RED                                                                                                                                                                               | PRIMARY | EOFPRIGP   |
| GREEN                                                                                                                                                                             | MEDIUM  | EOFMEDGP   |
| WHITE                                                                                                                                                                             | LARGE   | EOFLGGP    |
| Assorted: SM (20), PRI (20)<br>MED (10) LG (10)                                                                                                                                   |         | EOFGPAST-1 |
| <ul style="list-style-type: none"> <li>• Taper to match file design</li> <li>• Parallel when diameter meets 1.0mm</li> <li>• <b>Use with EdgeOne Fire™ Files:Pg 70</b></li> </ul> |         |            |

| EDGETAPER™<br>Gutta Percha                                                                                                                                                                                                                                                                      |    |            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----|------------|
| 29.5mm- 60 Pack                                                                                                                                                                                                                                                                                 |    | Item Codes |
| YELLOW                                                                                                                                                                                                                                                                                          | F1 | ETGPF1     |
| RED                                                                                                                                                                                                                                                                                             | F2 | ETGPF2     |
| BLUE                                                                                                                                                                                                                                                                                            | F3 | ETGPF3     |
| BLACK                                                                                                                                                                                                                                                                                           | F4 | ETGPF4     |
| YELLOW                                                                                                                                                                                                                                                                                          | F5 | ETGPF5     |
| Assorted: F1 (10), F2 (10), F3 (20), F4 (10), F5 (10)                                                                                                                                                                                                                                           |    | ETGPAST    |
| <ul style="list-style-type: none"> <li>• Taper to match file design</li> <li>• Parallel 16mm from tip</li> <li>• <b>Use with EdgeTaper Platinum™ Files: Pg 71</b></li> <li>• <b>Use with EdgeTaper™ Files: Pg 72</b></li> <li>• <b>Use with EdgeTaper Blaze Utopia™ Files: Pg 73</b></li> </ul> |    |            |

| EDGEONE R UTOPIA™<br>Gutta Percha                                                                                                                                                          |     |            |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|------------|
| 29.5mm- 60 Pack                                                                                                                                                                            |     | Item Codes |
| YELLOW                                                                                                                                                                                     | R20 | EORGP20-1  |
| RED                                                                                                                                                                                        | R25 | EORGP25-1  |
| GREEN                                                                                                                                                                                      | R35 | EORGP35-1  |
| WHITE                                                                                                                                                                                      | R45 | EORGP45-1  |
| <ul style="list-style-type: none"> <li>• Taper to match file design</li> <li>• Parallel when diameter meets 1.0mm</li> <li>• <b>Use with EdgeOne Blaze Utopia™ Files: Pg 73</b></li> </ul> |     |            |



| New Gutta Percha Item Code                                                          | File                                                                                | Size               | File Item Code (3 Pack) | Length |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------|-------------------------|--------|
|    |    | Small              | EOFSM21                 | 21mm   |
|                                                                                     |                                                                                     |                    | EOFSM25                 | 25mm   |
|                                                                                     |                                                                                     |                    | EOFSM31                 | 31mm   |
|    |    | Primary            | EOFPRI21                | 21mm   |
|                                                                                     |                                                                                     |                    | EOFPRI25                | 25mm   |
|                                                                                     |                                                                                     |                    | EOFPRI31                | 31mm   |
|   |   | Medium             | EOFMED21                | 21mm   |
|                                                                                     |                                                                                     |                    | EOFMED25                | 25mm   |
|                                                                                     |                                                                                     |                    | EOFMED31                | 31mm   |
|  |  | Large              | EOFLG21                 | 21mm   |
|                                                                                     |                                                                                     |                    | EOFLG25                 | 25mm   |
|                                                                                     |                                                                                     |                    | EOFLG31                 | 31mm   |
|  |  | Assorted<br>4 Pack | ASSTEOF21               | 21mm   |
|                                                                                     |                                                                                     |                    | ASSTEOF25               | 25mm   |
|                                                                                     |                                                                                     |                    | ASSTEOF31               | 31mm   |



## EDGE TAPER PLATINUM™

| New Gutta Percha Item Code                                                          | File                                                                                | Size            | File Item Code (6 Pack) | Length |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------------|-------------------------|--------|
|    |    | F1              | ETF121HF                | 21mm   |
|                                                                                     |                                                                                     |                 | ETF125HT                | 25mm   |
|                                                                                     |                                                                                     |                 | ETF131HT                | 31mm   |
|    |    | F2              | ETF221HF                | 21mm   |
|                                                                                     |                                                                                     |                 | ETF225HT                | 25mm   |
|                                                                                     |                                                                                     |                 | ETF231HT                | 31mm   |
|    |   | F3              | ETF321HF                | 21mm   |
|                                                                                     |                                                                                     |                 | ETF325HT                | 25mm   |
|                                                                                     |                                                                                     |                 | ETF331HT                | 31mm   |
|  |  | F4              | ETF421HF                | 21mm   |
|                                                                                     |                                                                                     |                 | ETF425HT                | 25mm   |
|                                                                                     |                                                                                     |                 | ETF431HT                | 31mm   |
|  |  | F5              | ETF521HF                | 21mm   |
|                                                                                     |                                                                                     |                 | ETF525HT                | 25mm   |
|                                                                                     |                                                                                     |                 | ETF531HT                | 31mm   |
|  |  | Assorted 6 Pack | ETP21MM                 | 21mm   |
|                                                                                     |                                                                                     |                 | ETP25MM                 | 25mm   |
|                                                                                     |                                                                                     |                 | ETP31MM                 | 31mm   |

# EDGE TAPER™

| New Gutta Percha Item Code                                                          | File                                                                                | Size            | File Item Code (6 Pack) | Length |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------------|-------------------------|--------|
|    |    | F1              | ETF121                  | 21mm   |
|                                                                                     |                                                                                     |                 | ETF125                  | 25mm   |
|                                                                                     |                                                                                     |                 | ETF131                  | 31mm   |
|    |    | F2              | ETF221                  | 21mm   |
|                                                                                     |                                                                                     |                 | ETF225                  | 25mm   |
|                                                                                     |                                                                                     |                 | ETF231                  | 31mm   |
|   |   | F3              | ETF321                  | 21mm   |
|                                                                                     |                                                                                     |                 | ETF325                  | 25mm   |
|                                                                                     |                                                                                     |                 | ETF331                  | 31mm   |
|  |  | F4              | ETF421                  | 21mm   |
|                                                                                     |                                                                                     |                 | ETF425                  | 25mm   |
|                                                                                     |                                                                                     |                 | ETF431                  | 31mm   |
|  |  | F5              | ETF521                  | 21mm   |
|                                                                                     |                                                                                     |                 | ETF525                  | 25mm   |
|                                                                                     |                                                                                     |                 | ETF531                  | 31mm   |
|  |  | Assorted 6 Pack | ETASST21                | 21mm   |
|                                                                                     |                                                                                     |                 | ETASST25                | 25mm   |
|                                                                                     |                                                                                     |                 | ETASST31                | 31mm   |

## EDGETAPER BLAZE UTOPIA™

| New Gutta Percha Item Code                                                        | File                                                                              | Size | File Item Code (6 Pack) | Length |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------|-------------------------|--------|
|  |  | F1   | ETBUF121SU              | 21mm   |
|                                                                                   |                                                                                   |      | ETBUF125SU              | 25mm   |
|                                                                                   |                                                                                   |      | ETBUF131SU              | 31mm   |
|  |  | F2   | ETBUF221SU              | 21mm   |
|                                                                                   |                                                                                   |      | ETBUF225SU              | 25mm   |
|                                                                                   |                                                                                   |      | ETBUF231SU              | 31mm   |
|  |  | F3   | ETBUF321SU              | 21mm   |
|                                                                                   |                                                                                   |      | ETBUF325SU              | 25mm   |
|                                                                                   |                                                                                   |      | ETBUF331SU              | 31mm   |

## EDGEONE BLAZE UTOPIA™

| New Gutta Percha Item Code                                                          | File                                                                                | Size | File Item Code (3 Pack) | Length |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|------|-------------------------|--------|
|  |  | R20  | EOBU2021SU              | 21mm   |
|                                                                                     |                                                                                     |      | EOBU2025SU              | 25mm   |
|                                                                                     |                                                                                     |      | EOBU2031SU              | 31mm   |
|  |  | R25  | EOBU2521SU              | 21mm   |
|                                                                                     |                                                                                     |      | EOBU2525SU              | 25mm   |
|                                                                                     |                                                                                     |      | EOBU2531SU              | 31mm   |
|  |  | R35  | EOBU3521SU              | 21mm   |
|                                                                                     |                                                                                     |      | EOBU3525SU              | 25mm   |
|                                                                                     |                                                                                     |      | EOBU3531SU              | 31mm   |
|  |  | R45  | EOBU4521SU              | 21mm   |
|                                                                                     |                                                                                     |      | EOBU4525SU              | 25mm   |
|                                                                                     |                                                                                     |      | EOBU4531SU              | 31mm   |

|                            |                |  |
|----------------------------|----------------|--|
| <b>EdgeCaOH™ - \$15.99</b> |                |  |
| 2.2G SYRINGE               | 00818365020003 |  |



Temporary root-canal filler for weeping canals that forms a hard-tissue barrier (apexification). Internal and external root resorption. Pre-mixed paste. Syringe tip for easy delivery of paste into canal. 10 Disposable tips included to prevent cross-contamination.

## EDGE TEMP™ PLUS

**\$7.46**

|     |             |           |
|-----|-------------|-----------|
| 40G | WHITE PASTE | ETEMPPLUS |
|-----|-------------|-----------|



Hydraulic temporary restorative material. Designed for fabricating temporary fillings and other temporary prosthetics until the permanent restoration can be delivered. Setting time: 30 minutes.

## EDGE BIO CERAMIC RETROFILL™ AND PERFORATION REPAIR MATERIAL

**\$157.88**

|      |             |              |
|------|-------------|--------------|
| 0.5G | WHITE PASTE | 818365020089 |
|------|-------------|--------------|



EdgeBioCeramic™ RetroFill is a unique bioceramic repair material that exhibits superior healing and handling properties that has been shown in studies to out perform MTA. In addition, EdgeBioCeramic RetroFill is premixed and ready to use unlike MTA. Simply syringe out a small amount of material and deliver it to the repair site.

# The Future of Root Repair



## Root Repair Material



EdgeUtopia™ RRM is the optimal solution for root repair, offering excellent clinical outcomes with its scientifically backed performance, versatility, and long-term efficacy. EdgeUtopia RRM is the ideal choice for all your root and pediatric repair needs.

This proven pre-mixed calcium silicate repair material has now been upgraded to feature a novel monodose delivery system to improve handling and reduce the risk of cross-contamination. Trust in EdgeUtopia™ RRM for tissue healing, bacteria prevention, and overall root canal repair success.



**Non Staining**



**Wide Range of Applications**



**Desired Biological Performance**



**Convenient Single Dose Delivery**



**Inhibits Bacterial Infiltration**



**Long-Lasting Results**



0.15g per capsule  
10 capsules per box

# EDGE UTOPIA RRM™

## Root Repair Material

### Indications

---



#### Pre-Mixed Ready-to-Use Monodose Capsule

Forget the hassle of your repair material drying out too soon. Experience effortless, reliable performance with every use. No manual mixing and no titration devices needed.



#### Reduced Cross-Contamination

Unlike other premixed, manually mixed or titrated bioceramics, EdgeUtopia RRM™ ensures hygienic, mess-free delivery with no risk of contamination between patients due to its single dose application and the elimination of extensive manipulation in the operatory.



#### Release of Ca<sup>2+</sup>

Proven track record with over 5 years of clinical cases and research studies\* demonstrating exceptional tissue healing and regeneration outcomes.



#### Wide Range of Applications

- Apexification
- Retrofilling
- Perforation Repairs
- Pulpotomy
- Pulp Capping
- External and Internal Resorption
- Deciduous Tooth Repairs



#### Cost-Effective Innovation

More product for a fraction of the price – helping you save money while delivering high-quality results. 3x the amount of our regular syringe formula per box.



0.15g per capsule  
10 capsules per box  
EURRM10MD  
**\$99.51/ 10 Pack**



**Pair with:**  
**Reusable EdgeUtopia™ RRM™**  
**Dispensing Gun**

EURRMDSG  
**\$18.35**

Available research and clinical case studies may refer to other private-label brands which are the same chemical formula as EdgeUtopia RRM™ Root Repair Material. For clarity, the brands EdgeUtopia RRM™ and Bio-C Repair are both manufactured by Angelus and are composed of equivalent chemical formulations.

Rev A 03/2025

|                                   |
|-----------------------------------|
| <b>EdgeFile® X7.04 - \$420.58</b> |
| 1 pack EdgeFile® X7 17.04 25mm    |
| 1 pack EdgeFile® X7 20.04 25mm    |
| 1 pack EdgeFile® X7 25.04 25mm    |
| 1 pack EdgeFile® X7 30.04 25mm    |
| 1 pack EdgeFile® X7 35.04 25mm    |
| 1 pack EdgeFile® X7 40.04 25mm    |

|                                |
|--------------------------------|
| 1 pack EdgeFile® X7 45.04 25mm |
| 3 packs EdgeGlide Path™ 25mm   |
| 4 pack Edge K-File™ #10, 25mm  |
| 4 pack Edge K-File™ #15, 25mm  |
| 1 bottle EdgeMix™              |
| 1 bottle EdgeEndo® 17% EDTA    |



|                                   |
|-----------------------------------|
| <b>EdgeFile® X7.06 - \$420.58</b> |
| 1 pack EdgeFile® X7 17.04 25mm    |
| 1 pack EdgeFile® X7 20.04 25mm    |
| 1 pack EdgeFile® X7 25.04 25mm    |
| 1 pack EdgeFile® X7 30.04 25mm    |
| 1 pack EdgeFile® X7 35.04 25mm    |
| 1 pack EdgeFile® X7 40.04 25mm    |

|                                |
|--------------------------------|
| 1 pack EdgeFile® X7 45.04 25mm |
| 3 packs EdgeGlide Path™ 25mm   |
| 4 pack Edge K-File™ #10, 25mm  |
| 4 pack Edge K-File™ #15, 25mm  |
| 1 bottle EdgeMix™              |
| 1 bottle EdgeEndo® 17% EDTA    |



|                                                  |
|--------------------------------------------------|
| <b>EdgeOne Fire™ - \$358.14</b>                  |
| 4 packs EdgeOne Fire™ Pri. 25mm                  |
| 1 pack EdgeOne Fire™ Small 25mm                  |
| 1 pack EdgeOne Fire™ Med. 25mm                   |
| 1 pack EdgeOne Fire™ Large 25mm                  |
| 3 packs EdgeOne GlidePath™ for EdgeOne Fire 25mm |

|                                |
|--------------------------------|
| 4 packs Edge K-File™ #10, 25mm |
| 4 packs Edge K-File™ #15, 25mm |
| 1 bottle EdgeMix™              |
| 1 bottle EdgeEndo® 17% EDTA    |
| 1 pack Edge K-File™ #15, 25mm  |



|                                                |
|------------------------------------------------|
| <b>EdgeTaper Platinum™ - \$338.13</b>          |
| 6 PACKS EdgeTaper Platinum™ Assorted Pack 25mm |
| 3 packs EdgeGlidePath™ with FireWire™ 25mm     |
| 4 packs Edge K-File™ #10, 25mm                 |

|                                |
|--------------------------------|
| 4 packs Edge K-File™ #15, 25mm |
| 1 EdgeMix™                     |
| 1 EdgeEndo® 17% EDTA           |



## Product Limits By Category

| Category                                                                                                                                     | Limit Qty |
|----------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| Power Equipment                                                                                                                              | Limit 2   |
| Power Accessories                                                                                                                            | Limit 10  |
| Power Equipment Tips & EdgePellets                                                                                                           | Limit 5   |
| Ancillary Products Group 1                                                                                                                   | Limit 5   |
| EdgeLube, EdgeGel, EdgeSeal, EdgeTemp, EdgeCaOH, Edgelce, EdgeMix, Sodium Hypochlorite, Chlorhexidine Gluconate, EDTA, EdgeUtopia BC Solvent |           |
| Handfiles                                                                                                                                    | No limit  |
| Files                                                                                                                                        | Limit 500 |
| Retreatment Files                                                                                                                            | Limit 100 |
| BioCeramic Sealers                                                                                                                           | Limit 10  |
| EdgeBurs 25 & 100 Packs                                                                                                                      | Limit 5   |
| EdgeBurs 5 & 10 Packs                                                                                                                        | Limit 10  |
| EdgeUltrasonic Tips                                                                                                                          | Limit 12  |
| Obturators - All Pack Quantities                                                                                                             | Limit 30  |
| Points- Paper Points, Gutta Percha Points, BioCeramic Points                                                                                 | Limit 30  |
| Supplies                                                                                                                                     | Limit 5   |
| Endo Rings, Edge Luer-Lock Syringes, Edge Irrigation Tips, Pre-Assembled Luer-Lock Syringes,                                                 |           |
| EdgeUtopia RRM                                                                                                                               | Limit 5   |
| EdgeUtopia RRM Dispensing Gun                                                                                                                | Limit 2   |

Limits are cumulative throughout your residency (see front cover for more information). Limit quantity is calculated per category, not by specific product lines, systems, or names.



**EDGEENDO®**

© 2025 EdgeEndo  
No copying without permission.  
Not responsible for typographical errors  
Rev. D 09/2025