



# New & Improved Fit Gutta Percha & Paper Point Comparison Chart



### Improved Navigation

Slimmer profile allows Gutta Percha to drop to length without hangups in the canal curvature. No more binding in canals due to taper.



### Seamless Integration

Tailored to work effortlessly with EdgeEndo file systems for streamlined endodontic procedures.



### Proprietary Stiffening Process

Provides improved handling, less tip crumpling, and improved displacement of sealer.



### Optimized Apical Fit with Coronal Venting

Designed for a tight apical fit with extra space in the coronal third to facilitate the upward venting of sealer thus minimizing the risk of apical extrusion.



## Improved Paper Points

- Highly absorbent
- Clear millimeter reference marks
- Specifically tapered to match corresponding finishing file
- Universal .04 and .06 available

### Gutta Percha















EdgeOne Platinum	2
EdgeOne Fire	3
EdgeOne Blaze Utopia	3
EdgeTaper	4
EdgeTaper Platinum	4
EdgeTaper Blaze Utopia	4
EdgeDiamond Utopia	4
Universal .04	5
EdgeSequel Sapphire .04	5
EdgeSequel Utopia .04	5
EdgeV-Taper HT	5
Universal .06	6

EdgeSequel Sapphire .06	6
EdgeSequel Utopia .06	6
EdgeV-Taper HT	6
DiaDent .04 & .06	8
<b>Paper Points</b>	
EdgeOne Platinum	9
EdgeOne Fire	10
EdgeOne Blaze Utopia	10
EdgeTaper	11
EdgeTaper Platinum	11
EdgeTaper Blaze Utopia	11
EdgeDiamond	11
Universal .04	12









EdgeSequel Sapphire .04	12
EdgeSequel Utopia .04	12
EdgeV-Taper HT	12
EdgeFile X7.04	12
EdgeEvolve .04	12
Universal .06	14
EdgeSequel Sapphire .06	14
EdgeSequel Utopia .06	14
EdgeV-Taper HT .06	14
EdgeFile X7.06	14
EdgeEvolve .06	14

\*Click to Jump to Page
















## For EdgeOne Platinum

Discontinued Classic EdgeOne Gutta Percha Points		Discontinued Slim EdgeOne Gutta Percha Points		Recommended Gutta Percha Points Replacement: EdgeEndo .06	
			EOPVERGP		SGP1506-1
	8806387801153_1		EOPSMGP		SGP2006-1
	8806387801160_1		EOPPRIGP		SGP2506-1
	8806387801184_1		EOPMEDGP		SGP3506-1
	8806387801207_1		EOPLGGP		SGP4506-1
				Available in Assorted Packs	

## For EdgeOne Fire

Discontinued EdgeOne Fire Gutta Percha Points		Recommended Gutta Percha Points Replacement: EdgeOne Fire	
	EOFSMGP		EOFSMGP-1
	EOFPRIGP		EOFPRIGP-1
	EOFMEDGP		EOFMEDGP-1
	EOFLGGP		EOFLGGP-1
		Available in Assorted Packs	

















## For EdgeTaper, EdgeTaper Platinum, EdgeTaper Blaze Utopia & EdgeDiamond Utopia

Discontinued Gutta Percha Points: Classic EdgeTaper		Discontinued Gutta Percha Points: Slim EdgeTaper		Recommended Gutta Percha Points Replacement: EdgeTaper	
	8806387803034		ETF1GP		ETGPF1-1
	8806387803041		ETF2GP		ETGPF2-1
	8806387803058		ETF3GP		ETGPF3-1
	8806387803126		ETF4GP		ETGPF4-1
	8806387803133		ETF5GP		ETGPF5-1
				Available in Assorted Packs	








## For Universal .04, EdgeSequel Sapphire, EdgeSequel Utopia and V-Taper HT

Discontinued Gutta Percha Points: Universal .04		Discontinued Gutta Percha Points: EdgeSequel Sapphire		Discontinued Gutta Percha Points: Edge V-File HT		Recommended Gutta Percha Points Replacement: EdgeEndo .04	
	EE1504GP		ESS1504GP				SGP1504-1
	EE2004GP		ESS2004GP		EVHTGP2006		SGP2004-1
	EE2504GP		ESS2504GP		EVHTGP2506		SGP2504-1
	EE3004GP		ESS3004GP		EVHTGP3006		SGP3004-1
	EE3504GP		ESS3504GP		EVHTGP3506		SGP3504-1
	EE4004GP		ESS4004GP		EVHTGP4006		SGP4004-1

## For Universal .04, EdgeSequel Sapphire, EdgeSequel Utopia and V-Taper HT


Discontinued Gutta Percha Points: Universal .04		Discontinued Gutta Percha Points: EdgeSequel Sapphire		Discontinued Gutta Percha Points: Edge V-File HT		Recommended Gutta Percha Points Replacement: EdgeEndo .04	
	EE4504GP		ESS4504GP		EVHTGP4506		SGP4504-1
	EE5004GP		ESS5004GP		EVHTGP5006		SGP5004-1
	EE5504GP		ESS5504GP				TBD
	EE6004GP		ESS6004GP				TBD
	EE7004GP		ESS7004GP				TBD
	EE8004GP		ESS8004GP				TBD
						Available in Assorted Packs	













# For Universal .06, EdgeSequel Sapphire & EdgeSequel Utopia

Discontinued Gutta Percha Points: EdgeEndo .06		Discontinued Gutta Percha Points: EdgeSequel Sapphire		Recommended Gutta Percha Points Replacement: EdgeEndo .06	
	EE1506GP		ESS1506GP		SGP1506-1
	EE2006GP		ESS2006GP		SGP2006-1
	EE2506GP		ESS2506GP		SGP2506-1
	EE3006GP		ESS3006GP		SGP3006-1
	EE3506GP		ESS3506GP		SGP3506-1
	EE4006GP		ESS4006GP		SGP4006-1
	EE4506GP		ESS4506GP		SGP4506-1
	EE5006GP		ESS5006GP		SGP5006-1
















Available in Assorted Packs

## For EdgeFile X3 and EdgeFile X7









Discontinued Gutta Percha Points: DiaDent .04		Recommended Gutta Percha Points Replacement: EdgeFile X7.04	
			X71704GP-1
	GTP2004GP		X72004GP-1
	GTP2504GP		X72504GP-1
	GTP3004GP		X73004GP-1
	GTP3504GP		X73504GP-1
	GTP4004GP		X74004GP-1
			X74504GP-1
		Available in Assorted Packs	

Discontinued Gutta Percha Points: DiaDent .06		Recommended Gutta Percha Points Replacement: EdgeFile X7.06	
			X71706GP-1
	GTP2006GP		X72006GP-1
	GTP2506GP		X72506GP-1
	GTP3006GP		X73006GP-1
	GTP3506GP		X73506GP-1
	GTP4006GP		X74006GP-1
			X74506GP-1
		Available in Assorted Packs	
















## For EdgeOne Platinum

Discontinued Paper Points: Classic EdgeOne		Discontinued Paper Points: Slim EdgeOne		Recommended Paper Points Replacement: EdgeEndo .06	
	8806387821922_1		EOPVERPP		SPP1506-1
	8806387821861_1		EOPSMPP		SPP2006-1
	8806387821663_1		EOPPRIPP		SPP2506-1
	8806387821687_1		EOPMEDPP		SPP3506-1
	8806387821946_1		EOPLGPP		SPP4506-1
				Available in Assorted Packs	





























## For EdgeOne for EdgeOne Fire & EdgeOne Blaze Utopia

Discontinued Paper Points: EdgeOne Fire		Recommended Paper Points Replacement: EdgeOne Fire	
	EOFSMPP		EOFSMPP-1
	EOFPRIPP		EOFPRIPP-1
	EOFMEDPP		EOFMEDPP-1
	EOF LGPP		EOF LGPP-1
		Available in Assorted Packs	

















## For EdgeTaper, EdgeTaper Platinum, EdgeTaper Blaze Utopia & EdgeDiamond Utopia

Discontinued Paper Points: Classic EdgeTaper		Discontinued Paper Points: Slim EdgeTaper		Recommended Paper Points Replacement: EdgeTaper	
	8806387821977		ETF1PP		ETPPF1-1
	8806387821984		ETF2PP		ETPPF2-1
	8806387821991		ETF3PP		ETPPF3-1
	8806387822004		ETF4PP		ETPPF4-1
	8806387822011		ETF5PP		ETPPF5-1
				Available in Assorted Packs	






























## For Universal .04, EdgeSequel Sapphire, EdgeSequel Utopia, EdgeV-Taper, EdgeFile X7 & EdgeEvolve

Discontinued Paper Points: Universal .04		Discontinued Paper Points: EdgeSequel .04		Discontinued Paper Points: V-File HT		Discontinued Paper Points: DiaDent .04		Recommended Paper Points Replacement: EdgeEndo .04	
	EE1504PP		ESS1504PP						SPP1504-1
	EE2004PP		ESS2004PP		EVHTPP2006		GTP2004PP		SPP2004-1
	EE2504PP		ESS2504PP		EVHTPP2506		GTP2504PP		SPP2504-1
	EE3004PP		ESS3004PP		EVHTPP3006		GTP3004PP		SPP3004-1
	EE3504PP		ESS3504PP		EVHTPP3506		GTP3504PP		SPP3504-1
	EE4004PP		ESS4004PP		EVHTPP4006		GTP4004PP		SPP4004-1
									Available in Assorted Packs

## For Universal .04, EdgeSequel Sapphire, EdgeSequel Utopia, EdgeV-Taper, EdgeFile X7 & EdgeEvolve

Discontinued Paper Points: Universal .04		Discontinued Paper Points: EdgeSequel .04		Discontinued Paper Points: V-File HT		Recommended Paper Points Replacement: EdgeEndo .04	
	EE4504PP		ESS4504PP		EVHTPP4506		SPP4504-1
	EE5004PP		ESS5004PP		EVHTPP5006		SPP5004-1
	EE5504PP		ESS5504PP				TBD
	EE6004PP		ESS6004PP				TBD
	EE7004PP		ESS7004PP				TBD
	EE8004PP		ESS8004PP				TBD
						Available in Assorted Packs	

# For Universal .06, EdgeSequel Sapphire, EdgeSequel Utopia, EdgeV-Taper, EdgeFile X7 & EdgeEvolve

Discontinued Paper Points: Universal .06		Discontinued Paper Points: EdgeSequel .06		Discontinued Paper Points: DiaDent .06		Recommended Paper Points Replacement: EdgeEndo .06	
	EE1506PP		ESS1506PP				SPP1506-1
	EE2006PP		ESS2006PP		GTP2006PP		SPP2006-1
	EE2506PP		ESS2506PP		GTP2506PP		SPP2506-1
	EE3006PP		ESS3006PP		GTP3006PP		SPP3006-1
	EE3506PP		ESS3506PP		GTP3506PP		SPP3506-1
	EE4006PP		ESS4006PP		GTP4006PP		SPP4006-1
	EE4506PP		ESS4506PP				SPP4506-1
	EE5006PP		ESS5006PP				SPP5006-1
						Available in Assorted Packs	