


EDGE X7  
UTOPIA™

Utopia Alternative to EdgeFile® X7



**\$38.80**  
Pre-Sterilized 6 Pack

- Premium version of EdgeFile® X7
- More Efficient Cutting
- Constant .04 or .06 Taper
- 21, 25, 29mm Lengths
- Assorted Pack Available

EDGE DIAMOND  
UTOPIA™

Utopia Alternative to ProTaper Ultimate™




**\$38.80**  
Pre-Sterilized 6 Pack

- Rotary or Forward Reciprocation
- Continuously Changing Diamond Shaped Cross-section
- Minimally Invasive Design
- Progressive and Regressive Variable Taper
- 21mm & 25mm lengths
- Assorted Pack Available

**Our Next Big....Little Thing!**

EDGE TAPER BLAZE  
UTOPIA™

Utopia Alternative to ProTaper Gold® & EdgeTaper Platinum™




**\$38.80**  
Pre-Sterilized 6 Pack

- Premium Version Of EdgeTaper Platinum™
- Electropolished Blade For Smoother Navigation
- Improved Cutting Efficiency
- Variable Taper System
- 21, 25, 31mm Lengths
- Assorted Pack Available

EDGE SEQUEL  
UTOPIA™

Utopia Alternative to Vortex® and EndoSequence® & EdgeSequel Sapphire



**\$38.80**  
Pre-Sterilized 6 Pack

- Premium Version Of EdgeSequel Sapphire™
- Shorter Handle For Improved Posterior Access
- Constant .02, .04 Or .06 Taper
- #25/.08 19mm Also Available
- 21, 25, 31mm Lengths
- Assorted Pack Available

EDGE ONE BLAZE  
UTOPIA™


Utopia Alternative to WaveOne® Gold & EdgeOne Fire



**\$27.95**  
Pre-Sterilized 3 Pack

- Premium Version Of EdgeOne Fire™
- Our Greatest Single File Shaping System Yet
- Resilient To Unwinding
- Variable Taper System
- 21, 25, 31mm Lengths
- Assorted 4 Pack and R-Glide 3 Pack Available

EDGE FILE® X7



#1 Fan Favorite File System

**\$38.80**  
6 Pack

HEAT-TREATED FIREWIRE™ NiTi

- Universal file system
- 2x more cyclic fatigue resistant than other file systems
- Tips 17-45, 04 & 06 constant taper
- 21, 25, & 29mm
- Available Assorted Pack

EDGE V-TAPER HT™




VTaper® 2H & Endosequence® VT Scout™ Alternative

**\$38.80**  
6 Pack

HEAT-TREATED FIREWIRE™ NiTi

- Identical to SS White V-Taper® 2H
- Minimally Invasive
- Expansive size portfolio available
- 21, 25, 31mm

EDGE TAPER PLATINUM™



ProTaper® Gold Alternative

**\$38.80**  
6 Pack

HEAT-TREATED FIREWIRE™ NiTi

- 2x more resistant to cyclic fatigue\* compared to ProTaper® Gold.
- SX, S1, S2, F1, F2, F3, F4, F5
- 21, 25, & 31mm
- Available Assorted Pack

EDGE TAPER™



ProTaper® Alternative

**\$38.80**  
6 Pack

- 6x more resistant than ProTaper®
- SX, S1, S2, F1, F2, F3, F4, F5
- 21, 25, & 31mm
- Available Assorted Pack

EDGE SEQUEL SAPPHIRE™



Vortex® Blue and EndoSequence® CM Alternative

**\$27.95**  
4 Pack

HEAT-TREATED FIREWIRE™ NiTi

- 2x more resistant to cyclic fatigue\* compared to Vortex® Blue, 8x more resistant than Vortex® and EndoSequence®
- Constant taper
- Tips 04 15-90 taper, 06 15-50 taper
- 21, 25, & 31mm
- Available Assorted Pack

EDGE ONE PLATINUM™



WaveOne® Gold Alternative

**\$27.95**  
3 Pack

HEAT-TREATED FIREWIRE™ NiTi

- Reverse Reciprocating
- Triangular Cross-section
- Verifile, Small, Primary, Medium & Large
- 21, 25, & 31mm
- Available Assorted Pack

EDGE ONE FIRE™



WaveOne® Gold Alternative

**\$27.95**  
3 Pack

HEAT-TREATED FIREWIRE™ NiTi

- 5X more resistant to cyclic fatigue\* compared to WaveOne® Gold
- Reverse Reciprocating
- Parallelogram Cross-section
- Small, Primary, Medium & Large
- 21, 25, & 31mm
- Assorted 4 Pack and EdgeOne GlidePath 3 Pack Available

\*Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph: <https://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> ProTaper®, ProTaper Gold®, Vortex®, Vortex Blue®, ProTaper Ultimate™ and WaveOne® Gold are registered trademarks of Dentsply Tulsa Dental Specialties. EndoSequence®, Endosequence® CM, EndoSequence VT Scout™ are registered trademarks of Brasseler USA. V-Taper® 2H is a registered trademark of SS White.

ACCESS • SHAPING • IRRIGATION • OBTURATION • SEALING

EDGE SURGICAL TIP™



\$53.85  
Per Tip

EDGE TIP™



\$53.85  
Per Tip



Take Advantage of  
Our Newly Lowered  
Prices on Select Burs!



EDGE K-FILE™

\$5.60  
6 pack

- Sizes: 6-80
- Alternative to Tulsa Dental K-file®
- 21, 25, & 31mm Lengths



EDGE C-FILE™

\$6.45  
6 pack

- Sizes: 6, 8, 10, 12, 15
- 21 & 25mm

Also Available: EdgeHeadstoms - \$5.60 • EdgeK-File NiTi - \$8.25

EDGEAPEX HP™

\$1179.00

- Pre-programmed with top selling EdgeEndo files
- 10 fully customizable Memory Program slots
- Integrated Apex Locator during instrumentation
- Option to use Apex locator with handfile
- 360° rotation of contra angle



EDGE ICE™

Pulp Vitality  
Refrigerant Spray

\$17.15



EDGE MIX™

EDTA-based formula  
enhanced with  
Chlorhexidine for added  
bacterial properties

\$64.65



EDGE UTOPIA™  
BC SOLVENT

Dissolves Bioceramic  
Sealers and Calcium  
Hydroxide.

\$107.85

Also Available: EdgeLube, EDTA, Sodium Hypochlorite, Chlorhexidine



EDGE CORE™

Pink Gutta Percha Core  
Thermal Obturator

\$38.80  
5 pack & 1 verifier



EDGE FILL™

Plastic Carrier  
Thermal Obturator

\$38.80  
6 Pack

Also Available: Universal sized EdgeCore and EdgeFill 20pks Available - \$105.67

EDGEFLOW™

\$749.00

- Simple preset heating with the press of a button
- Uniform and smooth gutta-percha output without bubbles
- Dual batteries enable filling ~100 root canals
- Re-usable injection needles available in various sizes and lengths
- Use with EdgePellets™ - \$30.20




NEW!



EdgeEndo® Gutta Percha Points - \$13.61

Optimized Apical Fit with Coronal Venting


- Tight Apical Fit
- Extra Space In The Coronal Third To Facilitate The Upward Venting Of Sealer
- Minimizes The Risk Of Apical Extrusion



EdgeEndo® Paper Points

Universal Constant Taper

\$8.65



EdgeEndo® Gutta Points

Universal Constant Taper

\$14.05

EDGE BIO CERAMIC™  
SEALER



\$140.15

2g syringe +  
15 reduced waste tips  
#1 Fan Favorite sealer with  
zero shrinkage

EDGE BIO CERAMIC™  
THERMAL FLOW™



\$140.15

1.5g syringe +  
15 reduced waste tips  
Low viscosity ideal for  
warm obturation

EDGE UTOPIA™  
BioCeramic Sealer



\$78.65

2g syringe + 20 tips  
High pH of 12.5 for optimal  
antibacterial properties

EDGE UTOPIA™ RRM™



\$102.50

0.15g per capsule - 10 capsules per box  
Pre-Mixed Monodose Root Repair Material  
Dispensing Gun Sold Separately - \$29.10

EDGE TEMP™ PLUS



Hydraulic Temporary Restorative  
Alternative to Cavit™ by 3M

\$11.85

EDGEPACK™

\$589.00

- High precision temperature control
- 4 hour battery life with each charge
- Variety of Tapers Available in Tip Sizes 35-60
- Improved penetrability into narrow root canals
- Ensures smooth cutting of gutta percha
- Warm up within 0.2 seconds auto shut off after 10 seconds



Cavit™ is a registered trademark of 3M. Tulsa Dental K-File® is a registered trademark of DENTSPLY Tulsa Dental Specialties

Rev. D 09/2025



Product  
Evaluation  
Request

Fast & Free Shipping On All Orders  
Hassle Free Returns • Straightforward Pricing  
• Exceptional Customer Service



Full  
Product  
Portfolio



EDGEENDO®



Scan For  
Digital Copy