



EDGEENDO®

2026 RESIDENT PRICING PRODUCT CATALOG

The Endodontic Graduate Program from EdgeEndo® is a special purchase program for first and second year residents. We're excited to provide you this opportunity to work with our products at a special discounted price. You may use this opportunity during residency and until December 31st of your graduation year.

Participants must create an account if they do not already have one. All products must be ordered by phone 1-855-985-3636. Limits are shown for each product category and are available to Mix/Match.

[Click For More Information](#)

Created by an Endodontist for Endodontists

Established in 2012 by American Endodontist Charles Goodis, EdgeEndo® was founded with the vision to make root canal treatments safe, easy, and cost-effective overcoming the shortfalls of manufacturers offering outdated technology with a high price tag.

The solution – Innovation combined with high-quality and cost-effectiveness to support the growth of endodontic practices by saving chair time and delivering outstanding patient outcomes. In just ten years, EdgeEndo® revolutionized the endodontic industry in the USA becoming a leader in endodontic solutions!

Safe – Healthy Patients - Happy Dentists. EdgeEndo® solutions are developed with the patient at the center. Our FireWire heat treatment processes provides excellent resistance to cyclic fatigue and flexibility to carefully follow the canal, preserving canal anatomy and dentin and allowing for faster procedures. All for increased safety and productivity.

Easy – Changing habits is always a challenge for any dental professional. Our solutions are designed for ease of use and are compatible with most endodontic techniques. EdgeEndo® continues expanding its product offering to cover the entire endodontic procedure; diagnosis, instrumentation, obturation and restoration. No need for multiple sources and complexity. We make it easy for the dentists.

Cost-Effective – EdgeEndo® has proven that quality does not need to come with a premium price tag. Our goal is to make root canals accessible to anyone, without compromises.

EdgeEndo® remains steadfast in its mission to provide game-changing products and solutions

Innovation - EdgeEndo® is the inventor of FireWire and FireWire Blaze, a unique combination of heat treatment and cryogenic applications allowing NiTi instruments to deliver incredible flexibility and resistance to cyclic fatigue.

Quality – Thousands of clinicians trust EdgeEndo® and our mission to deliver a seamless and consistent experience with all our products. The control of manufacturing processes is our obsession, and our instruments are all manufactured in the USA, Switzerland, and France under strict quality controls.

**Revolutionize endodontic procedures with EdgeEndo®:
the partner dental professionals have been searching for!**




Table of Contents

EdgePower Suite

EdgeApex HP	4
EdgeFlow	6
EdgeApex	9
EdgePack	10

Diagnosis & Isolation

EdgeIce.	13
EdgeDam	13

Supplies

EndoRing	14
--------------------	----

Initial Access

EdgeTip.	15
EdgeBurs - Carbide	16
EdgeGates.	19
EdgeBroaches	20
EdgeHedstroms	21
EdgeTaper Retreat.	21
BC Solvent.	22

Glidepath

EdgeFind.	23
EdgeGlidepath	23
EdgeOneGlidepath	24
EdgeOne Blaze Utopia - RGlide.	24

Access Refinement

EdgeSurgical Tip	25
EdgeBurs - Diamond	26

Working Length Determination

EdgeK-File NiTi	30
EdgeK-File	31
EdgeC-File.	32

File Systems

Rotary • Constant Taper

EdgeX7 Utopia	33
EdgeFile X7	35
EdgeFile X3	37
EdgeSequel Utopia	39

EdgeSequel Sapphire.	42
EdgeEvolve	44

Rotary • Variable Taper

EdgeDiamond Utopia.	47
EdgeV-Taper HT	49
EdgeTaper Blaze Utopia	51
EdgeTaper Platinum.	53
EdgeTaper	55

Reciprocating • Variable Taper

EdgeOne Blaze Utopia	57
EdgeOne Fire	59

Reciprocating • Constant Taper

EdgeOne Platinum	61
EdgeFile X1	63

Irrigation

EdgeLube	65
EdgeGel	65
Irrigation Tips & Syringes.	66
Clensers	67

Sealers

EdgeBioCeramic Sealer.	68
EdgeUtopia BioCeramic Sealer	68
EdgeBioCeramic Thermal Flow	68
EdgeSeal.	69

Obturation

EdgeCore	69
Edge Fill	69
Points.	70

Intracanal Medicament & Repair

EdgeCaOH.	75
EdgeTemp Plus.	75
EdgeBioCeramic Retrofill & Repair.	75
EdgeUtopia RRM.	76

Starter Bundles 78

Category Limit Quantities. 79

EDGEAPEX HP™



- Built in Apex Locator with option to use as a standalone Apex locator with handfile
- Pre-programmed with top selling EdgeEndo files PLUS 10 fully customizable Memory Program slots
- Efficient brushless motor, low noise, long service life
- Cordless portable endo motor with combined length determination
- 360° rotation of contra angle
- Adopt real-time feedback technology and dynamic torque control, effectively preventing file separation
- Wireless Charging
- Torque Control
- Auto On/Off, Torque Reverse, Apical Reverse, Apical Slow Down
- Visual torque on file in use
- Max Torque 5.0 Ncm
- Change Dominant Hand
- Range: 100 - 2500 RPM

\$825.30

EDGEAPEX HP™

\$825.30

1 UNIT	EAPXHP
--------	--------

Accessories included in package:

- Measuring Wire
- Lip Hook (x2)
- File Clip (x4)
- Touch Probe (x2)
- Power Adapter
- Protective Silicon Cover (x2)
- Spray Nozzle
- O-Rings (x2)
- Disposable Sleeves (1 pack)

Warranty: 2 years for handpiece motor, charging base, and power adapter. 1 year for contra angle attachment. 6 months for accessories.



Accessories

MEASURING WIRE - \$13.51
1 PACK
EAPXMW



CONTRA ANGLE - \$150.99
1 PACK
EAPXCA



TOUCH PROBE - \$6.72
1 PACK
EAPXTP



FILE CLIP - \$45.26
5 PACK
EAPXFC



LIP HOOK - \$11.27
5 PACK
EAPXLH



CHARGING BASE - \$55.06
1 PACK
EAPXCB



EDGEFLOW™



- Wireless and electric
- Excellent gutta-percha fluidity, uniform and smooth gutta-percha output without bubbles.
- Precise and intuitive display of remaining gutta-percha
- Auto heating stop after timeout; Auto shutdown after timeout
- Charging base with battery
- Temperature monitoring design
- Dual batteries deliver longer endurance, enabling filling about 100 root canals.
- Re-usable Injection needles available in various sizes and lengths
- Preset temperatures: 100°C, 120°C, 150°C, 180°C, 200°C
- Injection Speeds: Slow, Medium, Fast

\$524.30

EDGEFLOW™

\$524.30

1 UNIT	EEPTO
--------	-------

Accessories included in package:

- Main Unit
- Charging Stand
- Power Adapter
- Battery (x2)
- Heat Insulation Hood
- Wrench of Needle
- Gutta-Percha Ring
- Instruction Manual
- Qualified Certification
- Packing List
- Injection Needles (x4)
 - E23G-M (x1)
 - E23G-L (x1)
 - E25G-L (x2)

Warranty: 1 year for device.



Accessories

OBTURATOR WRENCH - \$11.27
1 PACK
EEPTOW



EDGEFLOW™ BATTERY - \$22.58
1 PACK
EEPTOB



EDGEFLOW™ HEAT INSULATION HOOD - \$13.51
3 PACK
EEPTOHIH3



EdgeFlow™ Re-Useable Needle Tips

NEEDLE TIP E20G-S - \$25.62
1 PACK

EE20GSIN



NEEDLE TIP E23G-L - \$25.62
1 PACK

EE23GLIN



NEEDLE TIP E20G-M - \$25.62
1 PACK

EE20GMIN



NEEDLE TIP E23G-M - \$25.62
1 PACK

EE25GMIN



NEEDLE TIP E20G-L - \$25.62
1 PACK

EE20GLIN



NEEDLE TIP E23G-L - \$25.62
1 PACK

EE25GLIN



NEEDLE TIP E23G-M - \$25.62
1 PACK

EE23GMIN



EDGEPELLETS™

100 PACK - \$24.16

00810184450092



EDGEAPEX™

\$377.30

1 UNIT	EEALU
--------	-------



Accessories included in package:

- Base
- Measuring Wire
- File Clip (x4)
- Touch Probe (x2)
- Pulp Tester Probe (x2)
- Adapter
- Tester
- Enlarged display of apical foramen
- Instruction card for apical stop
- Apex Locator troubleshooting card

Warranty: 1 year for device. 6 months for accessories.

Apex Locator Accessories

MEASURING WIRE - \$13.51
1 PACK
EEALUMW



TESTER - \$11.27
1 PACK
EEALUT



LIP HOOK - \$11.27
5 PACK
EEALULH5



ENDO RULER - \$11.27
1 PACK
EEALUER



TOUCH PROBE - \$8.23
1 PACK
EEALUMF



FILE CLIP - \$45.26
5 PACK
EEALUFC5



EDGE PACK™



- Wireless Gutta Percha Obturation Device
- Ergonomic design and ambidextrous accessible
- Easy to use
- 0.35mm pre-bend tip
- Improved penetrability into narrow root canals
- High precision temperature control
- Smooth cutting of GP from high-efficiency conversion system
- Compact size for easy cleaning and storage
- Visual and auditory prompts for safe operation
- Preset temperatures: 150°C, 180°C, 200°C, 230°C
- Time related features: Warm up time within 0.2 seconds, Auto shut off after 10 seconds,
- 4 hour battery life with each charge

\$412.30

EDGEPACK™

\$412.30

1 UNIT	EEDPU
--------	-------

Accessories included in package:

- Charging Stand
- Power Adapter
- WP4004 Tip
- WP4504 Tip
- WP5506 Tip
- WP5508 Tip
- Heat Insulation Sleeve of Tips (x4)
- Instruction Manual and Qualified Certification

Warranty: 1 year for device and tips.



Accessories

EDGEPACK CHARGING STAND - \$45.26	
1 UNIT	EEDPUCS



EdgePack Tips

EDGEPACK 3504 TIP - \$37.70

1 PACK

EEDPUT3504



EDGEPACK 5508 TIP - \$37.70

1 PACK

EEDPUT5508



EDGEPACK 4004 TIP - \$37.70

1 PACK

EEDPUT4004



EDGEPACK 5508L TIP - \$37.70

1 PACK

EEDPUT5508L



EDGEPACK 4504 TIP - \$37.70

1 PACK

EEDPUT4504



EDGEPACK 5510 TIP - \$37.70

1 PACK


EEDPUT5510



EDGEPACK 5004 TIP - \$37.70

1 PACK

EEDPUT5004



EDGEPACK 6012 TIP - \$37.70

1 PACK

EEDPUT6012



EDGEPACK 5506 TIP - \$37.70

1 PACK

EEDPUT5506




EdgeIce™ - \$13.72	
6 OZ. SPRAY BOTTLE	818365020065

A quick way to check cold sensitivity and pulp vitality. Significantly colder than ethyl chloride or ice. EdgeIce™ features a directional nozzle for easy placement on micro applicator or cotton pellet.



EDGE DAM™

Latex 5x5	52 Pack - \$9.00
Medium	EDBLM5552
Heavy	EDBLH5552
Latex 6x6	36 Pack- \$9.00
Medium	EDBLM6636
Heavy	EDBLH6636



Latex Free 5x5	50 Pack - \$53.55
Medium	EDVLFM5550
Heavy	EDVLFH5550
Latex Free 6x6	50 Pack - \$71.68
Medium	EDVLFM6650
Heavy	EDVLFH6650



ENDORING®

EndoRing® II Kit with REDUCED GLARE Metal Ruler and e-Foam® Rotary HD Inserts \$43.08

1 KIT	ER2GSEE-S
-------	-----------

1 ring assembly & 6 silver HD foam inserts



EndoRing® II Kit with Metal Ruler \$36.16

1 KIT	ER2BEE-S
-------	----------

1 ring assembly, 3 blue foam inserts & 3 yellow foam inserts



Foam Refills 48 Pack \$36.16

Silver Rotary HD e-foam® Insert Refill	ERFSHDEE-S
--	------------

Blue e-foam® Insert Refill	ERFBEE-S
----------------------------	----------

Yellow e-foam® Insert Refill	ERFYEE-S
------------------------------	----------

Combo e-foam® Insert Refill (Yellow & blue) -	ERFEE-S
---	---------



EdgeTips™ 1-5 are diamond-coated for more efficiency. We don't use zirconium nitride. EdgeTips™ 6-8 are a special titanium alloy and their special geometries allow precision work in areas of limited access.

- Find calcified or hidden canals.
- Remove intracanal obstructions.
- Remove dental restorations and cores.
- Dislodge posts and broken instruments.
- Remove remove pulp stones.
- Unique shaft designs for better access and visibility.
- Length and diameters designed to progress down the root.

EdgeTip #1 - Diamond-coated ball tip, length 20mm, taper 5%.

ETBD



EdgeTip #2 - Diamond-coated steel tip 76µm, length 18mm, taper 5%

ET18D



EdgeTip #3 - Diamond-coated retreatment tip, 30µm, length 20mm, taper 5%

ET20D



EdgeTip #4 - Long Diamond-coated retreatment tip, 30µm, 40mm, taper 4%

ET40D4



EdgeTip #6 - Short Titanium- Niobium tip, length 15mm taper 4%

ED25S



EdgeTip #7 - Titanium- Niobium tip, length 20mm, taper 3%

ET25



EdgeTip #8 - Titanium- Niobium tip, 25mm, taper 3%

ET25L



OPERATIVE CARBIDES



ROUND



US NO.	ISO NO.	1/4	1/2	1	2	3	4
FG (19MM)	500 314 001 001	005	006	008	010	012	014
	10 PACK	\$13.72 FG1QTR10	\$13.72 FG1HLF10	\$13.72 FG110	\$13.72 FG210	\$13.72 FG310	\$13.72 FG410
FG SURG. (25MM)	500 316 001 001		006		010		014
	10 PACK	\$13.72 FG1HLF-SURG10		\$13.72 FG2SURG10		\$13.72 FG4SURG10	
RA (22MM)	500 204 001 001		006		010	012	014
	10 PACK		\$13.72 RA1HLF10		\$13.72 RA210	\$13.72 RA310	\$13.72 RA410



ROUND

US NO.	ISO NO.	5	6	7	8	10
FG (19MM)	500 314 001 001	016	018		023	
	10 PACK	\$13.72 FG510	\$13.72 FG610		\$13.72 FG810	
FG SURG. (25MM)	500 316 001 001		018		023	
	10 PACK		\$13.72 FG6SURG10		\$13.72 FG8SURG10	
RA (22MM)	500 204 001 001		018	021	023	027
	10 PACK		\$13.72 RA610	\$13.72 RA710	\$13.72 RA810	\$13.72 RA1010



INVERTED CONE



US NO.	ISO NO.	3/4	1
FG (19MM)	500 314 010 001	008	010
	10 PACK	\$13.72 FG3410	\$13.72 FG3510



AMALGAM PREP



US NO.	ISO NO.	1/4
FG (19MM)	500 314 234 006	009
	10 PACK	\$13.72 FG24510
FG SHORT (16MM)	500 313 234 006	009
	10 PACK	\$13.72 FG245SS10



PEAR

US NO.	ISO NO.	330	331	332
FG (19MM)	500 314 237 001	008	010	012
	10 PACK	\$13.72 FG33010	\$13.72 FG33110	\$13.72 FG33210
FG SHORT (16MM)	500 313 237 001	008		
	10 PACK	\$13.72 FG330SS10		



PEAR LONG HEAD

US NO.	ISO NO.	331L
FG (19MM)	500 314 238 006	010
	10 PACK	\$13.72 FG331L10



FLAT FISSURE CROSS CUT

US NO.	ISO NO.	556	557	558
FG (19MM)	500 314 107 007	008	010	012
	10 PACK	\$13.72 FG55610	\$13.72 FG55710	\$13.72 FG55810
FG SURG (25MM)	500 316 107 007		010	012
	10 PACK		\$13.72 FG557SURG10	\$13.72 FG558SURG10
FG SHORT (16MM)	500 313 107 007		010	
	10 PACK		\$13.72 FG557SS10	



FLAT FISSURE
CROSS CUT LONG HEAD

US NO.	ISO NO.	557L
FG (19MM)	500 314 110 007	010
	10 PACK	\$13.72 FG557L10



DOMED FISSURE PLAIN
CUT

US NO.	ISO NO.	1157	1158
FG (19MM)	500 314 137 006	010	012
	10 PACK	\$13.72 FG115710	\$13.72 FG115810



DOMED FISSURE CROSS
CUT

US NO.	ISO NO.	1557	1558
FG (19MM)	500 314 137 007	010	012
	10 PACK	\$13.72 FG155710	\$13.72 FG155810

SPECIALTY



151-Z

FIG NO.	151-Z
DIAMETER (1/10MM)	016
H/LENGTH (MM)	11.0
OVERALL LENGTH (MM)	23.0
5 PACK	\$19.76 151ZB5



ENDO ACCESS

FIG NO.	152
DIAMETER (1/10MM)	016
H/LENGTH (MM)	9.0
OVERALL LENGTH (MM)	23
5 PACK	\$19.76 152B5



LINDEMANN
CROSS CUT

US NO.	162
H/LENGTH (MM)	10.0
FG SURG. (25MM)	016
5 PACK	\$19.76 FGSURG162B5



D LINDEMANN
CROSS CUT

US NO.	162A
H/LENGTH (MM)	9.2
FG SURG. (25MM)	016
5 PACK	\$19.76 FGSURG162AB5



LINDEMANN

US NO.	161
H/LENGTH (MM)	10.0
FG SURG. (25MM)	016
5 PACK	\$19.76 FGSURG161B5

6 Pack	\$7.65	28 mm	32 mm
#1		EG0128	EG0132
#2		EG0228	EG0232
#3		EG0328	EG0332
#4		EG0428	EG0432
#5		EG0528	EG0532
#6		EG0628	EG0632



ASSORTED PACK		\$7.65
28MM	EGASST28	
32MM	EGASST32	



6 Pack \$4.47			
Size	Color	21mm	25mm
0	 PURPLE	EB4XF021	EB4XF025
1	 WHITE	EB3XF121	EB3XF125
2	 YELLOW	EB2XF221	EB2XF225
3	 RED	EB1XF321	EB1XF325
4	 BLUE	EB0XF421	EB0XF425
5	 GREEN	EBM0521	EBM0525
6	 BLACK	EBC0621	EBC0625

Assorted Pack 10 Pack \$4.47		
Size	21mm	25mm
0-6	EBASST0621	EBASST0625
INCLUDED: #0 (X1), #1 (X1), #2 (X2), #3 (X2), #4 (X2), #5 (X1), #6 (X1)		
2-4	EBASST2421	EBASST2425
INCLUDED: #2 (X3), #4 (X4), #5 (X3)		

EDGEHEDSTROMS™

Stainless Steel



6 Pack				\$3.36
Size	Color	21mm	25mm	31mm
15	☐ WHITE	HED1521	HED1525	HED1531
20	☐ YELLOW	HED2021	HED2025	HED2031
25	☐ RED	HED2521	HED2525	HED2531
30	☐ BLUE	HED3021	HED3025	HED3031
35	☐ GREEN	HED3521	HED3525	HED3531
40	☐ BLACK	HED4021	HED4025	HED4031
45	☐ WHITE	HED4521	HED4525	HED4531
50	☐ YELLOW	HED5021	HED5025	HED5031

EDGETAPER RETREAT™

Rotary Retreatment Files



EdgeTaper Retreat™ files are designed for the different needs of unfilling the coronal third, the middle third and the apical third before canal reshaping.

16 mm 6 Pack			\$21.34
D1	☐ WHITE	ETRD116	

18 mm 6 Pack			\$21.34
D2	☐ WHITE	ETRD218	

22 mm 6 Pack			\$21.34
D3	☐ WHITE	ETRD322	

Assorted 6 Pack			\$21.34
D1-D3	☐ WHITE	ETRASSTD123	

EDGE U UTOPIA™ BC SOLVENT

**EdgeUtopia™ BC Solvent is not just another solution,
it's *THE* solution for dissolving Bioceramic Sealers**

Dissolves Bioceramic Sealers and Calcium Hydroxide.

EdgeUtopia™ BC Solvent is engineered to dissolve the toughest bioceramic sealers and calcium hydroxide, outperforming other traditional solutions on the market. No more struggle with difficult to remove materials.

Tough on Sealer. Gentle on Teeth.


While other products fall short on dissolving tough materials without damaging sensitive structures, EdgeUtopia BC Solvent works precisely where it's needed most, while maintaining the integrity of the surrounding tooth tissues.

EdgeUtopia™ BC Solvent - Preserving What Matters.

\$86.28

240mL Bottle

10818207024937



EDGE FIND™

Non Heat-Treated Endonol™ NiTi

6 Pack		\$19.58		
Size	Color	21mm	25mm	31mm
P1	 PURPLE	EFNDP121	EFNDP131	EFNDP125
P2	 WHITE	EFNDP221	EFNDP231	EFNDP225
P3	 YELLOW	EFNDP321	EFNDP331	EFNDP325



Assorted 3 Pack	\$19.58	21mm	25mm	31mm
P1, P2, P3		EFNDASST21	EFNDASST25	EFNDASST31

EDGE GLIDE PATH™

Heat-Treated FireWire™ NiTi

Rotary GlidePath



3 pack		\$15.37
Size	Color	
21mm	 WHITE	EGP0221
25mm	 WHITE	EGP0225
31mm	 WHITE	EGP0231

The EdgeGlidePath™ heat treated FireWire NiTi file features the flexibility and strength that EdgeEndo® is known for, replacing as many as 8 instruments to create a glidepath faster, with far less effort and with a 1mm maximum flute diameter. Replaces #15, 20, 25, 30, 35 hand files in a step-back technique and a #1, 2 & 3 Gates Glidden. Used exactly like Tulsa Dental ProGlider®* at substantially lower cost.

*ProGlider® is a registered trademark of Dentsply Tulsa Dental Specialties.

EDGEONE GLIDEPATH™

For EdgeOne Fire™



Reciprocating GlidePath

3 pack		\$15.37
Size	Color	
21mm	☐ WHITE	EOFGP21
25mm	☐ WHITE	EOFGP25
31mm	☐ WHITE	EOFGP31

EDGEONE GLIDEPATH™

For EdgeOne Platinum™



Reciprocating GlidePath

3 pack		\$15.37
Size	Color	
21mm	☐ WHITE	EOGP21
25mm	☐ WHITE	EOGP25
31mm	☐ WHITE	EOGP31

EDGEONE BLAZE UTOPIA™

R-Glide



Reciprocating GlidePath

3 pack		\$15.37
Size	Color	
21mm	☐ WHITE	EOBURG1421SU
25mm	☐ WHITE	EOBURG1425SU
31mm	☐ WHITE	EOBURG1431SU

*ProGlider® is a registered trademark of Dentsply Tulsa Dental Specialties.

EDGE SURGICAL TIP™

\$37.70

EdgeSurgical Tips™ are designed for use in endodontic microsurgery. Made with a diamond coating, these instruments cut dentin faster without need for repeated ultrasonics. The dentin surface is left smooth yet microscopically rough which results in better adaptation of filling materials.

- Diamond coating for efficient cutting of dentin.
- Irrigation strategically placed at the tip for control.
- Tip angles and lengths designed for access to difficult areas.
- We don't use zirconium nitride.

EdgeSurgical Tip #1

Diamond-coated universal tip 30µm, length 3mm, taper 9%

AS3D1



EdgeSurgical Tip #2

Diamond-coated universal tip 30µm, taper 9%

AS3D2



EdgeSurgical Tip #3

Left-oriented retro surgery tip, diamond-coated 30µm, length 5mm, taper 7%

P15LD



EdgeSurgical Tip #4

Left-oriented tip, diamond-coated 30µm, length 3mm, taper 10%

ASLD



EdgeSurgical Tip #5

Right-oriented tip, diamond-coated 30µm, length 6mm, taper 7%

P15RD



EdgeSurgical Tip #6

Right-oriented tip, diamond-coated 30µm, length 3mm, taper 10%

ASRD



DIAMONDS



001 | 801

Round

Head Diameter 1/10mm	010	012	014	016	018	021	023	025
L1: Head Length L1mm	1.0	1.2	1.4	1.6	1.8	2.1	2.3	2.5
L2: OAL mm	19.0	19.0	19.0	19.0	19.0	19.0	19.0	19.0
ISO								
806 314 001 514...								025 \$34.44 001025F
806 314 001 524...	010 \$34.44 001010M	012 \$34.44 001012M	014 \$34.44 001014M	016 \$34.44 001016M	018 \$34.44 001018M		023 \$34.44 001023M	
806 314 001 534...	010 \$34.44 001010C	012 \$34.44 001012C	014 \$34.44 001014C	016 \$34.44 001016C	018 \$34.44 001018C	021 \$34.44 001021C	023 \$34.44 001023C	
806 314 001 544...					018 \$34.44 001018XC		023 \$34.44 001023XC	



032 | 813

Diabolo

Head Diameter 1/10mm	014
L1: Head Length L1mm	1.5
L2: OAL mm	19.0
ISO	
806 314 032 534...	014 \$34.44 032014C



544 | 845KR Short

Round Edge Taper

Head Diameter 1/10mm	019	026
L1: Head Length L1mm	4.0	4.0
L2: OAL mm	18.0	18.0
ISO	Short	Short
806 313 544 514...		026 \$34.44 544026SF
806 313 544 534...	019 \$34.44 544019SC	



546 | 847KR

Round Edge Taper



Head Diameter 1/10mm	017
L1: Head Length L1mm	8.0
L2: OAL mm	20.0
ISO	Short
806 314 546 514...	017 \$34.44 546017SF
806 314 546 524...	017 \$34.44 546017SM
806 314 546 534...	017 \$34.44 546017SC
806 314 546 544...	017 \$34.44 546017SXC



546 | 847KR

Round Edge Taper



Head Diameter 1/10mm	016
L1: Head Length L1mm	8.0
L2: OAL mm	22.0
ISO	
806 313 546 534...	016 \$34.44 546016C



197 | 855 Short

Round End Taper



Head Diameter 1/10mm	017	024
L1: Head Length L1mm	6.0	7.0
L2: OAL mm	18.5	18.5
ISO	Short	Short
806 314 197 544...	017 \$34.44 197017SXC	024 \$34.44 197024SXC



197 | 855 XShort

Round End Taper



Head Diameter 1/10mm	016
L1: Head Length L1mm	6.0
L2: OAL mm	16.5
ISO	XShort
806 313 197 534...	016 \$34.44 197016SHC



198 | 856

Round End Taper

							
Head Diameter 1/10mm	012	014	015	016	018	019	021
L1: Head Length L1mm	8.0	8.0	8.0	8.0	8.0	8.0	8.0
L2: OAL mm	22.0	22.0	22.0	22.0	22.0	22.0	22.0
ISO							
806 316 198 514...		 014 \$34.44 198014F		 016 \$34.44 198016F	 018 \$34.44 198018F		
806 316 198 524...	 012 \$23.97 198012M	 014 \$34.44 198014M		 016 \$23.97 198016M	 018 \$23.97 198018M		
806 316 198 534...	 012 \$34.44 198012C	 014 \$34.44 198014C	 015 \$34.44 198015C	 016 \$34.44 198016C	 018 \$34.44 198018C		 021 \$34.44 198021C
806 316 198 544...	 012 \$34.44 198012XC	 014 \$34.44 198014XC	 015 \$34.44 198015XC	 016 \$34.44 198016XC	 018 \$34.44 198018XC	 019 \$23.97 198019XC	 021 \$34.44 198021XC



198 | 856

Round End Taper

		
Head Diameter 1/10mm	022	025
L1: Head Length L1mm	8.0	8.0
L2: OAL mm	22.0	22.0
ISO		
806 315 198 534...	 022 \$34.44 (468-0698)	 025 \$34.44 (468-0700)
806 315 198 544...	 022 \$34.44 (468-0699)	 025 \$34.44 (468-0701)



198 | 856 Short

Round End Taper

	
Head Diameter 1/10mm	023
L1: Head Length L1mm	8.5
L2: OAL mm	20.0
ISO	Short
806 314 198 534...	 023 \$34.44 198023SC



198 | 856L Long Round End Taper

Head Diameter 1/10mm	016	018	020
L1: Head Length L1mm	9.0	9.0	9.0
L2: OAL mm	22.0	22.0	22.0
ISO			
806 315 198L 514...	016 \$34.44 198L016F		
806 315 198L 534...	016 \$34.44 198L016C	018 \$34.44 198L018C	020 \$34.44 198L020C
806 315 198L 544...	016 \$34.44 198L016XC	018 \$34.44 198L018XC	



199 | 850 Round End Taper

Head Diameter 1/10mm	010	012	014	016	018
L1: Head Length L1mm	10.0	10.0	10.0	10.0	10.0
L2: OAL mm	24.0	24.0	24.0	24.0	24.0
ISO					
806 316 199 524...					018 \$34.44 199018M
806 316 199 534...	010 \$34.44 199010C	012 \$34.44 199012C	014 \$34.44 199014C	016 \$34.44 199016C	018 \$25.01 199018C
806 316 199 544...				016 \$34.44 199016XC	018 \$34.44 199018XC



199 | 850 Short Round End Taper

Head Diameter 1/10mm	014	016
L1: Head Length L1mm	10.0	10.0
L2: OAL mm	22.0	22.0
ISO	Short	Short
806 315 199 514...	014 \$34.44 199014SXC	016 \$34.44 199016SXC



263 | 369 Bud Slender

Head Diameter 1/10mm	025
L1: Head Length L1mm	5.5
L2: OAL mm	20.0
ISO	
806 314 263 534...	025 \$34.44 263025C

EDGEK-FILE™ NITI



Non Heat-Treated NiTi

Hand File

Alternative to Tulsa Dental K-File®

6 Pack				\$5.10
Size	Color	21mm	25mm	31mm
15	 WHITE	NITIKF1521	NITIKF1525	NITIKF1531
20	 YELLOW	NITIKF2021	NITIKF2025	NITIKF2031
25	 RED	NITIKF2521	NITIKF2525	NITIKF2531
30	 BLUE	NITIKF3021	NITIKF3025	NITIKF3031
35	 GREEN	NITIKF3521	NITIKF3525	NITIKF3531
40	 BLACK	NITIKF4021	NITIKF4025	NITIKF4031
45	 WHITE	NITIKF4521	NITIKF4525	NITIKF4531
50	 YELLOW	NITIKF5021	NITIKF5025	NITIKF5031
55	 RED	NITIKF5521	NITIKF5525	NITIKF5531
60	 BLUE	NITIKF6021	NITIKF6025	NITIKF6031

Tulsa Dental K-File® is a registered trademark of Dentsply Tulsa Dental Specialties



6 Pack					\$5.10
Size	Color	21mm	25mm	31mm	
6	 PINK	KF621	KF625	KF631	
8	 GREY	KF821	KF825	KF831	
10	 PURPLE	KF1021	KF1025	KF1031	
15	 WHITE	KF1521	KF1525	KF1531	
20	 YELLOW	KF2021	KF2025	KF2031	
25	 RED	KF2521	KF2525	KF2531	
30	 BLUE	KF3021	KF3025	KF3031	
35	 GREEN	KF3521	KF3525	KF3531	
40	 BLACK	KF4021	KF4025	KF4031	
45	 WHITE	KF4521	KF4525	KF4531	
50	 YELLOW	KF5021	KF5025	KF5031	
55	 RED	KF5521	KF5525	KF5531	
60	 BLUE	KF6021	KF6025	KF6031	
70	 GREEN	KF7021	KF7025	KF7031	
80	 BLACK	KF8021	KF8025	KF8031	

EDGE C-FILE™





Stainless Steel

Hand File

6 Pack		\$3.87	
Size	Color	21mm	25mm
6	 PINK	CF621	CF625
8	 GREY	CF821	CF825
10	 PURPLE	CF1021	CF1025
12	 BROWN	CF1221	CF1225
15	 WHITE	CF1521	CF1525

Enhanced Technology. Remarkably Beautiful.

- Enhanced version of EdgeFile® X7
- Features proprietary FireWire Blaze Heat Treatment for better cutting efficiency
- Constant Taper File
- Pre-Sterilized



- Parabolic Shapped Cross-Section
- Available in 3 lengths: 21, 25 & 29 mm
- File Size: 04 & 06
- Maximum flute diameter 1 mm
- Speed: 500 RPM
- Torque: 4.0 Ncm



High Cutting Efficiency



Pre-Bendable



High Strength

EDGE X7



UTOPIA™



Pre-Sterilized 6 Pack					Points		Obturators	
Size	Color	21mm	25mm	29mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
17.04	PURPLE	X7U170421SE	X7U170425SE	X7U170429SE	SPP1504-1	X71704GP-1	N/A	N/A
20.04	YELLOW	X7U170621SE	X7U170625SE	X7U170629SE	SPP2004-1	X72004GP-1	X72006C	X72006O
25.04	RED	X7U200421SE	X7U200425SE	X7U200429SE	SPP2504-1	X72504GP-1	X72506C	X72506O
30.04	BLUE	X7U200621SE	X7U200625SE	X7U200629SE	SPP3004-1	X73004GP-1	X73006C	X73006O
35.04	GREEN	X7U250421SE	X7U250425SE	X7U250429SE	SPP3504-1	X73504GP-1	X73506C	X73506O
40.04	BLACK	X7U250621SE	X7U250625SE	X7U250629SE	SPP4004-1	X74004GP-1	X74006C	X74006O
45.04	WHITE	X7U300421SE	X7U300425SE	X7U30049SE	SPP4504-1	X74504GP-1	X74506C	X74506O
17.06	PURPLE	X7U300621SE	X7U300625SE	X7U300629SE	SPP1506-1	X71706GP-1	N/A	N/A
20.06	YELLOW	X7U350421SE	X7U350425SE	X7U350429SE	SPP2006-1	X72006GP-1	X72006C	X72006O
25.06	RED	X7U350621SE	X7U350625SE	X7U350629SE	SPP2506-1	X72506GP-1	X72506C	X72506O
30.06	BLUE	X7U400421SE	X7U400425SE	X7U400429SE	SPP3006-1	X73006GP-1	X73006C	X73006O
35.06	GREEN	X7U400621SE	X7U400625SE	X7U40069SE	SPP3506-1	X73506GP-1	X73506C	X73506O
40.06	BLACK	X7U450421SE	X7U450425SE	X7U450429SE	SPP4006-1	X74006GP-1	X74006C	X74006O
45.06	WHITE	X7U450621SE	X7U450625SE	X7U450629SE	SPP4506-1	X74506GP-1	X74506C	X74506O

Pair EdgeDiamond Utopia™ with:



BioCeramic Sealer 2g

818365020072



Assorted Packs .04 - \$21.34/ 6 Pack

(17, 20, 25, 25, 30, 35)     

Length	Files	Paper Points	Gutta Percha
21MM	X7UASST0421SE	SPP04AST-1	X7GP04AST-1
25MM	X7UASST0425SE		

Assorted Packs .06 - \$21.34/ 6 Pack

(17, 20, 25, 25, 30, 35)     

Length	Files	Paper Points	Gutta Percha
21MM	X7UASST0621SE	SPP06AST-1	X7GP06AST-1
25MM	X7UASST0625SE		


EDGEFILE[®] X7



\$21.34
6 Pack

Unmatched Strength and Flexibility


- Universal File System
- #1 Fan Favorite
- FireWire NiTi Heat Treatment provides super elasticity with no shape memory
- Constant Tapered Rotary Files
- Non-Sterile
- Tighter pitch and flutes which deflect more torsional forces along the length of the file
- Non-cutting tip with flutes extending to 0.25mm—offering precision and safety in one design



- Parabolic Shapped Cross-Section
- Available in 3 lengths: 21, 25 & 29 mm
- Constant Taper: .04 and .06
- Rotary files
- Speed: 450 RPM
- Torque: 2.9 NCm



Capable of
90° Curves



Pre-Bendable



High
Strength

EDGEFILE® X7



6 Pack					Points		Obturators	
Size	Color	21mm	25mm	29mm	Gutta Percha	Paper Points	EdgeCore	EdgeFill
17.04	WHITE	X7170421RF	X7170425RF	X7170429RF	SPP1504-1	X71704GP-1	N/A	N/A
20.04	YELLOW	X7170621RF	X7170625RF	X7170629RF	SPP2004-1	X72004GP-1	X72006C	X72006O
25.04	RED	X7200421RF	X7200425RF	X7200429RF	SPP2504-1	X72504GP-1	X72506C	X72506O
30.04	BLUE	X7200621RF	X7200625RF	X7200629RF	SPP3004-1	X73004GP-1	X73006C	X73006O
35.04	GREEN	X7250421RF	X7250425RF	X7250429RF	SPP3504-1	X73504GP-1	X73506C	X73506O
40.04	BLACK	X7250621RF	X7250625RF	X7250629RF	SPP4004-1	X74004GP-1	X74006C	X74006O
45.04	WHITE	X7300421RF	X7300425RF	X7300429RF	SPP4504-1	X74504GP-1	X74506C	X74506O
17.06	WHITE	X7300621RF	X7300625RF	X7300629RF	SPP1506-1	X71706GP-1	N/A	N/A
20.06	YELLOW	X7350421RF	X7350425RF	X7350429RF	SPP2006-1	X72006GP-1	X72006C	X72006O
25.06	RED	X7350621RF	X7350625RF	X7350629RF	SPP2506-1	X72506GP-1	X72506C	X72506O
30.06	BLUE	X7400421RF	X7400425RF	X7400429RF	SPP3006-1	X73006GP-1	X73006C	X73006O
35.06	GREEN	X7400621RF	X7400625RF	X7400629RF	SPP3506-1	X73506GP-1	X73506C	X73506O
40.06	BLACK	X7450421RF	X7450425RF	X7450429RF	SPP4006-1	X74006GP-1	X74006C	X74006O
45.06	WHITE	X7450621RF	X7450625RF	X7450629RF	SPP4506-1	X74506GP-1	X74506C	X74506O

Pair EdgeFile® X7 with:

EDGEBIOCERAMIC™
SEALER

818365020072



Assorted Pack .04 - \$21.34 6 Pack

(17, 20, 25, 30, 35, 40)

Length	Files	Paper Points	Gutta Percha
21mm	ASSTX70421	SPP04AST-1	X7GP04AST-1
25mm	ASSTX70425		
29mm	ASSTX70429		

Assorted Pack .06 - \$21.34 6 Pack

(17, 20, 25, 30, 35, 40)

Length	Files	Paper Points	Gutta Percha
21mm	ASSTX70621	SPP06AST-1	X7GP06AST-1
25mm	ASSTX70625		
29mm	ASSTX70629		



EDGEFILE[®] X3



\$21.34
6 Pack

Engineered for Strength with Flexible Precision

- Designed with a Safe-Core[™] Parabolic Cross- Section
- .04 or .06 constant taper file until flute diameter reaches 1.0 then taper is 0
- Safe unwinding feature
- Super elasticity. No Shape Memory
- FireWire NiTi Heat Treated
- Variable Tapered Rotary Files
- Features its own dedicated Glidepath system
- Non-Sterile
- Alternative to ProTaper[®] Next¹ systems



- Parabolic shaped cross-section
- Available in 3 lengths: 21, 25 & 29 mm
- Speed 300-500 RPM
- Torque 2.94 Ncm



Highly Flexible




Pre-Bendable



High Strength

1) ProTaper[®] Next is a registered trademark of Dentsply Sirona.

EDGEFILE® X3



6 Pack					Points		Obturators	
Size/Taper	Color	21mm	25mm	29mm	Paper Points	Gutta Points	EdgeCore	EdgeFill
N1 17.06	PURPLE	X3170621RF	X3170625RF	X3170629RF	SPP1506-1	X71706GP-1	N/A	N/A
N2 17.04	WHITE	X3170421RF	X3170425RF	X3170429RF	SPP1504-1	X71704GP-1	N/A	N/A
C1 20.06	YELLOW	X3200621RF	X3200625RF	X3200629RF	SPP2006-1	X72006GP-1	X72006C	X72006O
C2 25.06	RED	X3250621RF	X3250625RF	X3250629RF	SPP2506-1	X72506GP-1	X72506C	X72506O
C3 30.06	BLUE	X3300621RF	X3300625RF	X3300629RF	SPP3006-1	X73006GP-1	X73006C	X73006O
C4 40.06	BLACK	X3400621RF	X3400625RF	X3400629RF	SPP4006-1	X74006GP-1	X74006C	X74006O

Pair EdgeFile® X3 with:

EDGEBIOCERAMIC™
SEALER


818365020072



EdgeVerify™ for EdgeFile® X3 6 Pack	
Size	Item Code
20	X32006V
25	X32506V
30	X33006V
40	X34006V
50	X35006V

Enhanced Technology. Remarkably Beautiful.

- Features proprietary FireWire Blaze Blue Heat Treatment
- Improved cutting efficiency compared to EdgeSequel Sapphire™
- 12mm handle length for better posterior access
- Constant Taper File
- Pre-Sterilized
- Alternative to EndoSequence® and Vortex®



- Triangular Cross Section
- Available in 3 lengths: 21, 25 & 31mm w/ 19mm orifice opener
- Speed: 800-1000 RPM
- Torque: 1.5 Ncm



Flexible



High
Cutting
Efficiency



Electropolished
Finish



High
Strength


*19mm orifice opener available in #25/.08

Vortex® is a registered trademark of Dentsply Tulsa Dental Specialties. EndoSequence® is a registered trademark of Brasseler USA

EDGESEQUEL



UTOPIA™



19mm Pre-Sterilized 6 Pack								
Size/Taper	Color							
25.08	RED	ES250819SE						
Pre-Sterilized 6 Pack					Points		Obturators	
Size/Taper	Color	21mm	25mm	31mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
10.02	PURPLE	ES100221SE	ES100225SE	N/A	N/A	N/A	N/A	N/A
15.02	WHITE	ES150221SE	ES150225SE	N/A	N/A	N/A	N/A	N/A
20.02	YELLOW	ES200221SE	ES200225SE	N/A	N/A	N/A	N/A	N/A
30.02	BLUE	ES300221SE	ES300425SE	N/A	N/A	N/A	N/A	N/A
10.04	PURPLE	ES100421SE	ES100425SE	N/A	N/A	N/A	N/A	N/A
15.04	WHITE	ES150421SE	ES150425SE	ES150431SE	SPP1504-1	SGP1504-1	N/A	N/A
20.04	YELLOW	ES200421SE	ES200421SE	N/A	SPP2004-1	SGP2004-1	ES2006C	ES2006O
25.04	RED	ES250421SE	ES250425SE	ES250431SE	SPP2504-1	SGP2504-1	ES2506C	ES2506O
30.04	BLUE	ES300421SE	ES300425SE	ES300431SE	SPP3004-1	SGP3004-1	ES3006C	ES3006O
35.04	GREEN	ES350421SE	ES350425SE	N/A	SPP3504-1	SGP3504-1	ES3506C	ES3506O
40.04	BLACK	ES400421SE	ES400425SE	ES400431SE	SPP4004-1	SGP4004-1	ES4006C	ES4006O
50.04	YELLOW	ES500421SE	ES500425SE	ES500431SE	SPP5004-1	SGP5004-1	ES5006C	ES5006O
15.06	WHITE	ES150621SE	ES150625SE	N/A	SPP1506-1	SGP1506-1	N/A	N/A
25.06	RED	ES250621SE	ES250625SE	ES250631SE	SPP2506-1	SGP2506-1	ES2506C	ES2506O
30.06	BLUE	ES300621SE	ES300625SE	ES300631SE	SPP3006-1	SGP3006-1	ES3006C	ES3006O
35.06	GREEN	ES350621SE	ES350625SE	ES350631SE	SPP3506-1	SGP3506-1	ES3506C	ES3506O
40.06	BLACK	ES400621SE	ES400625SE	N/A	SPP4006-1	SGP4006-1	ES4006C	ES4006O

Assorted Pack - \$15.37/ 3 Pack (15.04, 25.04, 30.04)   	
21mm	ASST15300421SE
25mm	ASST15300425SE
31mm	ASST15300431SE

Assorted Pack - \$15.37/ 3 Pack (15.04, 25.04, 25.06)   	
21mm	ASST15250621SE
25mm	ASST15250625SE
31mm	ASST15250631SE

Assorted Pack - \$21.34/ 6 Pack (15.04, 20.04, 25.04, 30.04, 35.04 & 25.08 19mm)      	
21mm	V2ASST0421SE
25mm	V2ASST0425SE
31mm	V2ASST0431SE

Pair EdgeSequel Utopia™ with:

EDGE  UTOPIA™
BioCeramic Sealer 2g
818365020072




EDGESEQUEL SAPPHIRE™



\$15.37
4 Pack

Bodybuilder Strength. Gymnast Flexibility.

- Excellent Cyclic Fatigue Compared To Other Leading Files*
- FireWire NiTi Heat Treated
- Constant Tapered Files
- Non-Sterile
- Alternative to Vortex®, VortexBlue®, and EndoSequence®



Cyclic Fatigue Testing - Average Time to Failure in Seconds

EDGESEQUEL
SAPPHIRE™

633Seconds

Vortex Blue®

299.9 Seconds

EndoSequence®

80.9 Seconds

*Vortex®, Vortex® Blue, are registered trademarks of Dentsply Tulsa Dental Specialties. EndoSequence® is a registered trademark of Brassler USA *Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph:<http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit www.edgeendo.com/matter-of-temperature



4 Pack					Points		Obturators	
Size	Color	21mm	25mm	31mm	Paper Points	Gutta Points	EdgeCore	EdgeFill
25.08	RED	ES250821	ES250825	ES250831	SPP1504-1	SGP1504-1	N/A	N/A
15.04	WHITE	ES150421	ES150425	ES150431	SPP2004-1	SGP2004-1	ES2006C	ES2006O
20.04	YELLOW	ES150621	ES150625	ES150631	SPP2504-1	SGP2504-1	ES2506C	ES2506O
25.04	RED	ES200421	ES200425	ES200431	SPP3004-1	SGP3004-1	ES3006C	ES3006O
30.04	BLUE	ES200621	ES200625	ES200631	SPP3504-1	SGP3504-1	ES3506C	ES3506O
35.04	GREEN	ES250421	ES250425	ES250431	SPP4004-1	SGP4004-1	ES4006C	ES4006O
40.04	BLACK	ES250621	ES250625	ES250631	SPP4504-1	SGP4504-1	ES4506C	ES4506O
45.04	WHITE	ES300421	ES300425	ES300431	SPP5004-1	SGP5004-1	ES5006C	ES5006O
50.04	YELLOW	ES300621	ES300625	ES300631	ESS5504GP	ESS5504PP	ES5506C	ES5506O
55.04	RED	ES350421	ES350425	ES350431	ESS6004GP	ESS6004PP	ES6006C	ES6006O
60.04	BLUE	ES350621	ES350625	ES350631	ESS7004GP	ESS7004PP	ES7006C	ES7006O
70.04	GREEN	ES400421	ES400425	ES400431	ESS8004GP	ESS8004PP	ES8006C	ES8006O
80.04	BLACK	ES400621	ES400625	ES400631	N/A	N/A	ES9006C	ES9006O
90.04	WHITE	ES450421	ES450425	ES450431	SPP1506-1	SGP1506-1	N/A	N/A
15.06	WHITE	ES450621	ES450625	ES450631	SPP2006-1	SGP2006-1	ES2006C	ES2006O
20.06	YELLOW	ES500421	ES500425	ES500431	SPP2506-1	SGP2506-1	ES2506C	ES2506O
25.06	RED	ES500621	ES500625	ES500631	SPP3006-1	SGP3006-1	ES3006C	ES3006O
30.06	BLUE	ES550421	ES550425	ES550431	SPP3506-1	SGP3506-1	ES3506C	ES3506O
35.06	GREEN	ES600421	ES600425	ES600431	SPP4006-1	SGP4006-1	ES4006C	ES4006O
40.06	BLACK	ES700421	ES700425	ES700431	SPP4506-1	SGP4506-1	ES4506C	ES4506O
45.06	WHITE	ES800421	ES800425	ES800431	SPP5006-1	SGP5006-1	ES5006C	ES5006O
50.06	YELLOW	ES900421	ES900425	ES900431	SPP5006-1	SGP5006-1	ES5006C	ES5006O

Pair EdgeSequel Sapphire[™] with:

EDGEBIOCERAMIC[™]
SEALER

818365020072



Assorted Packs .04 - \$15.37 4 Pack	
Sizes: (15, 20, 25, 30)	
21mm	ASST15300421
25mm	ASST15300425
31mm	ASST15300431
Sizes: (35, 40, 45, 50)	
21mm	ASST35500421
25mm	ASST35500425
31mm	ASST35500431

Assorted Packs .06 - \$15.37 4 Pack	
Sizes: (15, 20, 25, 30)	
21mm	ASST15300621
25mm	ASST15300625
31mm	ASST15300631
Sizes: (35, 40, 45, 50)	
21mm	ASST35500621
25mm	ASST35500625
31mm	ASST35500631

EDGE EVOLVE™



\$21.34
6 Pack

The Evolution of the Revolution

- Minimal dentin with a conservative flute and taper. All EdgeEvolve files, even our largest 40.10, only have a max diameter of 1mm. FireWire NiTi Heat Treated
- Extremely Flexible Constant Taper File with a non-cutting tip
- Available in .04, .06, .08, .10 Tapers
- Alternative to GTX®¹




- Triangular shaped cross-section
- Available in 3 lengths: 21, 25 & 29mm
- Speed: 300-500 RPM
- Torque: 2.94 Ncm



Extremely Flexible




Pre-Bendable



High Strength

1) GTX® is a registered trademark of Dentsply Tulsa Dental Specialties.



6 Pack					Points		Obturators	
Size/Taper	Color	21mm	25mm	29mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
17.04	WHITE	EE170421	EE170425	EE170429	SPP1504-1	SGP1504-1	N/A	N/A
17.06	WHITE	EE170621	EE170625	EE170629	SPP1506-1	SGP1506-1	N/A	N/A
17.08	WHITE	EE170821	EE170825	EE170829	N/A	N/A	N/A	N/A
17.10	WHITE	EE171021	EE171025	EE171029	N/A	N/A	N/A	N/A
20.04	YELLOW	EE200421	EE200425	EE200429	SPP2004-1	SGP2004-1	EE2006C	EE2006O
20.06	YELLOW	EE200621	EE200625	EE200629	SPP2006-1	SGP2006-1	EE2006C	EE2006O
20.08	YELLOW	EE200821	EE200825	EE200829	N/A	N/A	EE2006C	EE2006O
20.10	YELLOW	EE201021	EE201025	EE201029	N/A	N/A	EE2006C	EE2006O
25.04	RED	EE250421	EE250425	EE250429	SPP2504-1	SGP2504-1	EE2506C	EE2506O
25.06	RED	EE250621	EE250625	EE250629	SPP2506-1	SGP2506-1	EE2506C	EE2506O
25.08	RED	EE250821	EE250825	EE250829	N/A	N/A	EE2506C	EE2506O
25.10	RED	EE251021	EE251025	EE251029	N/A	N/A	EE2506C	EE2506O
30.04	BLUE	EE300421	EE300425	EE300429	SPP3004-1	SGP3004-1	EE3006C	EE3006O
30.06	BLUE	EE300621	EE300625	EE300629	SPP3006-1	SGP3006-1	EE3006C	EE3006O
30.08	BLUE	EE300821	EE300825	EE300829	N/A	N/A	EE3006C	EE3006O
30.10	BLUE	EE301021	EE301025	EE301029	N/A	N/A	EE3006C	EE3006O
35.04	GREEN	EE350421	EE350425	EE350429	SPP3504-1	SGP3504-1	EE3506C	EE3506O
35.06	GREEN	EE350621	EE350625	EE350629	SPP3506-1	SGP3506-1	EE3506C	EE3506O
35.08	GREEN	EE350821	EE350825	EE350829	N/A	N/A	EE3506C	EE3506O
35.10	GREEN	EE351021	EE351025	EE351029	N/A	N/A	EE3506C	EE3506O
40.04	BLACK	EE400421	EE400425	EE400429	SPP4004-1	SGP4004-1	EE4006C	EE4006O
40.06	BLACK	EE400621	EE400625	EE400629	SPP4006-1	SGP4006-1	EE4006C	EE4006O
40.08	BLACK	EE400821	EE400825	EE400829	N/A	N/A	EE4006C	EE4006O
40.10	BLACK	EE401021	EE401025	EE401029	N/A	N/A	EE4006C	EE4006O

Assorted Pack 20 4 Pack - \$21.34 (20.04, 20.06, 20.08, 20.10) 	
21mm	EEASST2021
25mm	EEASST2025
29mm	EEASST2029

Assorted Pack 25 4 Pack - \$21.34 (25.04, 25.06, 25.08, 25.10) 	
21mm	EEASST2521
25mm	EEASST2525
29mm	EEASST2529

Assorted Pack 30 4 Pack - \$21.34 (30.04, 30.06, 30.08, 30.10) 	
21mm	EEASST3021
25mm	EEASST3025
29mm	EEASST3029

Assorted Pack 40 4 Pack - \$21.34 (40.04, 40.06, 40.08, 40.10) 	
21mm	EEASST4021
25mm	EEASST4025
29mm	EEASST4029

Pair EdgeEvolve™ with:

EDGE BIO CERAMIC™
SEALER


818365020072



\$21.34
6 Pack

Elegance in Innovation, Strength in Versatility

- Minimally Invasive Progressive/Regressive Variable Tapered File System
- Ability to perform in **both rotary or forward reciprocation** settings
- Rotary Motion: Delivers maximum cutting efficiency and debris removal
- Forward Reciprocation Motion: Maximizes safety while still allowing for efficient cutting
- Unique geometry, innovative cutting angles and optimized FireWire™ Blaze Heat-Treatment** allows files to navigate curvatures with ease while maintaining dentin preservation in the coronal third
- Alternative to ProTaper Ultimate™* & Minimally Invasive Alternative to ProTaper™*



- Diamond Cross-Section
- Available in 21mm and 25mm lengths
- Progressive/Regressive Variable Taper File
- Pre-Sterilized



Incredibly Strong



Very Versatile



Minimally Invasive

*ProTaper Ultimate™ and ProTaper™ are registered trademarks of Dentsply Sirona. **Slider is not heat-treated Best Product Rating Dental Product Shopper Vol 19 No 2

EDGE DIAMOND



UTOPIA™



Best Product Rating Dental Product Shopper Vol 19 No 2

21mm Pre-Sterilized 6 Pack					EdgeTaper™ Points		EdgeTaper™ Obturators	
	Size	Taper	Color	Item Code	Paper Points	Gutta Points	EdgeCore	EdgeFill
SLIDER*	16	02v	PURPLE	EDU1602V21	N/A	N/A	N/A	N/A
SHAPER	20	04v	WHITE	EDU2004V21	N/A	N/A	N/A	N/A
F1	20	07v	YELLOW	EDU2007V21	ETPPF1-1	ETGPF1-1	ET2006C	ET2006O
F2	25	07v	RED	EDU2507V21	ETPPF2-1	ETGPF2-1	ET2506C	ET2506O
F3	30	07v	BLUE	EDU3007V21	ETPPF3-1	ETGPF3-1	ET3006C	ET3006O
FX	35	07v	GREEN	EDU3507V21	ETPPF4-1	ETGPF4-1	ET4006C	ET4006O

25mm Pre-Sterilized 6 Pack					EdgeTaper™ Points		EdgeTaper™ Obturators	
	Size	Taper	Color	Item Code	Paper Points	Gutta Points	EdgeCore	EdgeFill
SLIDER*	16	02v	PURPLE	EDU1602V25	N/A	N/A	N/A	N/A
SHAPER	20	04v	WHITE	EDU2004V25	N/A	N/A	N/A	N/A
F1	20	07v	YELLOW	EDU2007V25	ETPPF1-1	ETGPF1-1	ET2006C	ET2006O
F2	25	07v	RED	EDU2507V25	ETPPF2-1	ETGPF2-1	ET2506C	ET2506O
F3	30	07v	BLUE	EDU3007V25	ETPPF3-1	ETGPF3-1	ET3006C	ET3006O
FX	35	07v	GREEN	EDU3507V25	ETPPF4-1	ETGPF4-1	ET4006C	ET4006O

*Slider is not Heat Treated

Assorted Pack - \$21.34/ 6 Pack	
21mm	EDUASST21V
25mm	EDUASST25V

Pair EdgeDiamond Utopia™ with:

EDGE UTOPIA™

BioCeramic Sealer 2g

818365020072






\$21.34
6 Pack

Strong. Advanced. Flexible. Efficient.

- Designed with a Safe-Core™ Parabolic Cross- Section
- The Safe-Core™ Parabolic Cross- Section combines the attributes of being a highly efficient and flexible instrument while being resistant to fracture, minimizing your stress and resulting in a more reliable procedure.
- FireWire NiTi Heat Treated
- Variable Tapered Rotary Files
- Features its own dedicated Glidepath system
- Non-Sterile
- Alternative to DC Taper®¹



- Parabolic shaped cross-section
- Available in 3 lengths: 21, 25 & 31 mm
- Speed 300-400 RPM
- Torque - 4.0- 5.2 Ncm



Highly
Flexible



Pre-Bendable



High
Strength

1) DC Taper® is a registered trademark of SS White Dental®.

Rotary • Variable Taper • File System

EDGE V-TAPER^{HT}



Orifice Opener - 6 Pack		
Size/Taper	Color	17mm
25.09	GOLD	VFGP250917

Glidepaths 6 Pack					
Size/Taper	Color	21mm	25mm	29mm	31mm
13.03 GP	PINK	VFGP130321	VFGP130325	N/A	VFGP130331
14.03 GP	WHITE	VFGP140321	VFGP140325	VFGP140329	N/A
17.04 GP	WHITE	VFGP170421	VFGP170425	N/A	VFGP170431

6 Pack					Points		Obturators	
Size/Taper	Color	21mm	25mm	31mm	Paper Points	Gutta Points	EdgeCore	EdgeFill
20.06	YELLOW	VF200621	VF200625	VF200631	EVHTPP2006	EVHTGP2006	VF2006C	VF2006O
25.06	RED	VF250521	VF250525	VF250531	EVHTPP2506	EVHTGP2506	VF2506C	VF2506O
30.06	BLUE	VF220721	VF220725	VF220731	EVHTPP3006	EVHTGP3006	VF3006C	VF3006O
35.06	GREEN	VF250821	VF250825	VF250831	EVHTPP3506	EVHTGP3506	VF3506C	VF3506O
40.06	BLACK	VF300621	VF300625	VF300631	EVHTPP4006	EVHTGP4006	VF4006C	VF4006O
45.06	WHITE	VF300821	VF300825	VF300831	EVHTPP4506	EVHTGP4506	VF4506C	VF4506O
50.06	YELLOW	VF350621	VF350625	VF350631	EVHTPP5006	EVHTGP5006	VF5006C	VF5006O
22.07	YELLOW	VF400621	VF400625	VF400631	EVHTPP2207	EVHTGP2207	VF2006C	VF2006O
25.08	RED	VF450621	VF450625	VF450631	EVHTPP2508	EVHTGP2508	VF2506C	VF2506O
30.08	BLUE	VF500621	VF500625	VF500631	EVHTPP3008	EVHTGP3008	VF3006C	VF3006O

Pair EdgeV-Taper HTTM with:

EDGE BIO CERAMICTM
SEALER

818365020072



**EdgeVerifyTM for
EdgeV-Taper HTTM 6 Pack**

Size	Item Code
20	VF2006V
25	VF2506V
30	VF3006V
35	VF3506V
40	VF4006V
45	VF4506V
50	VF5006V


EDGETAPER BLAZE UTOPIA™



Enhanced Technology. Remarkably Beautiful.

\$21.34
6 Pack

- Features proprietary FireWire Blaze Heat Treatment
- Electropolished for better cutting efficiency
- Shorter handle length for better posterior access
- Non-cutting tip easily glides to the working length
- Pre-Sterilized
- Alternative to ProTaper Gold®
- Optimized version of EdgeTaper Platinum™



- Triangular Shaped Cross-Section
- Available in 3 lengths: 21, 25 & 31 mm
- Speed: 300 RPM
- Torque: 5.1 Ncm (SX&S1), 2.5 Ncm (S2&F1), 3.1 Ncm (F2&F3)

*ProTaper® Gold is a registered trademark of Dentsply Tulsa Dental Specialties



High
Cutting
Efficiency



High
Strength



Flexible

Rotary • Variable Taper • File System

EDGETAPER BLAZE



UTOPIA™



Pre-Sterilized 6 Pack		
Size	Color	19mm
SX	 YELLOW	ETBUSX19SU

Pre-Sterilized 6 Pack					EdgeTaper™ Points		EdgeTaper™ Obturators	
Size	Color	21mm	25mm	31mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
S1	 PURPLE	ETBUS121SU	ETBUS125SU	ETBUS131SU	N/A	N/A	N/A	N/A
S2	 WHITE	ETBUS221SU	ETBUS225SU	ETBUS231SU	N/A	N/A	N/A	N/A
F1	 YELLOW	ETBUF121SU	ETBUF125SU	ETBUF131SU	ETPPF1-1	ETGPF1-1	ET2006C	ET2006O
F2	 RED	ETBUF221SU	ETBUF225SU	ETBUF231SU	ETPPF2-1	ETGPF2-1	ET2506C	ET2506O
F3	 BLUE	ETBUF321SU	ETBUF325SU	ETBUF331SU	ETPPF3-1	ETGPF3-1	ET3006C	ET3006O

Assorted Packs - \$21.34/ 6 Pack (SX, S1, S2, F1, F2, F3)     	
21mm	ETBUA21SU
25mm	ETBUA25SU
31mm	ETBUA31SU

Pair EdgeTaper Blaze Utopia™ with:



BioCeramic Sealer 2g

818365020072



EDGE TAPER


PLATINUM™



\$21.34
6 Pack

Edge-ified Sequence with Enhanced Strength and Flexibility

- 2x The Cyclic Fatigue as ProTaper Gold®
- 6x The Cyclic Fatigue as ProTaper®
- FireWire NiTi Heat Treated
- Variable Tapered Files
- Non-Sterile



- Convex Triangular Cross- Section
- Available in 3 lengths: 21, 25 & 31 mm
- Speed: 300-400 RPM
- Torque: 4.0-5.2 Ncm



Highly Flexible




Pre-Bendable



High Strength

Rotary • Variable Taper • File System

EDGE TAPER PLATINUM™



6 Pack		
Size	Color	19mm
SX	 YELLOW	ETSX19HT

6 Pack					Points		Obturators	
Size	Color	21mm	25mm	31mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
S1	 PURPLE	ETS121HT	ETS125HT	ETS131HT	N/A			
S2	 WHITE	ETS221HT	ETS225HT	ETS231HT	N/A			
F1	 YELLOW	ETS121HT	ETF125HT	ETF131HT	ETPPF1-1	ETGPF1-1	ET2006C	ET2006O
F2	 RED	ETF221HT	ETF225HT	ETF231HT	ETPPF2-1	ETGPF2-1	ET2506C	ET2506O
F3	 BLUE	ETF321HT	ETF325HT	ETF331HT	ETPPF3-1	ETGPF3-1	ET3006C	ET3006O
F4	 BLACK	ETF421HT	ETF425HT	ETF431HT	ETPPF4-1	ETGPF4-1	ET4006C	ET4006O
F5	 YELLOW	ETF521HT	ETF525HT	ETF531HT	ETPPF5-1	ETGPF5-1	ET5006C	ET5006O

Assorted Packs - \$21.34				
Files      			Paper Points	Gutta Percha
6 Pack: (SX, S1, S2, F1, F2, F3)			    	
21mm	ETP21MM		ETPPAST	ETGPAST
25mm	ETP25MM			
31mm	ETP31MM			

Pair EdgeTaper Platinum™ with:

EDGE BIO CERAMIC™
SEALER

818365020072




EDGETAPER™



\$21.34
6 Pack

Edge-ified Sequence with Enhanced Strength

- 6x The Cyclic Fatigue as ProTaper®¹
- Variable Tapered Files
- Non-Sterile
- Non-Heat Treated



- Convex Triangular Cross- Section
- Available in 3 lengths: 21, 25 & 31 mm
- Speed: 300-400 RPM
- Torque: 4.0-5.2 Ncm



Pre-Bendable



High Strength

Cyclic Fatigue Testing - Average Time to Failure in Seconds

EDGETAPER
PLATINUM™

677 Seconds

ProTaper Gold®

312 Seconds

EDGETAPER™

127 Seconds

ProTaper®

119 Seconds

1) ProTaper® is a registered trademark of Dentsply Tulsa Dental Specialties. *Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph:<http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit www.edgeendo.com/matter-of-temperature

Rotary • Variable Taper • File System

EDGE TAPER™

Non-Heat Treated



6 Pack		
Size	Color	19mm
SX	 YELLOW	ETXS19

6 Pack					Points		Obturators	
Size	Color	21mm	25mm	31mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
S1	 PURPLE	ETS121	ETS125	ETS131	N/A			
S2	 WHITE	ETS221	ETS225	ETS231	N/A			
F1	 YELLOW	ETF121	ETF125	ETF131	ETPPF1-1	ETGPF1-1	ET2006C	ET2006O
F2	 RED	ETF221	ETF225	ETF231	ETPPF2-1	ETGPF2-1	ET2506C	ET2506O
F3	 BLUE	ETF321	ETF325	ETF331	ETPPF3-1	ETGPF3-1	ET3006C	ET3006O
F4	 BLACK	ETF421	ETF425	ETF431	ETPPF4-1	ETGPF4-1	ET4006C	ET4006O
F5	 YELLOW	ETF521	ETF525	ETF531	ETPPF5-1	ETGPF5-1	ET5006C	ET5006O

EDGE BIO CERAMIC™
SEALER



Assorted Packs - \$21.34			
Files      		Paper Points	Gutta Percha
6 Pack: (SX, S1, S2, F1, F2, F3)		    	
21mm	ETASST21	ETPPAST	ETGPAST
25mm	ETASST25		
31mm	ETASST31		


EDGEONE BLAZE UTOPIA™



Enhanced Technology. Remarkably Beautiful.

\$15.37
3 Pack

- Premier Single File Shaping Reciprocating System
- Features proprietary FireWire Blaze Heat Treatment
- Electropolished for better cutting efficiency
- Shorter handle length for better posterior access
- Variable Taper File
- Pre-Sterilized
- Alternative to WaveOne Gold® and EdgeOne Fire™



- S-Shaped Cross-Section - For Increased Debris Removal
- Available in 3 lengths: 21, 25 & 31 mm
- CCW 170°/ CW 60°



Flexible



High
Strength



High
Cutting
Efficiency

VOL. 18 NO. 4 Dental Product Shopper. "WaveOne® Gold is a registered trademarks of Dentsply Tulsa Dental Specialties.

EDGEONE BLAZE



UTOPIA™



VOL. 18 NO. 4 Dental Product Shopper.



Pre-Sterilized 3 Pack					Points		Obturators	
Size	Color	21mm	25mm	31mm	Paper Points	Gutta Points	EdgeCore	EdgeFill
R-GLIDE®	☐ WHITE	EOBURG1421SU	EOBURG1425SU	EOBURG1431SU	N/A		N/A	
SM R20	☒ YELLOW	EOBU2021SU	EOBU2025SU	EOBU2031SU	SPP2006-1	EORGP20-1	EOPSM20C	EOPSM20O
PRI R25	☒ RED	EOBU2521SU	EOBU2525SU	EOBU2531SU	SPP2506-1	EORGP25-1	EOPPRI25C	EOPPRI25O
MED R35	☒ GREEN	EOBU3521SU	EOBU3525SU	EOBU3531SU	SPP3506-1	EORGP35-1	EOPMED35C	EOPMED35O
LG R45	☐ WHITE	EOBU4521SU	EOBU4525SU	EOBU4531SU	SPP4506-1	EORGP45-1	EOPLG45C	EOPLG45O

Assorted Pack* - \$20.71/ 4 Pack
(R20, R25, R35, R45) ☒ ☒ ☒ ☐

21MM	EOBUA21SU
25MM	EOBUA25SU
31MM	EOBUA31SU

Pair EdgeOne Blaze Utopia™ with:

*R-Glide Sold Separately

EDGE UTOPIA™

BioCeramic Sealer 2g

818365020072




EDGEONE FIRE™



\$15.37
3 Pack

The Standard in Reciprocating Technology

- 5x the Cyclic Fatigue as WaveOne® Gold
- FireWire NiTi Heat Treated
- Variable Tapered Reverse Reciprocating Files
- Non-Sterile
- Alternative to WaveOne® Gold



- Parallelogram shaped cross-section
- Available in 3 lengths: 21, 25 & 31 mm
- Designed to shape canals in reverse-reciprocating motion CCW 150°/CW 60° or WaveOne® Gold¹ settings



Highly Flexible



Pre-Bendable



High Strength

Cyclic Fatigue Testing - Average Time to Failure in Seconds

EDGEONE
FIRE™

2844 Seconds

WaveOne® Gold

486 Seconds

WaveOne® Gold is a registered trademarks of Dentsply Tulsa Dental Specialties. *Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph: <http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit www.edgeendo.com/matter-of-temperature

EDGEONE FIRE™



3 Pack				
Size	Color	21mm	25mm	31mm
GLIDEPATH	☐ WHITE	EOFGP21	EOFGP25	EOFGP31

3 Pack					Points		Obturators	
Size	Color	21mm	25mm	31mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
SM 20	☒ YELLOW	EOFSM21	EOFSM25	EOFSM31	EOFSMPP-1	EOFSMGP-1	EOPSM20C	EOPSM20O
PRI 25	☒ RED	EOFPRI21	EOFPRI25	EOFPRI31	EOFPRIPP-1	EOFPRIGP-1	EOPPRI25C	EOPPRI25O
MED 35	☒ GREEN	EOFMED21	EOFMED25	EOFMED31	EOFMEDPP-1	EOFMEDGP-1	EOPMED35C	EOPMED35O
LG 45	☐ WHITE	EOFLG21	EOFLG25	EOFLG31	EOFLGGPP-1	EOFLGGP-1	EOPLG45C	EOPLG45O

Pair EdgeOne Fire™ with:

EDGE BIO CERAMIC™ SEALER

818365020072



Assorted Packs - \$20.30			
Files ☒ ☒ ☒ ☐		Paper Points	Gutta Percha
4 Pack (SM, PRI, MED, LG)		☒ ☒ ☒ ☐	
21mm	ASSTEOF21	EOFPAST-1	EOFGPAST-1
25mm	ASSTEOF25		
31mm	ASSTEOF31		

Constant Taper. Reliable Reciprocation.

- Constant taper reciprocating system
- Verifile dictates sequence
- Verifile immediately following 10 and 15 handfiles
- Non-Sterile
- Alternative to WaveOne® Gold¹



- Triangular shaped cross-section
- Available in 3 lengths: 21, 25 & 31 mm
- Designed to shape canals in reverse-reciprocating motion
- CCW 150°/CW 60° or WaveOne® Gold¹ settings



Highly Flexible



Pre-Bendable



High Strength

¹WaveOne® Gold is a registered trademarks of Dentsply Tulsa Dental Specialties.



3 Pack				
Size	Color	21mm	25mm	31mm
GLIDEPATH	☐ WHITE	EOGP21	EOGP25	EOGP31

3 Pack					Points		Obturators	
Size	Color	21mm	25mm	31mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
VERIFILE	☐ WHITE	EOPVER21	EOPVER25	EOPVER31	N/A	N/A	N/A	N/A
SMALL	☒ YELLOW	EOPSM21	EOPSM25	EOPSM31	SPP2006-1	X72006GP-1	EOPSM20C	EOPSM20O
PRIMARY	☒ RED	EOPPRI21	EOPPRI25	EOPPRI31	SPP2506-1	X72506GP-1	EOPPRI25C	EOPPRI25O
MEDIUM	☒ GREEN	EOPMED21	EOPMED25	EOPMED31	SPP3506-1	X73506GP-1	EOPMED35C	EOPMED35O
LARGE	☐ WHITE	EOPLG21	EOPLG25	EOPLG31	SPP4506-1	X74506GP-1	EOPLG45C	EOPLG45O

Pair EdgeOne Platinum[™] with:

EDGEBIOCERAMIC[™]
SEALER



818365020072



Assorted Pack - \$25.22			
Files ☐ ☒ ☒ ☒ ☐		Paper Points	Gutta Percha
5 Pack (VER, SM, PRI, MED, LG)		☒ ☒ ☒ ☐	
21mm	ASSTEOP21	EOFPPAST-1	EOFGPAST-1
25mm	ASSTEOP25		
31mm	ASSTEOP31		

Constant Taper. Consistent Reciprocation.

- Tighter pitch and flutes which deflect more torsional forces along the length of file due to the parabolic cross-section
- All X1 files (20.06 small, 25.06 primary, 40.06 large) have a constant 6% taper until the flute diameter is 1.0 mm then it goes to 0 taper
- FireWire NiTi Heat Treated
- Short handle 11mm
- Non-Sterile
- Alternative to WaveOne®¹



- Parabolic Shaped Cross-Section
- Non - Cutting Tip
- Available in 3 lengths: 21, 25 & 31 mm
- Constant Tapered Reciprocating file system should use the same motor settings as WaveOne®¹



Highly Flexible



Pre-Bendable



High Strength

¹) WaveOne® is a registered trademark of Dentsply Sirona.

EDGEFILE® X1



3 Pack					Points		Obturators	
Size/Taper	Color	21mm	25mm	29mm	Gutta Percha	Paper Points	EdgeCore	EdgeFill
20.06	 YELLOW	X1200621RF	X1200625RF	X1200629RF	SPP2006-1	X72006GP-1	X72006C	X72006O
25.06	 RED	X1250621RF	X1250625RF	X1250629RF	SPP2506-1	X72506GP-1	X72506C	X72506O
40.06	 BLACK	X1400621RF	X1400625RF	X1400629RF	SPP4006-1	X74006GP-1	X74006C	X74006O

Pair EdgeFile® X1 with:

EDGE BIO CERAMIC™
SEALER

818365020072



Assorted Pack - \$21.34/ 3 Pack
(20, 25, 40)



21mm	ASSTX121
25mm	ASSTX125
29mm	ASSTX129

EdgeLube™ - \$9.52 (17% EDTA SOLUTION)	
100ML BOTTLE	8806387892083



Presence of calcification in the canal. Narrow canals on posterior teeth. Smear layer removal.

Removal of inorganic substances by Chelation (EDTA-Ca²⁺). Easier instrumentation by softening dentin canal wall. Produces a canal wall free of smear layer.

EdgeGel™ - \$11.24 (19% EDTA GEL)	
5.4ML BOTTLE	8806387892076



Root canal preparation with files and reamers. Presence of calcification in the canal. Narrow canals in posterior teeth. Smear layer removal.

Removal of inorganic substances by chelation (EDTA-Ca²⁺). Easier filling with lubricating action. Produces a canal wall free of smear layer.

Irrigation Tips & Syringes

Edge® Irrigation Needle Tips - \$8.56 100 Pack		
23ga. Blue		10818365020260
27ga. Yellow		10818365020277
30ga. Purple		10818365020284



Edge® Irrigation Probe Tips - \$19.68 100 Pack		
23ga. Blue		10818365020291
27ga. Yellow		10818365020307
30ga. Purple		10818365020314



Edge® Pre-Assembled Irrigation Syringe with 1" Notched Tip - \$18.04 100 Pack		
23ga. Blue		00818365021062
27ga. Yellow		00818365021055
30ga. Purple		00818365021079

Edge® Luer-Lock Syringe - \$10.28 100 Pack		
3CC		10818365020383

2% Chlorhexidine Gluconate - \$25.80	
480ML BOTTLE	818365020041



17% EDTA - \$25.80	
480ML BOTTLE	818365020058



5% Sodium Hypochlorite- \$28.40	
480ML BOTTLE	818365020034



EDGEMIX™

\$51.72	
16OZ BOTTLE	18365020416



EDTA-based formula enhanced with Chlorhexidine for added anti-bacterial properties. Easy to use. 2-3 mL per canal. In instances where calcifications or a lot of debris is present, the higher end dosage should be utilized.

- Clears 27% more dentin tubules than EDTA
- Kills 99.99% bacteria in 10 seconds
- Does not form a precipitant when mixed with NaOCl

EDGE BIO CERAMIC™ SEALERS

EDGE BIO CERAMIC™ SEALER

2g syringe + 15 disposable tips

Cold Obturation

EdgeBioCeramic Sealer™ trusted in over **50 million cases worldwide** since 2008, EdgeBioCeramic Sealer™ offers proven, reliable performance. Easy-to-use, premixed, and resin-free, it's hydrophilic with zero shrinkage. The reduced waste tip preserves material by 62%, maximizing efficiency and cost savings.

Features:

- Radiopacity: High
- Backed by over 17 years of clinical use and extensive studies
- Ideal viscosity for single cone hydraulic condensation
- Composition: Tricalcium Silicate, Dicalcium Silicate, Zirconium Oxide, Tantalum Oxide

\$112.12

Sealer Sku: 818365020072
Disposable Tips Sku: EEBCT15

\$15.44/ 15 Tips



EDGE UTOPIA™ BioCeramic Sealer

2g syringe + 20 disposable tips

Warm & Cold Obturation


EdgeUtopia™ Sealer with its superior flow makes it **optimal for both warm or cold obturation techniques**. Its high pH of 12.5 offers powerful antibacterial properties, keeping the canal infection-free. Designed for moist environments, its hydrophilic nature ensures a secure seal in wet conditions.

Features:

- Radiopacity: Higher
- Our most cost effective Bioceramic sealer
- Ideal for both cold and warm obturation techniques
- Composition: Tricalcium silicate, dicalcium silicate, tricalcium aluminate, calcium oxide, zirconium oxide, silicon dioxide, polyethylene glycol and iron oxide

Sealer Sku: EUBCSEAL2G

\$62.92



EDGE BIO CERAMIC™ THERMALFLOW™

1.5g syringe + 15 disposable tips

Warm Obturation

EdgeBioCeramic™ ThermalFlow™ is an **optimized version** of EdgeBioCeramic Sealer™ with the same resin-free, hydrophilic, and zero shrinkage benefits, but with **lower viscosity for enhanced warm obturation techniques**. The reduced waste tip preserves material by 62%, maximizing efficiency and cost savings.

Features:

- Radiopacity: Highest
- Ideal for Warm Obturation with enhanced flowability
- Lowest viscosity and highest radiopacity
- Composition: Tricalcium Silicate, Dicalcium Silicate, Zirconium Oxide, Tantalum Oxide

\$112.12

Sealer Sku: 818365020232
Disposable Tips Sku: EEBCT15

\$15.44/ 15 Tips



EDGESEAL™



13.5g Dual Syringe

Resin-based root canal sealer designed to work with various obturation methods. Provides excellent biocompatibility. Easy to mix paste. Hermetic sealing ability. Non-staining to teeth. Insoluble in tissue fluids. Good radiopacity.

Working Time: 35 minutes at 23°C (73°F)

Setting Time: 45 minutes at 37°C (99°F)

Film Thickness: 3.3μ

\$25.88

818365020010

EDGECORE™

Pink Gutta Percha Core
Thermal Obturator



5 Pack & 1 Verifier

EDGEFILL™

Plastic Carrier
Thermal Obturator



6 Pack

20 Pack Obturators - \$59.97

Size	EdgeCore	EdgeFill
20	EE2020C	EE2020O
25	EE2520C	EE2520O
30	EE3020C	EE3020O
35	EE3520C	EE3520O
40	EE4020C	EE4020O
45	EE4520C	EE4520O
50	EE5020C	EE5020O

NEW EdgeEndo® Gutta Percha: Optimized Performance, Seamless Results Optimized For Use with Our Most Popular NiTi File Systems!



Seamless Integration

Tailored to work effortlessly with EdgeEndo file systems for streamlined endodontic procedures.



Proprietary Stiffening Process

Provides improved handling, less tip crumpling, and improved displacement of sealer.



Optimized Apical Fit with Coronal Venting

Designed for a tight apical fit with extra space in the coronal third to facilitate the upward venting of sealer thus minimizing the risk of apical extrusion. Additionally, this design helps prevent gutta percha from “hanging up” in longer canals promoting more efficient obturation.

Gutta Percha Quick Reference

These are the suggested gutta percha points for EdgeOne Fire, all EdgeTaper Families and EdgeOne Blaze Utopia NiTi Systems.


EDGEONE FIRE™ Gutta Percha		
29.5mm- 60 Pack		Item Codes
YELLOW	SMALL	EOFSGMP
RED	PRIMARY	EOFPRIGP
GREEN	MEDIUM	EOFMEDGP
WHITE	LARGE	EOFLGGP
Assorted: SM (20), PRI (20) MED (10) LG (10)		EOFGPAST-1
• Taper to match file design • Parallel when diameter meets 1.0mm • Use with EdgeOne Fire™ Files:Pg 70		

EDGETAPER™ Gutta Percha		
29.5mm- 60 Pack		Item Codes
YELLOW	F1	ETGPF1
RED	F2	ETGPF2
BLUE	F3	ETGPF3
BLACK	F4	ETGPF4
YELLOW	F5	ETGPF5
Assorted: F1 (10), F2 (10), F3 (20), F4 (10), F5 (10)		ETGPAST
• Taper to match file design • Parallel 16mm from tip • Use with EdgeTaper Platinum™ Files: Pg 71 • Use with EdgeTaper™ Files: Pg 72 • Use with EdgeTaper Blaze Utopia™ Files: Pg 73		

EDGEONE R UTOPIA™ Gutta Percha		
29.5mm- 60 Pack		Item Codes
YELLOW	R20	EORGP20-1
RED	R25	EORGP25-1
GREEN	R35	EORGP35-1
WHITE	R45	EORGP45-1
• Taper to match file design • Parallel when diameter meets 1.0mm • Use with EdgeOne Blaze Utopia™ Files: Pg 73		





\$9.84/ 60 Pack

New Gutta Percha Item Code	File	Size	File Item Code (3 Pack)	Length
		Small	EOFSM21	21mm
			EOFSM25	25mm
			EOFSM31	31mm
		Primary	EOFPRI21	21mm
			EOFPRI25	25mm
			EOFPRI31	31mm
		Medium	EOFMED21	21mm
			EOFMED25	25mm
			EOFMED31	31mm
		Large	EOFLG21	21mm
			EOFLG25	25mm
			EOFLG31	31mm
		Assorted 4 Pack	ASSTEOF21	21mm
			ASSTEOF25	25mm
			ASSTEOF31	31mm

EDGE TAPER PLATINUM™

\$9.84/ 60 Pack



New Gutta Percha Item Code	File	Size	File Item Code (6 Pack)	Length
		F1	ETF121HF	21mm
			ETF125HT	25mm
			ETF131HT	31mm
		F2	ETF221HF	21mm
			ETF225HT	25mm
			ETF231HT	31mm
		F3	ETF321HF	21mm
			ETF325HT	25mm
			ETF331HT	31mm
		F4	ETF421HF	21mm
			ETF425HT	25mm
			ETF431HT	31mm
		F5	ETF521HF	21mm
			ETF525HT	25mm
			ETF531HT	31mm
		Assorted 6 Pack	ETP21MM	21mm
			ETP25MM	25mm
			ETP31MM	31mm

EDGE TAPER™

\$9.84/ 60 Pack

New Gutta Percha Item Code	File	Size	File Item Code (6 Pack)	Length
		F1	ETF121	21mm
			ETF125	25mm
			ETF131	31mm
		F2	ETF221	21mm
			ETF225	25mm
			ETF231	31mm
		F3	ETF321	21mm
			ETF325	25mm
			ETF331	31mm
		F4	ETF421	21mm
			ETF425	25mm
			ETF431	31mm
		F5	ETF521	21mm
			ETF525	25mm
			ETF531	31mm
		Assorted 6 Pack	ETASST21	21mm
			ETASST25	25mm
			ETASST31	31mm

EDGETAPER BLAZE UTOPIA™
\$9.84/ 60 Pack

New Gutta Percha Item Code	File	Size	File Item Code (6 Pack)	Length
		F1	ETBUF121SU	21mm
			ETBUF125SU	25mm
			ETBUF131SU	31mm
		F2	ETBUF221SU	21mm
			ETBUF225SU	25mm
			ETBUF231SU	31mm
		F3	ETBUF321SU	21mm
			ETBUF325SU	25mm
			ETBUF331SU	31mm

EDGEONE BLAZE UTOPIA™
\$9.84/ 60 Pack

New Gutta Percha Item Code	File	Size	File Item Code (3 Pack)	Length
		R20	EOBU2021SU	21mm
			EOBU2025SU	25mm
			EOBU2031SU	31mm
		R25	EOBU2521SU	21mm
			EOBU2525SU	25mm
			EOBU2531SU	31mm
		R35	EOBU3521SU	21mm
			EOBU3525SU	25mm
			EOBU3531SU	31mm
		R45	EOBU4521SU	21mm
			EOBU4525SU	25mm
			EOBU4531SU	31mm

EdgeCaOH™ - \$17.24		
2.2G SYRINGE	818365020003	



Temporary root-canal filler for weeping canals that forms a hard-tissue barrier (apexification). Internal and external root resorption. Pre-mixed paste. Syringe tip for easy delivery of paste into canal. 10 Disposable tips included to prevent cross-contamination.

EDGE TEMP™ PLUS

\$9.48

40G	WHITE PASTE	ETEMPPLUS
-----	-------------	-----------



Hydraulic temporary restorative material. Designed for fabricating temporary fillings and other temporary prosthetics until the permanent restoration can be delivered. Setting time: 30 minutes.

EDGE BIO CERAMIC RETROFILL™ AND PERFORATION REPAIR MATERIAL

\$163.92

0.5G	WHITE PASTE	818365020089
------	-------------	--------------



EdgeBioCeramic™ RetroFill is a unique bioceramic repair material that exhibits superior healing and handling properties that has been shown in studies to out perform MTA. In addition, EdgeBioCeramic RetroFill is premixed and ready to use unlike MTA. Simply syringe out a small amount of material and deliver it to the repair site.

The Future of Root Repair



Root Repair Material



EdgeUtopia™ RRM is the optimal solution for root repair, offering excellent clinical outcomes with its scientifically backed performance, versatility, and long-term efficacy. EdgeUtopia RRM is the ideal choice for all your root and pediatric repair needs.

This proven pre-mixed calcium silicate repair material has now been upgraded to feature a novel monodose delivery system to improve handling and reduce the risk of cross-contamination. Trust in EdgeUtopia™ RRM for tissue healing, bacteria prevention, and overall root canal repair success.



Non Staining



Wide Range of Applications



Desired Biological Performance



Convenient Single Dose Delivery



Inhibits Bacterial Infiltration



Long-Lasting Results



0.15g per capsule
10 capsules per box

EDGE UTOPIA RRM™

Root Repair Material

Indications



Pre-Mixed Ready-to-Use Monodose Capsule

Forget the hassle of your repair material drying out too soon. Experience effortless, reliable performance with every use. No manual mixing and no titration devices needed.



Reduced Cross-Contamination

Unlike other premixed, manually mixed or titrated bioceramics, EdgeUtopia RRM™ ensures hygienic, mess-free delivery with no risk of contamination between patients due to its single dose application and the elimination of extensive manipulation in the operatory.



Release of Ca²⁺

Proven track record with over 5 years of clinical cases and research studies* demonstrating exceptional tissue healing and regeneration outcomes.



Wide Range of Applications

- Apexification
- Retrofilling
- Perforation Repairs
- Pulpotomy
- Pulp Capping
- External and Internal Resorption
- Deciduous Tooth Repairs



Cost-Effective Innovation

More product for a fraction of the price – helping you save money while delivering high-quality results. 3x the amount of our regular syringe formula per box.



0.15g per capsule
10 capsules per box
EURRM10MD
\$71.75/ 10 Pack



Pair with:
Reusable EdgeUtopia™ RRM™
Dispensing Gun

EURRMDSG
\$23.28

Available research and clinical case studies may refer to other private-label brands which are the same chemical formula as EdgeUtopia RRM™ Root Repair Material. For clarity, the brands EdgeUtopia RRM™ and Bio-C Repair are both manufactured by Angelus and are composed of equivalent chemical formulations.

Rev A 03/2025

EdgeFile® X7.04	
1 pack	EdgeFile® X7 17.04 25mm
1 pack	EdgeFile® X7 20.04 25mm
1 pack	EdgeFile® X7 25.04 25mm
1 pack	EdgeFile® X7 30.04 25mm
1 pack	EdgeFile® X7 35.04 25mm
1 pack	EdgeFile® X7 40.04 25mm

1 pack	EdgeFile® X7 45.04 25mm
3 packs	EdgeGlide Path™ 25mm
4 pack	Edge K-File™ #10, 25mm
4 pack	Edge K-File™ #15, 25mm
1 bottle	EdgeMix™
1 bottle	EdgeEndo® 17% EDTA



EdgeFile® X7.06	
1 pack	EdgeFile® X7 17.04 25mm
1 pack	EdgeFile® X7 20.04 25mm
1 pack	EdgeFile® X7 25.04 25mm
1 pack	EdgeFile® X7 30.04 25mm
1 pack	EdgeFile® X7 35.04 25mm
1 pack	EdgeFile® X7 40.04 25mm

1 pack	EdgeFile® X7 45.04 25mm
3 packs	EdgeGlide Path™ 25mm
4 pack	Edge K-File™ #10, 25mm
4 pack	Edge K-File™ #15, 25mm
1 bottle	EdgeMix™
1 bottle	EdgeEndo® 17% EDTA



EdgeOne Fire™	
4 packs	EdgeOne Fire™ Pri. 25mm
1 pack	EdgeOne Fire™ Small 25mm
1 pack	EdgeOne Fire™ Med. 25mm
1 pack	EdgeOne Fire™ Large 25mm
3 packs	EdgeOne GlidePath™ for EdgeOne Fire 25mm

4 packs	Edge K-File™ #10, 25mm
4 packs	Edge K-File™ #15, 25mm
1 bottle	EdgeMix™
1 bottle	EdgeEndo® 17% EDTA
1 pack	Edge K-File™ #15, 25mm



EdgeTaper Platinum™	
6 PACKS	EdgeTaper Platinum™ Assorted Pack 25mm
3 packs	EdgeGlidePath™ with FireWire™ 25mm
4 packs	Edge K-File™ #10, 25mm

4 packs	Edge K-File™ #15, 25mm
1	EdgeMix™
1	EdgeEndo® 17% EDTA



Product Limits By Category

Category	Limit Qty
Power Equipment	Limit 2
Power Accessories	Limit 10
Power Equipment Tips & EdgePellets	Limit 5
Ancillary Products Group 1	Limit 5
EdgeLube, EdgeGel, EdgeSeal, EdgeTemp, EdgeCaOH, EdgeGelce, EdgeMix, Sodium Hypochlorite, Chlorhexidine Gluconate, EDTA, EdgeUtopia BC Solvent	
Handfiles	No limit
Files	Limit 500
Retreatment Files	Limit 100
BioCeramic Sealers	Limit 10
EdgeBurs 25 & 100 Packs	Limit 5
EdgeBurs 5 & 10 Packs	Limit 10
EdgeUltrasonic Tips	Limit 12
Obturators - All Pack Quantities	Limit 30
Points- Paper Points, Gutta Percha Points, BioCeramic Points	Limit 30
Supplies	Limit 5
Endo Rings, Edge Luer-Lock Syringes, Edge Irrigation Tips, Pre-Assembled Luer-Lock Syringes,	
EdgeUtopia RRM	Limit 5
EdgeUtopia RRM Dispensing Gun	Limit 2

Price Breakdown:

20%	Carbide & Diamond Burs, Sealer, Solution, Syringes, Supplies
30%	Dam, Power, Ultrasonic Tips, Gutta Percha & Paper Points
40%	Obturators, Verifiers, Handfiles
45%	Files and Gates

Limits are cumulative throughout your residency (see front cover for more information). Limit quantity is calculated per category, not by specific product lines, systems, or names.



EDGEENDO®

© 2026 Edge Endo, a division of HSG-S Corp.
No copying without permission.
Not responsible for typographical errors
Rev. A 01/2026