

EDGEENDO[®]

PRODUCT CATALOG

2026

Created by an Endodontist for Endodontists

Established in 2012 by American Endodontist Charles Goodis, EdgeEndo® was founded with the vision to make root canal treatments safe, easy, and cost-effective overcoming the shortfalls of manufacturers offering outdated technology with a high price tag.

The Solution – Innovation combined with high-quality and cost-effectiveness to support the growth of endodontic practices by saving chair time and delivering outstanding patient outcomes. In just ten years, EdgeEndo revolutionized the endodontic industry in the USA becoming a leader in endodontic solutions!

Safe – Healthy Patients - Happy Dentists. EdgeEndo solutions are developed with the patient at the center. Our various treatment processes provide excellent resistance to cyclic fatigue and flexibility to carefully follow the canal, preserving canal anatomy and dentin and allowing for faster procedures. All for increased safety and productivity.

Easy – Changing habits is always a challenge for any dental professional. Our solutions are designed for ease of use and are compatible with most endodontic techniques. EdgeEndo continues expanding its product offering to cover the entire endodontic procedure; diagnosis, instrumentation, obturation and restoration. No need for multiple sources and complexity. We make it easy for the dentists.

Cost-Effective – EdgeEndo has proven that quality does not need to come with a premium price tag. Our goal is to make root canals accessible to anyone, without compromises.

EdgeEndo® remains steadfast in its mission to provide game-changing products and solutions.

Innovation - EdgeEndo is the inventor of FireWire and FireWire Blaze, two unique combinations of heat treatment and cryogenic applications allowing NiTi instruments to deliver incredible flexibility and resistance to cyclic fatigue.

Quality – Thousands of clinicians trust EdgeEndo and our mission to deliver a seamless and consistent experience with all our products. The control of manufacturing processes is our obsession, and our instruments are all manufactured under strict quality controls.

**Revolutionize endodontic procedures with EdgeEndo:
the partner dental professionals have been searching for!**




Table of Contents

EdgePower Suite

EdgeApex HP	4
EdgeFlow	6
EdgeApex	9
EdgePack	10

Diagnosis & Isolation

EdgeIcse	13
EdgeDam	13

Supplies

EndoRing	14
--------------------	----

Initial Access

EdgeTip.	15
EdgeBurs - Carbide	16
EdgeGates.	19
EdgeBroaches	19
EdgeHedstroms	20
EdgeTaper Retreat.	20
BC Solvent.	21

Glidepath

EdgeFind.	22
EdgeOneGlidepath	22
EdgeOne Blaze Utopia - RGlid.	23

Access Refinement

EdgeSurgical Tip	24
EdgeBurs - Diamond	25

Working Length Determination

EdgeK-File NiTi	29
EdgeUtopia K-File/EdgeK-File.	30
EdgeUtopia C-File/EdgeC-File	31

File Systems

Rotary • Constant Taper

EdgeOnyx	32
EdgeX7 Utopia	34
EdgeFile X7	36

EdgeFile X3	38
EdgeSequel Utopia	40
EdgeSequel Sapphire.	43
EdgeEvolve	45

Rotary • Variable Taper

EdgeDiamond Utopia.	48
EdgeV-Taper HT	50
EdgeTaper Blaze Utopia	52
EdgeTaper Platinum.	54
EdgeTaper	56

Reciprocating • Variable Taper

EdgeOne Blaze Utopia	58
EdgeOne Fire	60

Reciprocating • Constant Taper

EdgeOne Platinum	62
EdgeFile X1	64

Irrigation

EdgeLube	66
EdgeGel	66
Irrigation Tips & Syringes.	67

Sealers

EdgeBioCeramic Sealer.	69
EdgeUtopia BioCeramic Sealer	69
EdgeBioCeramic Thermal Flow	69
EdgeSeal	70

Obturation

EdgeCore	70
Edge Fill	70
Points.	71

Intracanal Medicament & Repair

EdgeCaOH.	72
EdgeTemp Plus.	72
EdgeBioCeramic Retrofill & Repair	72
EdgeUtopia RRM.	78

EDGEAPEX HP™



- Built in Apex Locator with option to use as a standalone apex locator with handfile
- Pre-programmed with top selling EdgeEndo files PLUS 10 fully customizable memory program slots
- Efficient brushless motor, low noise, long service life
- Cordless portable endo motor with combined length determination
- 360° rotation of contra angle
- Adopt real-time feedback technology and dynamic torque control, effectively preventing file separation
- Wireless charging
- Torque control
- Auto On/Off, Torque Reverse, Apical Reverse, Apical Slow Down
- Visual torque on file in use
- Max Torque 5.0 Ncm
- Change dominant hand
- Range: 100 - 2500 RPM

\$1,179.00

EDGEAPEX HP™

\$1179.00

1 UNIT	EAPEXHP
--------	---------

Accessories included in package:

- Measuring Wire
- Lip Hook (2x)
- File Clip (4x)
- Touch Probe (2x)
- Power Adapter
- Protective Silicon Cover (2x)
- Spray Nozzle
- O-Rings (2x)
- Disposable Sleeves (1 pack)

Warranty: 2 years for handpiece motor, charging base, and power adapter. 1 year for contra angle attachment. 6 months for accessories.



Accessories

USB-C EDGEAPEX HP™ MEASURING WIRE 1 PACK
EAPEXNMW \$19.30



CONTRA ANGLE 1 PACK
EAPEXCA \$215.70



TOUCH PROBE 1 PACK
EAPEXTP \$9.60



FILE CLIP 5 PACK
EAPEXFC \$64.65



LIP HOOK 5 PACK
EAPEXLH \$16.10



CHARGING BASE 1 PACK
EAPEXCB \$78.65



PROTECTIVE SILICONE COVER 3 PACK
EAPEXSC \$11.75



EDGEFLOW™



- Wireless and electric
- Excellent gutta percha fluidity, uniform and smooth gutta percha output without bubbles
- Precise and intuitive display of remaining gutta-percha
- Auto heating stop after timeout; auto shutdown after timeout
- Charging base with battery
- Temperature monitoring design
- Dual batteries deliver longer endurance, enabling filling about 100 root canals
- Reusable injection needles available in various sizes and lengths
- Preset temperatures: 100°C, 120°C, 150°C, 180°C, 200°C
- Injection speeds: slow, medium, fast

\$749.00

EDGEFLOW™

\$749.00

1 UNIT	EEPTO
--------	-------

Accessories included in package:

- Main Unit
- Charging Stand
- Power Adapter
- Battery (2x)
- Heat Insulation Hood
- Wrench of Needle
- Gutta-Percha Ring
- Instruction Manual
- Qualified Certification
- Packing List
- Injection Needles (4x)
 - E23G-M (1x)
 - E23G-L (1x)
 - E25G-L (2x)

Warranty: 1 year for device.



Accessories

CHARGING BASE 1 UNIT
EEPTOCH \$64.65



OBTURATOR WRENCH 1 PACK
EEPTOW \$16.10



EDGEFLOW™ BATTERY 1 PACK
EEPTOB \$32.25



EDGEFLOW™ HEAT INSULATION HOOD 3 PACK
EEPTOHIH3 \$19.30



EdgeFlow™ Reusable Needle Tips

NEEDLE TIP E20G-S 1 PACK 20 GAUGE; 22mm LENGTH	
EE20GSIN	\$36.60



NEEDLE TIP E23G-L 1 PACK 23 GAUGE; 28mm LENGTH	
EE23GLIN	\$36.60



NEEDLE TIP E20G-M 1 PACK 20 GAUGE; 24mm LENGTH	
EE20GMIN	\$36.60



NEEDLE TIP E25G-M 1 PACK 25 GAUGE; 24mm LENGTH	
EE25GMIN	\$36.60



NEEDLE TIP E20G-L 1 PACK 20 GAUGE; 28mm LENGTH	
EE20GLIN	\$36.60



NEEDLE TIP E25G-L 1 PACK 25 GAUGE; 28mm LENGTH	
EE25GLIN	\$36.60



NEEDLE TIP E23G-M 1 PACK 23 GAUGE; 24mm LENGTH (x1)	
EE23GMIN	\$36.60



Note:
Tip item name refers to gauge and tip length. (e.g. Needle Tip E20G - S is a 20 gauge needle with a small tip length)
Small = 22mm • Medium = 24mm • Large = 28mm

EDGE PELLETS™

100 PACK	
00810184450092	\$30.20



EDGEAPEX™

\$539.00

1 UNIT	EEALU
--------	-------

Features

- Apex Locator with full color LED display
- Built in Electronic Pulp Vitality Tester
- Omnidirectional base

Accessories included in package:

- Base
- Measuring Wire
- File Clip (4x)
- Touch Probe (2x)
- Pulp Tester Probe (2x)
- Adapter
- Tester
- Enlarged display of apical foramen
- Instruction card for apical stop
- Apex Locator troubleshooting card

Warranty: 1 year for device. 6 months for accessories.



Apex Locator Accessories

MEASURING WIRE 1 PACK	
EEALUMW	\$19.30



TESTER 1 PACK	
EEALUT	\$16.10



LIP HOOK 5 PACK	
EEALULH5	\$16.10



ENDO RULER 1 PACK	
EEALUER	\$16.10



MEASURING FILE 1 PACK	
EEALUMF	\$11.75



FILE CLIP 5 PACK	
EEALUFC5	\$64.65



EDGEPACK™



- Wireless gutta percha obturation device
- Ergonomic design and ambidextrous accessible
- Easy-to-use
- Pre-bent tip in a variety of tip sizes and tapers
- Improved penetrability into narrow root canals
- High precision temperature control
- Smooth cutting of gp from high-efficiency conversion system
- Compact size for easy cleaning and storage
- Visual and auditory prompts for safe operation
- Preset temperatures: 150°C, 180°C, 200°C, 230°C
- Time related features: warm up time within 0.2 seconds, auto shut off after 10 seconds
- 4-Hour battery life with each charge

\$589.00

EDGEPACK™

\$589.00

1 UNIT	EEDPU
--------	-------

Accessories included in package:

- Charging Stand
- Power Adapter
- WP4004 Tip
- WP4504 Tip
- WP5506 Tip
- WP5508 Tip
- Heat Insulation Sleeve of Tips (4x)
- Instruction Manual and Qualified Certification

Warranty: 1 year for device



Accessories

EDGEPACK CHARGING STAND		
1 UNIT	EEDPUCS	\$64.65



EdgePack Tips

EDGEPACK 3504 TIP 1 PACK	
EEDPUT3504	\$53.85



EDGEPACK 5508 TIP 1 PACK	
EEDPUT5508	\$53.85



EDGEPACK 4004 TIP 1 PACK	
EEDPUT4004	\$53.85



EDGEPACK 5508L TIP 1 PACK	
EEDPUT5508L	\$53.85



EDGEPACK 4504 TIP 1 PACK	
EEDPUT4504	\$53.85



EDGEPACK 5510 TIP 1 PACK	
EEDPUT5510	\$53.85



EDGEPACK 5004 TIP 1 PACK	
EEDPUT5004	\$53.85



EDGEPACK 6012 TIP 1 PACK	
EEDPUT6012	\$53.85



EDGEPACK 5506 TIP 1 PACK	
EEDPUT5506	\$53.85



Note:

Tip numbering of EdgePack Tips refer to tip size and taper. (e.g. 3504 tip is a tip size 35 with an .04 taper)

EdgeIce™		
6 OZ. SPRAY BOTTLE	818365020065	\$17.15



A quick way to check cold sensitivity and pulp vitality. Significantly colder than ethyl chloride or ice. EdgeIce™ features a directional nozzle for easy placement on micro applicator or cotton pellet.

EDGE DAM™

Available in Latex or Latex-Free options

Latex 5x5	52 Pack	\$12.85
Medium	EDBLM5552	
Heavy	EDBLH5552	
Latex 6x6	36 Pack	\$12.85
Medium	EDBLM6636	
Heavy	EDBLH6636	



Latex-Free 5x5	50 Pack	\$76.50
Medium	EDVLFM5550	
Heavy	EDVLFH5550	
Latex-Free 6x6	50 Pack	\$102.40
Medium	EDVLFM6650	
Heavy	EDVLFH6650	



ENDORING[™]

EndoRing[®] II Kit with REDUCED GLARE Metal Ruler and e-Foam[®] Rotary HD Inserts

1 KIT	ER2GSEE-S	\$53.85
-------	-----------	----------------

1 ring assembly & 6 silver HD foam inserts



EndoRing[®] II Kit with Metal Ruler

1 KIT	ER2BEE-S	\$45.20
-------	----------	----------------



1 ring assembly, 3 blue foam inserts & 3 yellow foam inserts



Foam Refills - 48 Pack

\$43.05

Silver Rotary HD e-foam [®] Insert Refill	ERFSHDEE-S
Blue e-foam [®] Insert Refill	ERFBEE-S
Yellow e-foam [®] Insert Refill	ERFYEE-S
Combo e-foam [®] Insert Refill (24 Yellow & 24 Blue)	ERFEE-S



EdgeTips™ 1-5 are diamond-coated for more efficiency. We don't use zirconium nitride. EdgeTips™ 6-8 are a special titanium alloy and their special geometries allow precision work in areas of limited access.

- Find calcified or hidden canals
- Remove intracanal obstructions
- Remove dental restorations and cores
- Dislodge posts and broken instruments
- Remove remove pulp stones
- Unique shaft designs for better access and visibility
- Length and diameters designed to progress down the root

EdgeTip #1 - ETBD
Diamond-coated ball tip, length 20mm, taper 5%

Used for searching for canals and locating calcified canals

EdgeTip #2 - ET18D
Diamond-coated steel tip 76µm, length 18mm, taper 5%

Used for finishing access cavity, removing overhangs, calcifications, and filling materials

EdgeTip #3 - ET20D
Diamond-coated retreatment tip, 30µm, length 20mm, taper 5%

Used in 1st coronal 3rd to remove very hard materials by brushing walls

EdgeTip #4 & #5 - ET40D4
Long Diamond-coated retreatment tip, 30µm, 40mm, taper 4%

Retreatment middle 3rd

EdgeTip #6 - ED25S
Short Titanium- Niobium tip, length 15mm taper 4%

Retreatment in coronal 3rd and isthmuses

EdgeTip #7 - ET25
Titanium- Niobium tip, length 20mm, taper 3%

Retreatment middle 3rd and used for broken instrument retrieval

EdgeTip #8 - ET25L
Titanium- Niobium tip, 25mm, taper 3%

Retreatment apical 3rd and used in long straight canals





OPERATIVE CARBIDES




ROUND



US NO.	ISO NO.	1/4	1/2	1	2	3	4
FG (19mm)	500 314 001 001	005	006	008	010	012	014
	10 PACK	\$17.15 FG1QTR10	\$17.15 FG1HLF10	\$17.15 FG110	\$17.15 FG210	\$17.15 FG310	\$17.15 FG410
FG SURG. (25mm)	500 316 001 001		006		010		014
	10 PACK		\$17.15 FG1HLFSURG10		\$17.15 FG2SURG10		\$17.15 FG4SURG10
RA (22mm)	500 204 001 001		006		010	012	014
	10 PACK		\$17.15 RA1HLF10		\$17.15 RA210	\$17.15 RA310	\$17.15 RA410

ROUND



US NO.	ISO NO.	5	6	7	8	10
FG (19mm)	500 314 001 001	016	018		023	
	10 PACK	\$17.15 FG510	\$17.15 FG610		\$17.15 FG810	
FG SURG. (25mm)	500 316 001 001		018		023	
	10 PACK		\$17.15 FG6SURG10		\$17.15 FG8SURG10	
RA (22mm)	500 204 001 001		018	021	023	027
	10 PACK		\$17.15 RA610	\$17.15 RA710	\$17.15 RA810	\$17.15 RA1010



INVERTED CONE




US NO.	ISO NO.	34	35
FG (19mm)	500 314 010 001	008	010
	10 PACK	\$17.15 FG3410	\$17.15 FG3510



AMALGAM PREP



US NO.	ISO NO.	245
FG (19mm)	500 314 234 006	009
	10 PACK	\$17.15 FG24510
FG SHORT (16mm)	500 313 234 006	009
	10 PACK	\$17.15 FG245SS10




PEAR

US NO.	ISO NO.	330	331	332
FG (19mm)	500 314 237 001	008	010	012
	10 PACK	\$17.15 FG33010	\$17.15 FG33110	\$17.15 FG33210
FG SHORT (16mm)	500 313 237 001	008		
	10 PACK	\$16.63 FG330SS10		



PEAR LONG HEAD

US NO.	ISO NO.	331L
FG (19mm)	500 314 238 006	010
	10 PACK	\$17.15 FG331L10



FLAT FISSURE CROSS CUT

US NO.	ISO NO.	556	557	558
FG (19mm)	500 314 107 007	008	010	012
	10 PACK	\$17.15 FG55610	\$17.15 FG55710	\$17.15 FG55810
FG SURG (25mm)	500 316 107 007		010	012
	10 PACK		\$17.15 FG557SURG10	\$17.15 FG558SURG10
FG SHORT (16mm)	500 313 107 007		010	
	10 PACK		\$17.15 FG557SS10	



FLAT FISSURE CROSS CUT LONG HEAD

US NO.	ISO NO.	557L
FG (19mm)	500 314 110 007	010
	10 PACK	\$17.15 FG557L10



DOMED FISSURE PLAIN CUT

US NO.	ISO NO.	1157	1158
FG (19mm)	500 314 137 006	010	012
	10 PACK	\$17.15 FG115710	\$17.15 FG115810



DOMED FISSURE CROSS CUT

US NO.	ISO NO.	1557	1558
FG (19mm)	500 314 137 007	010	012
	10 PACK	\$17.15 FG155710	\$17.15 FG155810

SPECIALTY



151-Z

FIG NO.	151-Z
DIAMETER (1/10mm)	016
H/LENGTH (mm)	11.0
OVERALL LENGTH (mm)	23.0
5 PACK	\$24.70 151ZB5



ENDO ACCESS

FIG NO.	152
DIAMETER (1/10mm)	016
H/LENGTH (mm)	9.0
OVERALL LENGTH (mm)	23
5 PACK	\$24.70 152B5



LINDEMANN
CROSS CUT

US NO.	162
H/LENGTH (mm)	10.0
FG SURG. (25mm)	016
5 PACK	\$24.70 FGSURG162B5



D LINDEMANN
CROSS CUT

US NO.	162A
H/LENGTH (mm)	9.2
FG SURG. (25mm)	016
5 PACK	\$24.70 FGSURG162AB5



LINDEMANN

US NO.	161
H/LENGTH (mm)	10.0
FG SURG. (25mm)	016
5 PACK	\$24.70 FGSURG161B5


EDGE GATES™

Stainless Steel



6 Pack		\$13.90 each	
	28mm	32mm	
#1	EG0128	EG0132	
#2	EG0228	EG0232	
#3	EG0328	EG0332	
#4	EG0428	EG0432	
#5	EG0528	EG0532	
#6	EG0628	EG0632	
ASSORTED PACK		\$13.90 each	
6 Pack	EGASST28	EGASST32	

EDGE BROACHES™



10 Pack		\$7.45 each	
Size	Color	21mm	25mm
0	 PURPLE	EB4XF021	EB4XF025
1	 WHITE	EB3XF121	EB3XF125
2	 YELLOW	EB2XF221	EB2XF225
3	 RED	EB1XF321	EB1XF325
4	 BLUE	EB0XF421	EB0XF425
5	 GREEN	EBM0521	EBM0525
6	 BLACK	EBC0621	EBC0625
Assorted Pack 10 Pack		\$7.45 each	
Size		21mm	25mm
INCLUDED: #0 (X1), #1 (X1), #2 (X2), #3 (X2), #4 (X2), #5 (X1), #6 (X1)		EBASST0621	EBASST0625
INCLUDED: #2 (X3), #4 (X4), #5 (X3)		EBASST2421	EBASST2425

EDGEHEDSTROMS™

Stainless Steel



6 Pack		\$5.60 each		
Size	Color	21mm	25mm	31mm
15	<input type="checkbox"/> WHITE	HED1521	HED1525	HED1531
20	<input checked="" type="checkbox"/> YELLOW	HED2021	HED2025	HED2031
25	<input checked="" type="checkbox"/> RED	HED2521	HED2525	HED2531
30	<input checked="" type="checkbox"/> BLUE	HED3021	HED3025	HED3031
35	<input checked="" type="checkbox"/> GREEN	HED3521	HED3525	HED3531
40	<input checked="" type="checkbox"/> BLACK	HED4021	HED4025	HED4031
45	<input type="checkbox"/> WHITE	HED4521	HED4525	HED4531
50	<input checked="" type="checkbox"/> YELLOW	HED5021	HED5025	HED5031

EDGETAPER RETREAT™

Rotary Retreatment Files



EdgeTaper Retreat™ files are designed for the different needs of unfilling the coronal third, the middle third and the apical third before canal reshaping.

6 Pack		\$38.80 each	
Size	Color	Length	SKU
D1	<input type="checkbox"/> White	16mm	ETRD116
D2	<input type="checkbox"/> White	18mm	ETRD218
D3	<input type="checkbox"/> White	22mm	ETRD322
D1-D3	<input type="checkbox"/> White	D1 - 16mm (X2) D2 - 18mm (X2) D3 - 22mm (X2)	ETRASSTD123

EDGE U UTOPIA™ BC SOLVENT

**EdgeUtopia™ BC Solvent is not just another solution,
it's *THE* solution for dissolving Bioceramic Sealers**

Dissolves Bioceramic Sealers and Calcium Hydroxide.

EdgeUtopia™ BC Solvent is engineered to dissolve the toughest bioceramic sealers and calcium hydroxide, outperforming other traditional solutions on the market. No more struggle with difficult to remove materials.

Tough on Sealer. Gentle on Teeth.


While other products fall short on dissolving tough materials without damaging sensitive structures, EdgeUtopia BC Solvent works precisely where it's needed most, while maintaining the integrity of the surrounding tooth tissues.

EdgeUtopia™ BC Solvent - Preserving What Matters.

\$107.85




240mL Bottle

10818207024937



EDGEFIND™

Non Heat-Treated Endonol™ NiTi


6 Pack		\$35.60 each		
Size	Color	21mm	25mm	31mm
P1	 PURPLE	EFNDP121	EFNDP131	EFNDP125
P2	 WHITE	EFNDP221	EFNDP231	EFNDP225
P3	 YELLOW	EFNDP321	EFNDP331	EFNDP325
Assorted 6 Pack		\$35.60 each		
P1 (X2), P2 (X2), P3 (X2)		EFNDASST21	EFNDASST25	EFNDASST31



EDGEGLIDEPATH™

Heat-Treated FireWire™ NiTi

Rotary GlidePath

4 Pack		\$27.95 each		
Color	21mm	25mm	31mm	
 WHITE	EGP0221	EGP0225	EGP0231	



The EdgeGlidePath™ Heat Treated FireWire NiTi file features the flexibility and strength that EdgeEndo® is known for, replacing as many as 8 instruments to create a glidepath faster, with far less effort and with a 1mm maximum flute diameter. Replaces #15, 20, 25, 30, 35 hand files in a step-back technique and a #1, 2 & 3 Gates Glidden. Used exactly like Tulsa Dental ProGlider® at substantially lower cost.

*ProGlider® is a registered trademark of Dentsply Tulsa Dental Specialties.

EDGEONE GLIDEPATH™

For EdgeOne Fire™ and EdgeOne Platinum™

Reciprocating GlidePath

3 pack			\$27.95 each
Color	21mm	25mm	31mm
<input type="checkbox"/> WHITE	EOFGP21	EOFGP25	EOFGP31



EDGEONE BLAZE UTOPIA™

R-Glide

Reciprocating GlidePath

3 pack			\$27.95 each
Color	21mm	25mm	31mm
<input type="checkbox"/> WHITE	EOBURG1421SU	EOBURG1425SU	EOBURG1431SU



EDGE SURGICAL TIP™ \$53.85 each

EdgeSurgical Tips™ are designed for use in endodontic microsurgery. Made with a diamond coating, these instruments cut dentin faster without the need for repeated ultrasonics. The dentin surface is left smooth yet microscopically rough which results in better adaptation of filling materials.

- Diamond coating for efficient cutting of dentin
- Irrigation strategically placed at the tip for control
- Tip angles and lengths designed for access to difficult areas
- We don't use zirconium nitride

EdgeSurgical Tip #1 AS3D1 Diamond-coated universal tip 30µm, length 3mm, taper 9%
Intended for apical surgery of anterior teeth. Tip should be used at the lowest possible effective power.



EdgeSurgical Tip #2 AS3D2 Diamond-coated universal tip 30µm, taper 9%
For apical surgery of anterior teeth.



EdgeSurgical Tip #3 P15LD Left-oriented retro surgery tip, diamond-coated 30µm, length 5mm, taper 7%
Double angle 75 degrees for mandibular left and maxillary right premolars.



EdgeSurgical Tip #4 ASLD Left-oriented tip, diamond-coated 30µm, length 3mm, taper 10%
Double angle 110 degrees for mandibular right and maxillary left posterior teeth. Tip should be used with very light pressure.



EdgeSurgical Tip #5 P15RD Right-oriented tip, diamond-coated 30µm, length 6mm, taper 7%
Double angle 75 degrees for mandibular right and maxillary left premolars.



EdgeSurgical Tip #6 ASRD Right-oriented tip, diamond-coated 30µm, length 3mm, taper 10%
Double angle 110 degrees for mandibular right and maxillary left posterior teeth. Tip should be used with very light pressure.



DIAMONDS



001 | 801

Round



Head Diameter 1/10mm	010	012	014	016	018	021	023	025
L1: Head Length L1mm	1.0	1.2	1.4	1.6	1.8	2.1	2.3	2.5
L2: OAL mm	19.0	19.0	19.0	19.0	19.0	19.0	19.0	19.0
ISO								
806 314 001 514...								025 \$43.05 001025F
806 314 001 524...	010 \$43.05 001010M	012 \$43.05 001012M	014 \$43.05 001014M	016 \$43.05 001016M	018 \$43.05 001018M		023 \$43.05 001023M	
806 314 001 534...	010 \$43.05 001010C	012 \$43.05 001012C	014 \$43.05 001014C	016 \$43.05 001016C	018 \$43.05 001018C	021 \$43.05 001021C	023 \$43.05 001023C	
806 314 001 544...					018 \$43.05 001018XC		023 \$43.05 001023XC	



032 | 813

Diablo



Head Diameter 1/10mm	014
L1: Head Length L1mm	1.5
L2: OAL mm	19.0
ISO	
806 314 032 534...	014 \$43.05 032014C



544 | 845KR Short

Round Edge Taper



Head Diameter 1/10mm	019	026
L1: Head Length L1mm	4.0	4.0
L2: OAL mm	18.0	18.0
ISO	Short	Short
806 313 544 514...		026 \$43.05 544026SF
806 313 544 534...	019 \$43.05 544019SC	



546 | 847KR

Round Edge Taper

Head Diameter 1/10mm	017
L1: Head Length L1mm	8.0
L2: OAL mm	20.0
ISO	Short
806 314 546 514...	017 \$43.05 546017SF
806 314 546 524...	017 \$43.05 546017SM
806 314 546 534...	017 \$43.05 546017SC
806 314 546 544...	017 \$43.05 546017SXC



546 | 847KR

Round Edge Taper

Head Diameter 1/10mm	016
L1: Head Length L1mm	8.0
L2: OAL mm	22.0
ISO	
806 313 546 534...	016 \$43.05 546016C



197 | 855 Short

Round End Taper

Head Diameter 1/10mm	017	024
L1: Head Length L1mm	6.0	7.0
L2: OAL mm	18.5	18.5
ISO	Short	Short
806 314 197 544...	017 \$43.05 197017SXC	024 \$43.05 197024SXC



197 | 855 XShort

Round End Taper

Head Diameter 1/10mm	016
L1: Head Length L1mm	6.0
L2: OAL mm	16.5
ISO	XShort
806 313 197 534...	016 \$43.05 197016SHC





198 | 856

Round End Taper

Head Diameter 1/10mm	012	014	015	016	018	019	021
L1: Head Length L1mm	8.0	8.0	8.0	8.0	8.0	8.0	8.0
L2: OAL mm	22.0	22.0	22.0	22.0	22.0	22.0	22.0
ISO							
806 316 198 514...		■ 014 \$43.05 198014F		■ 016 \$43.05 198016F	■ 018 \$43.05 198018F		
806 316 198 524...	■ 012 \$23.97 198012M	■ 014 \$43.05 198014M		■ 016 \$43.05 198016M	■ 018 \$43.05 198018M		
806 316 198 534...	■ 012 \$43.057 198012C	■ 014 \$43.05 198014C	■ 015 \$43.05 198015C	■ 016 \$43.05 198016C	■ 018 \$43.05 198018C		■ 021 \$43.05 198021C
806 316 198 544...	■ 012 \$43.05 198012XC	■ 014 \$43.05 198014XC	■ 015 \$43.05 198015XC	■ 016 \$43.05 198016XC	■ 018 \$43.05 198018XC	■ 019 \$43.05 198019XC	■ 021 \$43.05 198021XC



198 | 856

Round End Taper

Head Diameter 1/10mm	022	025
L1: Head Length L1mm	8.0	8.0
L2: OAL mm	22.0	22.0
ISO		
806 315 198 534...	■ 022 \$43.05 198022C	■ 025 \$43.05 198025C
806 315 198 544...	■ 022 \$43.05 198022XC	■ 025 \$43.05 198025XC




198 | 856 Short

Round End Taper

Head Diameter 1/10mm	023
L1: Head Length L1mm	8.5
L2: OAL mm	20.0
ISO	Short
806 314 198 534...	■ 023 \$43.05 198023SC



198 | 856L Long
Round End Taper



Head Diameter 1/10mm	016	018	020
L1: Head Length L1mm	9.0	9.0	9.0
L2: OAL mm	22.0	22.0	22.0
ISO			
806 315 198L 514...	■ 016 \$43.05 198L016F		
806 315 198L 534...	■ 016 \$43.05 198L016C	■ 018 \$43.05 198L018C	■ 020 \$43.05 198L020C
806 315 198L 544...	■ 016 \$43.05 198L016XC	■ 018 \$43.05 198L018XC	




199 | 850
Round End Taper



Head Diameter 1/10mm	010	012	014	016	018
L1: Head Length L1mm	10.0	10.0	10.0	10.0	10.0
L2: OAL mm	24.0	24.0	24.0	24.0	24.0
ISO					
806 316 199 524...					■ 018 \$43.05 199018M
806 316 199 534...	■ 010 \$43.05 199010C	■ 012 \$43.05 199012C	■ 014 \$43.05 199014C	■ 016 \$43.05 199016C	018 \$25.80 199018C
806 316 199 544...				■ 016 \$43.05 199016XC	■ 018 \$43.05 199018XC



199 | 850 Short
Round End Taper



Head Diameter 1/10mm	014	016
L1: Head Length L1mm	10.0	10.0
L2: OAL mm	22.0	22.0
ISO	Short	Short
806 315 199 514...	■ 014 \$43.05 199014SXC	■ 016 \$43.05 199016SXC



263 | 369
Bud Slender













Head Diameter 1/10mm	025
L1: Head Length L1mm	5.5
L2: OAL mm	20.0
ISO	
806 314 263 534...	■ 025 \$43.05 263025C

EDGEK-FILE™ NITI

Non Heat-Treated NiTi
Alternative to Tulsa Dental K-File®



\$8.50 each

6 Pack		Non-Sterile		
Size	Color	21mm	25mm	31mm
15	 WHITE	NITIKF1521	NITIKF1525	NITIKF1531
20	 YELLOW	NITIKF2021	NITIKF2025	NITIKF2031
25	 RED	NITIKF2521	NITIKF2525	NITIKF2531
30	 BLUE	NITIKF3021	NITIKF3025	NITIKF3031
35	 GREEN	NITIKF3521	NITIKF3525	NITIKF3531
40	 BLACK	NITIKF4021	NITIKF4025	NITIKF4031
45	 WHITE	NITIKF4521	NITIKF4525	NITIKF4531
50	 YELLOW	NITIKF5021	NITIKF5025	NITIKF5031
55	 RED	NITIKF5521	NITIKF5525	NITIKF5531
60	 BLUE	NITIKF6021	NITIKF6025	NITIKF6031

Tulsa Dental K-File® is a registered trademark of Dentsply Tulsa Dental Specialties



\$6.50 each

\$5.60 each

EdgeUtopia K-File™/ 6pk				Sterile Bliser Pack			EdgeK-File™/ 6pk			Non-Sterile
Alternative to Lexicon® Pre-Sterilized K-Files and Roydent™ K-Files						Alternative to Tulsa Dental K-File®				
Stiffer with more efficient cutting compared to EdgeK-Files™										
Size	Color	21mm	25mm	21mm	25mm	31mm				
6	PINK	EUKF621	EUKF625	KF621	KF625	KF631				
8	GRAY	EUKF821	EUKF825	KF821	KF825	KF831				
10	PURPLE	EUKF1021	EUKF1025	KF1021	KF1025	KF1031				
15	WHITE	EUKF1521	EUKF1525	KF1521	KF1525	KF1531				
20	YELLOW	EUKF2021	EUKF2025	KF2021	KF2025	KF2031				
25	RED	EUKF2521	EUKF2525	KF2521	KF2525	KF2531				
30	BLUE	EUKF3021	EUKF3025	KF3021	KF3025	KF3031				
35	GREEN	EUKF3521	EUKF3525	KF3521	KF3525	KF3531				
40	BLACK	EUKF4021	EUKF4025	KF4021	KF4025	KF4031				
45	WHITE	EUKF4521	EUKF4525	KF4521	KF4525	KF4531				
50	YELLOW	EUKF5021	EUKF5025	KF5021	KF5025	KF5031				
55	RED	EUKF5521	EUKF5525	KF5521	KF5525	KF5531				
60	BLUE	EUKF6021	EUKF6025	KF6021	KF6025	KF6031				
70	GREEN	EUKF7021	EUKF7025	KF7021	KF7025	KF7031				
80	BLACK	EUKF8021	EUKF8025	KF8021	KF8025	KF8031				
Assorted 6 Pack			\$6.50 each							
#10 - #35		EUKF103521	EUKF103525							

Lexicon® is a registered trademark of Dentsply Sirona Inc. Roydent™ is a registered trademark of Roydent Inc. Tulsa Dental K-File® is a registered trademark of Dentsply Tulsa Dental Specialties





C-FILE



Stainless Steel

\$7.25 each

\$6.45 each

EdgeUtopia C-File™/ 6pk				Sterile Bliser Pack		EdgeC-File™/ 6pk		Non-Sterile	
Alternative to Lexicon® Pre-Sterilized C-Files and Roydent™ C-Files						Alternative to Tulsa Dental C-File®			
Stiffer with more efficient cutting compared to EdgeC-Files™									
Size	Color	21mm	25mm	21mm	25mm	21mm	25mm	21mm	25mm
6	 PINK	EUCF621	EUCF625	CF621	CF625	CF621	CF625	CF621	CF625
8	 GRAY	EUCF821	EUCF825	CF821	CF825	CF821	CF825	CF821	CF825
10	 PURPLE	EUCF1021	EUCF1025	CF1021	CF1025	CF1021	CF1025	CF1021	CF1025
15	 WHITE	EUCF1521	EUCF1525	CF1521	CF1525	CF1521	CF1525	CF1521	CF1525

Lexicon® is a registered trademark of Dentsply Sirona Inc. Roydent™ is a registered trademark of Roydent Inc. Tulsa Dental C-File® is a registered trademark of Dentsply Tulsa Dental Specialties.

EDGEONYX™

\$28.80*
6 Pack


Balanced Performance

EdgeOnyx™ is the newest evolution in EdgeEndo's NiTi portfolio. Engineered over 10+ years to provide balance between flexibility, cutting efficiency and control.

- Cutting efficiency
- Flexibility
- Clinical control

Featuring:

- New safety tip
- Anti-wobble shank
- Proprietary heat treatment for pre-bending without sacrificing flute fidelity
- Minimally invasive, < 0.8 max flute diameter
- Non-Sterile Blister Pack



- Parabolic cross section
- Available in 2 Lengths: 21 & 25mm w/ 15mm orifice opener**
- Speed: 300-500 RPM
- Torque: 2.5 Ncm



Precise Control



Pre-Bendable










Minimally Invasive

*Not Available For Discounts **15mm orifice opener available in #09.09

EDGEONYX™



X-FIND Orifice Opener - 6 Pack		
Size	Color	15mm only
09.09	 YELLOW	EO090915HTX

6 Pack				Points	
Size	Color	21mm	25mm	Paper Points*	Gutta Percha
X-FOLLOW Glidepath Options					
15.04	<input type="checkbox"/> WHITE	EO150421HTX	EO150425HTX	N/A	N/A
17.04	 PURPLE	EO170421HTX	EO170425HTX	SPP1504-1	X71704GP-1
X-FINISH Shaping Options					
20.04	 YELLOW	EO200421HTX	EO200425HTX	SPP1504-1	X72004GP-1
25.04	 RED	EO250421HTX	EO250425HTX	SPP2004-1	X72504GP-1
30.04	 BLUE	EO300421HTX	EO300425HTX	SPP2504-1	X73004GP-1
35.04	 GREEN	EO350421HTX	EO350425HTX	SPP3004-1	X73504GP-1
40.04	 BLACK	EO400421HTX	EO400425HTX	SPP3504-1	X74004GP-1
45.04	<input type="checkbox"/> WHITE	EO450421HTX	EO450425HTX	SPP4004-1	X74504GP-1

*It Is Recommended To Drop Down One Tip Size For Paper Points






Pair EdgeOnyx™ with:



BioCeramic Sealer 2g

818365020072



Assorted Packs - \$28.80/ 6 Pack		Points	
X-Find (09/.09) 15mm, X-Follow (15/.04), X-Finish (20-30/.04)		Paper Points	Gutta Percha
 <input type="checkbox"/>    			
21mm	EOASST21HTX	SPP04AST-1	X7GP04AST-1
25mm	EOASST25HTX		

26-0705 | 2026-26-03 | Rev 0



EDGE X7 UTOPIA™



\$38.80
6 Pack

Enhanced Technology. Remarkably Beautiful.

- Enhanced version of EdgeFile® X7
- Features proprietary FireWire™ Blaze heat treatment for better cutting efficiency
- Constant taper file
- Sterile, ready-to-use



- Parabolic shaped cross-section
- Available in 3 lengths: 21, 25 & 29mm
- File Size: 04 & 06
- Maximum flute diameter 1mm
- Speed: 500 RPM
- Torque: 4.0 Ncm



High Cutting Efficiency



Pre-Bendable



High Strength

EDGE X7



UTOPIA™



Sterile 6 Pack					Points		Universal Obturators	
Size	Color	21mm	25mm	29mm	Standard .04/.06 Paper Points	X7 Gutta Percha	EdgeCore	EdgeFill
17.04	PURPLE	X7U170421SE	X7U170425SE	X7U170429SE	SPP1504-1	X71704GP-1	N/A	N/A
20.04	YELLOW	X7U200421SE	X7U200425SE	X7U200429SE	SPP2004-1	X72004GP-1	X72006C	X72006O
25.04	RED	X7U250421SE	X7U250425SE	X7U250429SE	SPP2504-1	X72504GP-1	X72506C	X72506O
30.04	BLUE	X7U300421SE	X7U300425SE	X7U300429SE	SPP3004-1	X73004GP-1	X73006C	X73006O
35.04	GREEN	X7U350421SE	X7U350425SE	X7U350429SE	SPP3504-1	X73504GP-1	X73506C	X73506O
40.04	BLACK	X7U400421SE	X7U400425SE	X7U400429SE	SPP4004-1	X74004GP-1	X74006C	X74006O
45.04	WHITE	X7U450421SE	X7U450425SE	X7U450429SE	SPP4504-1	X74504GP-1	X74506C	X74506O
17.06	PURPLE	X7U170621SE	X7U170625SE	X7U170629SE	SPP1506-1	X71706GP-1	N/A	N/A
20.06	YELLOW	X7U200621SE	X7U200625SE	X7U200629SE	SPP2006-1	X72006GP-1	X72006C	X72006O
25.06	RED	X7U250621SE	X7U250625SE	X7U250629SE	SPP2506-1	X72506GP-1	X72506C	X72506O
30.06	BLUE	X7U300621SE	X7U300625SE	X7U300629SE	SPP3006-1	X73006GP-1	X73006C	X73006O
35.06	GREEN	X7U350621SE	X7U350625SE	X7U350629SE	SPP3506-1	X73506GP-1	X73506C	X73506O
40.06	BLACK	X7U400621SE	X7U400625SE	X7U400629SE	SPP4006-1	X74006GP-1	X74006C	X74006O
45.06	WHITE	X7U450621SE	X7U450625SE	X7U450629SE	SPP4506-1	X74506GP-1	X74506C	X74506O
Assorted Packs .04 \$38.80/ 6 Pack					Paper Points	Gutta Percha	See page 70 for 20pk obturators	
17, 20, 25, 25, 30, 35	X7UASST0421SE	X7UASST0425SE	N/A	SPP04AST-1	X7GP04AST-1			
Assorted Packs .06 \$38.80/ 6 Pack					Paper Points	Gutta Percha	See page 70 for 20pk obturators	
17, 20, 25, 25, 30, 35	X7UASST0621SE	X7UASST0621SE	N/A	SPP06AST-1	X7GP06AST-1			

*Stepping down one size on paper points from size listed may be desired.

Pair EdgeDiamond Utopia™ with:

EDGE BIO CERAMIC™

SEALER

818365020072




EDGEFILE[®] X7




\$38.80
6 Pack

Unmatched Strength and Flexibility

- Universal file system
- #1 Fan Favorite
- FireWire™ NiTi heat treatment provides super elasticity with no shape memory
- Constant tapered rotary files
- Sterile, ready-to-use blister packaging*
- Tighter pitch and flutes which deflect more torsional forces along the length of the file
- Non-cutting tip with flutes extending to 0.25mm—offering precision and safety in one design



- Parabolic shaped cross-section
- Available in 3 lengths: 21, 25 & 29mm
- Constant taper: .04 and .06
- Rotary files
- Speed: 450 RPM
- Torque: 2.9 NCm



Capable of 90° Curves



Pre-Bendable



High Strength

* Product Packaging Upgrade: X7 Files will be transitioning to sterile, ready-to-use blister packs. A new [REF] will appear on the back top left of the blister next to the [REF] symbol which has been updated from ending "RF" to "HTS" representing "Heat-Treated Sterile" or "ASST" for assorted packs. Invoices and packing slips will continue to display former REF ending "RF" as featured in table on the following page. Searching either number when ordering will return the correct product.

EDGEFILE® X7



Sterile 6 Pack					Points		Universal Obturators	
Size	Color	21mm	25mm	29mm	Standard .04/.06 Paper Points	X7 Gutta Percha	EdgeCore	EdgeFill
17.04	WHITE	X7170421RF	X7170425RF	X7170429RF	SPP1504-1	X71704GP-1	N/A	N/A
20.04	YELLOW	X7200421RF	X7200425RF	X7200429RF	SPP2004-1	X72004GP-1	X72006C	X72006O
25.04	RED	X7250421RF	X7250425RF	X7250429RF	SPP2504-1	X72504GP-1	X72506C	X72506O
30.04	BLUE	X7300421RF	X7300425RF	X7300429RF	SPP3004-1	X73004GP-1	X73006C	X73006O
35.04	GREEN	X7350421RF	X7350425RF	X7350429RF	SPP3504-1	X73504GP-1	X73506C	X73506O
40.04	BLACK	X7400421RF	X7400425RF	X7400429RF	SPP4004-1	X74004GP-1	X74006C	X74006O
45.04	WHITE	X7450421RF	X7450425RF	X7450429RF	SPP4504-1	X74504GP-1	X74506C	X74506O
17.06	WHITE	X7170621RF	X7170625RF	X7170629RF	SPP1506-1	X71706GP-1	N/A	N/A
20.06	YELLOW	X7200621RF	X7200625RF	X7200629RF	SPP2006-1	X72006GP-1	X72006C	X72006O
25.06	RED	X7250621RF	X7250625RF	X7250629RF	SPP2506-1	X72506GP-1	X72506C	X72506O
30.06	BLUE	X7300621RF	X7300625RF	X7300629RF	SPP3006-1	X73006GP-1	X73006C	X73006O
35.06	GREEN	X7350621RF	X7350625RF	X7350629RF	SPP3506-1	X73506GP-1	X73506C	X73506O
40.06	BLACK	X7400621RF	X7400625RF	X7400629RF	SPP4006-1	X74006GP-1	X74006C	X74006O
45.06	WHITE	X7450621RF	X7450625RF	X7450629RF	SPP4506-1	X74506GP-1	X74506C	X74506O
Assorted Pack .04 <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> \$21.34/ 6 Pack					Paper Points	Gutta Percha	See page 70 for 20pk obturators	
17, 20, 25, 30, 35, 40		ASSTX70421	ASSTX70425	ASSTX70429	SPP04AST-1	X7GP04AST-1		
Assorted Pack .06 <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> \$21.34/ 6 Pack					Paper Points	Gutta Percha	See page 70 for 20pk obturators	
17, 20, 25, 30, 35, 40		ASSTX70621	ASSTX70625	ASSTX70629	SPP06AST-1	X7GP06AST-1		

*Stepping down one size on paper points from size listed may be desired.

Pair EdgeFile® X7 with:

EDGE BIOCERAMIC™
SEALER





EDGEFILE[®] X3



\$38.80
6 Pack

Engineered for Strength with Flexible Precision

- Designed with a Safe-Core™ parabolic cross-section
- .04 or .06 constant taper file until flute diameter reaches 1.0 Then taper is 0
- Safe unwinding feature
- Super elasticity. No shape memory
- FireWire NiTi™ heat treated
- Variable tapered rotary files
- Features its own dedicated glidepath system
- Sterile, ready-to-use blister packaging*
- Alternative to ProTaper® Next¹ systems



- Parabolic shaped cross-section
- Available in 3 lengths: 21, 25 & 29mm
- Speed 300-500 RPM
- Torque 2.94 Ncm



Highly Flexible




Pre-Bendable



High Strength

1) ProTaper® Next is a registered trademark of Dentsply Sirona. * Product Packaging Upgrade: X7 Files will be transitioning to sterile, ready-to-use blister packs. A new [REF] will appear on the back top left of the blister next to the [REF] symbol which has been updated from ending "RF" to "HTS" representing "Heat-Treated Sterile" or "ASST" for assorted packs. Invoices and packing slips will continue to display former REF ending "RF" as featured in table on the following page. Searching either number when ordering will return the correct product.

EDGEFILE® X3



Sterile 6 Pack					Points		Universal Obturators	
Size/Taper	Color	21mm	25mm	29mm	Standard .06 Paper Points*	X7 Gutta Percha	EdgeCore	EdgeFill
N1 17.06	PURPLE	X3170621RF	X3170625RF	X3170629RF	SPP1506-1	X71706GP-1	N/A	N/A
N2 17.04	WHITE	X3170421RF	X3170425RF	X3170429RF	SPP1504-1	X71704GP-1	N/A	N/A
C1 20.06	YELLOW	X3200621RF	X3200625RF	X3200629RF	SPP2006-1	X72006GP-1	X72006C	X72006O
C2 25.06	RED	X3250621RF	X3250625RF	X3250629RF	SPP2506-1	X72506GP-1	X72506C	X72506O
C3 30.06	BLUE	X3300621RF	X3300625RF	X3300629RF	SPP3006-1	X73006GP-1	X73006C	X73006O
C4 40.06	BLACK	X3400621RF	X3400625RF	X3400629RF	SPP4006-1	X74006GP-1	X74006C	X74006O

*Stepping down one size on paper points from size listed may be desired.

Pair EdgeFile® X3 with:

EDGEBIOCERAMIC™
SEALER

818365020072



EdgeVerify™ for EdgeFile® X3 6 Pack	
Size	Item Code
20	X32006V
25	X32506V
30	X33006V
40	X34006V
50	X35006V


EDGESEQUEL UTOPIA™



\$38.80
6 Pack

Enhanced Technology. Remarkably Beautiful.

- Features proprietary FireWire™ Blaze Blue heat treatment
- Improved cutting efficiency compared to EdgeSequel Sapphire™
- 12mm handle length for better posterior access
- Constant taper file
- Sterile, ready-to-use
- Alternative to EndoSequence® and Vortex®



- Triangular cross section
- Available in 3 lengths: 21, 25 & 31 mm w/ 19mm orifice opener
- Speed: 800-1000 RPM
- Torque: 1.5 Ncm



Flexible



High Cutting Efficiency



Electropolished Finish




High Strength

*19mm orifice opener available in #25/.08



EDGESEQUEL



UTOPIA™



Sterile 6 Pack								
Size/Taper	Color	19mm						
25.08	RED	ES250819SE						
Sterile 6 Pack				Points		Obturators		
Size/Taper	Color	21mm	25mm	31mm	Standard .04/.06 Points	Gutta Percha	EdgeCore	EdgeFill
10.02	PURPLE	ES100221SE	ES100225SE	N/A	N/A	N/A	N/A	N/A
15.02	WHITE	ES150221SE	ES150225SE	N/A	N/A	N/A	N/A	N/A
20.02	YELLOW	ES200221SE	ES200225SE	N/A	N/A	N/A	N/A	N/A
30.02	BLUE	ES300221SE	ES300225SE	N/A	N/A	N/A	N/A	N/A
10.04	PURPLE	ES100421SE	ES100425SE	N/A	N/A	N/A	N/A	N/A
15.04	WHITE	ES150421SE	ES150425SE	ES150431SE	SPP1504-1	SGP1504-1	N/A	N/A
20.04	YELLOW	ES200421SE	ES200425SE	N/A	SPP2004-1	SGP2004-1	ES2006C	ES2006O
25.04	RED	ES250421SE	ES250425SE	ES250431SE	SPP2504-1	SGP2504-1	ES2506C	ES2506O
30.04	BLUE	ES300421SE	ES300425SE	ES300431SE	SPP3004-1	SGP3004-1	ES3006C	ES3006O
35.04	GREEN	ES350421SE	ES350425SE	N/A	SPP3504-1	SGP3504-1	ES3506C	ES3506O
40.04	BLACK	ES400421SE	ES400425SE	ES400431SE	SPP4004-1	SGP4004-1	ES4006C	ES4006O
50.04	YELLOW	ES500421SE	ES500425SE	ES500431SE	SPP5004-1	SGP5004-1	ES5006C	ES5006O
15.06	WHITE	ES150621SE	ES150625SE	N/A	SPP1506-1	SGP1506-1	N/A	N/A
25.06	RED	ES250621SE	ES250625SE	ES250631SE	SPP2506-1	SGP2506-1	ES2506C	ES2506O
30.06	BLUE	ES300621SE	ES300625SE	ES300631SE	SPP3006-1	SGP3006-1	ES3006C	ES3006O
35.06	GREEN	ES350621SE	ES350625SE	ES350631SE	SPP3506-1	SGP3506-1	ES3506C	ES3506O
40.06	BLACK	ES400621SE	ES400625SE	N/A	SPP4006-1	SGP4006-1	ES4006C	ES4006O

Assorted Packs	21mm	25mm	31mm	Paper Points	Gutta Percha
15.04, 20.04, 25.04, 30.04, 35.04 & 25.08 19mm  \$38.80/ 6 pack	V2ASST0421SE	V2ASST0425SE	V2ASST0431SE	SPP04AST-1 \$8.65	SGP04AST-1 \$14.05
15.04, 25.04, 30.04  \$27.95/ 3 pack	ASST15300421SE	ASST15300425SE	ASST15300431SE		
15.04, 25.04, 25.06  \$27.95/ 3 pack	ASST15250621SE	ASST15250625SE	ASST15250631SE		

Pair EdgeSequel Utopia™ with:

EDGE  UTOPIA™

BioCeramic Sealer 2g

818365020072




EDGESEQUEL SAPPHIRE™



\$27.95
4 Pack

Bodybuilder Strength. Gymnast Flexibility.

- Excellent cyclic fatigue compared to other leading files*
- FireWire™ NiTi heat treated
- Constant tapered files
- Sterile, ready-to-use blister packaging**
- Alternative to Vortex®, VortexBlue®, and EndoSequence®



- Triangular shaped cross-section
- Available in 3 lengths: 21, 25 & 31 mm
- Speed: 500 RPM
- Torque: 4.0-5.2 Ncm



Highly Flexible




Pre-Bendable



High Strength

Cyclic Fatigue Testing - Average Time to Failure in Seconds



*Vortex®, Vortex® Blue, are registered trademarks of Dentsply Tulsa Dental Specialties. EndoSequence® is a registered trademark of Brassler USA. **Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph: <http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit www.edgeendo.com/matter-of-temperature. **Product Packaging Upgrade: EdgeSequel Sapphire™ Files will be transitioning to sterile, ready-to-use blister packs. A new [REF] will appear on the back top left of the blister next to the [REF] symbol which has been updated from ending "RF" to "HTS" representing "Heat-Treated Sterile" or "ASST" for assorted packs. Invoices and packing slips will continue to display former REF ending "RF" as featured in table on the following page. Searching either number when ordering will return the correct product.



4 Pack					Points		EdgeSequel Obturators	
Size	Color	21mm	25mm	31mm	Standard .04/.06 Paper Points*	Gutta Points	EdgeCore	EdgeFill
25.08	RED	ES250821	ES250825	ES250831	N/A	N/A	ES2506C	ES2506O
15.04	WHITE	ES150421	ES150425	ES150431	SPP1504-1	SGP1504-1	N/A	N/A
20.04	YELLOW	ES200421	ES200425	ES200431	SPP2004-1	SGP2004-1	ES2006C	ES2006O
25.04	RED	ES250421	ES250425	ES250431	SPP2504-1	SGP2504-1	ES2506C	ES2506O
30.04	BLUE	ES300421	ES300425	ES300431	SPP3004-1	SGP3004-1	ES3006C	ES3006O
35.04	GREEN	ES350421	ES350425	ES350431	SPP3504-1	SGP3504-1	ES3506C	ES3506O
40.04	BLACK	ES400421	ES400425	ES400431	SPP4004-1	SGP4004-1	ES4006C	ES4006O
45.04	WHITE	ES450421	ES450425	ES450431	SPP4504-1	SGP4504-1	ES4506C	ES4506O
50.04	YELLOW	ES500421	ES500425	ES500431	SPP5004-1	SGP5004-1	ES5006C	ES5006O
55.04	RED	ES550421	ES550425	ES550431	ESS5504PP	ESS5504GP	ES5506C	ES5506O
60.04	BLUE	ES600421	ES600425	ES600431	ESS6004PP	ESS6004GP	ES6006C	ES6006O
70.04	GREEN	ES700421	ES700425	ES700431	ESS7004PP	ESS7004GP	ES7006C	ES7006O
80.04	BLACK	ES800421	ES800425	ES800431	ESS8004PP	ESS8004GP	ES8006C	ES8006O
90.04	WHITE	ES900421	ES900425	ES900431	N/A	N/A	ES8006C	ES8006O
15.06	WHITE	ES150621	ES150625	ES150631	SPP1506-1	SGP1506-1	ES2006C	ES2006O
20.06	YELLOW	ES200621	ES200625	ES200631	SPP2006-1	SGP2006-1	ES2006C	ES2006O
25.06	RED	ES250621	ES250625	ES250631	SPP2506-1	SGP2506-1	ES2506C	ES2506O
30.06	BLUE	ES300621	ES300625	ES300631	SPP3006-1	SGP3006-1	ES3006C	ES3006O
35.06	GREEN	ES350621	ES350625	ES350631	SPP3506-1	SGP3506-1	ES3506C	ES3506O
40.06	BLACK	ES400621	ES400625	ES400631	SPP4006-1	SGP4006-1	ES4006C	ES4006O
45.06	WHITE	ES450621	ES450625	ES450631	SPP4506-1	SGP4506-1	ES4506C	ES4506O
50.06	YELLOW	ES500621	ES500625	ES500631	SPP5006-1	SGP5006-1	ES5006C	ES5006O

*Stepping down one size on paper points from size listed may be desired

Assorted Packs .04	21mm	25mm	31mm
15, 20, 25, 30	\$27.95/ 4 pack	ASST15300421	ASST15300431
35, 40, 45, 50	\$27.95/ 4 pack	ASST35500421	ASST35500431
Assorted Packs .06	21mm	25mm	31mm
15, 20, 25, 30	\$27.95/ 4 pack	ASST15300621	ASST15300631
35, 40, 45, 50	\$27.95/ 4 pack	ASST35500621	ASST35500631
Assorted Paper Points .04 (20, 25, 30, 35, 40, 45)	\$8.65	SPP04AST-1	
Assorted Gutta Percha .04 (20, 25, 30, 35, 40, 45)	\$14.05	SGP04AST-1	
Assorted Paper Points .06 (20, 25, 30, 35, 40, 45)	\$8.65	SPP06AST-1	
Assorted Gutta Percha .06 (20, 25, 30, 35, 40, 45)	\$14.05	SGP06AST-1	

EDGE EVOLVE™

The Evolution of the Revolution



\$38.80
6 Pack

- Minimal dentin with a conservative flute and taper. All EdgeEvolve files, even our largest 40.10, only have a max diameter of 1mm. FireWire™ NiTi Heat Treated
- Extremely flexible constant taper file with a non-cutting tip
- Available in .04, .06, .08, .10 tapers
- Alternative to GTX®¹
- Non-sterile



- Triangular shaped cross-section
- Available in 3 lengths: 21, 25 & 29mm
- Speed: 300-500 RPM
- Torque: 2.94 Ncm



Extremely Flexible




Pre-Bendable







High Strength

1) GTX® is a registered trademark of Dentsply Tulsa Dental Specialties.

EDGE EVOLVE™



6 Pack					Points		Obturators	
Size/Taper	Color	21mm	25mm	29mm	Standard .04/.06 Paper Points	Gutta Percha	EdgeCore	EdgeFill
17.04	WHITE	EE170421	EE170425	EE170429	SPP1504-1	SGP1504-1	N/A	N/A
17.06	WHITE	EE170621	EE170625	EE170629	SPP1506-1	SGP1506-1		
17.08	WHITE	EE170821	EE170825	EE170829	N/A	N/A		
17.10	WHITE	EE171021	EE171025	EE171029	N/A	N/A		
20.04	YELLOW	EE200421	EE200425	EE200429	SPP2004-1	SGP2004-1	EE2006C	EE2006O
20.06	YELLOW	EE200621	EE200625	EE200629	SPP2006-1	SGP2006-1		
20.08	YELLOW	EE200821	EE200825	EE200829	N/A	N/A		
20.10	YELLOW	EE201021	EE201025	EE201029	N/A	N/A		
25.04	RED	EE250421	EE250425	EE250429	SPP2504-1	SGP2504-1	EE2506C	EE2506O
25.06	RED	EE250621	EE250625	EE250629	SPP2506-1	SGP2506-1		
25.08	RED	EE250821	EE250825	EE250829	N/A	N/A		
25.10	RED	EE251021	EE251025	EE251029	N/A	N/A		
30.04	BLUE	EE300421	EE300425	EE300429	SPP3004-1	SGP3004-1	EE3006C	EE3006O
30.06	BLUE	EE300621	EE300625	EE300629	SPP3006-1	SGP3006-1		
30.08	BLUE	EE300821	EE300825	EE300829	N/A	N/A		
30.10	BLUE	EE301021	EE301025	EE301029	N/A	N/A		
35.04	GREEN	EE350421	EE350425	EE350429	SPP3504-1	SGP3504-1	EE3506C	EE3506O
35.06	GREEN	EE350621	EE350625	EE350629	SPP3506-1	SGP3506-1		
35.08	GREEN	EE350821	EE350825	EE350829	N/A	N/A		
35.10	GREEN	EE351021	EE351025	EE351029	N/A	N/A		
40.04	BLACK	EE400421	EE400425	EE400429	SPP4004-1	SGP4004-1	EE4006C	EE4006O
40.06	BLACK	EE400621	EE400625	EE400629	SPP4006-1	SGP4006-1		
40.08	BLACK	EE400821	EE400825	EE400829	N/A	N/A		
40.10	BLACK	EE401021	EE401025	EE401029	N/A	N/A		

Assorted Packs			21mm	25mm	29mm
20 (20.04, 20.06, 20.08, 20.10)  \$38.80/ 4 pack			EEASST2021	EEASST2025	EEASST2029
30 (30.04, 30.06, 30.08, 30.10)  \$38.80/ 4 pack			EEASST3021	EEASST3025	EEASST3029
25 (25.04, 25.06, 25.08, 25.10)  \$38.80/ 4 pack			EEASST2521	EEASST2525	EEASST2529
40 (40.04, 40.06, 40.08, 40.10)  \$38.80/ 4 pack			EEASST4021	EEASST4025	EEASST4029

Pair EdgeEvolve™ with:

EDGE BIO CERAMIC™
SEALER

818365020072






Elegance in Innovation, Strength in Versatility

- Minimally Invasive Progressive/Regressive Variable Tapered File System
- Ability to perform in **both rotary or forward reciprocation** settings
- Rotary Motion: Delivers maximum cutting efficiency and debris removal
- Forward Reciprocation Motion: Maximizes safety while still allowing for efficient cutting
- Unique geometry, innovative cutting angles and optimized FireWire™ Blaze Heat-Treatment** allows files to navigate curvatures with ease while maintaining dentin preservation in the coronal third
- Alternative to ProTaper Ultimate™* & Minimally Invasive Alternative to ProTaper™*

\$38.80
6 Pack



- Diamond cross-section
- Available in 21mm and 25mm lengths
- Progressive/Regressive Variable Taper File
- Sterile, ready-to-use

Non-continuous rotation: (+340° / -20°)
Recommended speed: 350 - 500 rpm
Torque: 3.5 Ncm



Incredibly Strong



Very Versatile



Minimally Invasive

*ProTaper Ultimate™ and ProTaper™ are registered trademarks of Dentsply Sirona. **Slider is not heat-treated Best Product Rating Dental Product Shopper Vol 19 No 2

EDGE DIAMOND



Slider 16.02v (not heat-treated)

Shaper 20.04v

F1 20.07v

F2 25.07v

F3 30.07v

FX 35.07v

Sterile 6 Pack				EdgeTaper™ Points		EdgeTaper™ Obturators		
Size/Taper		Color	21mm	25mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
Slider*	16.02	PURPLE	EDU1602V21	EDU1602V25	N/A	N/A	N/A	N/A
Shaper	20.04	WHITE	EDU2004V21	EDU2004V25	N/A	N/A	N/A	N/A
F1	20.07	YELLOW	EDU2007V21	EDU2007V25	ETPPF1-1	ETGPF1-1	ET2006C	ET2006O
F2	25.07	RED	EDU2507V21	EDU2507V25	ETPPF2-1	ETGPF2-1	ET2506C	ET2506O
F3	30.07	BLUE	EDU3007V21	EDU3007V25	ETPPF3-1	ETGPF3-1	ET3006C	ET3006O
FX	35.07	GREEN	EDU3507V21	EDU3507V25	ETPPF4-1	ETGPF4-1	ET4006C	ET4006O
Assorted 6 Pack				ETPPAST	ETGPAST			
Slider*, Shaper, F1, F2, F3, FX		EDUASST21V	EDUASST25V					

*Slider is not Heat Treated

Pair EdgeDiamond Utopia™ with:



BioCeramic Sealer 2g

818365020072




EDGE V-TAPER^{HT}



\$38.80
6 Pack

Strong. Advanced. Flexible. Efficient.

- Designed with a Safe-Core™ parabolic cross-section
- The Safe-Core™ parabolic cross-section combines the attributes of being a highly efficient and flexible instrument while being resistant to fracture, minimizing your stress and resulting in a more reliable procedure.
- FireWire™ NiTi heat treated
- Variable tapered rotary files
- Features its own dedicated glidepath system
- Non-sterile
- Alternative to DC Taper®¹



- Parabolic shaped cross-section
- Available in 3 lengths: 21, 25 & 31mm
- Speed 300-400 RPM
- Torque - 4.0- 5.2 Ncm



Highly Flexible



Pre-Bendable



High Strength

1) DC Taper® is a registered trademark of SS White Dental®.


EDGETAPER BLAZE UTOPIA™



Enhanced Technology. Remarkably Beautiful.

\$38.80
6 Pack

- Features proprietary FireWire™ Blaze heat treatment
- Electropolished for better cutting efficiency
- Shorter handle length for better posterior access
- Non-cutting tip easily glides to the working length
- Sterile, ready-to-use
- Alternative to ProTaper Gold®
- Optimized version of EdgeTaper Platinum™



- Triangular shaped cross-section
- Available in 3 lengths: 21, 25 & 31mm
- Speed: 300 RPM
- Torque: 5.1 Ncm (SX & S1), 2.5 Ncm (S2 & F1), 3.1 Ncm (F2 & F3)



High
Cutting
Efficiency



High
Strength



Flexible

*ProTaper® Gold is a registered trademark of Dentsply Tulsa Dental Specialties


EDGETAPER BLAZE UTOPIA™



Enhanced Technology. Remarkably Beautiful.

\$38.80
6 Pack

- Features proprietary FireWire™ Blaze heat treatment
- Electropolished for better cutting efficiency
- Shorter handle length for better posterior access
- Non-cutting tip easily glides to the working length
- Sterile, ready-to-use
- Alternative to ProTaper Gold®
- Optimized version of EdgeTaper Platinum™



- Triangular shaped cross-section
- Available in 3 lengths: 21, 25 & 31mm
- Speed: 300 RPM
- Torque: 5.1 Ncm (SX & S1), 2.5 Ncm (S2 & F1), 3.1 Ncm (F2 & F3)



High
Cutting
Efficiency



High
Strength



Flexible

*ProTaper® Gold is a registered trademark of Dentsply Tulsa Dental Specialties

EDGETAPER BLAZE



Sterile 6 Pack		
Size	Color	19mm
SX	YELLOW	ETBUSX19SU

Sterile 6 Pack					EdgeTaper™ Points		EdgeTaper™ Obturators	
Size	Color	21mm	25mm	31mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
S1	PURPLE	ETBUS121SU	ETBUS125SU	ETBUS131SU	N/A	N/A	N/A	N/A
S2	WHITE	ETBUS221SU	ETBUS225SU	ETBUS231SU	N/A	N/A	N/A	N/A
F1	YELLOW	ETBUF121SU	ETBUF125SU	ETBUF131SU	ETPPF1-1	ETGPF1-1	ET2006C	ET2006O
F2	RED	ETBUF221SU	ETBUF225SU	ETBUF231SU	ETPPF2-1	ETGPF2-1	ET2506C	ET2506O
F3	BLUE	ETBUF321SU	ETBUF325SU	ETBUF331SU	ETPPF3-1	ETGPF3-1	ET3006C	ET3006O
Assorted Pack 6 Pack					Paper Points	Gutta Percha		
(SX, S1, S2, F1, F2, F3)		ETBUA21SU	ETBUA25SU	ETBUA31SU	ETPPAST	ETGPAST		

Pair EdgeTaper Blaze Utopia™ with:



BioCeramic Sealer 2g

818365020072





EDGE TAPER PLATINUM™



\$38.80
6 Pack

Edge-ified Sequence with Enhanced Strength and Flexibility

- 2x the cyclic fatigue as ProTaper Gold®
- 6x the cyclic fatigue as ProTaper®
- FireWire™ NiTi heat treated
- Variable tapered files
- Non-sterile



- Convex triangular cross-section
- Available in 3 lengths: 21, 25 & 31mm
- Speed: 300-400 RPM
- Torque: 4.0-5.2 Ncm



Highly Flexible



Pre-Bendable

















High Strength

*ProTaper® Gold and ProTaper are registered trademarks of Dentsply Tulsa Dental Specialties

EDGETAPER PLATINUM™



6 Pack		
Size	Color	19mm
SX	 YELLOW	ETSX19HT

6 Pack					EdgeTaper™ Points		EdgeTaper™ Obturators		
Size	Color	21mm	25mm	31mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill	
S1	 PURPLE	ETS121HT	ETS125HT	ETS131HT	N/A				
S2	 WHITE	ETS221HT	ETS225HT	ETS231HT	N/A				
F1	 YELLOW	ETS121HT	ETF125HT	ETF131HT	ETPPF1-1	ETGPF1-1	ET2006C	ET2006O	
F2	 RED	ETF221HT	ETF225HT	ETF231HT	ETPPF2-1	ETGPF2-1	ET2506C	ET2506O	
F3	 BLUE	ETF321HT	ETF325HT	ETF331HT	ETPPF3-1	ETGPF3-1	ET3006C	ET3006O	
F4	 BLACK	ETF421HT	ETF425HT	ETF431HT	ETPPF4-1	ETGPF4-1	ET4006C	ET4006O	
F5	 YELLOW	ETF521HT	ETF525HT	ETF531HT	ETPPF5-1	ETGPF5-1	ET5006C	ET5006O	
Assorted 6 pack		     	\$38.80			Paper Points	Gutta Percha		
(SX, S1, S2, F1, F2, F3)		ETP21MM	ETP25MM	ETP31MM	ETPPAST	ETGPAST			

Pair EdgeTaper Platinum™ with:

EDGE BIOCERAMIC™
SEALER

818365020072





EDGETAPER™



\$38.80
6 Pack

Edge-ified Sequence with Enhanced Strength

- 6x The Cyclic Fatigue as ProTaper®¹
- Variable Tapered Files
- Non-Sterile
- Non-Heat Treated



- Convex triangular cross-section
- Available in 3 lengths: 21, 25 & 31mm
- Speed: 300-400 RPM
- Torque: 4.0-5.2 Ncm




Pre-Bendable



High Strength

Cyclic Fatigue Testing - Average Time to Failure in Seconds



EDGE TAPER™

Non-Heat Treated



6 Pack		
Size	Color	19mm
SX	YELLOW	ETSX19

6 Pack					EdgeTaper™ Points		EdgeTaper™ Obturators	
Size	Color	21mm	25mm	31mm	Paper Points	Gutta Percha	EdgeCore	EdgeFill
S1	PURPLE	ETS121	ETS125	ETS131	N/A			
S2	WHITE	ETS221	ETS225	ETS231	N/A			
F1	YELLOW	ETF121	ETF125	ETF131	ETPPF1-1	ETGPF1-1	ET2006C	ET2006O
F2	RED	ETF221	ETF225	ETF231	ETPPF2-1	ETGPF2-1	ET2506C	ET2506O
F3	BLUE	ETF321	ETF325	ETF331	ETPPF3-1	ETGPF3-1	ET3006C	ET3006O
F4	BLACK	ETF421	ETF425	ETF431	ETPPF4-1	ETGPF4-1	ET4006C	ET4006O
F5	YELLOW	ETF521	ETF525	ETF531	ETPPF5-1	ETGPF5-1	ET5006C	ET5006O
Assorted 6 pack		\$38.80			Paper Points	Gutta Percha		
(SX, S1, S2, F1, F2, F3)		ETASST21	ETASST25	ETASST31	ETPPAST	ETGPAST		

Pair EdgeTaper™ with:

EDGE BIO CERAMIC™
SEALER

818365020072




EDGEONE BLAZE UTOPIA™

Enhanced Technology. Remarkably Beautiful.



\$27.95
3 Pack

- Premier single file shaping reciprocating system
- Features proprietary firewire blaze heat treatment
- Electropolished for better cutting efficiency
- Shorter handle length for better posterior access
- Variable taper file
- Sterile, ready-to-use
- Alternative to WaveOne Gold® and EdgeOne Fire™



- S-shaped cross-section - for increased debris removal
- Available in 3 lengths: 21, 25 & 31mm
- CCW 170*/ CW 60*



Flexible




High Strength



High Cutting Efficiency

VOL. 18 NO. 4 Dental Product Shopper. *WaveOne® Gold is a registered trademarks of Dentsply Tulsa Dental Specialties.

EDGEONE BLAZE



Sterile 3 Pack					Points		EdgeOne Obturators	
Size	Color	21mm	25mm	31mm	Standard .06 Paper Points	EdgeOne-R Gutta Points	EdgeCore	EdgeFill
R-GLIDE	<input type="checkbox"/> WHITE	EOBURG1421SU	EOBURG1425SU	EOBURG1431SU	N/A		N/A	
R20 SM	<input checked="" type="checkbox"/> YELLOW	EOBU2021SU	EOBU2025SU	EOBU2031SU	SPP2006-1	EORGP20-1	EOPSM20C	EOPSM200
R25 PRI	<input checked="" type="checkbox"/> RED	EOBU2521SU	EOBU2525SU	EOBU2531SU	SPP2506-1	EORGP25-1	EOPPRI25C	EOPPRI250
R35 MED	<input checked="" type="checkbox"/> GREEN	EOBU3521SU	EOBU3525SU	EOBU3531SU	SPP3506-1	EORGP35-1	EOPMED35C	EOPMED350
R45 LG	<input type="checkbox"/> WHITE	EOBU4521SU	EOBU4525SU	EOBU4531SU	SPP4506-1	EORGP45-1	EOPLG45C	EOPLG450
Assorted Pack* <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>					\$37.65/ 4 Pack			
(R20, R25, R35, R45)		EOBUA21SU	EOBUA25SU	EOBUA31SU				

*R-Glide Sold Separately

Pair EdgeOne Blaze Utopia™ with:

EDGE  **UTOPIA™**
BioCeramic Sealer 2g

818365020072




EDGEONE FIRE™



\$27.95
3 Pack

The Standard in Reciprocating Technology

- 5x the cyclic fatigue as WaveOne® Gold
- FireWire™ NiTi heat treated
- Variable tapered reverse reciprocating files
- Non-sterile
- Alternative to WaveOne® Gold



- Parallelogram shaped cross-section
- Available in 3 lengths: 21, 25 & 31mm
- Designed to shape canals in reverse-reciprocating motion
CCW 150*/CW 60* or WaveOne® Gold¹ settings



Highly Flexible



Pre-Bendable




High Strength

Cyclic Fatigue Testing - Average Time to Failure in Seconds

EDGEONE
FIRE™



WaveOne® Gold



WaveOne® Gold is a registered trademarks of Dentsply Tulsa Dental Specialties. *Based on resistance cyclic testing as shown on the FireWire NiTi Strength Graph:<http://edgeendo.com/comparative-study-of-cyclic-fatigue-resistance/> and price comparisons vs retail price. For more information, please visit www.edgeendo.com/matter-of-temperature



3 Pack				
Size	Color	21mm	25mm	31mm
GLIDEPATH	<input type="checkbox"/> WHITE	EOFGP21	EOFGP25	EOFGP31

3 Pack					Points		EdgeOne Obturators	
Size	Color	21mm	25mm	31mm	EdgeOne Fire™ Paper Points	Gutta Percha	EdgeCore	EdgeFill
SM (20)	<input checked="" type="checkbox"/> YELLOW	EOFSM21	EOFSM25	EOFSM31	EOFSMPP-1	EOFSMGP-1	EOPSM20C	EOPSM20O
PRI (25)	<input checked="" type="checkbox"/> RED	EOFPRI21	EOFPRI25	EOFPRI31	EOFPRIPP-1	EOFPRIGP-1	EOPPRI25C	EOPPRI25O
MED (35)	<input checked="" type="checkbox"/> GREEN	EOFMED21	EOFMED25	EOFMED31	EOFMEDPP-1	EOFMEDGP-1	EOPMED35C	EOPMED35O
LG (45)	<input type="checkbox"/> WHITE	EOFLG21	EOFLG25	EOFLG31	EOFLGGPP-1	EOFLGGP-1	EOPLG45C	EOPLG45O
Assorted 4 pack <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>								
					\$36.50			
(SM, PRI, MED, LG)		ASSTEOF21	ASSTEOF25	ASSTEOF31	EOFPAST-1	EOFGPAST-1		

Pair EdgeOne Fire™ with:

EDGE BIOCERAMIC™
SEALER



EDGEONE PLATINUM™



\$27.95
3 Pack

Constant Taper. Reliable Reciprocation.

- Constant taper reciprocating system
- Verifile dictates sequence
- Verifile immediately following 10 and 15 handfiles
- Non-Sterile
- Alternative to WaveOne® Gold¹



- Triangular shaped cross-section
- Available in 3 lengths: 21, 25 & 31mm
- Designed to shape canals in reverse-reciprocating motion
- CCW 150*/CW 60* or WaveOne® Gold¹ settings



Highly Flexible




Pre-Bendable



High Strength

¹WaveOne® Gold is a registered trademarks of Dentsply Tulsa Dental Specialties.



3 Pack				
Size	Color	21mm	25mm	31mm
GLIDEPATH	<input type="checkbox"/> WHITE	EOGP21	EOGP25	EOGP31

3 Pack					Points		EdgeOne Obturators	
Size	Color	21mm	25mm	31mm	Standard .06 PaperPoints	Gutta Percha	EdgeCore	EdgeFill
VERIFILE	<input type="checkbox"/> WHITE	EOPVER21	EOPVER25	EOPVER31	SPP1506-1	SGP1506-1	N/A	N/A
SMALL	<input checked="" type="checkbox"/> YELLOW	EOPSM21	EOPSM25	EOPSM31	SPP2006-1	SGP2006-1	EOPSM20C	EOPSM20O
PRIMARY	<input checked="" type="checkbox"/> RED	EOPPRI21	EOPPRI25	EOPPRI31	SPP2506-1	SGP2506-1	EOPPRI25C	EOPPRI25O
MEDIUM	<input checked="" type="checkbox"/> GREEN	EOPMED21	EOPMED25	EOPMED31	SPP3506-1	SGP3506-1	EOPMED35C	EOPMED35O
LARGE	<input type="checkbox"/> WHITE	EOPLG21	EOPLG25	EOPLG31	SPP4506-1	SGP4506-1	EOPLG45C	EOPLG45O
Assorted 5 pack <input type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>					\$45.85	Paper Points	Gutta Percha	
(VER, SM, PRI, MED, LG)		ASSTEOP21	ASSTEOP25	ASSTEOP31	SPP06AST-1	SGP06AST-1		

Pair EdgeOne Platinum™ with:

EDGEBIOCERAMIC™
SEALER
818365020072





EDGEFILE® X1



\$38.80
3 Pack

Constant Taper. Consistent Reciprocation.

- Tighter pitch and flutes which deflect more torsional forces along the length of file due to the parabolic cross-section
- All X1 files (20.06 small, 25.06 primary, 40.06 large) have a constant 6% taper until the flute diameter is 1.0 mm then it goes to 0 taper
- FireWire NiTi™ heat treated
- Short handle 11mm
- Non-sterile
- Alternative to WaveOne®¹



- Parabolic shaped cross-section
- Non - cutting tip
- Available in 3 lengths: 21, 25 & 31mm
- Constant tapered reciprocating file system should use the same motor settings as WaveOne®¹



Highly Flexible



Pre-Bendable



High Strength

¹) WaveOne® is a registered trademark of Densply Sirona.

EDGEFILE® X1



3 Pack					Points		Obturator	
Size/Taper	Color	21mm	25mm	29mm	Standard .06 Paper Points	Gutta Percha	EdgeCore	EdgeFill
20.06	■ YELLOW	X1200621RF	X1200625RF	X1200629RF	SPP2006-1	SGP2006-1	X72006C	X72006O
25.06	■ RED	X1250621RF	X1250625RF	X1250629RF	SPP2506-1	SGP2506-1	X72506C	X72506O
40.06	■ BLACK	X1400621RF	X1400625RF	X1400629RF	SPP4006-1	SGP4006-1	X74006C	X74006O

Pair EdgeFile® X1 with:

EDGE BIOCERAMIC™
SEALER
818365020072



Assorted Pack - \$38.80/ 3 Pack (20, 25, 40)	
21mm	ASSTX121 ■ ■ ■
25mm	ASSTX125
29mm	ASSTX129

EdgeLube™ (17% EDTA SOLUTION) \$11.90 each	
100ML BOTTLE	8806387892083



EdgeLube™ 17% Solution comes in a 100mL bottle accompanied by a 1cc syringe and three flexible disposable tips.

- Ideal for Challenging Canal Conditions
- Effective in cases with:
 - Canal calcification
 - Narrow posterior canals
 - Where smear layer removal is required
- Utilizes chelation with EDTA (Ca^{2+}) to remove inorganic debris, facilitating easier instrumentation by softening dentin walls
- Delivers exceptionally clean dentinal surfaces following smear layer removal for optimal canal preparation

Composition: EDTA, Ammonia, Water

EdgeGel™ (19% EDTA GEL) \$14.05 each	
5.4ML BOTTLE	8806387892076







EdgeGel™ comes in two 7g (5.4mL) syringes with a pouch of 10 disposable tips




- Enhanced Root Canal Preparation
- Formulated to support instrumentation with files and reamers, especially in challenging cases such as calcified canals and narrow posterior canals.
- Features chelating action with EDTA (Ca^{2+}) to effectively remove inorganic debris and the smear layer
- Provides a lubricating effect for smoother instrumentation and easier filling, while ensuring exceptionally clean canal walls for optimal treatment outcomes

Composition: EDTA, PEG




Irrigation Tips & Syringes




Edge® Irrigation Needle Tips		\$10.70 each
23ga. Blue		10818365020260
27ga. Yellow		10818365020277
30ga. Purple		10818365020284



Edge® Irrigation Probe Tips		\$24.60 each
23ga. Blue		10818365020291
27ga. Yellow		10818365020307
30ga. Purple		10818365020314



Edge® Pre-Assembled Irrigation Syringe with 1\"		\$22.55 each
23ga. Blue		00818365021062
27ga. Yellow		00818365021055
30ga. Purple		00818365021079

Edge® Luer-Lock Syringe		
3CC	 \$12.85/ 100 Pack	10818365020383
6CC	 \$12.85/ 50 Pack	10818365020390
12CC	 \$16.10 / 100 Pack	10818365020406

2% Chlorhexidine Gluconate		\$32.25 each
480ML BOTTLE	818365020041	



17% EDTA		\$32.25 each
480ML BOTTLE	818365020058	



5% Sodium Hypochlorite		\$35.50 each
480ML BOTTLE	818365020034	



EDGE MIX™

EDGE MIX		\$64.65 each
16OZ BOTTLE	818365020416	



EDTA-based formula enhanced with Chlorhexidine for added anti-bacterial properties. Easy to use. 2-3ML per canal. In instances where calcifications or a lot of debris is present, the higher end dosage should be utilized.

- Clears 27% more dentin tubules than EDTA
- Kills 99.99% bacteria in 10 seconds
- Does not form a precipitant when mixed with NaOCl

EDGE BIO CERAMIC™ SEALERS

EDGE BIO CERAMIC™ SEALER

2g syringe + 15 disposable tips

Cold Obturation

EdgeBioCeramic Sealer™ trusted in over **50 million cases worldwide** since 2008, EdgeBioCeramic Sealer™ offers proven, reliable performance. Easy-to-use, premixed, and resin-free, it's hydrophilic with zero shrinkage. The reduced waste tip preserves material by 62%, maximizing efficiency and cost savings.

\$140.15 each

Features:

- Radiopacity: High
- Backed by over 17 years of clinical use and extensive studies
- Ideal viscosity for single cone hydraulic condensation
- Composition: Tricalcium Silicate, Dicalcium Silicate, Zirconium Oxide, Tantalum Oxide

Sealer Sku: 818365020072
Disposable Tips Sku: EEBCT15

\$19.30/ 15 tips



EDGE UTOPIA™ BioCeramic Sealer

2g syringe + 20 disposable tips

Warm & Cold Obturation


EdgeUtopia™ Sealer with its superior flow makes it **optimal for both warm or cold obturation techniques**. Its high pH of 12.5 offers powerful antibacterial properties, keeping the canal infection-free. Designed for moist environments, its hydrophilic nature ensures a secure seal in wet conditions.

\$78.65 each

Features:

- Radiopacity: Higher
- Our most cost effective Bioceramic sealer
- Ideal for both cold and warm obturation techniques
- Composition: Tricalcium silicate, dicalcium silicate, tricalcium aluminate, calcium oxide, zirconium oxide, silicon dioxide, polyethylene glycol and iron oxide

Sealer Sku: EUBCSEAL2G



EDGE BIO CERAMIC™ THERMALFLOW™

1.5g syringe + 15 disposable tips

Warm Obturation

EdgeBioCeramic™ ThermalFlow™ is an **optimized version** of EdgeBioCeramic Sealer™ with the same resin-free, hydrophilic, and zero shrinkage benefits, but with **lower viscosity for enhanced warm obturation techniques**. The reduced waste tip preserves material by 62%, maximizing efficiency and cost savings.


\$140.15 each

Features:

- Radiopacity: Highest
- Ideal for Warm Obturation with enhanced flowability
- Lowest viscosity and highest radiopacity
- Composition: Tricalcium Silicate, Dicalcium Silicate, Zirconium Oxide, Tantalum Oxide

Sealer Sku: 818365020232
Disposable Tips Sku: EEBCT15

\$19.30/ 15 tips



EDGESEAL™



818365020010

\$32.35 each

Item Code: (468-0153) Resin-based root canal sealer designed to work with various obturation methods. Designed for optimal performance, this root canal sealing material offers excellent biocompatibility and a hermetic seal to ensure long-term success. The easy-to-mix paste allows for efficient handling, while its nonstaining formula preserves natural tooth color. EdgeSeal provides clear visibility on radiographs with its strong radiopacity.

Working Time: 35 minutes at 23°C (73°F)

Setting Time: 45 minutes at 37°C (99°F)

Film Thickness: 3.3µ

Composition: BASE: Epoxy oligomer resin, Ethylene glycol salicylate, Calcium phosphate, Bismuth subcarbonate, Zirconium oxide. CATALYST: Poly aminobenzoate, Tri ethanolamine, Calcium phosphate, Bismuth subcarbonate, Zirconium oxide, Calcium oxide

As Featured on File Product Page Tables

EDGECORE™

EDGEFILL™

Pink Gutta Percha Core
Thermal Obturator

\$37.50



5 Pack & 1 Verifier

Plastic Carrier
Thermal Obturator

\$37.50



6 Pack

20 Pack Obturators - \$99.95 each

Size	EdgeCore	EdgeFill
20	EE2020C	EE2020O
25	EE2520C	EE2520O
30	EE3020C	EE3020O
35	EE3520C	EE3520O
40	EE4020C	EE4020O
45	EE4520C	EE4520O
50	EE5020C	EE5020O

NEW EdgeEndo® Gutta Percha: Optimized Performance, Seamless Results Optimized For Use with Our Most Popular NiTi File Systems!



Seamless Integration

Tailored to work effortlessly with EdgeEndo file systems for streamlined endodontic procedures.



Proprietary Stiffening Process

Provides improved handling, less tip crumpling, and improved displacement of sealer.



Optimized Apical Fit with Coronal Venting

Designed for a tight apical fit with extra space in the coronal third to facilitate the upward venting of sealer thus minimizing the risk of apical extrusion. Additionally, this design helps prevent gutta percha from “hanging up” in longer canals promoting more efficient obturation.

Gutta Percha Quick Reference

These are the suggested gutta percha points for EdgeOne Fire, all EdgeTaper Families and EdgeOne Blaze Utopia NiTi Systems.

EDGEONE FIRE™ Gutta Percha			
29.5mm- 60 Pack		Item Codes	
	YELLOW	SMALL	EOFSMGP
	RED	PRIMARY	EOFPRIGP
	GREEN	MEDIUM	EOFMEDGP
	WHITE	LARGE	EOFLGGP
Assorted: SM (20), PRI (20) MED (10) LG (10)		EOFGPAST-1	
<ul style="list-style-type: none"> • Taper to match file design • Parallel when diameter meets 1.0mm • Use with EdgeOne Fire™ Files 			

EDGETAPER™ Gutta Percha			
29.5mm- 60 Pack		Item Codes	
	YELLOW	F1	ETGPF1
	RED	F2	ETGPF2
	BLUE	F3	ETGPF3
	BLACK	F4	ETGPF4
	YELLOW	F5	ETGPF5
Assorted: F1 (10), F2 (10), F3 (20), F4 (10), F5 (10)		ETGPAST	
<ul style="list-style-type: none"> • Taper to match file design • Parallel 16mm from tip • Use with EdgeTaper Platinum™ Files • Use with EdgeTaper™ Files • Use with EdgeTaper Blaze Utopia™ Files 			

EDGEONE R UTOPIA™ Gutta Percha			
29.5mm- 60 Pack		Item Codes	
	YELLOW	R20	EORGP20-1
	RED	R25	EORGP25-1
	GREEN	R35	EORGP35-1
	WHITE	R45	EORGP45-1
<ul style="list-style-type: none"> • Taper to match file design • Parallel when diameter meets 1.0mm • Use with EdgeOne Blaze Utopia™ Files 			

EDGECAOH™

EdgeCaOH™		\$20.90 each
2.2G SYRINGE	0081836502003	



Temporary root-canal filler for weeping canals that forms a hard-tissue barrier (apexification). Internal and external root resorption. Pre-mixed paste. Syringe tip for easy delivery of paste into canal. 10 Disposable tips included to prevent cross-contamination.

EDGETEMP™ PLUS

EdgeTemp™ Plus		\$11.85 each
40G	WHITE PASTE	ETEMPPLUS



Hydraulic temporary restorative material. Designed for fabricating temporary fillings and other temporary prosthetics until the permanent restoration can be delivered. Setting time: 30 minutes.

EDGEBIOCERAMIC™ RETROFILL™ AND PERFORATION REPAIR

EdgeBioCeramic Retrofill™ and Perforation Repair		\$204.90 each
0.5G	WHITE PASTE	818365020089



EdgeBioCeramic™ RetroFill is a unique bioceramic repair material that exhibits superior healing and handling properties that has been shown in studies to out perform MTA. In addition, EdgeBioCeramic RetroFill is premixed and ready to use unlike MTA. Simply syringe out a small amount of material and deliver it to the repair site.

The Future of Root Repair

EDGE  **UTOPIA RRM™**

Root Repair Material



EdgeUtopia™ RRM is the optimal solution for root repair, offering excellent clinical outcomes with its scientifically backed performance, versatility, and long-term efficacy. EdgeUtopia RRM is the ideal choice for all your root and pediatric repair needs.

This proven pre-mixed calcium silicate repair material has now been upgraded to feature a novel monodose delivery system to improve handling and reduce the risk of cross-contamination. Trust in EdgeUtopia™ RRM for tissue healing, bacteria prevention, and overall root canal repair success.



Non Staining



Wide Range of Applications



Desired Biological Performance



Convenient Single Dose Delivery



Inhibits Bacterial Infiltration



Long-Lasting Results



0.15g per capsule
10 capsules per box

EDGE UTOPIA RRM™

Root Repair Material

Indications



Pre-Mixed Ready-to-Use Monodose Capsule

Forget the hassle of your repair material drying out too soon. Experience effortless, reliable performance with every use. No manual mixing and no titration devices needed.



Reduced Cross-Contamination

Unlike other premixed, manually mixed or titrated bioceramics, EdgeUtopia RRM™ ensures hygienic, mess-free delivery with no risk of contamination between patients due to its single dose application and the elimination of extensive manipulation in the operatory.



Release of Ca²⁺

Proven track record with over 5 years of clinical cases and research studies* demonstrating exceptional tissue healing and regeneration outcomes.



Wide Range of Applications

- Apexification
- Retrofilling
- Perforation Repairs
- Pulpotomy
- Pulp Capping
- External and Internal Resorption
- Deciduous Tooth Repairs



Cost-Effective Innovation

More product for a fraction of the price – helping you save money while delivering high-quality results. 3x the amount of our regular syringe formula per box.



0.15g per capsule
10 capsules per box
EURRM10MD
\$102.50/ Per Box



Pair with:
Reusable EdgeUtopia™ RRM™
Dispensing Gun

EURRMDSG
\$29.10 each

Available research and clinical case studies may refer to other private-label brands which are the same chemical formula as EdgeUtopia RRM™ Root Repair Material. For clarity, the brands EdgeUtopia RRM™ and Bio-C Repair are both manufactured by Angelus and are composed of equivalent chemical formulations.



EDGEENDO®

© 2026 Edge Endo, a division of HSG-S Corp.
No copying without permission.

Not responsible for typographical errors.

California customers please use website when ordering for Prop 65 information.

26-1126 | 2026-04-08 | Rev 0